



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center – WO66-G609  
Silver Spring, MD 20993-0002

Barco NV  
% Mr. Lieven De Wandel  
Regulatory Affairs Officer  
35 President Kennedypark  
8500 Kortrijk  
BELGIUM

November 12, 2014

Re: K143157  
Trade/Device Name: Coronis Fusion 4MP LED  
Regulation Number: 21 CFR 892.2050  
Regulation Name: Picture archiving and communications system  
Regulatory Class: II  
Product Code: PGY  
Dated: October 21, 2014  
Received: November 3, 2014

Dear Mr. De Wandel:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Industry and Consumer Education at its toll-free number (800) 638 2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, “Misbranding by reference to premarket notification” (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH’s Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

A handwritten signature in black ink that reads "Robert A. Ochs". The signature is written in a cursive style. A faint, large "FDA" watermark is visible in the background behind the signature.

for

Janine M. Morris  
Director  
Division of Radiological Health  
Office of In Vitro Diagnostics  
and Radiological Health  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)  
K143157

Device Name  
Coronis Fusion 4MP LED (MDCC-4230)

Indications for Use (Describe)

The Coronis Fusion 4MP LED (MDCC-4230) Medical Flat Panel Display System is intended to be used as a tool in displaying and viewing digital images (excluding digital mammography) for review and analysis by trained medical practitioners.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

**PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON A SEPARATE PAGE IF NEEDED.**

### FOR FDA USE ONLY

Concurrence of Center for Devices and Radiological Health (CDRH) (Signature)

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*“An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number.”*

## 510(k) Summary (in accordance with 21 CFR 807.92)

1. Company	Barco N.V. Healthcare Division 35 President Kennedypark B-8500 Kortrijk BELGIUM																													
2. Contact person	Lieven De Wandel Regulatory Affairs Officer																													
3. Date of submission	October 21, 2014																													
4. Device information	Trade name/model: <b>Coronis Fusion 4MP LED (MDCC-4230)</b> Common name: 4MP color LCD display Classification name: System, image processing, Radiological Classification code: PGY Regulation number: 892.2050																													
5. Predicate device	<b>Coronis Fusion 4MP DL (MDCC-4130)</b> cleared under 510(K) <b>K111989</b>																													
6. Device description	The <b>Coronis Fusion 4MP LED (MDCC-4230)</b> is a high-resolution flat panel LCD display system for reviewing medical images. It consists of an LCD display ( <b>MDCC-4230</b> ), an optional high-resolution display controller board and QA software.  The display controller board is installed in a PACS workstation computer, connected to the display. The QA software helps to make and keep the displays DICOM compliant.  The display uses LED backlight technology.																													
7. Intended Use of the Device	The <b>Coronis Fusion 4MP LED (MDCC-4230) Medical Flat Panel Display System</b> is intended to be used as a tool in displaying and viewing digital images (excluding digital mammography) for review and analysis by trained medical practitioners.																													
8. Comparison of technological characteristics	<table border="1"> <thead> <tr> <th>Specification</th> <th>MDCC-4130</th> <th>MDCC-4230</th> </tr> </thead> <tbody> <tr> <td>Screen technology</td> <td>TFT AM LCD Dual Domain IPS Pro</td> <td>TFT AM LCD Dual Domain IPS Pro</td> </tr> <tr> <td>Active screen size (diagonal)</td> <td>756 mm (29.8")</td> <td>756 mm (29.8")</td> </tr> <tr> <td>Active screen size (H x V)</td> <td>641.28 x 400.8 mm (25.2 x 15.8")</td> <td>641.28 x 400.8 mm (25.2 x 15.8")</td> </tr> <tr> <td>Aspect ratio (H:V)</td> <td>16:10</td> <td>16:10</td> </tr> <tr> <td>Resolution</td> <td>4MP (2560 x 1600)</td> <td>4MP (2560 x 1600)</td> </tr> <tr> <td>Pixel pitch</td> <td>0.2505 mm</td> <td>0.2505 mm</td> </tr> <tr> <td>Color imaging</td> <td>Yes</td> <td>Yes</td> </tr> <tr> <td>Gray imaging</td> <td>No</td> <td>No</td> </tr> </tbody> </table>	Specification	MDCC-4130	MDCC-4230	Screen technology	TFT AM LCD Dual Domain IPS Pro	TFT AM LCD Dual Domain IPS Pro	Active screen size (diagonal)	756 mm (29.8")	756 mm (29.8")	Active screen size (H x V)	641.28 x 400.8 mm (25.2 x 15.8")	641.28 x 400.8 mm (25.2 x 15.8")	Aspect ratio (H:V)	16:10	16:10	Resolution	4MP (2560 x 1600)	4MP (2560 x 1600)	Pixel pitch	0.2505 mm	0.2505 mm	Color imaging	Yes	Yes	Gray imaging	No	No		
Specification	MDCC-4130	MDCC-4230																												
Screen technology	TFT AM LCD Dual Domain IPS Pro	TFT AM LCD Dual Domain IPS Pro																												
Active screen size (diagonal)	756 mm (29.8")	756 mm (29.8")																												
Active screen size (H x V)	641.28 x 400.8 mm (25.2 x 15.8")	641.28 x 400.8 mm (25.2 x 15.8")																												
Aspect ratio (H:V)	16:10	16:10																												
Resolution	4MP (2560 x 1600)	4MP (2560 x 1600)																												
Pixel pitch	0.2505 mm	0.2505 mm																												
Color imaging	Yes	Yes																												
Gray imaging	No	No																												

Viewing angle (H, V)	170°	178°
Uniform Luminance Technology (ULT)	Yes	Yes
Per Pixel Uniformity (PPU)	No	Yes
Ambient Light Compensation (ALC)	Yes	Yes
Backlight Output Stabilization (BLOS)	Yes	Yes
I-Guard	Yes	Yes
Maximum luminance	950 cd/m <sup>2</sup>	720 cd/m <sup>2</sup>
DICOM calibrated luminance (ULT off)	500 cd/m <sup>2</sup>	500 cd/m <sup>2</sup>
Contrast ratio (typical)	1100:1	1000:1
Response time (Tr + Tf)	20 ms	20 ms
Scanning frequency (H; V)	30-150 kHz; 15-80 Hz	30-150 kHz; 15-80 Hz
Video input signals	DVI-D Dual Link	DVI-D Dual Link / DisplayPort 1.1
USB ports	1 upstream (endpoint), 3 downstream	1 upstream (endpoint), 3 downstream
USB standard	1.0	2.0
Power requirements (nominal)	100-240V	100-240V
Power consumption (nominal)	135W	105W
Power save mode	Yes	Yes
Net weight with stand	28.1 kg	21.5 kg
Net weight w/o stand	21.6 kg	15 kg

<p>9. Performance testing</p>	<p>The bench tests mentioned below were performed to validate the device characteristics that differ from the predicate device:</p> <table border="1" data-bbox="391 1045 1430 1581"> <thead> <tr> <th>Modification to device</th> <th>Test performed</th> </tr> </thead> <tbody> <tr> <td>LED backlight instead of CCFL</td> <td>DICOM calibration and Luminance Uniformity tests</td> </tr> <tr> <td>Different platform (including firmware)</td> <td>Firmware tests</td> </tr> <tr> <td>Additional DisplayPort video input</td> <td>Firmware tests</td> </tr> <tr> <td>Uniformity correction: Per pixel uniformity (vs per zone unif. on predicate device)</td> <td>Luminance Uniformity tests</td> </tr> <tr> <td>Other material of front filter</td> <td>Impact test in IEC 60601-1 tests</td> </tr> <tr> <td>Other material for sheet metal parts</td> <td>Shock and Vibration tests in Environmental test report</td> </tr> </tbody> </table> <p>Additional tests performed: Electrical Safety test (IEC 60601-1), EMC test (IEC 60601-1-2)</p> <p>The tests showed that the device has similar or superior characteristics compared to the predicate device and did not reveal new issues of safety and performance.</p> <p>Animal or clinical testing have not been performed.</p>	Modification to device	Test performed	LED backlight instead of CCFL	DICOM calibration and Luminance Uniformity tests	Different platform (including firmware)	Firmware tests	Additional DisplayPort video input	Firmware tests	Uniformity correction: Per pixel uniformity (vs per zone unif. on predicate device)	Luminance Uniformity tests	Other material of front filter	Impact test in IEC 60601-1 tests	Other material for sheet metal parts	Shock and Vibration tests in Environmental test report
Modification to device	Test performed														
LED backlight instead of CCFL	DICOM calibration and Luminance Uniformity tests														
Different platform (including firmware)	Firmware tests														
Additional DisplayPort video input	Firmware tests														
Uniformity correction: Per pixel uniformity (vs per zone unif. on predicate device)	Luminance Uniformity tests														
Other material of front filter	Impact test in IEC 60601-1 tests														
Other material for sheet metal parts	Shock and Vibration tests in Environmental test report														
<p>10. Conclusion</p>	<p>The <b>Coronis Fusion 4MP LED (MDCC-4230)</b> was found to be substantially equivalent to</p>														

the predicate device, due to the following reasons:

- a) Device and predicate device have the same intended use
- b) The technological characteristics differences from the predicate device do not affect safety or effectiveness
- c) Bench testing showed that the device has similar or superior characteristics compared to the predicate device and did not reveal new issues of safety and performance.