



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center – WO66-G609  
Silver Spring, MD 20993-0002

ROCHE DIAGNOSTICS  
KELLI TURNER  
REGULATORY AFFAIRS PRINCIPAL  
9115 HAGUE ROAD, PO BOX 50416  
INDIANAPOLIS IN 46250-0416

December 19, 2014

Re: K143342  
Trade/Device Name: Elecsys T-Uptake CalSet  
Regulation Number: 21 CFR 862.1150  
Regulation Name: Calibrator  
Regulatory Class: II  
Product Code: JIT  
Dated: November 20, 2014  
Received: November 21, 2014

Dear Ms. Kelli Turner:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulations (21 CFR Parts 801 and 809), please contact the Division of Industry and Consumer Education at its toll-free number (800) 638 2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, “Misbranding by reference to premarket notification” (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH’s Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

  
**Katherine Serrano -S**

For: Courtney H. Lias, Ph.D.  
Director  
Division of Chemistry and Toxicology Devices  
Office of In Vitro Diagnostics  
and Radiological Health  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)

k143342

Device Name

Elecsys T-Uptake CalSet

Indications for Use (Describe)

T-Uptake CalSet is used for calibrating the quantitative Elecsys T-Uptake assay on the Elecsys and cobas e immunoassay analyzers.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

## 510(k) Summary, Elecsys T-Uptake CalSet

---

**Introduction** According to the requirements of 21 CFR 807.92, the following information provides sufficient detail to understand the basis for a determination of substantial equivalence.

The purpose of this premarket notification is to obtain FDA review and clearance for the Elecsys T-Uptake CalSet.

---

**Submitter name, address, contact** Roche Diagnostics  
9115 Hague Road, P.O. Box 50416  
Indianapolis, IN 46250-0416

Contact Person: Kelli Turner  
Phone: 317-521-4515  
Fax: 317-521-2324  
Email: kelli.turner@roche.com

Secondary Contact: Michael Leuther  
Phone: 317-521-3930  
Fax: 317-521-2324  
Email: michael.leuther@roche.com

Date Prepared: December 18, 2014

---

**Device Name** Proprietary name: Elecsys T-Uptake CalSet  
Common name: T-Uptake CalSet  
Classification: 21 CFR 862.1150, Calibrator, Secondary  
Product Code: JIT

---

**Establishment Registration** For the T-Uptake CalSet, the establishment registration number for Roche Diagnostics GmbH in Mannheim, Germany, is 9610126 and for Penzberg, Germany, is 9610529. The establishment registration number for Roche Diagnostics United States is 1823260.

---

*Continued on next page*

## 510(k) Summary, Elecsys T-Uptake CalSet, *continued*

**Classification** FDA has classified the product as a Class II device.

<b>Product Name</b>	<b>Panel</b>	<b>Product Code</b>	<b>Classification Name</b>	<b>Regulation Citation</b>
Elecsys T-Uptake CalSet	Clinical Chemistry	JIT	Calibrator, Secondary	21 CFR 862.1150

**Device Description**

Elecsys T-Uptake CalSet:

- T-Uptake CalSet is a lyophilized buffer/protein/TBG matrix with added L-thyroxine. The Elecsys T-Uptake assay is traceable to a clinically defined human serum panel with a mean TBI (thyroxine-binding index) of 1.

**Intended use**

Elecsys T-Uptake CalSet:

- T-Uptake CalSet is used for calibrating the quantitative Elecsys T-Uptake assay on the Elecsys and **cobas e** immunoassay analyzers.

**Predicate device**

The Elecsys T-Uptake CalSet is substantially equivalent to the predicate device Elecsys T-Uptake CalSet (k961488).

**Substantial Equivalence Comparison**

The following tables compare the Elecsys T-Uptake CalSet with the predicate device.

*Continued on next page*

## 510(k) Summary, Elecsys T-Uptake CalSet, *continued*

**Comparison Table** The table below compares Elecsys T-Uptake CalSet with the predicate device, Elecsys T-Uptake CalSet (k961488).

The change in the new product was in the format, going from liquid to lyophilized.

**Comparison Table**

Characteristic	Elecsys T-Uptake CalSet (Candidate)	Elecsys T-Uptake CalSet (k961488)
Intended Use	T-Uptake CalSet is used for calibrating the quantitative Elecsys T-Uptake assay on the Elecsys and <b>cobas e</b> immunoassay analyzers.	Same.
Analyte	L-Thyroxine, syntheticThyroxine binding globuline (TBG), bovine	Same
Matrix	Buffer / protein (bovine serum albumin) / TBG matrix	Same
Levels	Two	Same
Target Ranges	Cal 1: 0.5 TBI conc. Cal 2: 1.3 TBI conc.	Cal 1: Same. Cal 2: Same.
Format	Lyophilized	Liquid, ready-for-use.
Stability	<p><u>Unopened:</u></p> <ul style="list-style-type: none"> <li>Store at 2-8°C until expiration date</li> </ul> <p><u>Opened:</u></p> <ul style="list-style-type: none"> <li>2-8°C: 84 days</li> <li>-20°C: 84 days (freeze only once)</li> </ul> <p>Elecsys 2010 and <b>cobas e</b> 411 analyzers at 20-25°C:</p> <ul style="list-style-type: none"> <li>up to 5 hours</li> </ul> <p>MODULAR ANALYTICS E170, <b>cobas e</b> 601 and <b>cobas e</b> 602:</p> <ul style="list-style-type: none"> <li>use only once</li> </ul>	<p><u>Unopened:</u></p> <ul style="list-style-type: none"> <li>Store at 2-8°C until expiration date</li> </ul> <p><u>Opened:</u></p> <ul style="list-style-type: none"> <li>2-8°C: 24 hours</li> <li>20-25°C: use only once</li> <li>-20°C: 3 month</li> </ul>

*Continued on next page*

## 510(k) Summary, Elecsys T-Uptake CalSet, *continued*

Table *continued*

Characteristic	Elecsys T-Uptake CalSet (Candidate)	Elecsys T-Uptake CalSet (k961488)
Handling	<p>Carefully dissolve the contents of one bottle by adding exactly 1.0 mL of distilled or deionized water and allow to stand closed for 15 minutes to reconstitute. Mix carefully, avoiding the foam formation.</p> <p>Transfer the reconstituted calibrators into the supplied empty labeled snap-cap bottles .</p> <p>Elecsys 2010 and <b>cobas e 411</b> analyzers: The calibrators should only be left on the analyzers during calibration at 20-25 °C. After use, close the bottles as soon as possible and store upright at 2-8 °C.</p> <p>Due to possible evaporation effects, not more than 5 calibration procedures per bottle set should be performed. If necessary, freeze in aliquots; see section on</p> <p>MODULAR ANALYTICS E170, <b>cobas e 601</b> and <b>cobas e 602</b> analyzers.</p> <p>MODULAR ANALYTICS E170, <b>cobas e 601</b> and <b>cobas e 602</b> analyzers:</p> <p>Unless the entire volume is necessary for calibration on the analyzers, transfer aliquots of the reconstituted calibrators into empty snap-cap bottles (CalSet Vials). Attach the supplied labels to these additional bottles. Store the aliquots at 2-8 °C or -20 °C for later use.</p> <p>Perform <b>only one</b> calibration procedure per aliquot. .</p>	<p>The calibrators are supplied ready for use in bottles compatible with the system.</p> <p>Elecsys 2010 and <b>cobas e 411</b> analyzers: The calibrators should only be left on the analyzers during calibration at 20-25°C. After use, close the bottles as soon as possible and store at 2-8°C. Due to possible evaporation effects, not more than 5 calibration procedures per bottle set should be performed.</p> <p>MODULAR ANALYTICS E170, <b>cobas e 601</b>, <b>cobas e 602</b> analyzers: Unless the entire volume is necessary for calibration on the analyzers, transfer aliquots of the ready for use calibrators into empty snap-cap bottles (CalSet Vials). Attach the supplied labels to these additional bottles. Store the aliquots at 2-8°C for later use. Perform only one calibration procedure per aliquot.</p>
Traceability	Standardized using a clinically defined human serum panel with a mean TBI of 1.0.	Same

## 510(k) Summary, Elecsys T-Uptake CalSet, *continued*

---

**Traceability** The T-Uptake CalSet has been standardized against a clinically defined human serum panel with a mean TBI of 1.

---

**Evaluations Summary** The Elecsys T-Uptake CalSet was evaluated for value assignment, stability and reconstitution.

---

**CalSet Value Assignment** Value assignment testing was conducted and passed pre-defined acceptance criteria. The target values for the two levels of the T-Uptake CalSet kit are chosen to obtain the best fit with the Reference Standard Curve. For each Elecsys T-Uptake CalSet lot manufactured, the calibrators are run in duplicate on at least three (3) Elecsys 2010 analyzers and at least three (3) MODULAR ANALYTICS E170 analyzers with all T-Uptake reagent lots available. The assigned value of each calibrator is defined as the median value obtained over at least six (6) runs on at least three (3) analyzers) of the respective calibrator.

Measurement values for PreciControl Universal (Level 1 & 2) are read off from the calibration curves generated. The pre-defined acceptance criteria for PreciControl Universal have to be met to release the Assigned Values for T-Uptake CalSet.

---

*Continued on next page*

## 510(k) Summary, Elecsys T-Uptake CalSet, *continued*

---

### **Stability Studies**

Three studies were performed in order to verify the stability claims for the T-Uptake CalSet. Stability studies after reconstitution and an accelerated stability study were completed on the **cobas e 411**. Additionally, a real-time stability study is planned.

#### **Study 1. Stability at 2-8°C, -20°C, in open vial and freeze/thaw cycles (after reconstitution):**

The on-test and reference materials were tested in duplicate. The on-test material was reconstituted and stored in closed vials for 85 days at 2 to 8°C, for 699 days at -15 to -25°C and for 6 hours at 20 to 25°C in open vial. In addition, the stability of the T-Uptake CalSet for one (1) freeze/thaw cycles was evaluated.

The on-test signal recovery was calculated as percent of the reference value. The reference material was a freshly reconstituted set of T-Uptake CalSet.

T-Uptake CalSet was evaluated in duplicate on the **cobas e 411**. The acceptance criterion was 95 to 105 % signal recovery of the reference material value.

#### **Study 2. Accelerated Stability:**

The on-test material was stored lyophilized (as supplied to the user) at 35°C for 3 weeks. The reference material was a freshly reconstituted set of T-Uptake CalSet (stored at 2 to 8°C). After 3 weeks, the on-test and reference materials were tested in duplicate. The on-test recovery was calculated as a percent of the reference value.

One T-Uptake CalSet lot was evaluated in duplicate on the **cobas e 411** analyzer.

The acceptance criterion was 95 to 105 % recovery of the reference material value.

---

*Continued on next page*

## 510(k) Summary, Elecsys T-Uptake CalSet, *continued*

---

### Stability Studies, *continued*

#### **Study 3. Real-Time Stability:**

In addition, real-time stability is being evaluated as follows:

In the real-time stability study, the T-Uptake CalSet test material is stored at 2-8°C. The CalSets are tested in duplicate at specified intervals over the shelf life of the device up to the planned shelf life plus one month (25 months).

Real-time stability is calculated based on the recovery of signal of stressed calibrator (stored at 2-8°C) vs. unstressed calibrator (stored at -20°C). At the specified intervals over the shelf life, the mean value of the stressed calibrator was calculated as percent recovery of the unstressed value (each tested in duplicates at the same time point).

The acceptance criterion for T-Uptake Calibrator 1 and 2 is recovery of 95-105 % of the reference value.

---

### Reconstitution Study

Reconstitution time for the lyophilized T-Uptake CalSet was tested. Two sets of T-Uptake CalSet were reconstituted, one for 15 minutes and the other for 30 minutes. Signal recovery after 30 minutes reconstitution was compared to the signal value after 15 minutes.

T-Uptake CalSet was evaluated in duplicate on the **cobas e 411** analyzer as a reference.

The acceptance criterion was 95 to 105 % signal recovery of the reference material value.

---

### Conclusion

We trust that the data and information provided in this Premarket Notification (510(k)) will support a determination of substantial equivalence for the Elecsys T-Uptake CalSet. The data supports a safe effective device which performs as well as or better than the predicate device.

---