



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center – WO66-G609  
Silver Spring, MD 20993-0002

June 9, 2016

Spinal Elements, Incorporated  
Julie Lamothe, Ph.D., M.B.A.  
Regulatory Affairs and Quality Assurance Director  
3115 Melrose Drive, Suite 200  
Carlsbad, CA 92010

Re: K152011

Trade/Device Name: Lucent® Intervertebral Body Fusion Device  
Regulation Number: 21 CFR 888.3080  
Regulation Name: Intervertebral body fusion device  
Regulatory Class: Class II  
Product Code: MAX  
Dated: May 6, 2016  
Received: May 9, 2016

Dear Dr. Lamothe:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set

forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to

<http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

**Mark N. Melkerson -S**

Mark N. Melkerson  
Director  
Division of Orthopedic Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)

K152011

Device Name

Lucent® Intervertebral Body Fusion Device

### Indications for Use (Describe)

Lucent intervertebral body fusion devices are intended for spinal fusion procedures at one or two contiguous levels (L2-S1) in skeletally mature patients with degenerative disc disease (DDD). DDD is defined as back pain of discogenic origin with degeneration of the disc confirmed by history and radiographic studies. DDD patients may also have up to Grade 1 spondylolisthesis or retrolisthesis at the involved levels. These patients may have had a previous non-fusion spinal surgery at the involved spinal level(s).

This device is intended to be used with supplemental spinal fixation systems that have been cleared for use in the lumbosacral spine (i.e., posterior pedicle screw and rod systems, anterior plate systems, and anterior screw and rod systems).

This device is intended to be used with autogenous bone graft. Patients must have undergone a regimen of at least six (6) months non-operative treatment prior to being treated with this device.

Type of Use (Select one or both, as applicable)

☒ Prescription Use (Part 21 CFR 801 Subpart D)

☐ Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

#### **\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
PRAStaff@fda.hhs.gov

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

*Spinal Elements, Inc.*  
*Premarket Notification –Lucent®XP Intervertebral Body Fusion Device*

**510(k) Summary**  
**Lucent® XP and Lucent XP Ti-Bond®**

**510(k) Number** \_\_K152011\_\_

**Manufacturer Identification**

**Submitted by:**

Spinal Elements, Inc.  
3115 Melrose Dr., Suite 200  
Carlsbad, CA 92010  
P. 760-607-0121  
F. 760-607-0125

**Contact Information:**

Julie Lamothe  
Director Regulatory Affairs  
Spinal Elements, Inc.  
3115 Melrose Dr., Suite 200  
Carlsbad, CA 92010  
760-607-1816  
jlamothe@spinalelements.com

**Date Prepared:**

July 17<sup>th</sup>, 2015

**Proprietary Name**

Lucent® Intervertebral Body Fusion Device

**Device Classification**

21 CFR 888.3080 (Appliance, Fixation Spinal  
Intervertebral Body)

**Proposed Regulatory Class**

Class II

**Device Product Code**

MAX

**Purpose of this 510(k)**

This 510(k) seeks clearance for line additions to the Lucent, Lucent Ti-Bond and Lucent Lateral previously cleared for use under K071724 and K150061, K110632, K122967, respectively.

**Device Description**

The Lucent® XP device is an intervertebral body fusion device for use in lumbar spinal surgery. It may also be referred to as an interbody device or interbody cage. The Lucent® XP device is comprised of two PEEK endplates with teeth for engaging the vertebral body endplates on the outside of the device and an internal titanium mechanism for expanding the PEEK endplates of the device. The device is generally box-shaped with various holes throughout its design to allow for the placement of autograft. The exterior

of the device has “teeth” or other generally sharp engagement members on the superior and inferior surfaces to help prevent the device from migrating once it is surgically positioned. The device is made from PEEK-Optima® grade LT1 conforming to ASTM F2026 with or without a coating of commercially pure titanium conforming to ASTM F1580. The internal mechanism is made from titanium alloy conforming to ASTM F136.

### **Indications for Use**

Lucent intervertebral body fusion devices are intended for spinal fusion procedures at one or two contiguous levels (L2-S1) in skeletally mature patients with degenerative disc disease (DDD). DDD is defined as back pain of discogenic origin with degeneration of the disc confirmed by history and radiographic studies. DDD patients may also have up to Grade 1 spondylolisthesis or retrolisthesis at the involved levels. These patients may have had a previous non-fusion spinal surgery at the involved spinal level(s).

This device is intended to be used with supplemental spinal fixation systems that have been cleared for use in the lumbosacral spine (i.e., posterior pedicle screw and rod systems, anterior plate systems, and anterior screw and rod systems).

This device is intended to be used with autogenous bone graft. Patients must have undergone a regimen of at least six (6) months non-operative treatment prior to being treated with this device.

### **Substantial Equivalence**

Lucent® XP and Lucent XP Ti-Bond® is identical in indications for use, manufacturing method, raw material and operating principles to the predicate devices cleared in:

K071724 – Primary Predicate

K093242 – Reference Device

K110632 – Additional Predicate

K122967 – Additional Predicate

K150061 – Additional Predicate

K123231 – Additional Predicate

### **Performance Data**

Performance testing included:

- Static Compression Testing per ASTM F 2077-14
- Dynamic Compression Testing per ASTM F 2077-14
- Subsidence Testing per ASTM F 2267-11
- Static & Dynamic Compression Shear Testing per ASTM F2077-14

All data indicates that the device will perform as intended.

Based on the indications for use, technological characteristics, performance testing, and comparison to predicate devices, the subject device has been shown to be substantially equivalent to the aforementioned predicate devices cleared by FDA for commercial distribution in the United States.