



Food and Drug Administration
10903 New Hampshire Avenue
Document Control Center – WO66-G609
Silver Spring, MD 20993-0002

July 7, 2016

Medline Industries, Inc.
Ms. Jennifer Mason
Senior Regulatory Affairs Specialist
1 Medline Place
Mundelein, IL 60060

Re: K152564
Trade/Device Name: Gemini Sterilization Wrap
Regulation Number: 21 CFR 880.6850
Regulation Name: Sterilization wrap
Regulatory Class: Class II
Product Code: FRG
Dated: May 31, 2016
Received: June 1, 2016

Dear Ms. Jennifer Mason:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA).

You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical

device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to

<http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

Tejashri Purohit-Sheth, M.D.

Tejashri Purohit-Sheth, M.D.
Clinical Deputy Director
DAGRID/ODE/CDRH FOR

Erin I. Keith, M.S.
Director
Division of Anesthesiology, General Hospital,
Respiratory, Infection Control and Dental Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

K152564

Device Name

Gemini Sterilization Wrap

Indications for Use (Describe)

Gemini Sterilization Wraps are intended to be used to enclose another medical device that is to be sterilized by a health care provider. It is intended to allow sterilization of the enclosed medical device and also to maintain sterility of such content until used. Gemini Sterilization Wrap is validated for use in steam or STERRAD sterilization processes in the following sterilization modes and cycles:

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

TABLE 1: Pre-vacuum steam at 270°F/132°C for 4 minutes with the following dry times by weight

Gemini Wrap Model	Gemini Wrap Weight	Dry Time
GEM11XX	Lightweight	20 Minutes
GEM21XX	Regular Weight	20 Minutes
GEM31XX	Medium Weight	20 Minutes
GEM41XX	Heavyweight	30 Minutes
GEM51XX	Super Heavyweight	30 Minutes

TABLE 2: Gravity steam at 250°F/121°C for 30 minutes with the following dry times by weight

Gemini Wrap Model	Gemini Wrap Weight	Dry Time
GEM11XX	Lightweight	20 Minutes
GEM21XX	Regular Weight	20 Minutes
GEM31XX	Medium Weight	20 Minutes
GEM41XX	Heavyweight	30 Minutes
GEM51XX	Super Heavyweight	30 Minutes

TABLE 3: STERRAD® Sterilization

Gemini Wrap Model	Gemini Wrap Weight	STERRAD® Cycles
GEM11XX	Lightweight	STERRAD® 50, 200S and 100NX™ DUO Cycles
GEM21XX	Regular Weight	STERRAD® 50, 200S and 100NX™ DUO Cycles
GEM31XX	Medium Weight	STERRAD® 50, 200S and 100NX™ DUO Cycles
GEM41XX	Heavyweight	STERRAD® 50, 200S and 100NX™ DUO Cycles
GEM51XX	Super Heavyweight	STERRAD® 50, 200S and 100NX™ DUO Cycles

The Gemini Sterilization Wrap intended loads with pre-vacuum steam sterilization, gravity steam sterilization and the Advanced Sterilization Products STERRAD® Sterilization systems are provided in Table 4.

The Gemini Sterilization Wrap Recommendations for use with Pre-vacuum steam cycles, Gravity steam cycles, and the Advanced Sterilization Products STERRAD® Sterilization systems are provided in Table 5.

TABLE 4: Validated Pre-Vacuum Steam, Gravity Steam, and Advanced Sterilization Products STERRAD® 50, 200S and STERRAD® 100 NX™ DUO Cycles

Cycle	Maximum recommended Chamber Load		Intended Load
Pre-Vacuum Steam & Gravity Steam Cycles	Lightweight	6 lbs.	2 huck towels (17 in. x 29 in.), 3 fluid-resistant drapes (108 in. x 72 in.)
	Regular Weight	9 lbs.	16 huck towels (17 in. x 29 in.), 1 fluid resistant table cover (90 in. x 60 in.). 5 lbs. of metal mass.
	Medium Weight	13 lbs.	4 tray liners (20 in. x 25 in.) stacked 23 in. x 11 in. x 3.5 in. tray containing 11 lbs. of metal mass.

	Heavyweight	17 lbs.	4 tray liners (20 in. x 25 in.) stacked 23 in. x 11 in. x 3.5 in. tray containing 15 lbs. of metal mass.
	Super Heavyweight	25 lbs.	4 tray liners (20 in. x 25 in.) stacked 23in. x 11in. x 3.5 in. tray containing 20 lbs. of metal mass.
STERRAD® 50	10.7 lbs.		<p>Reusable metal and non-metal medical devices, including up to 10 lumens of the following dimensions per chamber load:</p> <ul style="list-style-type: none"> • An inside diameter of 1 mm or larger and a length of 125 mm or shorter of single-channel stainless steel lumens. • An inside diameter of 2 mm or larger and a length of 250 mm or shorter of single-channel stainless steel lumens. • An inside diameter of 3 mm or larger and length of 400 mm or shorter of single-channel stainless steel lumens • An inside diameter of 6 mm or larger and a length of 310 mm or shorter of single-channel Teflon-Polyethylene lumens. • Refer to the STERRAD® 50 Sterilizer User’s Guide for complete instructions on load(s) and cycle(s), including chamber loading instructions (i.e., 10 lumens per load).
STERRAD® 200S	10.7 lbs.		<p>Reusable metal and non-metal medical devices, including up to 12 lumens of the following dimensions per chamber load:</p> <ul style="list-style-type: none"> • An inside diameter of 1 mm or larger and a length of 125 mm or shorter of single-channel stainless steel lumens. • An inside diameter of 2 mm or larger and a length of 250 mm or shorter of single-channel stainless steel lumens. • An inside diameter of 3 mm or larger and length of 400 mm or shorter of single-channel stainless steel lumens • An inside diameter of 6 mm or larger and a length of 310 mm or shorter of single-channel Teflon-Polyethylene lumens. <p>Refer to the STERRAD® 200 Sterilizer User’s Guide for complete instructions on load(s) and cycle(s), including chamber loading instructions (i.e., 36.48 lbs. per load).</p>
STERRAD® 100NX™ DUO cycle	10.7 lbs.		<p>One or two single-channel Flexible Endoscope with accessory devices that are normally connected to it, with or without a silicone mat. The flexible endoscope may contain:</p> <ul style="list-style-type: none"> • A single-channel Teflon®/Polyethylene lumen with an inside diameter of 1 mm or larger and a length of 850 mm or shorter. • Accessory devices that are normally connected to a flexible endoscope during use. • Flexible endoscopes without lumens. <p>Refer to the STERRAD® 100NX™ Sterilizer User’s Guide for</p>

		complete instructions on load(s) and cycle(s), including chamber loading instructions (i.e., 13.2 lbs. per load).
--	--	---

TABLE 5: Wrap Model Recommendations for Use with Pre-vacuum steam cycles, Gravity steam cycles, and the Advanced Sterilization Products STERRAD® Sterilization Systems

Gemini Wrap Weight	Gemini Wrap Model	Intended Load	Maximum Recommended Wrapped Package Content ²	
			Pre-Vacuum, Gravity, and EO	ASP STERRAD® 50, 200S and 100NX™ DUO Cycles
Light Weight	GEM11XX	Light weight package (for example: standard linen packs)	6 lbs.	10.7 lbs.
Regular Weight	GEM21XX	Light to moderate weight package (for example: general use medical instruments)	9 lbs.	10.7 lbs.
Medium Weight	GEM31XX	Moderate to heavy weight package (for example: general use medical instruments)	13 lbs.	10.7 lbs.
Heavy Weight	GEM41XX	Heavy weight package (for example: general use medical instruments)	17 lbs.	10.7 lbs.
Super Heavy Weight	GEM51XX	Very heavy weight package (for example: general use medical instruments)	25 lbs.	10.7 lbs.

¹ Individual results may differ due to factors such as variations in handling practices, wrapping techniques, and folding methods. Results may also differ due to the use of irregularly shaped contents, which may put added stress on the wrap. Each healthcare facility should determine for itself which wrap model is most appropriate for each intended use.

² It is recommended to not exceed the maximum wrapped package content weights indicated for each wrap model. Furthermore, it is recommended to not exceed the number, weight, and size of individual content types that were validated for the Gemini Sterilization Wraps (i.e.: the weight of the metal mass).



Medline Industries, Inc.
One Medline Place
Mundelein, IL 60060

K152564

510(k) SUMMARY

[AS REQUIRED BY 21CFR807.92(c)]

Submitter / 510(k) Sponsor

Medline Industries, Inc.
1 Medline Place
Mundelein, IL 60060

Registration Number: 1417592

Contact Person

Jennifer Mason
Senior Regulatory Affairs Specialist
Phone: 847-643-3652
Email: jamason@medline.com

Summary Preparation Date

July 7, 2016

Type of 510(k) Submission

Traditional

Device Name / Classification

Name of Device: Gemini Sterilization Wrap
Proprietary Name: Gemini Sterilization Wrap
Common Name: Sterilization Wrap
Classification Name: Wrap, Sterilization
Product Code: FRG
Classification Panel: General Hospital
Regulation: 21 CFR 880.6850

Predicate Device

Gemini Sterilization Wrap
K113353

Device Description

There have been no changes to the Gemini Sterilization Wrap. The materials of construction, colorants, sizes and product specifications have not changed and are identical to what was cleared under K113353.



Medline Industries, Inc.
One Medline Place
Mundelein, IL 60060

As previously described in K113353, Gemini Sterilization Wrap is offered to the marketplace as bulk packages of single ply sheets for use by customers in accordance with standard hospital practices which require that two sheets are used each time a medical device or collection of medical devices are wrapped.

Gemini Sterilization Wrap items are square or rectangular sheets of fabric produced using a five-layer SSMMS (spunbond-spunbond-meltblown-meltblown-spunbond) process.

The standard blue wrap fabric is made of polypropylene with the addition of less than 2% of phthalocyanine blue pigmentation and less than 0.35% titanium dioxide white pigmentation. The wrap allows for aseptic opening of the sterilized package.

Gemini Sterilization Wrap is available in sizes ranging from 12”x12” to 54”x90” across the five different material weights/models.

Indications for Use

Gemini Sterilization Wrap is intended to be used to enclose another medical device that is to be sterilized by a health care provider. It is intended to allow sterilization of the enclosed medical device and also to maintain sterility of such content until used. Gemini Sterilization Wrap is validated for use in steam or STERRAD® sterilization processes in the following sterilization modes and cycles:

TABLE 1: Pre-vacuum steam at 270°F/132°C for 4 minutes with the following dry times by weight

Gemini Wrap Model	Gemini Wrap Weight	Dry Time
GEM11XX	Lightweight	20 Minutes
GEM21XX	Regular Weight	20 Minutes
GEM31XX	Medium Weight	20 Minutes
GEM41XX	Heavyweight	30 Minutes
GEM51XX	Super Heavyweight	30 Minutes

TABLE 2: Gravity steam at 250°F/121°C for 30 minutes with the following dry times by weight

Gemini Wrap Model	Gemini Wrap Weight	Dry Time
GEM11XX	Lightweight	20 Minutes
GEM21XX	Regular Weight	20 Minutes
GEM31XX	Medium Weight	20 Minutes
GEM41XX	Heavyweight	30 Minutes
GEM51XX	Super Heavyweight	30 Minutes



Medline Industries, Inc.
 One Medline Place
 Mundelein, IL 60060

TABLE 3: STERRAD® Sterilization

Gemini Wrap Model	Gemini Wrap Weight	STERRAD® Cycles
GEM11XX	Lightweight	STERRAD® 50, 200 S and 100NX™ DUO Cycles
GEM21XX	Regular Weight	STERRAD® 50, 200 S and 100NX™ DUO Cycles
GEM31XX	Medium Weight	STERRAD® 50, 200 S and 100NX™ DUO Cycles
GEM41XX	Heavyweight	STERRAD® 50, 200 S and 100NX™ DUO Cycles
GEM51XX	Super Heavyweight	STERRAD® 50, 200 S and 100NX™ DUO Cycles

The Gemini Sterilization Wrap intended loads with pre-vacuum steam sterilization, gravity steam sterilization and the Advanced Sterilization Products STERRAD® Sterilization systems are provided in Table 4.

The Gemini Sterilization Wrap Recommendations for use with Pre-vacuum steam cycles, Gravity steam cycles, and the Advanced Sterilization Products STERRAD® Sterilization systems are provided in Table 5.

TABLE 4: Validated Pre-Vacuum Steam, Gravity Steam, and Advanced Sterilization Products STERRAD® 50, 200S, 100NX™ DUO Cycles

Cycle	Maximum recommended Chamber Load		Intended Load
Pre-Vacuum Steam & Gravity Steam Cycles	Lightweight	6 lbs.	2 huck towels (17 in. x 29 in.), 3 fluid-resistant drapes (108 in. x 72 in.)
	Regular Weight	9 lbs.	16 huck towels (17 in. x 29 in.), 1 fluid resistant table cover (90 in. x 60 in.). 5 lbs. of metal mass.
	Medium Weight	13 lbs.	4 tray liners (20 in. x 25 in.) stacked 23 in. x 11 in. x 3.5 in. tray containing 11 lbs. of metal mass.
	Heavyweight	17 lbs.	4 tray liners (20 in. x 25 in.) stacked 23 in. x 11 in. x 3.5 in. tray containing 15 lbs. of metal mass.
	Super Heavyweight	25 lbs.	4 tray liners (20 in. x 25 in.) stacked 23in. x 11in. x 3.5 in. tray containing 20 lbs. of metal mass.
STERRAD® 50	10.7 lbs.		Reusable metal and non-metal medical devices, including up to 10 lumens of the following dimensions per chamber load: <ul style="list-style-type: none"> • An inside diameter of 1mm or larger and a length of 125 mm or shorter of single-channel stainless steel lumens. • An inside diameter of 2 mm or larger and a length of 250 mm or shorter of single-channel stainless steel lumens.



Medline Industries, Inc.
 One Medline Place
 Mundelein, IL 60060

		<ul style="list-style-type: none"> • An inside diameter of 3mm or larger and length of 400 mm or shorter of single-channel stainless steel lumens • An inside diameter of 6 mm or larger and a length of 310 mm or shorter of single-channel Teflon-Polyethylene lumens. • Refer to the STERRAD® 50 Sterilizer User’s Guide for complete instructions on load(s) and cycle(s), including chamber loading instructions (i.e., 10 lumens per load).
STERRAD® 200 S	10.7 lbs.	<p>Reusable metal and non-metal medical devices, including up to 12 lumens of the following dimensions per chamber load:</p> <ul style="list-style-type: none"> • An inside diameter of 1mm or larger and a length of 125 mm or shorter of single-channel stainless steel lumens. • An inside diameter of 2 mm or larger and a length of 250 mm or shorter of single-channel stainless steel lumens. • An inside diameter of 3mm or larger and length of 400 mm or shorter of single-channel stainless steel lumens • An inside diameter of 6 mm or larger and a length of 310 mm or shorter of single-channel Teflon-Polyethylene lumens. <p>Refer to the STERRAD® 200 Sterilizer User’s Guide for complete instructions on load(s) and cycle(s), including chamber loading instructions (i.e., 36.48 lbs. per load).</p>
STERRAD® 100NX™ DUO cycle	10.7 lbs.	<p>One or two single-channel Flexible Endoscope with accessory devices that are normally connected to it, with or without a silicone mat. The flexible endoscope may contain:</p> <ul style="list-style-type: none"> • A single-channel Teflon®/Polyethylene lumen with an inside diameter of 1 mm or larger and a length of 850 mm or shorter. • Accessory devices that are normally connected to a flexible endoscope during use. • Flexible endoscopes without lumens. <p>Refer to the STERRAD® 100NX™ Sterilizer User’s Guide for complete instructions on load(s) and cycle(s), including chamber loading instructions (i.e., 13.2 lbs. per load).</p>



Medline Industries, Inc.
One Medline Place
Mundelein, IL 60060

TABLE 5: Wrap Model Recommendations for Use with Pre-vacuum steam cycles, Gravity steam cycles, and the Advanced Sterilization Products STERRAD® Sterilization Systems

Gemini Wrap Weight	Gemini Wrap Model	Intended Load	Maximum Recommended Wrapped Package Content ²	
			Pre-Vacuum, Gravity, and EO	ASP STERRAD® 50, 200 S and 100NX™ DUO Cycles
Light Weight	GEM11XX	Light weight package (for example: standard linen packs)	6 lbs.	10.7 lbs.
Regular Weight	GEM21XX	Light to moderate weight package (for example: general use medical instruments)	9 lbs.	10.7 lbs.
Medium Weight	GEM31XX	Moderate to heavy weight package (for example: general use medical instruments)	13 lbs.	10.7 lbs.
Heavy Weight	GEM41XX	Heavy weight package (for example: general use medical instruments)	17 lbs.	10.7 lbs.
Super Heavy Weight	GEM51XX	Very heavy weight package (for example: general use medical instruments)	25 lbs.	10.7 lbs.

¹ Individual results may differ due to factors such as variations in handling practices, wrapping techniques, and folding methods. Results may also differ due to the use of irregularly shaped contents, which may put added stress on the wrap. Each healthcare facility should determine for itself which wrap model is most appropriate for each intended use.

² It is recommended to not exceed the maximum wrapped package content weights indicated for each wrap model. Furthermore, it is recommended to not exceed the number, weight, and size of individual content types that were validated for the Gemini Sterilization Wraps (i.e.: the weight of the metal mass).

Maintenance of Sterility (Shelf Life)

The Gemini Sterilization Wrap has a shelf life of two years following steam sterilization and 180 days following STERRAD® sterilization.



Medline Industries, Inc.
One Medline Place
Mundelein, IL 60060

Validated Sterilization Cycles

Gemini Sterilization Wrap is validated for use in steam or STERRAD® sterilization processes in the following sterilization modes and cycles:

TABLE 6: Pre-vacuum steam at 270°F/132°C for 4 minutes with the following dry times by weight (cleared under K113353)

Gemini Wrap Model	Gemini Wrap Weight	Dry Time
GEM11XX	Lightweight	20 Minutes
GEM21XX	Regular Weight	20 Minutes
GEM31XX	Medium Weight	20 Minutes
GEM41XX	Heavyweight	30 Minutes
GEM51XX	Super Heavyweight	30 Minutes

TABLE 7: Gravity steam at 250°F/121°C for 30 minutes with the following dry times by weight (cleared under K113353)

Gemini Wrap Model	Gemini Wrap Weight	Dry Time
GEM11XX	Lightweight	20 Minutes
GEM21XX	Regular Weight	20 Minutes
GEM31XX	Medium Weight	20 Minutes
GEM41XX	Heavyweight	30 Minutes
GEM51XX	Super Heavyweight	30 Minutes

Advanced Sterilization Products STERRAD® Sterilization System

- STERRAD® 50, 200 S and 100NX™ DUO – Subject of this 510(k)



Medline Industries, Inc.
 One Medline Place
 Mundelein, IL 60060

Summary of Technological Characteristics

TABLE 8: COMPARISON OF PROPOSED AND PREDICATE DEVICES

Device Characteristic	Proposed Device	Predicate Device	Comparison Analysis
Product Name	Gemini Sterilization Wrap	Gemini Sterilization Wrap	Same
510(k) Reference		K113353	
Product Owner	Medline Industries, Inc.	Medline Industries, Inc.	Same
Product Code	FRG	FRG	Same
Intended Use	Gemini Sterilization Wrap is intended to be used to enclose another medical device that is to be sterilized by a health care provider. It is intended to allow sterilization of the enclosed medical device and also to maintain sterility of such content.	Gemini Sterilization Wrap is intended to be used to enclose another medical device that is to be sterilized by a health care provider. It is intended to allow sterilization of the enclosed medical device and also to maintain sterility of such content.	Same
Regulation Number	21 CFR 880.6850	21 CFR 880.6850	Same
Design Features	Square or rectangular sheets mfd. by spunbond-meltblown process	Square or rectangular sheets mfd. by spunbond-meltblown process	Same
Design Configurations	12 in. x 12 in. to 54 in. x 90 in.	12 in. x 12 in. to 54 in. x 90 in.	Same
Materials	Polypropylene with phthalocyanine blue and titanium dioxide	Polypropylene with phthalocyanine blue and titanium dioxide	Same
Wrapping Technique	Sequential/simultaneous double wrapping	Sequential/simultaneous double wrapping	Same
Prescription vs. OTC	OTC	OTC	Same
Sterilization	Pre-Vacuum Steam Cycle: 4 Minutes Exposure at 270°F/132°C with minimum 20 minutes dry time	Pre-Vacuum Steam Cycle: 4 Minutes Exposure at 270°F/132°C with minimum 20 minutes dry time	Same
	Gravity Steam Cycle: 30 Minutes Exposure at 250°F/121°C with minimum 20 minutes dry time	Gravity Steam Cycle: 30 Minutes Exposure at 250°F/121°C with minimum 20 minutes dry time	Same
	STERRAD® 50, 200 S and 100NX™ DUO cycle	N/A	Different



Medline Industries, Inc.
 One Medline Place
 Mundelein, IL 60060

	STERRAD® NX™, (Standard Cycle, Advanced Cycle)	STERRAD® NX™, (Standard Cycle, Advanced Cycle)	Same
	STERRAD® 100S, Standard Cycle	STERRAD® 100S, Standard Cycle	Same
	STERRAD® 100NX™, (Standard Cycle, Flex Cycle, EXPRESS Cycle)	STERRAD® 100NX™, (Standard Cycle, Flex Cycle, EXPRESS Cycle)	Same
Maintenance of Sterility (Shelf Life)	Steam – 2 years	Steam – 180 days	Different
	STERRAD® – 180 days	STERRAD® – 180 days	Same
Single Use vs. Reusable	Single Use	Single Use	Same

Discussion of Similarities and Differences

The proposed Gemini Sterilization Wrap is substantially equivalent in intended use, materials, device features / specifications and function in comparison to the predicate (K113353) Gemini Sterilization Wrap. As noted in the table above, these characteristic comparisons render the devices identical. The difference between subject/predicate devices is related to:

- Maintenance of Sterility (Shelf Life): Extending the maintenance of sterility from 180 days to two years for steam sterilization based on real time stability studies. The completed stability studies were identical to the in-process stability studies described within K113353.
- Addition of STERRAD® Systems: STERRAD® 50 Cycle, STERRAD® 200 S Cycle and STERRAD® 100NX™ Duo Cycle.

This product modification involving the subject and predicate devices does not significantly alter the Gemini Sterilization Wrap or raise questions regarding safety or effectiveness.

Information within this Premarket Notification demonstrates that there are no significant differences in technological characteristics between Medline’s Gemini Sterilization Wrap and the cited predicate device K113353.



Medline Industries, Inc.
 One Medline Place
 Mundelein, IL 60060

Summary of Non-Clinical Testing

TABLE 9: Summary of Performance Testing

Study		Performance Results
Maintenance of Sterility (Package Integrity)	Steam Pre-Vacuum & Steam Gravity Cycles (Real Time Event Related Two Years Shelf Life)	Passed
	STERRAD® (Real Time Event Related 180 Days Shelf Life)	Passed
Pre-Vacuum Steam and Gravity Sterilant Penetration		Passed
STERRAD® 50 Cycle, STERRAD® 200 S Cycle and STERRAD® 100NX™ DUO Cycle) Sterilant Penetration		Passed
Post Sterilization Biocompatibility Testing (Primary Skin Irritation Testing - ISO 10993-10)		Passed

The following performance testing was conducted on the lightest weight and the heaviest weight Gemini Sterilization Wrap both pre-sterilization and after exposure to a STERRAD® 100NX™ DUO cycle.

Air permeability – per ASTM D737-04 R2012

Basis weight – per ASTM D3776

Material burst strength – per ASTM D3786

Resistance to water penetration – per AATCC127

Tensile strength and elongation – per ASTM D5034

Tear strength – per ASTM D5587

The safety and effectiveness of Medline’s Gemini Sterilization Wrap is adequately supported by the substantial equivalence information, materials information, and Design Control activities referenced within this Premarket Notification.

Summary of Clinical Testing

Not applicable.

Basis for Determination of Substantial Equivalence/Conclusion

Based on the intended use, technological characteristics, performance data, and nonclinical tests performed, the subject Gemini Sterilization Wrap are substantially equivalent to, and are as safe and as effective as, the legally marketed predicate device, Gemini Sterilization Wrap, cleared under K113353 under regulation 21 CFR 880.6850, product code FRG.