



Food and Drug Administration
10903 New Hampshire Avenue
Document Control Center – WO66-G609
Silver Spring, MD 20993-0002

G.N.S neoLaser Ltd.
Mr. Gil Shapira
Chief Executive Officer
7 Ha'Eshel Str. 38900
Caesarea ISRAEL, P.O. Box 3203

December 8, 2015

Re: K152722
Trade/Device Name: neoV980 & neoV1470 Diode Lasers
Regulation Number: 21 CFR 878.4810
Regulation Name: Laser surgical instrument for use in general and
plastic surgery and in dermatology
Regulatory Class: Class II
Product Code: GEX
Dated: July 15, 2015
Received: September 21, 2015

Dear Mr. Shapira:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you; however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the [Federal Register](#).

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set

forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to

<http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

Joshua C. Nipper -S

For Binita S. Ashar, M.D., M.B.A., F.A.C.S.
Director
Division of Surgical Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

K152722

Device Name

neoV980 & neoV1470 Diode Lasers

Indications for Use (Describe)

The neoV980 & neoV1470 Diode Lasers, (and their delivery accessories used to deliver optical energy) are indicated for use in surgical applications requiring the vaporization, incision, excision, ablation, cutting and hemostasis, or coagulation of soft tissue in conjunction with endoscopic equipment for medical specialties including: Urology, Thoracic Surgery, Plastic Surgery and Dermatology, General Surgery, Ophthalmology, Orthopedics, Podiatry, Arthroscopy, Spinal Surgery, Gynecology, Pulmonary Surgery, Neurosurgery, Gastroenterology, Head/neck/ENT and Radiology, Oral Surgery and Dental procedures, Endovascular coagulation, and endovenous occlusion of the greatest saphenous vein in patients with superficial vein reflux.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

	neoV FDA 510(k) Submission K152722 Section 5 – 510(k) Summary	Revision #2
Ha'Eshel St. 7, 38900, Caesarea ISRAEL	SOP Number: DHF-001-R-025-5	Effective Date: July 15, 2015
PO Box 3203, Phone:+972 4 677 9919, Fax:+972 4 8591505	Page 10 Section 5	

510(k) SUMMARY

Title: neoV980 & neoV1470 Diode Lasers

Submitter: G.N.S neoLaser Ltd.
Ha'Eshel St. 7, 38900, Caesarea ISRAEL, PO Box 3203

Contact: Gil Shapira, CEO
G.N.S neoLaser Ltd.
Ha'Eshel St. 7, 38900, Caesarea ISRAEL, PO Box 3203
Phone: +972 52 2246965
Fax: +972 4 859 1505
Email: shapirag@neo-laser.com

Date
Prepared: July 15, 2015

Device Trade
Name: neoV980 & neoV1470 Diode Lasers

Common
Name: Laser surgical instrument for use in general surgery and dermatology

Classification
Name: Instrument, surgical, powered laser
GEX
21 CFR 878.4810

Predicate
Devices: Quanta Diode Laser Family (K100558), Ceralas 1470nm Diode Laser (K112253) and Ceralas 980nm Diode Laser (K112324)

Device
Description: neoV980 & neoV1470, members of the neoV Diode Lasers family, cleared under K133006, are medical grade, solid-state, infrared Diode lasers, designed to deliver continuous or pulsed, infrared laser energy at wavelengths of 980nm and 1470nm respectively, with power level ranging from 10 to 20 Watts respectively. The lasers are controlled via a high-resolution color touch screen. The touch screen display includes a user interface allowing selection of continuous, repeat pulse, or single pulse modes of operation as well as repetition rates, aiming beam settings, password key protection, and standby/ready mode selection. The units have an emergency shut off button on

	neoV980 & neoV1470 Diode Lasers 510(k) Submission 510(k) Summary	Revision #2	
Ha'Eshel St. 7, 38900, Caesarea ISRAEL	SOP Number: DHF-001-R-025-	Effective Date: July 15, 2015	Page 11
PO Box 3203, Phone:+972 4 677 9919, Fax:+972 4 8591505			Section 5

the front of the unit.

The Laser System: The laser system consists of an optical block which contains the laser diode, mirrors, lens, and aiming beam diode, an air cooling system, and electronics which include the color touch screen control panel. The unit utilizes an external low voltage power supply, as well as an external wired foot switch for laser activation.

The Delivery System: The delivery system consists of either sterile fibers (not provided with the system) or non-sterile fibers and hand pieces. Safety goggles, non sterile fibers, handpieces and a safety sign are provided with the unit.

**Intended
Use:**

The neoV980 & neoV1470 Diode lasers are indicated for use in surgical applications requiring the vaporization, incision, excision, ablation, cutting and hemostasis, or coagulation of soft tissue in conjunction with endoscopic equipment for medical specialties including: Urology, Thoracic Surgery, Plastic Surgery and Dermatology, General Surgery, Ophthalmology, Orthopedics, Podiatry, Arthroscopy, Spinal Surgery, Gynecology, Pulmonary Surgery, Neurosurgery, Gastroenterology, Head/neck/ENT and Radiology, Oral Surgery and Dental procedures, endovascular coagulation and endovenous occlusion of the greatest saphenous vein in patients with superficial vein reflux.

Substantial

Equivalence: The neoV980 & neoV1470 Diode Lasers share the same intended use, technical and performance characteristics of the predicate devices Quanta System QUANTA Diode Laser Family (K100558), the Ceralas 1470nm Diode Laser (K112253) and the Ceralas 980nm Diode Laser (K112324).

Conclusion: From a design and clinical perspective, the predicates and subject laser devices, have the same technological characteristics and share the same intended use. Accordingly, the safety and effectiveness of the neoV980 & neoV1470, is based upon a determination of the substantial equivalence to the predicate devices.

	<p>neoV980 & neoV1470 Diode Lasers 510(k) Submission</p> <p>510(k) Summary</p>		<p>Revision #2</p>
<p>Ha'Eshel St. 7, 38900, Caesarea ISRAEL</p>	<p>SOP Number: DHF-001-R-025-</p>	<p>Effective Date: July 15, 2015</p>	<p>Page 12 Section 5</p>
<p>PO Box 3203, Phone:+972 4 677 9919, Fax:+972 4 8591505</p>			

Non clinical
Performance

Data: The neoV Diode Laser family including neoV980 & neoV1470 has been tested for compliance to IEC 60601-2-22: 2007 (Third Edition), Medical electrical equipment Part 2: Particular requirements for the safety of diagnostic and therapeutic laser equipment, IEC 60601-1: 2005 (Third Edition) Medical electrical equipment Part 1: General requirements for basic safety and essential performance and tested for compliance with all functional requirements, IEC 60825-1:2007 Safety of Laser Products – Part 1: Equipment classification, requirements and user's guide, EN 60601-1-2: 2007, Medical electrical equipment Part 1-2: General requirements for safety - Collateral standard: Electromagnetic compatibility – Requirements and tests.

Clinical
Performance

Data: None