



Food and Drug Administration
10903 New Hampshire Avenue
Document Control Center – WO66-G609
Silver Spring, MD 20993-0002

Centers For Disease Control And Prevention
Yon Yu, Pharm. D.
Associate Director For Regulatory Affairs
1600 Clifton Road, NE MS-C18
Atlanta, GA 30329-4027

December 1, 2015

Re: K153148

Trade/Device Name: CDC Human Influenza Virus Real- Time RT-PCR Diagnostic Panel,
Influenza A/H5 Subtyping Kit (VER 3)

Regulation Number: 21 CFR 866.3980

Regulation Name: Respiratory viral panel multiplex nucleic acid assay

Regulatory Class: II

Product Code: OZE, NSU, NXD

Dated: October 29, 2015

Received: October 30, 2015

Dear Dr. Yu:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulations (21 CFR Parts 801 and 809), please contact the Division of Industry and Consumer Education at its toll-free number (800) 638 2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, “Misbranding by reference to premarket notification” (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH’s Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

Tamara V. Feldblyum -S for

Uwe Scherf, M.Sc., Ph.D.
Director
Division of Microbiology Devices
Office of In Vitro Diagnostics
and Radiological Health
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

K153148

Device Name

CDC Human Influenza Virus Real-Time RT-PCR Diagnostic Panel, Influenza A/H5 Subtyping Kit (VER 3)

Indications for Use (Describe)

The Influenza A/H5 Subtyping Kit contains reagents and controls of the CDC Human Influenza Virus Real-Time RT-PCR Diagnostic Panel and is intended for use in real-time RT-PCR (rRT-PCR) assays on an Applied Biosystems (ABI) 7500 Fast Dx Real-Time PCR instrument in conjunction with clinical and epidemiological information:

- For the presumptive identification of virus in patients who may be infected with influenza A subtype A/H5 (Asian lineage) from viral RNA in human respiratory specimens and viral culture in conjunction with clinical and epidemiological risk factors;
- To provide epidemiologic information for surveillance of circulating influenza viruses.

Performance characteristics for influenza were established during a season when seasonal influenza viruses A/H1 and A/H3 were the predominant influenza A viruses in circulation and during a season when the A/H1pdm09 influenza virus was the predominant influenza A virus in circulation. Performance characteristics may vary with other emerging influenza A viruses.

Testing with the influenza H5a and H5b primer and probe sets should not be performed unless the patient meets the most current U.S. Department of Health and Human Services (DHHS) clinical and epidemiologic criteria for testing suspect A/H5 specimens. The definitive identification of influenza A/H5 (Asian lineage) either directly from patient specimens or from virus cultures requires additional laboratory testing, along with clinical and epidemiological assessment in consultation with national influenza surveillance experts.

Negative results do not preclude influenza virus infection and should not be used as the sole basis for treatment or other patient management decisions. Conversely, positive results do not rule out bacterial infection or co-infection with other viruses. The agent detected may not be the definite cause of disease.

If infection with a novel influenza A virus is suspected based on current clinical and epidemiological screening criteria recommended by public health authorities, specimens should be collected with appropriate infection control precautions for novel virulent influenza viruses and sent to a state or local health department for testing. Viral culture should not be attempted unless a BSL 3+ facility is available to receive and culture specimens.

All users, analysts, and any person reporting results from use of this device should be trained to perform and interpret the results from this procedure by a competent instructor prior to use. CDC Influenza Division will limit the distribution of this device to only those users who have successfully completed a training course provided by CDC instructors or designees.

Type of Use (Select one or both, as applicable)

☒ Prescription Use (Part 21 CFR 801 Subpart D)

☐ Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

9. **510(k) Summary**

I. GENERAL INFORMATION

Submitter:

Centers for Disease Control and Prevention
1600 Clifton Road, NE
Atlanta, GA 30333

Contact Person:

CDC Yon Yu, Pharm. D.
Associate Director for Regulatory Affairs
Office of the Director
National Center for Emerging and Zoonotic Infectious Diseases
Centers for Disease Control and Prevention
1600 Clifton Road, MS-C18
Atlanta, GA 30329-4027
Phone: 404-639-3046
Fax: 404-639-1275
Email: fk8@cdc.gov

Date Prepared: October 22, 2015

II. DEVICE INFORMATION

| | |
|--|---|
| Proprietary Name: | CDC Human Influenza Virus Real-Time RT-PCR Diagnostic Panel, Influenza A/H5 Subtyping Kit (VER 3) |
| Common Name: | Influenza A/H5 Subtyping Kit |
| Regulation Section: | 866.3980-Respiratory viral panel multiplex nucleic acid assay |
| Subsequent Regulation Sections: | 862.2570-Instrumentation for clinical multiplex test systems 866.3332-Reagents for detection of specific novel influenza A viruses |
| Device Classification: | Class II |
| Product Code: | OZE |
| Subsequent Product Codes: | NSU, NXD |
| Panel: | Microbiology |

III. PREDICATE DEVICE

CDC Human Influenza Virus Real-Time RT-PCR Diagnostic Panel, Influenza A/H5 Subtyping Kit (K141859)

IV. DEVICE DESCRIPTION

The CDC Human Influenza Real-Time RT-PCR Diagnostic Panel Influenza A/H5 Subtyping Kit is a real-time RT-PCR (rRT-PCR) assay that utilizes the Applied Biosystems® (ABI) 7500 Fast Dx Real-time PCR system. The Influenza A/H5 Subtyping Kit consists of oligonucleotide primers, fluorescently labeled hydrolysis probes, and controls which are used in rRT-PCR assays for the in vitro qualitative detection and characterization of influenza virus RNA in respiratory specimens from patients presenting with influenza-like illness (ILI). The oligonucleotide primers and probes for detection of influenza A viruses (InfA) were selected from highly conserved regions of the matrix (M) protein. Oligonucleotide primers and probes for characterization and differentiation of avian influenza A/H5 (Asian lineage) viruses (H5a and H5b) were selected from highly conserved regions of their HA genes. The Influenza A/H5 Subtyping kit also contains primers and probes to detect the human RNase P gene (RP) in control samples and clinical specimens.

V. INTENDED USE

The Influenza A/H5 Subtyping Kit contains reagents and controls of the CDC Human Influenza Virus Real-Time RT-PCR Diagnostic Panel and is intended for use in real-time RT-PCR (rRT-PCR) assays on an Applied Biosystems (ABI) 7500 Fast Dx Real-Time PCR instrument in conjunction with clinical and epidemiological information:

- For the presumptive identification of virus in patients who may be infected with influenza A subtype A/H5 (Asian lineage) from viral RNA in human respiratory specimens and viral culture in conjunction with clinical and epidemiological risk factors;
- To provide epidemiologic information for surveillance of circulating influenza viruses.

Performance characteristics for influenza were established during a season when seasonal influenza viruses A/H1 and A/H3 were the predominant influenza A viruses in circulation and during a season when the A/H1pdm09 influenza virus was the predominant influenza A virus in circulation. Performance characteristics may vary with other emerging influenza A viruses.

Testing with the influenza H5a and H5b primer and probe sets should not be performed unless the patient meets the most current U.S. Department of Health and Human Services (DHHS) clinical and epidemiologic criteria for testing suspect A/H5 specimens. The definitive identification of influenza A/H5 (Asian lineage) either directly from patient specimens or from virus cultures requires additional laboratory testing, along with clinical and epidemiological assessment in consultation with national influenza surveillance experts.

Negative results do not preclude influenza virus infection and should not be used as the sole basis for treatment or other patient management decisions. Conversely, positive results do not rule out bacterial infection or co-infection with other viruses. The agent detected may not be the definite cause of disease.

If infection with a novel influenza A virus is suspected based on current clinical and epidemiological screening criteria recommended by public health authorities, specimens should be collected with appropriate infection control precautions for novel virulent influenza viruses and sent to a state or

local health department for testing. Viral culture should not be attempted unless a BSL 3+ facility is available to receive and culture specimens.

All users, analysts, and any person reporting results from use of this device should be trained to perform and interpret the results from this procedure by a competent instructor prior to use. CDC Influenza Division will limit the distribution of this device to only those users who have successfully completed a training course provided by CDC instructors or designees.

VI. TECHNOLOGICAL CHARACTERISTICS

As circulating Asian lineage highly pathogenic avian influenza (HPAI) viruses continue to evolve, the HA gene target of certain phylogenetic groups has acquired changes in the oligonucleotide primer and probe regions that require modification to ensure comprehensive detection of circulating influenza A/H5 clades. The predicate Influenza A/H5 Subtyping Kit (K141859) includes reagents for two assays, referred to as H5a and H5b, specifically detecting influenza A H5 Asian lineage viruses. The predicate H5a assay consists of two primers and two probes. The predicate H5b assay consists of two primers and one probe. The oligonucleotide primers and probes used in the modified H5a and H5b assays target the HA gene at the same locations and include a total of five primers and three probes for H5a and three primers and one probe for H5b. In addition to evaluating the modified oligonucleotide primers and probes, CDC has evaluated the ZENTM Double-Quenched Probe technology (manufactured by Integrated DNA Technologies) as an alternate fluorescent hydrolysis probe quencher chemistry. The Influenza A/H5 Subtyping Kit assays contain ZENTM double quenched probes (ZEN probes) that include an internal ZENTM quencher located nine nucleotides away from the 5' FAM reporter dye in addition to an Iowa Black[®] FQ quencher (IABkFQ) at the 3' end of the probe.

VII. SUBSTANTIAL EQUIVALENCE COMPARISON

The CDC Human Influenza Virus Real-Time RT-PCR Diagnostic Panel, Influenza A/H5 Subtyping Kit (K141859) will serve as the predicate for the intended change. See the table below for a detailed comparison.

Device Comparison

| | CDC Human Influenza Virus Real-Time RT-PCR Diagnostic Panel Diagnostic Panel, Influenza A/H5 Subtyping Kit (K141859) | CDC Human Influenza Virus Real-Time RT-PCR Diagnostic Panel Diagnostic Panel, Influenza A/H5 Subtyping Kit |
|---------------------|---|---|
| Intended Use | <p>The Influenza A/H5 Subtyping Kit contains reagents and controls of the CDC Human Influenza Virus Real-Time RT-PCR Diagnostic Panel and is intended for use in real-time RT-PCR (rRT-PCR) assays on an Applied Biosystems (ABI) 7500 Fast Dx Real-Time PCR instrument in conjunction with clinical and epidemiological information:</p> <ul style="list-style-type: none"> For the presumptive identification of virus in patients who may be infected with influenza A subtype A/H5 (Asian lineage) from viral RNA in human respiratory specimens and viral culture in conjunction with clinical and epidemiological risk factors; To provide epidemiologic information for surveillance of circulating influenza viruses. | Same |

| | | |
|---------------------------------|--|------|
| | <p>Performance characteristics for influenza were established during a season when seasonal influenza viruses A/H1 and A/H3 were the predominant influenza A viruses in circulation and during a season when the A/H1pdm09 influenza virus was the predominant influenza A virus in circulation. Performance characteristics may vary with other emerging influenza A viruses.</p> <p>Testing with the influenza H5a and H5b primer and probe sets should not be performed unless the patient meets the most current U.S. Department of Health and Human Services (DHHS) clinical and epidemiologic criteria for testing suspect A/H5 specimens. The definitive identification of influenza A/H5 (Asian lineage) either directly from patient specimens or from virus cultures requires additional laboratory testing, along with clinical and epidemiological assessment in consultation with national influenza surveillance experts.</p> <p>Negative results do not preclude influenza virus infection and should not be used as the sole basis for treatment or other patient management decisions. Conversely, positive results do not rule out bacterial infection or co-infection with other viruses. The agent detected may not be the definite cause of disease.</p> <p>If infection with a novel influenza A virus is suspected based on current clinical and epidemiological screening criteria recommended by public health authorities, specimens should be collected with appropriate infection control precautions for novel virulent influenza viruses and sent to state or local health department for testing. Viral culture should not be attempted unless a BSL 3+ facility is available to receive and culture specimens.</p> <div style="border: 1px solid black; padding: 5px;"> <p>All users, analysts, and any person reporting results from use of this device should be trained to perform and interpret the results from this procedure by a competent instructor prior to use. CDC Influenza Division will limit the distribution of this device to only those users who have successfully completed a training course provided by CDC instructors or designees.</p> </div> | |
| Organism Detected | Universal influenza A viruses (animal and human) and Influenza A/H5 Subtype (Asian lineage) viruses. | Same |
| Specimen Types | Human respiratory specimens and viral culture. | Same |
| Nucleic Acid Extraction | Yes | Same |
| Extraction Method | <ul style="list-style-type: none"> • QIAamp® DSP Viral RNA Mini Kit, Qiagen • MagNA Pure Compact –Nucleic Acid Isolation Kit I, Roche • MagNA Pure Compact – RNA Isolation Kit, Roche • MagNA Pure LC – Total Nucleic Acid Kit, Roche • Qiagen QIAcube – QIAamp® DSP Viral RNA Mini Kit, Qiagen • NucliSENS® easyMAG®, bioMérieux | Same |
| Enzyme Master Mix | <p>Invitrogen SuperScript™ III Platinum® One-Step Quantitative RT-PCR Kit (with or without ROX)</p> <p>OR</p> <p>Quanta BioSciences qScript™ One-Step qRT-PCR Kit, Low ROX</p> | Same |
| PCR Technology | Real-Time RT-PCR | Same |
| Required Instrumentation | Applied Biosystems 7500 Fast Dx Real-Time PCR Instrument with SDS software version 1.4 | Same |

| | | |
|---------------------------------|---|---|
| Probe Quenching Molecule | Black Hole Quencher Probe® (BHQ-1) | ZEN™ Double-Quenched Probe (InfA, H5a, H5b, and RP assays) OR Black Hole Quencher Probe® (InfA and RP assays) |
| Oligonucleotides | H5a assay-Targets a region of the HA gene H5b assay-Targets a region of the HA gene InfA assay-Targets a conserved region of the matrix gene in Influenza A viruses | Gene targets of the oligonucleotide assays are the same as the predicate; minor changes to the oligonucleotide sequences have been made |

VIII. ANALYTICAL PERFORMANCE EVALUATION

Analytical Sensitivity - Limit of Detection Study (LOD)

Analytical sensitivity of the Influenza A/H5 Subtyping Kit was demonstrated by determining the LOD using Quanta qScript™ and Invitrogen SuperScript™ enzyme kits. Characterized viruses of known 50% infectious dose titers (EID₅₀/mL) were extracted, and the RNA was serially diluted and tested (n=3 replicates) in order to determine an apparent endpoint range. The LOD for each primer and probe set was confirmed by testing extraction replicates (n=20) of the highest virus dilution where ≥95% of all replicates tested positive. Virus dilutions were prepared in virus transport medium containing human A549 cells to emulate clinical specimen matrix. The lowest concentration where the InfA and both H5a and H5b primer and probe sets demonstrate uniform detection was reported as the LOD. The results are summarized in the table below.

LOD Summary

| Influenza Virus Tested | Influenza Strain Designation | LOD (EID ₅₀ /mL) | |
|------------------------|---|-----------------------------|--------------------|
| | | Invitrogen SuperScript™ | Quanta qScript™ |
| A/H5N1 | A/Vietnam/1203/2004×A/Puerto Rico/8/34 reassortant (A/Vietnam/1203/2004 PR8- VN H5N1- PR8/CDC-RG) | 10 ^{3.8} | 10 ^{2.4} |
| | A/duck/Vietnam/NCVD-1544/2012 | 10 ^{3.1} | 10 ^{3.1} |
| A/H5N8 | A/gyrfalcon/Washington/41088-6/2014 | 10 ^{3.35} | 10 ^{3.35} |

Analytical Sensitivity - Inclusivity Testing

Inclusivity testing was conducted to demonstrate the capability of the oligonucleotide primers and probes in the Influenza A/H5 Subtyping Kit to detect strains of influenza A/H5 viruses (Asian lineage) representative of different geographic locations and phylogenetic clades at or near the established LOD. Inclusivity testing was performed with sixteen representative H5 viruses (Asian lineage). A virus of the phylogenetic clade 2.2.2.1 was unavailable for testing, therefore reactivity of the probes with A/Bangladesh/3222/2011 was performed *in silico*. The remaining fifteen viruses were grown to high titer, harvested, and serially diluted to near the LOD of the assays. The diluted influenza A/H5 viruses were extracted and tested in triplicate with the InfA, H5a, and H5b assays to demonstrate

reactivity. Inclusivity of the Influenza A/H5 Subtyping Kit was evaluated with both enzyme systems (i.e. Invitrogen SuperScript™ and Quanta qScript™) and one cleared extraction method.

The Influenza A/H5 Subtyping Kit was reactive with all H5 (Asian lineage) isolates that were tested and predicted to be reactive with the influenza A/Bangladesh/3222/2011 *in silico*. The inclusivity results are summarized in the table below.

Inclusivity Results of the Influenza A/H5 Subtyping Kit

| Influenza Virus Strain Identification | Subtype | Clade | EID ₅₀ /mL | Invitrogen SuperScript™ | | | Quanta qScript™ | | |
|---|---------|----------|-----------------------|-------------------------|---------|---------|-----------------|---------|---------|
| | | | | InfA | H5a | H5b | InfA | H5a | H5b |
| A/Northern pintail/Washington/40964/2014 | H5N2 | 2.3.4.4 | 10 ^{8.4} | 3/3 (+) | 3/3 (+) | 3/3 (+) | 3/3 (+) | 3/3 (+) | 3/3 (+) |
| A/common magpie/Hong Kong/645/2006 | H5N1 | 2.3.4 | 10 ^{8.2} | 3/3 (+) | 3/3 (+) | 3/3 (+) | 3/3 (+) | 3/3 (+) | 3/3 (+) |
| A/chicken/Bangladesh/11rs-1984-30/2011 | H5N1 | 2.3.4.2 | 10 ^{8.75} | 3/3 (+) | 3/3 (+) | 3/3 (+) | 3/3 (+) | 3/3 (+) | 3/3 (+) |
| A/chicken/Vietnam/NCVD-279/2009 | H5N1 | 2.3.4.3 | 10 ^{4.25} | 3/3 (+) | 3/3 (+) | 3/3 (+) | 3/3 (+) | 3/3 (+) | 3/3 (+) |
| A/Cambodia/W0526301/2012 | H5N1 | 1.1 | 10 ^{3.4} | 3/3 (+) | 3/3 (+) | 3/3 (+) | 3/3 (+) | 3/3 (+) | 3/3 (+) |
| A/Cambodia/R040505/2007 | H5N1 | 1.1 | 10 ^{3.5} | 3/3 (+) | 3/3 (+) | 3/3 (+) | 3/3 (+) | 3/3 (+) | 3/3 (+) |
| A/duck/Hunan/795/2002 | H5N1 | 2.1.1 | 10 ^{3.9} | 3/3 (+) | 3/3 (+) | 3/3 (+) | 3/3 (+) | 3/3 (+) | 3/3 (+) |
| A/Indonesia/NIHRD11771/2011 | H5N1 | 2.1.3.2a | 10 ^{4.4} | 3/3 (+) | 3/3 (+) | 3/3 (+) | 3/3 (+) | 3/3 (+) | 3/3 (+) |
| A/turkey/Turkey/1/2005 | H5N1 | 2.2.1 | 10 ^{4.1} | 3/3 (+) | 3/3 (+) | 3/3 (+) | 3/3 (+) | 3/3 (+) | 3/3 (+) |
| A/Egypt/1050/NAMRU3/2013 | H5N1 | 2.2.1 | 10 ^{3.4} | 3/3 (+) | 3/3 (+) | 3/3 (+) | 3/3 (+) | 3/3 (+) | 3/3 (+) |
| A/Bangladesh/3233/2011 (<i>in silico</i> analysis) | H5N1 | 2.2.2.1 | -- | + | + | + | + | + | + |
| A/common magpie/Hong Kong/5052/2007 | H5N1 | 2.3.2.1 | 10 ^{4.75} | 3/3 (+) | 3/3 (+) | 3/3 (+) | 3/3 (+) | 3/3 (+) | 3/3 (+) |
| A/chicken/Bangladesh/42010/2012 | H5N1 | 2.3.2.1a | 10 ^{8.2} | 3/3 (+) | 3/3 (+) | 3/3 (+) | 3/3 (+) | 3/3 (+) | 3/3 (+) |
| A/duck/Vietnam/NCVD-672/2011 | H5N1 | 2.3.2.1b | 10 ^{3.8} | 3/3 (+) | 3/3 (+) | 3/3 (+) | 3/3 (+) | 3/3 (+) | 3/3 (+) |
| A/goose/Guiyang/337/2006 | H5N1 | 4 | 10 ^{4.1} | 3/3 (+) | 3/3 (+) | 3/3 (+) | 3/3 (+) | 3/3 (+) | 3/3 (+) |
| A/chicken/Vietnam/NCVD-1088/2013 | H5N1 | 7.2 | 10 ^{3.9} | 3/3 (+) | 3/3 (+) | 3/3 (+) | 3/3 (+) | 3/3 (+) | 3/3 (+) |

Analytical Specificity –Cross-Reactivity

Cross-reactivity of the Influenza A/H5 Subtyping Kit was evaluated by testing influenza A viruses of different types and subtypes that include viruses representing diverse geographic locations and different sources. Samples were tested in triplicate using RNA extracted from high titer preparations

of viruses ($\geq 10^6$ EID₅₀/mL). Cross-reactivity testing of the Influenza A/H5 Subtyping Kit was evaluated with the Invitrogen Superscript™ enzyme system and one cleared extraction method. The results are summarized in the tables below.

Human Influenza Viruses for Cross-Reactivity Testing

| Virus Designation | Subtype | EID ₅₀ /mL | Invitrogen SuperScript™ | | |
|--------------------------|---------------------|-----------------------|-------------------------|-----|-----|
| | | | InfA | H5a | H5b |
| A/Brisbane/59/07 | A(H1N1) | 10 ^{8.4} | (+) 3/3 | - | - |
| A/Fujian Gulou/1896/2009 | | 10 ^{9.1} | (+) 3/3 | - | - |
| A/Perth/16/2009 | A(H3N2) | 10 ^{8.9} | (+) 3/3 | - | - |
| A/Texas/50/2012 | | 10 ^{9.2} | (+) 3/3 | - | - |
| A/California/07/09 | A(H1N1)pdm09 | 10 ^{8.4} | (+) 3/3 | - | - |
| A/Washington/24/2012 | | 10 ^{8.5} | (+) 3/3 | - | - |
| B/Brisbane/60/2008 | B(Victoria lineage) | 10 ^{9.3} | - | - | - |
| B/Montana/5/2012 | | 10 ^{8.4} | - | - | - |
| B/Wisconsin/01/2010 | B(Yamagata lineage) | 10 ^{9.2} | - | - | - |
| B/Massachusetts/02/2012 | | 10 ^{9.2} | - | - | - |

Animal Influenza Viruses for Cross-Reactivity Testing

| Species | Virus | Type | EID ₅₀ /mL | Invitrogen SuperScript™ | | |
|---------|--------------------------------------|---------|-----------------------|-------------------------|-----|-----|
| | | | | InfA | H5a | H5b |
| Swine | A/swine/Wisconsin/125/1997 | SW-H1N1 | 10 ^{7.9} | (+) 3/3 | - | - |
| Swine | A/Maryland/12/1991 | SW-H1N1 | 10 ^{9.0} | (+) 3/3 | - | - |
| Swine | A/Iowa/1/2006 | SW-H1N1 | 10 ^{9.0} | (+) 3/3 | - | - |
| Avian | A/chicken/Pennsylvania/298101-4/2004 | H2N2 | 10 ^{9.5} | (+) 3/3 | - | - |
| Canine | A/canine/Florida/43/2004 | H3N8 | 10 ^{7.2} | (+) 3/3 | - | - |
| Equine | A/equine/Ohio/01/2003 | EQ-H3N8 | 10 ^{7.5} | (+) 3/3 | - | - |
| Avian | A/chicken/Alabama/1975 | H4N8 | 10 ^{10.0} | (+) 3/3 | - | - |

| | | | | | | |
|-------|-----------------------------------|-------|-------------|------------|------------|------------|
| Avian | A/duck/Singapore-Q/F119-3/1997 | H5N3 | $10^{8.5}$ | (+) 3/3 | (+) 3/3 | (+) 3/3 |
| Avian | A/duck/Pennsylvania/1969 | H6N1 | $10^{9.2}$ | (+) 3/3 | - | - |
| Avian | A/chicken/California/32213-1/2000 | H6N2 | $10^{9.1}$ | (+) 3/3 | - | - |
| Avian | A/chicken/New York/13237-6/1998 | H6N8 | $10^{10.0}$ | (+) 3/3 | - | - |
| Avian | A/chicken/New Jersey/15906-9/1996 | H11N1 | $10^{6.5}$ | (+) 3/3 | - | - |
| Avian | A/duck/Memphis/546/1974 | H11N9 | $10^{9.8}$ | (+) 3/3 | - | - |
| Avian | A/Taiwan/2/2013 | H6N1 | $10^{10.2}$ | (+) 3/3 | - | - |
| Avian | A/Anhui/1/2013 | H7N9 | $10^{10.9}$ | (+) 3/3 | - | - |
| Avian | A/mallard/Netherlands/12/2000 | H7N3 | $10^{9.5}$ | (+) 3/3 | - | - |
| Avian | A/chicken/Arkansas/10/2008 | H7N3 | $10^{9.9}$ | (+) 3/3 | - | - |
| Avian | A/Bangladesh/0994/2011 | H9N2 | $10^{10.5}$ | (+) 3/3 | - | - |
| Swine | A/Indiana/21/2012 | H3N2v | $10^{9.4}$ | (+) 3/3 | - | - |
| Swine | A/Minnesota/11/2010 | H3N2v | $10^{9.2}$ | (+) 3/3 | - | - |

Exclusivity Testing

An exclusivity study was performed to demonstrate the specificity of each primer and probe set of the Influenza A/H5 Subtyping Kit when tested with common non-influenza human respiratory viruses, respiratory bacteria, and commensal organisms of the human respiratory tract. All organisms used in the study were propagated, titered, and characterized to confirm identity prior to testing. Nucleic acids were purified from thirty-five (35) non-influenza organisms (16 viruses, 18 bacteria, and 1 yeast) representing common respiratory pathogens or flora commonly present in specimens collected from the human nasopharynx region. High titer preparations of bacteria and yeast, generally greater than or equal to 10^6 cfu/mL, and non-influenza respiratory virus preparations at concentrations greater than 10^6 TCID₅₀/mL were tested (except in cases where production of high titer virus stock was not possible, e.g. parainfluenza virus type 2). The Influenza A/H5 Subtyping Kit was evaluated with the Invitrogen Superscript™ enzyme system and one cleared extraction method. The results are summarized in the table below.

Exclusivity results with Respiratory Pathogens and Flora

| Organism Tested | | | Invitrogen SuperScript™ | | |
|--|-------------|------------------------|-------------------------|-----|-----|
| Bacteria and Yeast | Strain | cfu/mL | InfA | H5a | H5b |
| <i>Bordetella pertussis</i> | A639 | 10 ^{8.3} | - | - | - |
| <i>Candida albicans</i> (yeast) | 2001-21-196 | 10 ^{8.8} | - | - | - |
| <i>Chlamydia pneumoniae</i> | TW183 | 40 IFU/mL ¹ | - | - | - |
| <i>Corynebacterium diphtheriae</i> | - | 10 ¹⁰ | - | - | - |
| <i>Escherichia coli</i> | K12 | 10 ^{9.6} | - | - | - |
| <i>Streptococcus pyogenes</i> | 7790-06 | 10 ^{7.5} | - | - | - |
| <i>Haemophilus influenzae</i> | M15709 | 10 ^{6.4} | - | - | - |
| <i>Lactobacillus plantarum</i> | - | 10 ^{8.8} | - | - | - |
| <i>Legionella pneumophila</i> | - | 10 ^{7.1} | - | - | - |
| <i>Moraxella catarrhalis</i> | M15757 | 10 ^{9.5} | - | - | - |
| <i>Mycobacterium tuberculosis</i> | H37Rv | 100 ng/uL ² | - | - | - |
| <i>Mycoplasma pneumoniae</i> | M129 | 10 ^{7.7} | - | - | - |
| <i>Neisseria elongata</i> | - | 10 ^{8.6} | - | - | - |
| <i>Neisseria meningitidis</i> | M2578 | 10 ^{7.9} | - | - | - |
| <i>Pseudomonas aeruginosa</i> | - | 10 ^{10.5} | - | - | - |
| <i>Staphylococcus aureus</i> | - | 10 ^{10.7} | - | - | - |
| <i>Staphylococcus epidermidis</i> | - | 10 ^{10.5} | - | - | - |
| <i>Streptococcus pneumoniae</i> | 249-06 | 10 ^{6.6} | - | - | - |
| <i>Streptococcus salivarius</i> | SS1672 | 10 ^{8.4} | - | - | - |
| Organism Tested | | | Invitrogen SuperScript™ | | |
| Viruses | Strain | TCID ₅₀ /mL | InfA | H5a | H5b |
| Human adenovirus type 1 | Ad.71 | 10 ^{9.2} | - | - | - |
| Human adenovirus type 7a | S-1058 | 10 ^{7.1} | - | - | - |
| Human parainfluenza virus type 1 | - | 3.0 ng/μL ² | - | - | - |
| Human parainfluenza virus type 2 | Greer | 10 ^{3.1} | - | - | - |
| Human parainfluenza virus type 3 | C-243 | 10 ^{7.9} | - | - | - |
| Respiratory syncytial virus | CH93-18b | 10 ^{6.8} | - | - | - |
| Human rhinovirus type A | 1A | 10 ^{5.8} | - | - | - |
| Enterovirus | Echo 6 | 10 ^{6.9} | - | - | - |
| Human coronavirus | 299E | 31.6ng/μL ² | - | - | - |
| Human coronavirus | OC43 | 50.4ng/μL ² | - | - | - |
| Herpes simplex virus | KOS | 5 X 10 ^{7.75} | - | - | - |
| Varicella-zoster virus | AV92-3:H | 5 X 10 ^{3.75} | - | - | - |
| Epstein-Barr virus | B95-8 | 1.7 ng/μL ² | - | - | - |
| Measles virus (<i>Paramyxoviridae</i>) | Edmonston | 5 X 10 ^{4.5} | - | - | - |
| Mumps virus | Enders | 5 X 10 ^{6.5} | - | - | - |
| Human Cytomegalovirus | AD-169 | 5 X 10 ^{6.25} | - | - | - |

¹Organism quantified by Infectious Forming Units (IFU)

²Organism nucleic acid quantified by spectrophotometry

IX. CLINICAL PERFORMANCE EVALUATION

The clinical performance of oligonucleotide primer and probe sets of the Influenza A/H5 Subtyping Kit were evaluated using contrived samples of grown virus added to an A549 cell suspension to simulate positive clinical samples. A total of fifty positive contrived samples in high, moderate, and low concentrations were evaluated. In addition, sixty-five specimens that tested negative for influenza A with the CDC Human Influenza rRT-PCR Diagnostic Panel that were obtained from a clinical study conducted during the 2011-2012 influenza season were evaluated. Testing was performed with both enzymes cleared for use with the kit. The results are summarized in the table below.

Clinical Performance Evaluation Results

| Enzyme Utilized | # of Positives ¹ | % Positive Agreement (95% CI) | # of Negatives ² | % Negative Agreement (95% CI) |
|-----------------------------|-----------------------------|-------------------------------|-----------------------------|-------------------------------|
| Quanta BioSciences qScript™ | 50/50 | 100 (92.9 – 100.00) | 65/65 | 100 (94.4 – 100.00) |
| Invitrogen SuperScript™ | 44/50 | 88.0 (76.2 - 94.4) | 65/65 | 100 (94.4 – 100.00) |

¹Proportion of contrived samples correctly identified as positive by both influenza A H5a and H5b primer and probe sets.

²Proportion of negative samples correctly identified versus the comparator.

X. CONCLUSION

The modification of the CDC Human Influenza Virus rRT-PCR Diagnostic Panel, Influenza A/H5 Subtyping kit to ensure comprehensive detection of circulating influenza A/H5 clades does not substantially change the device. Analytical and clinical data demonstrate that the performance of the device to detect Asian lineage influenza A H5 viruses is accomplished with high positive and negative percent agreement in a manner substantially equivalent to the predicate.

XI. REFERENCES

Lednicky JA, Villanueva JM, Burke SA, Shively R, Shaw MW, Daniels DE, Hamilton SB, Donis RO. 2010. Validation of a method for preparing influenza H5N1 simulated samples. J Virol Methods. 2010 Aug;167(2):125-31.