



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center - WO66-G609  
Silver Spring, MD 20993-0002

February 12, 2016

Acclarent, Inc.  
James Patrick Garvey  
Senior Manager, Regulatory Affairs  
1525-B O'Brien Drive  
Menlo Park, California 94025

Re: K153341

Trade/Device Name: Relieva Scout™ Multi-Sinus Dilation System  
Regulation Number: 21 CFR 874.4420  
Regulation Name: Ear, Nose, And Throat Manual Surgical Instrument  
Regulatory Class: Class I  
Product Code: LRC  
Dated: November 18, 2015  
Received: November 19, 2015

Dear Mr. Garvey:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

  
Eric A. Mann -S

for Malvina Eydelman, M.D.  
Director  
Division of Ophthalmic and Ear, Nose,  
and Throat Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)

K153341

Device Name

Relieva Scout™ Multi-Sinus Dilation System

Indications for Use (Describe)

For patients aged 18 and older, the Relieva Scout™ Multi-Sinus Dilation System is intended to provide a means to access the sinus space and to dilate the sinus ostia and spaces associated with the sphenoid, frontal, and maxillary sinus cavities for diagnostic and therapeutic procedures. In addition, the device is intended to illuminate within and transilluminate across nasal and sinus structures.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

**PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON A SEPARATE PAGE IF NEEDED.**

### FOR FDA USE ONLY

Concurrence of Center for Devices and Radiological Health (CDRH) (Signature)

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*“An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number.”*

---

**APPENDIX A: 510(k) SUMMARY****[807.92(a)(1)] Submitter Information**

**Sponsor/Submitter:** Acclarent, Inc.  
1525-B O'Brien Drive  
Menlo Park, California 94025

**Contact Person:** Patrick Garvey  
Sr. Manager, Regulatory Affairs  
Email: pgarvey@its.jnj.com  
Tel: 650-687-5888

**Date Summary Prepared:** February 12, 2016

**[807.92(a)(2)] Name of Device**

**Device Trade Name:** Relieva Scout™ Multi-Sinus Dilation System

**Common Name:** Sinus Balloon Catheter

**Device Classification:** Class I

**Regulation Number:** 21 CFR 874.4420

**Classification Name:** Ear, Nose, and Throat Manual Surgical Instrument (21 CFR 874.4420)

**Product Code:** LRC

**[807.92(a)(3)] Legally Marketed Devices**

**Predicate Devices:** RELIEVA SCOUT® Sinus Dilation System (K120280)  
XprESS™ Multi-Sinus Dilation Tool (K121174)  
RELIEVA® Spin Balloon Sinuplasty System (K111875)

**[807.92(a)(4)] Device Description**

**Device Description:** The Relieva Scout™ Multi-Sinus Dilation System is a sterile, single-use system that is an integrated device with a handle, rail, balloon catheter, and sinus illumination system with an illuminated ball tip. The packaged device contains the sinus balloon catheter an angle selection tool to enable angular bends for multiple sinus access (sphenoid, frontal, maxillary). The device may be used to access the sinus space and the balloon inflated with sterile water or sterile or saline to dilate the sinus ostia and infundibulum. The inflation device is provided separately. The device may be used in Operating Room and physician office settings.

**[807.92(a)(5)] Intended Use**

**Indications for Use:** For patients aged 18 and older, the Relieva Scout™ Multi-Sinus Dilation System is intended to provide a means to access the sinus space and to dilate

## Relieva Scout™ Multi-Sinus Dilation System

the sinus ostia and spaces associated with the sphenoid, frontal, and maxillary sinus cavities for diagnostic and therapeutic procedures. In addition, the device is intended to illuminate within and transilluminate across nasal and sinus structures.

Difference in Indications from Predicate Device

The difference in indications for use between the subject and predicate devices is supported is presented in Table 1.

**[807.92(a)(6)] Technical Characteristics****Technological Characteristics:**

The Relieva Scout™ Multi-Sinus Dilation System combines features of a rail-based balloon catheter with the tissue expansion effect of balloon dilation. The distal end of the device may be shaped with the angle selection tool to optimize sinus access. Light from an extendable integrated illumination system can be seen via transillumination.

See Table 1 for a comparison of the technological characteristics between the Relieva Scout and the predicate devices.



**Table 1: Comparison of Technological Characteristics between the Relieva Scout™ Multi-Sinus Dilation System and predicate devices.**

| <b>Attribute</b>       | <b><u>Primary Predicate Device:</u><br/>RELIEVA SCOUT® Sinus Dilation System (K120280)</b>  | <b><u>Secondary Predicate Device:</u><br/>XprESS™ Multi-Sinus Dilation Tool (K121174)</b>   | <b><u>Predicate Device:</u><br/>RELIEVA® Spin Sinus Dilation System (K111875)</b>   | <b><u>Subject Device:</u><br/>Relieva Scout™ Multi-Sinus (SMS) Dilation System</b>  |
|------------------------|---|---|---|---|
| Manufacturer           | Acclarent, Inc.   | Entellus Medical, Inc.  | Acclarent, Inc.   | Acclarent, Inc.   |
| Common Name            | Sinus Balloon Dilation System   | Sinus Balloon Dilation System   | Sinus Balloon Dilation System   | Sinus Balloon Dilation System   |
| Class                  | I   | I   | I   | I   |
| Product Code is LRC    | LRC   | LRC   | LRC   | LRC   |
| Classification Section | 21 CFR 874.4420   | 21 CFR 874.4420   | 21 CFR 874.4420   | 21 CFR 874.4420   |
| Indications for Use    | For patients aged 18 and older, the Relieva Scout™ Sinus Dilation System is intended to provide a means to access the frontal sinus space and to dilate the frontal recess, frontal sinus ostia and spaces within the frontal sinus cavity for diagnostic and therapeutic procedures. In addition, the device is intended to illuminate within and transilluminate across nasal and sinus structures. | To access and treat the frontal recesses, sphenoid sinus ostia and maxillary ostia/ethmoid infundibula in adults using a trans-nasal approach. The bony sinus outflow tracts are remodeled by balloon displacement of adjacent bone and paranasal sinus structures. | Intended to provide a means to access the sinus space and to dilate the sinus ostia and spaces associated with the paranasal sinus cavities for diagnostic and therapeutic procedures. For children aged 17 and under, the device is intended to dilate sinus ostia and spaces associated with the maxillary sinus for diagnostic and therapeutic procedures. | For patients aged 18 and older, the Relieva Scout™ Multi-Sinus Dilation System is intended to provide a means to access the sinus space and to dilate the sinus ostia and spaces associated with the sphenoid, frontal, and maxillary sinus cavities for diagnostic and therapeutic procedures. In addition, the device is intended to illuminate within and transilluminate across nasal and sinus structures. |
| Indicated for Children | No  | No  | Yes   | No  |



Traditional 510(k) Premarket Notification

Relieva Scout™ Multi-Sinus Dilation System

| <b>Attribute</b>   | <b><u>Primary Predicate Device:</u><br/>RELIEVA SCOUT® Sinus Dilation System (K120280)</b>   | <b><u>Secondary Predicate Device:</u><br/>XprESS™ Multi-Sinus Dilation Tool (K121174)</b>  | <b><u>Predicate Device:</u><br/>RELIEVA® Spin Sinus Dilation System (K111875)</b>   | <b><u>Subject Device:</u><br/>Relieva Scout™ Multi-Sinus (SMS) Dilation System</b>   |
|--|--|--|---|--|
| Single Patient Use   | Yes  | Yes  | Yes   | Yes  |
| Direct Patient Contact   | Yes  | Yes  | Yes   | Yes  |
| Labeled as Non-Pyrogenic?  | No   | No   | No  | No   |
| Technological Characteristics  | Combines features of a frontal ostium seeker with the tissue expansion effect of balloon dilation. The distal end of the device is permanently curved to optimize frontal ostium access. Light from the distal tip of the integrated sinus illumination system can be seen via transillumination. The device is connected to any standard light source via a light cable and an adapter. | Combines features of a curved suction tip and a frontal ostium seeker (access) with the tissue expansion effect of balloon dilation (treat). The distal end of the device is re-shapeable. | Combines a sinus balloon catheter and a sinus guide catheter to access the sinus space and dilate sinus ostia.<br><br>A Sinus Illumination System comes packaged with the device and is pre-loaded into the Spin Sinus Dilation System. | Combines a sinus balloon catheter with rail-based balloon guidance to access the sinus space and dilate the sinus ostia. The packaged device contains an angle selection tool to enable angular bends for multiple sinus access.<br><br>Light from the distal tip of the integrated sinus illumination system can be seen via transillumination. |
| Constructed of Materials Commonly Used in Patient Contacting Medical Devices | Yes  | Yes  | Yes   | Yes  |
| Balloon Diameter   | 6mm  | 5mm to 7mm   | 6mm   | 6mm  |
| Balloon Length   | 24mm   | 8-20mm   | 16mm  | 24mm   |



Traditional 510(k) Premarket Notification

Relieva Scout™ Multi-Sinus Dilation System

| <b>Attribute</b>   | <b><u>Primary Predicate Device:</u><br/>RELIEVA SCOUT® Sinus Dilation System (K120280)</b> | <b><u>Secondary Predicate Device:</u><br/>XprESS™ Multi-Sinus Dilation Tool (K121174)</b> | <b><u>Predicate Device:</u><br/>RELIEVA® Spin Sinus Dilation System (K111875)</b> | <b><u>Subject Device:</u><br/>Relieva Scout™ Multi-Sinus (SMS) Dilation System</b> |
|--|--|---|---|--|
| Maximum Inflation Pressure   | 12 ATM   | 12 ATM  | 12 ATM  | 12 ATM   |
| Single-Handed Use  | Yes  | Yes   | Yes   | Yes  |
| Balloon Slide Mechanism  | Yes  | Yes   | Yes   | Yes  |
| Rail-Based Design  | Yes, preset rail intended for frontal sinus access only                                    | Yes, malleable rail   | No  | Yes, malleable rail  |
| Handle and Slider System to Facilitate Advancement of a Balloon Catheter on the Distal End | Yes  | Yes   | Yes   | Yes  |
| Uses a Bending Tool to Achieve Target Angles for Access into Targeted Anatomy              | No   | Yes   | No  | Yes  |
| Guidewire Capability (Sinus Illumination System)   | Yes  | Yes   | Yes   | Yes  |
| Sinuses in which device is intended for use: [Sphenoid, Maxillary, Frontal]                | Frontal  | Sphenoid, Maxillary, Frontal  | Sphenoid, Maxillary, Frontal  | Sphenoid, Maxillary, Frontal   |
| EtO Sterilized   | Yes  | Yes   | Yes   | Yes  |



Traditional 510(k) Premarket Notification

Relieva Scout™ Multi-Sinus Dilation System

| <b>Attribute</b>        | <b><u>Primary Predicate Device:</u><br/>RELIEVA SCOUT® Sinus Dilation System (K120280)</b>                  | <b><u>Secondary Predicate Device:</u><br/>XprESS™ Multi-Sinus Dilation Tool (K121174)</b>                   | <b><u>Predicate Device:</u><br/>RELIEVA® Spin Sinus Dilation System (K111875)</b>                           | <b><u>Subject Device:</u><br/>Relieva Scout™ Multi-Sinus (SMS) Dilation System</b>                          |
|-------------------------|---|---|---|---|
| Packaging               | Thermoformed tray in pouch  | Backer card in pouch  | Backer card in pouch  | Thermoformed tray in pouch  |
| Principles of Operation | Manually operated device. Balloon inflated with sterile saline or water to mechanically dilate sinus ostia. | Manually operated device. Balloon inflated with sterile saline or water to mechanically dilate sinus ostia. | Manually operated device. Balloon inflated with sterile saline or water to mechanically dilate sinus ostia. | Manually operated device. Balloon inflated with sterile saline or water to mechanically dilate sinus ostia. |

**[807.92(b) (1)] Determination of Substantial Equivalence****Non-Clinical Performance Data:**

Bench testing met all acceptance criteria for attributes such as dimensional attributes, cycle fatigue, balloon burst, and bond separation. Testing in accordance with ISO 10993-1 AAMI ANSI ISO 10993-1:2009/(R) 2013 also showed that the Relieva Scout™ Multi-Sinus Dilation System is biocompatible.

The sterilization process was validated per AAMI/ANSI/ISO 11135-1:2007 and demonstrated a sterility assurance level of  $10^{-6}$ . The method used for sterilization validation was the overkill (half-cycle approach) in a fixed chamber. Ethylene oxide residuals were tested and met ISO 10993-7:2008 requirements. The subject device is not tested nor labeled as “non-pyrogenic”.

Packaging shelf life was established per ASTM F1980-07.

The performance data demonstrated that the device performs as intended.

**[807.92(b) (2)] Determination of Substantial Equivalence****Clinical Performance Data**

Clinical data was not necessary for the Relieva Scout™ Multi-Sinus Dilation System. The performance data demonstrated that the device performs as intended.

**[807.92(b) (3)] Conclusion****Conclusion from Non-Clinical and Clinical Tests**

The Relieva Scout™ Multi-Sinus Dilation System is substantially equivalent to the predicate devices.