



Food and Drug Administration
10903 New Hampshire Avenue
Document Control Center - WO66-G609
Silver Spring, MD 20993-0002

June 10, 2016

Tornier S.A.S.
Aymen Azaiez
Regulatory Affairs Specialist
161, Rue Lavoisier
38330 Montbonnot Saint Martin
FRANCE

Re: K160555

Trade/Device Name: BLUEPRINT Patient Specific Instrumentation
Regulation Number: 21 CFR 888.3660
Regulation Name: Shoulder Joint Metal/Polymer Semi-Constrained Cemented Prosthesis
Regulatory Class: Class II
Product Code: KWS
Dated: February 22, 2016
Received: February 29, 2016

Dear Aymen Azaiez:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in

the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to

<http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

Mark N. Melkerson -S

Mark N. Melkerson
Director
Division of Orthopedic Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)
K160555

Device Name
BLUEPRINT Patient Specific Instrumentation

Indications for Use (Describe)

The hardware

The Aequalis Glenoid Guides are patient-specific drill guides. They have been specially designed to assist in the intraoperative positioning of glenoid components used with total anatomic shoulder arthroplasty procedures using anatomic landmarks that are identifiable on patient-specific preoperative CT-scans.

Aequalis PerFORM Anatomic Glenoid Guide is used by surgeons to facilitate the placement of the Aequalis PerFORM glenoids.

The software

The BLUEPRINT 3D planning software is a medical device for surgeon composed of one software component. It is intended to be used as a pre-surgical planner for shoulder orthopedic surgery.

BLUEPRINT 3D planning software runs on standard personal and business computers running Microsoft Windows or Mac OS operating systems.

The software supports DICOM standard to import the CT-Scan (Computed Tomography) images of the patient. Only CT-Scan modality can be loaded with BLUEPRINT3D planning software.

BLUEPRINT 3D planning software allows surgeon to visualize, measure, reconstruct, and annotate anatomic data. It allows surgeon to design patient specific guides based on the presurgical plan.

This device is intended for use provided anatomic reference points necessary for positioning of the guide are present on the CT scan.

The software leads to the generation of a surgery report along with a 3D file of the patient-specific guide.

BLUEPRINT 3D planning software does not include any system to manufacture the guide.

BLUEPRINT 3D planning software is to be used for adult patients only and should not be used for Diagnostic purpose.

Type of Use (Select one or both, as applicable)

☒ Prescription Use (Part 21 CFR 801 Subpart D)

☐ Over-The-Counter Use (21 CFR 801 Subpart C)

PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON A SEPARATE PAGE IF NEEDED.

FOR FDA USE ONLY

Concurrence of Center for Devices and Radiological Health (CDRH) (Signature)

TORNIER

Implants Chirurgicaux

BLUEPRINT™ Patient Specific
Instrumentation K160555

Summary of Safety and Effectiveness information *Special 510(k) Premarket – BLUEPRINT™ Patient Specific Instrumentation*

1) Device name

Trade name: BLUEPRINT™ Patient Specific Instrumentation
510(k) Number: K160555
Common name: Patient specific instrumentation + 3D planning software
Classification name: Prosthesis, Shoulder, Semi-constrained, Metal/Polymer Cemented
 (§888.3660)

2) Submitter :

TORNIER SAS
 161 rue Lavoisier
 38330 Montbonnot Saint Martin- France
 Registration Number: 3000931034

3) Company contact :

Tornier SAS
 Mr Aymen AZAIEZ
 Regulatory Affairs Specialist
 161 rue Lavoisier
 38334 Montbonnot
 Tel: 00 33 4 76 61 35 00
 Fax: 00 33 4 76 61 35 59
 e-mail : aymen.azaiez@tornier.com

4) Classification

Device class: Class II
Classification panel: Orthopedic
Product code: KWS

5) Equivalent / Predicate device :

BLUEPRINT™ Patient Specific Instrumentation, Tornier SAS (K143374)
 Delta TT Acetabular System, Limacorporate S.p.A. (K112898)

TORNIER

Implants Chirurgicaux

*BLUEPRINT™ Patient Specific
Instrumentation K160555*

6) Device description :

BLUEPRINT™ Patient Specific Instrumentation is composed of two components: *Aequalis Glenoid Guides (hardware) and BLUEPRINT 3D planning software (software).*

BLUEPRINT™ Patient Specific Instrumentation which includes the *Aequalis Glenoid Guides* and *BLUEPRINT 3D planning software* is the responsibility of Tornier. Tornier is the legal manufacturer for the hardware and the software.

The hardware

The *Aequalis Glenoid Guides* are patient-specific instruments specially designed to facilitate the implantation of the *Aequalis PerFORM* shoulder prostheses and are exclusively reserved for this use.

The *Aequalis Glenoid Guides* are designed and manufactured based on a pre-operative plan generated only by the software *BLUEPRINT™ 3D planning software*.

The software

BluePrint 3D Planning software is composed of one software component connected to an Online Management System (OMS). The software installed on a computer is intended to be used by orthopedic surgeons, as a preoperative planning software for shoulder arthroplasty surgery (= total anatomic shoulder replacement).

It is intended to help to plan an operation by allowing surgeons to:

- position and select the glenoid implant,
- design a patient specific pin guide.

This submission seeks clearance for:

- **Hardware: a guide made of Titanium with an orientation hole which allows for controlling rotation of the commercially available implant, Aequalis™ PerFORM**
- **Software modified to:**
 - Integrate a guide made of titanium,
 - Add a glenoid sphere radius measurement.

These modifications do not affect the intended use or the indications for use of the device or alter the fundamental scientific technology of the device.

Intended Use

The hardware

The Aequalis Glenoid Guides are intended to be used as surgical instruments to assist in the intraoperative positioning of glenoid components used with total anatomic shoulder

TORNIER

Implants Chirurgicaux

BLUEPRINT™ Patient Specific
Instrumentation K160555

arthroplasty procedures using anatomic landmarks that are identifiable on patient-specific preoperative CT-scans.

The software

The *BLUEPRINT 3D planning software* is intended to be used as a medical software to assist in pre-operative surgical planning for shoulder surgery.

7) Materials :

The commercially available *Aequalis Glenoid Guides* are manufactured from medical grade polyamide 2200. This 510k seeks clearance for guides made of medical grade titanium (Ti6Al4V).

8) Indications :

The hardware

The *Aequalis Glenoid Guides* are patient-specific drill guides. They have been specially designed to assist in the intraoperative positioning of glenoid components used with total anatomic shoulder arthroplasty procedures using anatomic landmarks that are identifiable on patient-specific preoperative CT-scans.

Aequalis PerFORM Anatomic Glenoid Guide is used by surgeons to facilitate the placement of the *Aequalis PerFORM* glenoids.

The software

The *BLUEPRINT 3D planning software* is a medical device for surgeon composed of one software component. It is intended to be used as a pre-surgical planner for shoulder orthopedic surgery.

BLUEPRINT 3D planning software runs on standard personal and business computers running Microsoft Windows or Mac OS operating systems.

The software supports DICOM standard to import the CT-Scan (Computed Tomography) images of the patient. Only CT-Scan modality can be loaded with *BLUEPRINT 3D planning software*.

BLUEPRINT 3D planning software allows surgeon to visualize, measure, reconstruct, and annotate anatomic data. It allows surgeon to design patient specific guides based on the presurgical plan.

This device is intended for use provided anatomic reference points necessary for positioning of the guide are present on the CT scan.

The software leads to the generation of a surgery report along with a 3D file of the patient-specific guide.

BLUEPRINT 3D planning software does not include any system to manufacture the guide.

BLUEPRINT 3D planning software is to be used for adult patients only and should not be used for Diagnostic purpose.

TORNIER

Implants Chirurgicaux

*BLUEPRINT™ Patient Specific
Instrumentation K160555*

9) Summary of technological characteristics

Table 1 Main features comparison

Main features or system characteristics	<i>BLUEPRINT™ Patient Specific Instrumentation</i> SUBJECT DEVICE SYSTEM	<i>BLUEPRINT™ Patient Specific Instrumentation</i> (K143374) PRIMARY	Delta TT Acetabular System (K112898) REFERENCE
Material	Ti6Al4V	Polyamide 2200	Ti6Al4V
Standard	ISO 5832-3	USP Class VI compatible	ISO 5832-3 ASTM F 1472
Manufacturing	3D printing	3D printing	3D printing
Product Code	KWS	KWS	LPH, MBL
Surgical procedure	Total anatomic shoulder arthroplasty	Total anatomic shoulder arthroplasty	Hip replacement
Single-use	Yes	Yes	Yes
Sterile	No	No	Yes
Manufacturer	Tornier SAS	Tornier SAS	Limacorporate S.p.A.
Software	Assist surgeon in pre-operative surgical planning	Assist surgeon in pre-operative surgical planning	NA

Device comparison showed that the proposed device is substantially equivalent:

- in intended use and performance characteristics to the predicate device
- in material to the reference device

10) Non-clinical testing

BLUEPRINT™ Patient Specific Instrumentation was validated through studies using cadaver specimens. Non clinical testing was performed on *BLUEPRINT™ Patient Specific Instrumentation* to assess that no new safety or effectiveness questions were raised with this device.

Table 2: Non-clinical testing

Validation and / or Verification Method	Acceptance Criteria description	Verification and Validation Results
Patient Specific Guiding Wire test	Version angle error, inclination angle error, main (central) entry point error and rotation entry point error must be compliant with device specifications	Acceptable
Dimensional test of titanium guide	The titanium guide meets dimensional specifications	Acceptable

TORNIER

Implants Chirurgicaux

*BLUEPRINT™ Patient Specific
Instrumentation K160555*

This testing aimed to validate that the software measures in preoperative planning (implant positioning with a human cadaver) with the modified Blueprint Software generated a patient specific glenoid guide matching the patient anatomy.

11) Substantial equivalence conclusion

Based upon this comparative study, substantial equivalence of *BLUEPRINT™ Patient Specific Instrumentation* to the predicates can be demonstrated on the following grounds, according to the FDA's Guidelines for Substantial Equivalence Decision making Process:

- *BLUEPRINT™ Patient Specific Instrumentation (pending device)* is compared to the predicate device and the reference device.
- *BLUEPRINT™ Patient Specific Instrumentation (pending device)* has the same intended use as the predicate device: *BLUEPRINT™ Patient Specific Instrumentation (k143374)*,
- *BLUEPRINT 3D planning software (pending device)* is equivalent to the *BLUEPRINT 3D planning software (k143374)*.
- *BLUEPRINT 3D planning software (pending device)* user manual is similar in indications precautions, warnings, and instructions as the predicate device: *(k143374)*.
- Major technological characteristics are equivalent between *BLUEPRINT™ Patient Specific Instrumentation (pending device)* and the predicate device:
 - Equivalence of general features
 - Equivalent surgical procedures
 - Equivalent intended use, indications for use
- *BLUEPRINT™ Patient Specific Instrumentation (pending device)* is considered equivalent to the reference device *Delta TT Acetabular System (K112898)* in material.
- *BLUEPRINT™ Patient Specific Instrumentation (pending device)* differs from the predicate device *BLUEPRINT™ Patient Specific Instrumentation (K143374)* regarding two features: the orientation hole intended to help the surgeon to control the rotation of the commercially available implant, Aequalis™ PerFORM, and glenoid sphere radius measurement to facilitate choosing the implant Aequalis™ PerFORM. These differences do not, however, change the fundamental technology, principle of operation for the drill guide/software or intended use.

Therefore, in the light of the above information, the *BLUEPRINT™ Patient Specific Instrumentation* is considered equivalent to the predicate devices