



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center – WO66-G609  
Silver Spring, MD 20993-0002

Sumitomo Heavy Industries, Ltd.  
% Cynthia Nolte, Ph.D., RAC  
Senior Director, Regulatory Affairs  
ICON Clinical Research LCC  
62 Forest Street, Suite 300  
MARLBOROUGH MA 01752

August 17, 2016

Re: K160612

Trade/Device Name: Sumitomo Proton Therapy System  
Regulation Number: 21 CFR 892.5050  
Regulation Name: Medical charged-particle radiation therapy system  
Regulatory Class: II  
Product Code: LHN  
Dated: August 3, 2016  
Received: August 4, 2016

Dear Dr. Nolte:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Industry and Consumer Education at its toll-free number (800) 638 2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, “Misbranding by reference to premarket notification” (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH’s Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

A handwritten signature in blue ink that reads "Michael D. O'Hara". The signature is written over a large, light gray watermark of the FDA logo.

For

Robert Ochs, Ph.D.  
Director  
Division of Radiological Health  
Office of In Vitro Diagnostics  
and Radiological Health  
Center for Devices and Radiological Health

Enclosure

## 4. INDICATIONS FOR USE

DEPARTMENT OF HEALTH AND HUMAN SERVICES Food and Drug Administration <b>Indications for Use</b>	Form Approved: OMB No. 0910-0120 Expiration Date: January 31, 2017 <i>See PRA Statement below.</i>
---	--

510(k) Number (if known)

K160612

Device Name

Sumitomo Proton Therapy System

Indications for Use (Describe)

The Sumitomo Proton Therapy System is a medical device designed to produce and deliver a proton beam for treatment of patients with localized tumors and other conditions susceptible to treatment by radiation.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

**PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON A SEPARATE PAGE IF NEEDED.**

### FOR FDA USE ONLY

Concurrence of Center for Devices and Radiological Health (CDRH) (Signature)

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASTaff@fda.hhs.gov](mailto:PRASTaff@fda.hhs.gov)

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

**510(K) SUMMARY**  
**Sumitomo Heavy Industries, Ltd.**  
**Sumitomo Proton Therapy System**  
**(per 21CFR 807.92)**

**1. SUBMITTER/510(K) HOLDER**

Sumitomo Heavy Industries, Ltd.

5-2, Sobiraki-cho

Niihama, Ehime, 792-8588

Japan

Establishment Registration No: 3002674566

Contact: Katsuhiko Narasaki, Engineering Dept.

Contact Phone: +81-897-32-6370

Contact Email: [kth\\_narasaki@shi.co.jp](mailto:kth_narasaki@shi.co.jp)

Date Prepared: June 27, 2016

**2. DEVICE NAME**

Proprietary Name: Sumitomo Proton Therapy System

Common Name: Proton therapy system

Classification Name: Medical charged-particle radiation therapy system

Classification Regulation: 21 CFR 892.5050

Product Code: LHN

**3. PREDICATE DEVICE**

- Sumitomo Proton Therapy System, Sumitomo Heavy Industries, Ltd. (K130426)
- IBA Proton Therapy System- Proteus 235, Ion Beam Applications, S.A. (K082416)

**4. DEVICE DESCRIPTION**

The Sumitomo Proton Therapy System (PTS) is a large-scale medical electrical system that consists of an integrated system of medical electrical equipment and non-medical electrical components to provide proton beam radiation therapy. The Sumitomo PTS consists of a 230 MeV Cyclotron, an Energy Selection System, a Beam Transport System, and a Gantry Treatment System. The new Sumitomo PTS is a modification of the Sumitomo PTS cleared for marketing under K130426 in November 2013. The purpose of the current 510(k) premarket notification is to add the pencil beam scanning functionality to the cleared Sumitomo PTS. The pencil beam scanning function allows for precise treatment of complex cancers, such as tumors located near critical structures and tumors with special shapes, including multi-site targets, without the need for beam shaping devices such as collimators or compensators. The pencil beam scanning function also produces a radiation field directly from the beam scan path layer by layer, which shortens

the treatment duration.

## **5. INDICATION FOR USE/INTENDED USE**

The Sumitomo Proton Therapy System is a medical device designed to produce and deliver a proton beam for treatment of patients with localized tumors and other conditions susceptible to treatment by radiation.

## **6. SUMMARY OF TECHNOLOGICAL CHARACTERISTICS COMPARED TO THE PREDICATE DEVICE/S**

The Sumitomo PTS with pencil beam scanning and predicate Sumitomo PTS are identical except for the pencil beam scanning functionality. The major features of both the new device and the predicate Sumitomo PTS include:

- Cyclotron (230 MeV)
- Energy selection system (ESS)
- Beam transport system (BTS)
- Multipurpose Nozzle
- Respiratory Gating System
- Multi Leaf Collimator
- Online PET (Gamma ray detector)
- Compact Gantry
- Patient Positioning System

The operational characteristics of the new and predicate Sumitomo PTS are identical. The accelerator energies are in the same range (70 to 250 MeV). The beam transport systems are identical and deliver the accelerated protons to the gantries to be used for treatment in the same manner. The gantries are identical, and direct the radiation through treatment heads to patients in the same manner. The interface between the treatment planning software and the treatment heads, and the treatment planning software itself, are also identical.

The multi-purpose nozzle described in the current 510(k) premarket notification provides both wobbling and pencil beam scanning functions. The nozzle described in the original 510(k) for the predicate device was limited to the wobbling scanning function only. The user selects whether to use wobbling scanning or pencil beam scanning when planning for treatment.

In pencil beam scanning mode, proton beams delivered by the ESS, BTS, and GTS (Gantry Transport System) are irradiated to the isocenter without any scatter, focused by quadrupole magnets installed at the entrance of the nozzle. The beam is scanned laterally in two orthogonal directions continuously, controlled by a pair of scanning magnets.

## **7. SUMMARY OF NON-CLINICAL PERFORMANCE TESTING AS BASIS FOR SUBSTANTIAL EQUIVALENCE**

A series of factory and on-site safety and performance studies were completed to assess the performance of the new Sumitomo PTS with pencil beam scanning. Testing was performed and documented to ensure the design outputs of the system and its components have met design specifications in a complete and verifiable manner. Testing was performed at the unit level, subsystem level at the factory and on-site, and system level on-site. The extensive factory testing and on-site testing demonstrates that the Sumitomo PTS with pencil beam scanning met all performance requirement specifications for hardware, software, and safety requirements for a proton therapy system. The Sumitomo PTS with pencil beam scanning met the requirements of all design specifications.

Sumitomo conducted a usability assessment of the Sumitomo PTS with pencil beam scanning. Study participants completed two simulated clinical treatment scenarios as well as tasks identified as essential and critical. The results confirmed that users can safely and effectively perform the activities involved in pencil beam scanning.

## **8. SUMMARY OF CLINICAL TESTING AS BASIS FOR SUBSTANTIAL EQUIVALENCE**

No clinical testing was conducted in support of this 510(k) premarket notification.

## **9. SUMMARY OF OTHER INFORMATION**

No other information is submitted in support of this 510(k) premarket notification.

## **10. CONCLUSIONS DRAWN FROM NON-CLINICAL AND CLINICAL TESTS**

Based on the information provided in this 510(k), Sumitomo believes that the new Sumitomo PTS with pencil beam scanning is substantially equivalent to the previously cleared Sumitomo PTS. The Sumitomo PTS with pencil beam scanning has the same intended use, similar design, principle of operation, and technological characteristics as the predicate device. The only difference between the Sumitomo PTS with pencil beam scanning and the predicate Sumitomo PTS is the addition of the pencil beam scanning function to the multi-purpose nozzle. The addition of the pencil beam scanning function raises no new issues of safety or effectiveness. The Sumitomo PTS with pencil beam scanning met all performance requirements with no usability issues. Thus, the Sumitomo PTS with pencil beam scanning is substantially equivalent to the predicate Sumitomo PTS.

The new device raises no new issues of safety and effectiveness. The non-clinical safety and performance testing performed demonstrates that the new device met all test specifications and is suitable for its intended use. The usability evaluation confirmed that users can safely and effectively use the new Sumitomo PTS for pencil beam scanning.

**Table 5-1. Side-by-Side Comparison of New and Predicate Sumitomo Proton Therapy Systems with Reference Device**

<b>Characteristic</b>	<b>Sumitomo Proton Therapy System with Pencil Beam Scanning NEW</b>	<b>Sumitomo Proton Therapy System (K130426) PREDICATE</b>	<b>Proteus 235 Proton Therapy System (K082416) REFERENCE</b>
Manufacturer	Sumitomo Heavy Industries, Ltd.	Sumitomo Heavy Industries, Ltd.	IBA
Indication for Use	The Sumitomo Proton Therapy System is a medical device designed to produce and deliver a proton beam for treatment of patients with localized tumors and other conditions susceptible to treatment by radiation.	The Sumitomo Proton Therapy System is a medical device designed to produce and deliver a proton beam for treatment of patients with localized tumors and other conditions susceptible to treatment by radiation.	The PTS is a medical device designed to produce and deliver a proton beam for the treatment of patients with localized tumors and other conditions susceptible to treatment by radiation.
Accelerator Type	Cyclotron	Cyclotron	Cyclotron
Particle	Proton	Proton	Proton
Accelerator: Energy (MeV)	70 – 230 continuous	70 – 230 continuous	70 – 230 continuous
Depth of penetration	3 cm-32cm	3cm - 32cm	5cm - 28cm
Field Shaping System	Wobbling, Pencil Beam Scanning	Wobbling	Double and Single Scattering, Uniform Scanning, Pencil Beam Scanning
Dose Rate: (Gy/min)	For Wobbling > 2 Gy/min (110 to 230 MeV) > 0.5 Gy/min for 10 cm x 10 cm x 4 cm (SOBP) at 70 MeV For pencil beam scanning functionality > 2Gy/min/liter (>150 MeV) > 2Gy/min/(10cmx10cmx4cm) (110MeV) > 0.5 Gy/min/(10cmx10cmx4cm) (Max depth 5cm with snout degrader)	For Wobbling > 2 Gy/min (110 to 230 MeV) > 0.5 Gy/min for 10 cm x 10 cm x 4 cm (SOBP) at 70 MeV	> 2Gy/min.
Safety System	Hard-wired relay and programmable logic controller (PLC) based interlock systems	Hard-wired relay and programmable logic controller (PLC) based interlock systems	Hard-wired relay-based interlock system
Treatment Geometry	Isocentric gantry or fixed horizontal beam delivery	Isocentric gantry or fixed horizontal beam delivery	Isocentric gantry or fixed horizontal beam delivery
Proton Beam current stability	Stable (Cyclotron)	Stable (Cyclotron)	Stable (Cyclotron)
Proton beam time structure	CW beam	CW beam	CW (Continues Wave) beam
Respiratory gating for wobbling mode	Simple (Respiratory gating+ Wobbling timing ) (Passive ridge filter)	Simple (Respiratory gating+ Wobbling timing ) (Passive ridge filter)	Not simple (Range modulator + Respiratory gating + Wobbling timing)
Pencil beam scanning	Line scanning + spot scanning	N/A	Line scanning + spot scanning
Gantry size	Short axis gantry (Similar to Loma Linda gantry)	Short axis gantry (Similar to Loma Linda gantry)	Long axis gantry
Delivery proton dose distribution measurement	Off line phantom measurement + On line PET measurement	Off line phantom measurement + On line PET measurement	Off line phantom measurement
Accessibility to the patient in gantry in case of emergency	Quick access (Caterpillar floor)	Quick access (Caterpillar floor)	Need time to access patient (in case nozzle is lower position)

N/A: Not available