



June 29, 2016

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center – WO66-G609  
Silver Spring, MD 20993-0002

CEPHEID  
JIM KELLY, Ph.D.  
EXECUTIVE DIRECTOR, REGULATORY AFFAIRS  
904 CARIBBEAN DRIVE  
SUNNYVALE CA 94089-1189

Re: K160901  
Trade/Device Name: Xpert<sup>®</sup> Carba-R Assay  
Regulation Number: 21 CFR 866.1640  
Regulation Name: Antimicrobial susceptibility test powder  
Regulatory Class: II  
Product Code: POC, OOI  
Dated: March 31, 2016  
Received: April 1, 2016

Dear Dr. Kelly:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulations (21 CFR Parts 801 and 809), please contact the Division of Industry and Consumer Education at its toll-free number (800) 638 2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, “Misbranding by reference to premarket notification” (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH’s Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

**Ribhi Shawar -S**

For Uwe Scherf, M.Sc., Ph.D.  
Director  
Division of Microbiology Devices  
Office of *In Vitro* Diagnostics  
and Radiological Health  
Center for Devices and Radiological Health

Enclosure

---

---

510(k) Number (if known)

K160901

---

Device Name

Xpert Carba-R

---

Indications for Use (Describe)

The Xpert<sup>®</sup> Carba-R Assay, performed on the GeneXpert<sup>®</sup> Instrument Systems, is a qualitative *in vitro* diagnostic test designed for the detection and differentiation of the *bla*KPC, *bla*NDM, *bla*VIM, *bla*OXA-48, and *bla*IMP gene sequences associated with carbapenem-non-susceptibility. The test utilizes automated real-time polymerase chain reaction (PCR).

The Xpert Carba-R Assay is intended as an aid to infection control in the detection of carbapenem-non-susceptible bacteria that colonize patients in healthcare settings. The Xpert Carba-R Assay is not intended to guide or monitor treatment for carbapenem-non-susceptible bacterial infections. A negative Xpert Carba-R Assay result does not preclude the presence of other resistance mechanisms.

The Xpert Carba-R Assay is for use with the following sample types:

Rectal Swab Specimens

The assay is performed on rectal swab specimens from patients at risk for intestinal colonization with carbapenem-non-susceptible bacteria. Concomitant cultures are necessary to recover organisms for epidemiological typing, antimicrobial susceptibility testing, and for further confirmatory bacterial identification.

Pure Colonies

The assay is performed on carbapenem-non-susceptible pure colonies of *Enterobacteriaceae*, *Acinetobacter baumannii*, or *Pseudomonas aeruginosa*, when grown on blood agar or MacConkey agar. For testing pure colonies, the Xpert Carba-R Assay should be used in conjunction with other laboratory tests including phenotypic antimicrobial susceptibility testing.

---

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

---

---

**CONTINUE ON A SEPARATE PAGE IF NEEDED.**

---

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
PRAStaff@fda.hhs.gov

*“An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number.”*

## 510(k) Summary

As required by 21 CFR Section 807.92(c).

Submitted by: Cepheid  
904 Caribbean Drive  
Sunnyvale, CA 90489  
Phone number: (847) 228-3299  
Fax number: (847) 890-6589

Contact: Scott A. Campbell, PhD, MBA

Date of Preparation: June 28, 2016

Device:

Trade name: Xpert<sup>®</sup> Carba-R

Common name: Xpert Carba-R Assay

Type of Test: Qualitative nucleic acid amplification test of the *bla*<sub>KPC</sub>, *bla*<sub>NDM</sub>, *bla*<sub>VIM</sub>, *bla*<sub>OXA-48</sub>, and *bla*<sub>IMP</sub> gene sequences associated with carbapenem-non-susceptibility in gram-negative bacteria obtained from rectal swab specimens

Classification: II

Regulation number: 866.1640

Classification name: Antimicrobial susceptibility test powder

Product code: POC, OOI

Classification Advisory Panel: Microbiology (83)

Prescription Use: Yes

Predicate Device Assay: Cepheid Xpert<sup>®</sup> Carba-R  
[510(k) #K152614]

### Device Description:

The Xpert Carba-R Assay is an automated real-time polymerase chain reaction (PCR) *in vitro* diagnostic test for qualitative detection of the *bla*<sub>KPC</sub>, *bla*<sub>NDM</sub>, *bla*<sub>VIM</sub>, *bla*<sub>OXA-48</sub>, and *bla*<sub>IMP</sub> gene sequences from rectal swab specimens. The Xpert Carba-R Assay is intended as an aid for infection control for monitoring the spread of carbapenem-non-susceptible organisms in healthcare settings.

The Xpert Carba-R Assay is performed on the Cepheid GeneXpert<sup>®</sup> Instrument Systems (GeneXpert Dx, GeneXpert Infinity-48, GeneXpert Infinity-48s, and GeneXpert Infinity-80 systems). The GeneXpert Instrument System platform automates sample preparation, amplification and real-time detection.

The GeneXpert Instrument Systems require the use of single-use, disposable cartridges (the Xpert Carba-R cartridges) that hold the PCR reagents and host the PCR process. Because the cartridges are self-contained and specimens never come into contact with working parts of the instrument modules, cross-contamination between samples is minimized.

The Xpert Carba-R Assay cartridges contain reagents for the detection of *bla*<sub>KPC</sub>, *bla*<sub>NDM</sub>, *bla*<sub>VIM</sub>, *bla*<sub>OXA-48</sub>, and *bla*<sub>IMP</sub> gene sequences. A Sample Processing Control (SPC) and a Probe Check Control (PCC) are controls utilized by the GeneXpert Instrument System platform. The SPC is present to control for adequate processing of the target bacteria and to monitor the presence of inhibitors in the real-time PCR reaction to reduce the possibility of false negative results. The PCC verifies reagent rehydration, real-time PCR tube filling in the cartridge, probe integrity, and dye stability.

The single-use, multi-chambered fluidic cartridges are designed to complete sample preparation and real-time PCR for the detection of the *bla*<sub>KPC</sub>, *bla*<sub>NDM</sub>, *bla*<sub>VIM</sub>, *bla*<sub>OXA-48</sub>, and *bla*<sub>IMP</sub> gene sequences from rectal swab specimens in approximately 50 minutes. The GeneXpert Instrument Systems, comprised of the GeneXpert Dx Systems and the GeneXpert Infinity Systems, have 1 to 80 randomly accessible modules, depending upon the instrument, that are each capable of performing separate sample processing and real-time PCR and RT-PCR tests. Each module contains a syringe drive for dispensing fluids (i.e., the syringe drive activates the plunger that works in concert with the rotary valve in the cartridge to move fluids between chambers), an ultrasonic horn for lysing cells or spores, and a proprietary I-CORE<sup>®</sup> thermocycler for performing real-time PCR and RT-PCR and detection.

Rectal swab specimens are placed into a sample reagent. The sample is transferred to the sample chamber of the disposable fluidic cartridge (the Xpert Carba-R cartridge). The user initiates a test from the system user interface and places the cartridge into the GeneXpert instrument platform, which performs hands-off real-time, multiplex PCR for detection of the *bla*<sub>KPC</sub>, *bla*<sub>NDM</sub>, *bla*<sub>VIM</sub>, *bla*<sub>OXA-48</sub>, and *bla*<sub>IMP</sub> gene sequences. The results are automatically generated at the end of the process in a report that can be viewed and printed.

### Device Intended Use:

The Xpert<sup>®</sup> Carba-R Assay, performed on the GeneXpert<sup>®</sup> Instrument Systems, is a qualitative *in vitro* diagnostic test designed for the detection and differentiation of the *bla*<sub>KPC</sub>, *bla*<sub>NDM</sub>, *bla*<sub>VIM</sub>, *bla*<sub>OXA-48</sub>, and *bla*<sub>IMP</sub> gene sequences associated with carbapenem-non-susceptibility. The test utilizes automated real-time polymerase chain reaction (PCR).

The Xpert Carba-R Assay is intended as an aid to infection control in the detection of carbapenem-non-susceptible bacteria that colonize patients in healthcare settings. The Xpert Carba-R Assay is not intended to guide or monitor treatment for carbapenem-non-susceptible bacterial infections. A negative Xpert Carba-R Assay result does not preclude the presence of other resistance mechanisms.

The Xpert Carba R-Assay is for use with the following sample types:

#### Rectal Swab Specimens

The assay is performed on rectal swab specimens from patients at risk for intestinal colonization with carbapenem-non-susceptible bacteria. Concomitant cultures are necessary to recover organisms for epidemiological typing, antimicrobial susceptibility testing, and for further confirmatory bacterial identification.

#### Pure Colonies

The assay is performed on carbapenem-non-susceptible pure colonies of *Enterobacteriaceae*, *Acinetobacter baumannii*, or *Pseudomonas aeruginosa*, when grown on blood agar or MacConkey agar. For testing pure colonies, the Xpert Carba-R Assay should be used in conjunction with other laboratory tests including phenotypic antimicrobial susceptibility testing.

### Substantial Equivalence:

The Cepheid Xpert Carba-R Assay is substantially equivalent to the Xpert<sup>®</sup> Carba-R Assay, 510(k) #K152614. Both assays utilize the same GeneXpert cartridge and detect target gene sequences (*bla*<sub>KPC</sub>, *bla*<sub>NDM</sub>, *bla*<sub>VIM</sub>, *bla*<sub>OXA-48</sub>, and *bla*<sub>IMP</sub>) using real-time PCR amplification and fluorogenic target-specific hybridization detection. The performance of the Xpert Carba-R Assay for an expanded indication was determined in a multi-site clinical study in which the performance of the Xpert Carba-R Assay was evaluated relative to culture and reference DNA sequence analysis. The results of the study demonstrated that the performance of the Xpert Carba-R Assay is substantially equivalent to the predicate device. Please refer to K152614 for information on the performance of the Xpert Carba-R Assay with pure colonies of carbapenem non-susceptible *Enterobacteriaceae*, *Acinetobacter baumannii*, and *Pseudomonas aeruginosa*.

Table 5-1 shows the similarities and differences between the Xpert Carba-R Assay and the predicate device.

**Table 5-1: Comparison of Similarities and Differences of the Xpert Carba-R Assay with the Predicate Device**

Similarities		
Item	Device	Predicate Device
	Cepheid Xpert Carba-R Assay	Cepheid Xpert Carba-R Assay K152614
General Intended Use	<p>The Xpert<sup>®</sup> Carba-R Assay, performed on the GeneXpert<sup>®</sup> Instrument Systems, is a qualitative <i>in vitro</i> diagnostic test designed for the <b>detection and differentiation of the <i>bla</i><sub>KPC</sub>, <i>bla</i><sub>NDM</sub>, <i>bla</i><sub>VIM</sub>, <i>bla</i><sub>OXA-48</sub>, and <i>bla</i><sub>IMP</sub> gene sequences associated with carbapenem-non-susceptibility. The test utilizes automated real-time polymerase chain reaction (PCR).</b></p> <p>The Xpert Carba-R Assay is intended as <b>an aid to infection control</b> in the detection of carbapenem-non-susceptible bacteria that colonize patients <b>in healthcare settings. The Xpert Carba-R Assay is not intended to guide or monitor treatment for carbapenem-non-susceptible bacterial infections. A negative Xpert Carba-R Assay result does not preclude the presence of other resistance mechanisms.</b></p> <p>The Xpert Carba-R Assay is for use with the following sample types: Rectal Swab Specimens</p> <p>The assay is performed on rectal swab specimens from patients at risk for intestinal colonization with carbapenem-non-susceptible bacteria. Concomitant cultures are necessary to recover organisms for epidemiological typing, antimicrobial susceptibility testing, and for further confirmatory bacterial identification.</p> <p>Pure Colonies</p> <p>The assay is performed on carbapenem-non-susceptible pure colonies of <i>Enterobacteriaceae</i>, <i>Acinetobacter baumannii</i>, or <i>Pseudomonas aeruginosa</i>, when grown on blood agar or MacConkey agar. For testing pure colonies, the Xpert Carba-R Assay should be used in conjunction with other</p>	<p>The Xpert<sup>®</sup> Carba-R Assay, performed on the GeneXpert<sup>®</sup> Instrument Systems, is a qualitative <i>in vitro</i> diagnostic test for the <b>detection and differentiation of the <i>bla</i><sub>KPC</sub>, <i>bla</i><sub>NDM</sub>, <i>bla</i><sub>VIM</sub>, <i>bla</i><sub>OXA-48</sub>, and <i>bla</i><sub>IMP</sub> gene sequences associated with carbapenem-non-susceptible pure colonies of <i>Enterobacteriaceae</i>, <i>Acinetobacter baumannii</i>, or <i>Pseudomonas aeruginosa</i> grown on blood agar or MacConkey agar. The test utilizes automated real-time polymerase chain reaction (PCR).</b></p> <p><b>A negative Xpert Carba-R Assay result does not preclude the presence of other resistance mechanisms.</b> The Xpert Carba-R Assay <b>should be used in conjunction with other laboratory tests including phenotypic antimicrobial susceptibility testing.</b> The Xpert Carba-R Assay is intended <b>as an aid for infection control</b> in detecting and differentiating genetic markers of resistance to monitor the spread of carbapenem-non-susceptible organisms <b>in healthcare settings. The Xpert Carba-R Assay is not intended to guide or monitor treatment for carbapenem-non-susceptible bacterial infections.</b></p>

<b>Similarities</b>		
<b>Item</b>	<b>Device</b>	<b>Predicate Device</b>
	<b>Cepheid Xpert Carba-R Assay</b>	<b>Cepheid Xpert Carba-R Assay K152614</b>
	laboratory tests including phenotypic antimicrobial susceptibility testing.	
Type of test	Same	Qualitative
Technological Principles	Same	Fully-automated nucleic acid amplification (DNA); real-time PCR
Test Cartridge	Same	Disposable single-use, multi-chambered fluidic cartridge
Probes	Same	TaqMan <sup>®</sup> Probes
Controls	Same	Internal sample processing control (SPC) and probe check control (PCC)  External controls available
Instrument System	Same	GeneXpert Instrument System (includes GeneXpert Dx, Infinity-48, Infinity-48s, and Infinity-80)
Time to obtain test results	Same	Approximately 50 minutes to results
Interpretation of test results	Same	Diagnostic software of the GeneXpert Instrument System
Laboratory Users	Same	Operators in CLIA Moderate or High Complexity labs
<b>Differences</b>		
	<b>New Device</b>	<b>Predicate Device</b>
	<b>Cepheid Xpert Carba-R Assay</b>	<b>Cepheid Xpert Carba-R Assay K152614</b>
Sample Types	Bacterial isolates from culture, <b>rectal swab specimens</b>	Bacterial isolates from culture

The Xpert Carba-R Assay has the same general intended use as the predicate device and has the same technological characteristics as the predicate device. The differences between the Xpert Carba-R Assay and the predicate device do not raise questions of safety and effectiveness. The clinical study demonstrates that the Xpert Carba-R Assay is acceptable for its intended use with inexperienced laboratory users and is substantially equivalent to the predicate device described above.

Non-Clinical Studies:

**Analytical Sensitivity (Limit of Detection)**

The analytical sensitivity or Limit of Detection (LoD) of the Xpert Carba-R Assay was assessed using carbapenemase-producing organisms seeded into pooled negative human rectal swab matrix. The LoD was determined for two carbapenemase-producing bacteria for each gene analyte, i.e., the genes encoding KPC, NDM, VIM, OXA-48, and IMP. Bacteria were titrated by plate counts and spiked onto clean swabs. Swabs were placed into pooled negative rectal swab matrix and replicates of 20 were evaluated at a minimum of five different concentrations over four days. The LoD for each of the ten carbapenemase-producing organisms was estimated by probit analysis. The LoD is defined as the lowest concentration of target cells (CFU/swab) that can be reproducibly distinguished from negative samples with 95% confidence. The study was performed with two different lots of Xpert Carba-R reagents and the claimed LoD is the higher of the two determinations. The estimated LoDs were verified by preparing and testing 10 replicates from two independent dilutions of each bacterium at each estimated LoD.

The claimed LoD for each pair of carbapenemase-producing organism in rectal swab matrix are shown in Table 5-2.

**Table 5-2. LoD Estimates and Verification for Organisms Harboring Carbapenemase Genes using the Xpert Carba-R Assay in Rectal Swab Matrix**

Target Gene and Organism	LoD Estimates (Probit) CFU /swab		LoD Claim CFU/swab	Estimated LoD In Sample Reagent CFU/mL	Verification (Positives/20)
	Lot 1	Lot 2			
IMP-1 <i>Acinetobacter baumannii</i>	<b>174</b>	141	174	35	20/20
IMP-1 <i>Klebsiella pneumoniae</i>	303	<b>306</b>	306	61	20/20
VIM-1 <i>Klebsiella pneumoniae</i>	247	<b>305</b>	305	61	20/20
VIM-4 <i>Escherichia coli</i>	<b>815</b>	468	815	163	20/20
NDM-1 <i>Klebsiella pneumoniae</i> ATCC BAA-2146	117	<b>251</b>	251	50	20/20
NDM <i>Klebsiella pneumoniae</i>	<b>74</b>	57	74	15	19/20
KPC-3 <i>Klebsiella pneumoniae</i> NCTC 13438	<b>373</b>	292	373	75	20/20
KPC <i>Enterobacter cloacae</i>	<b>779</b>	537	779	156	20/20
OXA-48 <i>Enterobacter cloacae</i>	<b>154</b>	109	154	31	20/20
OXA-48 <i>Escherichia coli</i>	<b>104</b>	99	104	21	20/20

**Analytical Reactivity (Inclusivity)**

The analytical reactivity of the Xpert Carba-R Assay with rectal swab matrix was evaluated by testing a panel of 72 samples. This panel consisted of 11 *bla*<sub>KPC</sub> (KPC), 11 *bla*<sub>VIM</sub> (VIM), 13 *bla*<sub>NDM</sub> (NDM), 8 *bla*<sub>OXA-48</sub> (OXA-48), 5 *bla*<sub>NDM</sub>/*bla*<sub>OXA-181</sub> (NDM/OXA-181), 6 *bla*<sub>OXA-181</sub> (OXA-181), 17 *bla*<sub>IMP</sub> (IMP), and one *bla*<sub>KPC</sub>/*bla*<sub>VIM</sub> (KPC/VIM) well-characterized

bacterial strains. The strains tested in rectal swab matrix and their test concentrations are presented in Table 5-3.

For testing in rectal swab matrix, organisms were seeded into pooled negative rectal swab matrix. All bacterial strains were tested in triplicate at approximately 3x LoD for each specimen type. Xpert Carba-R Assay target genes were detected in 69 of 72 carbapenemase-producing bacterial strains although IMP-4 was detected only using a higher concentration (Table 5-3). Xpert Carba-R Assay target DNA sequences were not detected in three bacterial strains as shown in Table 5-3. In one of the three bacterial strains, the IMP-13 gene was not detected by the assay, although it was predicted to be detected by *in silico* analysis. In two of the other three bacterial strains, the IMP-7 and IMP-14 genes were not predicted to be detected by *in silico* analysis and were not detected by the assay. See Limitations in the package insert.

The variants detected, and predictions for detecting other subtypes of each resistance gene based on *in silico* analysis, are presented in Table 5-4.

**Table 5-3. Analytical Reactivity of the Xpert Carba-R Assay in Rectal Swab Matrix**

Strain ID	Organism	Resistance Marker with variant information	Concentration Tested in Rectal Swab Matrix (CFU/mL)
NCTC 13438	<i>Klebsiella pneumoniae</i>	KPC-3	153
31551	<i>Klebsiella pneumoniae</i>	KPC-4	50
ATCC BAA-1705	<i>Klebsiella pneumoniae</i>	KPC-2	130
PA-Col	<i>Pseudomonas aeruginosa</i>	KPC-2	250
KBM18	<i>Enterobacter aerogenes</i>	KPC-2	250
BM9	<i>Klebsiella pneumoniae</i>	KPC-3	330
PA3	<i>Klebsiella pneumoniae</i>	KPC-2	100
CGNC	<i>Serratia marcescens</i>	KPC-2	300
CFVL	<i>Enterobacter cloacae</i>	KPC-2	160
COL	<i>Escherichia coli</i>	KPC-2	147
GR-04/KP-69	<i>Klebsiella pneumoniae</i>	KPC-2, VIM	80
164-3	<i>Klebsiella oxytoca</i>	KPC	70
NCTC 13437	<i>Pseudomonas aeruginosa</i>	VIM-10	500
NCTC 13439	<i>Klebsiella pneumoniae</i>	VIM-1	130
NCTC 13440	<i>Klebsiella pneumoniae</i>	VIM-1	70
758	<i>Pseudomonas aeruginosa</i>	VIM	250
PA-87	<i>Klebsiella pneumoniae</i>	VIM	200
B92A	<i>Pseudomonas aeruginosa</i>	VIM	2000
Col1	<i>Pseudomonas aeruginosa</i>	VIM-2	500
BM19	<i>Serratia marcescens</i>	VIM-2	250
KOW7	<i>Escherichia coli</i>	VIM-4	250
DIH	<i>Klebsiella pneumoniae</i>	VIM-19	250
MSH2014-3	<i>Enterobacter cloacae</i>	VIM	500
NCTC 13443	<i>Klebsiella pneumoniae</i>	NDM-1	80

Strain ID	Organism	Resistance Marker with variant information	Concentration Tested in Rectal Swab Matrix (CFU/mL)
ATCC BAA-2146	<i>Klebsiella pneumoniae</i>	NDM-1	80
34262	<i>Klebsiella pneumoniae</i>	NDM	80
GEN	<i>Acinetobacter baumannii</i>	NDM-1	130
3047	<i>Enterobacter cloacae</i>	NDM-1	70
7892	<i>Proteus mirabilis</i>	NDM-1	30
CAN	<i>Salmonella spp.</i>	NDM-1	70
EGY	<i>Acinetobacter baumannii</i>	NDM-2	40
I5	<i>Escherichia coli</i>	NDM-4	30
405	<i>Escherichia coli</i>	NDM-5	30
CF-ABE	<i>Citrobacter freundii</i>	NDM	30
73999	<i>Pseudomonas aeruginosa</i>	NDM	50
39365	<i>Providencia rettgeri</i>	NDM-1	70
NCTC 13442	<i>Klebsiella pneumoniae</i>	OXA-48	40
OM11	<i>Klebsiella pneumoniae</i>	OXA-48	60
501	<i>Enterobacter cloacae</i>	OXA-48	80
DUW	<i>Klebsiella pneumoniae</i>	OXA-48	120
OM22	<i>Escherichia coli</i>	OXA-48	80
BOU	<i>Enterobacter cloacae</i>	OXA-48	80
TUR	<i>Enterobacter cloacae</i>	OXA-48	120
11670	<i>Escherichia coli</i>	OXA-48	100
166643	<i>Klebsiella pneumoniae</i>	OXA-181	20
42194	<i>Klebsiella pneumoniae</i>	OXA-181	20
MSH2014-64	<i>Klebsiella pneumoniae</i>	OXA-181	280
MSH2014-72	<i>Escherichia coli</i>	OXA-181	100
74	<i>Escherichia coli</i>	OXA-181	100
CDC0051	<i>Klebsiella ozaenae</i> <sup>a</sup>	OXA-181	250
B108A	<i>Klebsiella pneumoniae</i>	NDM, OXA-181	10
C10192-DISCS	<i>Enterobacter aerogenes</i>	NDM, OXA-181	10
KP-OMA3	<i>Klebsiella pneumoniae</i>	NDM, OXA-181	60
1300920	<i>Klebsiella pneumoniae</i>	NDM, OXA-181	15
MSH2014-69	<i>Klebsiella pneumoniae</i>	NDM, OXA-181	20
NCTC 13476	<i>Escherichia coli</i>	IMP-1	250
695	<i>Acinetobacter baumannii</i>	IMP-1	1720
2340	<i>Enterobacter cloacae</i>	IMP-1	250
IMPBMI	<i>Klebsiella pneumoniae</i>	IMP-1	100
Yonsei_1	<i>Acinetobacter baumannii</i>	IMP-1	1000
Yonsei_2	<i>Acinetobacter baumannii</i>	IMP-1	500
6852	<i>Klebsiella pneumoniae</i>	IMP-1	100
MKAM	<i>Pseudomonas aeruginosa</i>	IMP-1	500
70450-1	<i>Pseudomonas aeruginosa</i>	IMP-1	250
3994	<i>Pseudomonas spp.</i>	IMP-10	250

Strain ID	Organism	Resistance Marker with variant information	Concentration Tested in Rectal Swab Matrix (CFU/mL)
CDC0161	<i>Enterobacter aerogenes</i> <sup>a</sup>	IMP-4	5.00E+04
5344	<i>Pseudomonas aeruginosa</i>	IMP-2	60
3985	<i>Pseudomonas aeruginosa</i>	IMP-11	2000
4032	<i>Pseudomonas aeruginosa</i>	IMP-6	80
3424	<i>Pseudomonas aeruginosa</i>	IMP-7 <sup>b, c</sup>	1.00E+06
32443	<i>Klebsiella pneumoniae</i>	IMP-13 <sup>c</sup>	1.00E+06
92	<i>Pseudomonas aeruginosa</i>	IMP-14 <sup>b, c</sup>	1.00E+06

a. These organisms were not tested as bacterial isolates.

b. IMP-7 and IMP-14 genes (*Pseudomonas aeruginosa*) were not detected by the assay and were not predicted to be detected by *in silico* analysis (see Limitations).

c. IMP-13 gene (*Klebsiella pneumoniae*): although predicted to be detected by *in silico* analysis, the IMP-13 gene was not detected by the assay (see Limitations).

**Table 5-4. Summary of Variants Detected by Wet Testing or Predicted to be Detected Based on *In Silico* Analysis**

Marker (or Traditional Subgroup)	Wet testing			Not tested but predicted to be detected based on <i>in silico</i> analysis
	No. of Samples	Type(s) Detected	Type(s) not Detected	
KPC	12	KPC-2, 3, 4	--	KPC-5, 6, 7, 8, 9, 10, 11, 12, 13, 14, 15, 16
NDM	18	NDM-1, 2, 4, 5	--	NDM-3, 6, 7, 8, 9
VIM	12	VIM-1, 2, 4, 10, 19	--	VIM-5, 6, 7, 8, 9, 11, 12, 13, 14, 15, 16, 17, 18, 20, 23, 24, 25, 26, 27, 28, 29, 30, 31, 32, 33, 34, 35, 36, 37, 38
OXA-48	19	OXA-48, 181(OXA-48 variant)	--	OXA-162, 163, 204, 232, 244, 245, 247
IMP	17	IMP-1 (9 strains), IMP-2, 4, 6, 10, 11	IMP-7 <sup>a</sup> , 13 <sup>b</sup> , 14 <sup>a</sup>	IMP-3, 8, 9, 13 <sup>b</sup> , 19, 20, 21, 22, 24, 25, 27, 28, 30, 31, 33, 37, 40, 42

a. IMP-7 and IMP-14 genes (*Pseudomonas aeruginosa*) were not detected by the assay and were not predicted to be detected by *in silico* analysis (see Limitations).

b. IMP-13 gene (*Klebsiella pneumoniae*) was tested: although predicted to be detected by *in silico* analysis, the IMP-13 gene was not detected by the assay (see Limitations).

### Analytical Specificity (Cross-reactivity)

The analytical specificity of the Xpert Carba-R Assay was evaluated for bacterial isolates, organisms seeded into rectal swab matrix. For all three specimen types, a panel of 62 well-characterized bacterial strains of carbapenem-susceptible bacteria or bacteria with carbapenem non-susceptibility due to genes or mechanisms other than the Xpert Carba-R target genes (Table 5-5 and Table 5-6) and 24 commensal bacterial strains and other enteric

microorganisms were also evaluated in the study (Table 5-7). Human cells were also tested in rectal swab matrix (Table 5-8). Resistance mechanisms were determined by individual PCR assays, DNA sequence analysis, or Check-Points array version CT102.

For rectal swab matrix samples, 62 strains were tested at concentrations  $>1 \times 10^6$  CFU/mL with the exception of *Peptostreptococcus anaerobius* that was tested at  $5 \times 10^5$  CFU/mL. Viruses were tested at  $>1 \times 10^5$  TCID<sub>50</sub>/mL or greater than  $2.5 \times 10^7$  RNA copies/mL. A bladder cell line (human genomic DNA) was tested at  $1 \times 10^5$  cells/mL. Organisms were diluted into pooled negative rectal swab matrix and tested in triplicate. None of the 94 potentially cross-reactive organisms and nucleic acids tested was detected with the Xpert Carba-R Assay.

The Xpert Carba-R Assay did not cross react with any of the organisms tested (Table 5-6, Table 5-7, and Table 5-8). The analytical specificity of the assay was 100%.

**Table 5-5. Number of Carbapenem-susceptible and Non-susceptible Organisms for each Antibiotic**

	Ertapenem	Imipenem	Meropenem
<b>Susceptible</b>	19	30	24
<b>Intermediate</b>	0	8	4
<b>Resistant</b>	43	24	34

**Table 5-6. Cross-reactivity Panel**

Organism	Strain ID	Confirmed Resistance Mechanisms	Carbapenem susceptibility (S/I/R) <sup>a</sup>		
			ETP <sup>a</sup>	IMP <sup>a</sup>	MEM <sup>a</sup>
<i>Escherichia coli</i>	NCTC 13441	CTX-M (-1, -type 15 like); TEM	S	S	S
<i>Klebsiella pneumoniae</i>	NCTC 13465	CTX-M (25)	S	S	S
<i>Enterobacter aerogenes</i>	810	OmpC/OmpF deficient; TEM	R	R	R
<i>Citrobacter freundii</i>	1698	TEM (WT+164S)	S	S	S
<i>Enterobacter cloacae</i>	5557	AmpC (ACT/MIR)	R	R	R
<i>Klebsiella pneumoniae</i>	kpn5	CTX-M-2	R	S	R
<i>Klebsiella pneumoniae</i>	kpn12	TEM; SHV; CTX-M	R	R	R
<i>Escherichia coli</i>	eco1	TEM; CTX-M-2	R	R	R
<i>Escherichia coli</i>	eco2	CTX-M (2); TEM	R	S	S
<i>Enterobacter cloacae</i>	cor1	CTX-M (2); TEM	R	R	R
<i>Serratia marcescens</i>	hpp21	CTX-M (2); TEM	S	S	S
<i>Morganella morganii</i>	fer29	CTX-M (2); TEM	S	R	S
<i>Proteus mirabilis</i>	gut25	CTX-M (2); TEM	S	R	S
<i>Salmonella spp.</i>	3209	CTX-M (2); TEM	S	S	S

Organism	Strain ID	Confirmed Resistance Mechanisms	Carbapenem susceptibility (S/I/R) <sup>a</sup>		
			ETP <sup>a</sup>	IMP <sup>a</sup>	MEM <sup>a</sup>
<i>Shigella flexnerii</i>	3331	CTX-M (2); TEM	S	S	S
<i>Enterobacter cloacae</i>	PA_3	AmpC; CTX-M-15; TEM	S	S	S
<i>Klebsiella pneumoniae</i>	32189	SHV	S	S	S
<i>Klebsiella pneumoniae</i>	32443	CTX-M (1, -type 15 like); SHV	S	S	S
<i>Klebsiella pneumoniae</i>	32598	CTX-M (-1, -type 15 like); SHV; TEM	R	I	R
<i>Klebsiella pneumoniae</i>	33560	CTX-M (15); SHV-11; TEM-1	S	S	S
<i>Klebsiella pneumoniae</i>	33603	SHV-2	R	I	R
<i>Klebsiella pneumoniae</i>	33617	SHV-27	S	S	S
<i>Klebsiella pneumoniae</i>	33643	SHV (-5, -55); TEM	S	S	S
<i>Klebsiella pneumoniae</i>	34430	SHV; TEM; CTX-M-15	S	S	S
<i>Klebsiella pneumoniae</i>	34680	TEM; CTX-M-2	R	S	R
<i>Klebsiella pneumoniae</i>	34732	CTX-M (15); SHV; TEM	R	S	S
<i>Enterobacter cloacae</i>	PA_174	GX-/Culture+; SHV; TEM	S	S	S
<i>Enterobacter aerogenes</i>	STU 645	SHV (WT+238S+240K)	R	S	R
<i>Enterobacter aerogenes</i>	STU 669	SHV (WT+238S+240K)	R	R	R
<i>Escherichia coli</i>	C3015	AmpC (CMY II); TEM	R	R	R
<i>Enterobacter aerogenes</i>	RI_100	AmpC (DHA); SHV	R	R	R
<i>Klebsiella pneumoniae</i>	B4A	SHV (WT + 238S +240K)	R	R	R
<i>Klebsiella pneumoniae</i>	B13A	SHV (WT + 238S +240K)	R	S	S
<i>Enterobacter cloacae</i>	RI_474	AmpC (ACT/MIR)	R	I	I
<i>Enterobacter amnigenus</i>	B71	AmpC (ACT/MIR)	R	R	R
<i>Klebsiella pneumoniae</i>	DD82A	SHV (WT + 238S + 240K)	R	S	R
<i>Klebsiella pneumoniae</i>	B100	CTX-M (-1, type-15 like); SHV (WT+238S); TEM	R	S	R
<i>Enterobacter cloacae</i>	135B	TEM	S	S	S
<i>Klebsiella pneumoniae</i>	B157	SHV; TEM	R	R	R
<i>Escherichia coli</i>	T2914280	CTX-M (-1, -15); TEM	R	S	R
<i>Providencia stuartii</i>	DD188	TEM (104K + 164S)	R	I	I
<i>Enterobacter cloacae</i>	DD189	AmpC (ACT/MIR)	R	S	S

Organism	Strain ID	Confirmed Resistance Mechanisms	Carbapenem susceptibility (S/I/R) <sup>a</sup>		
			ETP <sup>a</sup>	IMP <sup>a</sup>	MEM <sup>a</sup>
<i>Escherichia coli</i>	B198B	CTX-M (-1, type -15 like); TEM	R	S	R
<i>Klebsiella pneumoniae</i>	T3019989-1	CTX-M (-1, type-15 like); SHV	R	I	R
<i>Klebsiella pneumoniae</i>	T3019989-2	CTX-M (-1, type-15 like); SHV	R	S	R
<i>Enterobacter cloacae</i>	ENC-THAI14	VEB-1, TEM	S	S	S
<i>Escherichia coli</i>	CB154006	CTX-M (9); TEM	R	I	I
<i>Enterobacter cloacae</i>	S35766	AmpC(ACT/MIR)	S	S	S
<i>Enterobacter cloacae</i>	X1856910	AmpC (ACT/MIR); TEM	R	I	I
<i>Klebsiella pneumoniae</i>	W3758164	CTX-M (-1, -15 like); SHV; TEM.	R	I	R
<i>Klebsiella pneumoniae</i>	X2135758	CTX-M (-1, -15 like); SHV	R	S	S
<i>Klebsiella pneumoniae</i>	W3809535	CTX-M (-1, -15 like); SHV	R	R	R
<i>Pseudomonas aeruginosa</i>	CDC0064	SPM	R	R	R
<i>Serratia marcescens</i>	CDC0099	SME	R	R	R
<i>Serratia marcescens</i>	CDC0121	SME	R	R	R
<i>Serratia marcescens</i>	CDC0122	SME	R	R	R
<i>Serratia marcescens</i>	CDC0123	SME	R	R	R
<i>Serratia marcescens</i>	CDC0124	SME	R	R	R
<i>Serratia marcescens</i>	CDC0130	SME	R	R	R
<i>Serratia marcescens</i>	CDC0131	SME	R	R	R
<i>Enterobacter cloacae</i> group	CDC0132	IMI	R	R	R
<i>Enterobacter cloacae</i> complex	CDC0164	IMI	R	R	R

a. S/I/R = Susceptible/Intermediate/Resistant, ETP = Ertapenem, IMP = Imipenem, MEM = Meropenem

**Table 5-7. Cross-reactivity Panel (Commensal and Other Enteric Microorganisms)**

Strain ID	Organism	Concentration Tested (CFU/mL unless otherwise specified)
ATCC 25922	<i>Escherichia coli</i>	2.67E+06
ATCC 29212	<i>Enterococcus faecalis</i>	3.15E+06
ATCC 700603	<i>Klebsiella pneumoniae</i>	5.20E+06
ATCC 35218	<i>Escherichia coli</i>	2.47E+06
ATCC 25923	<i>Staphylococcus aureus</i>	4.53E+06
ATCC 27853	<i>Pseudomonas aeruginosa</i>	3.17E+06

Strain ID	Organism	Concentration Tested (CFU/mL unless otherwise specified)
ATCC 9689	<i>Clostridium difficile</i>	1.80E+07
ATCC 700621	<i>Enterobacter cloacae</i>	8.95E+06
ATCC 9756	<i>Enterococcus faecium</i>	6.54E+06
ATCC 13182	<i>Klebsiella oxytoca</i>	4.76E+06
ATCC BAA-747	<i>Acinetobacter baumannii</i>	2.27E+06
ATCC 33128	<i>Citrobacter freundii</i>	2.01E+06
ATCC 49948	<i>Morganella morganii</i>	8.19E+06
ATCC 51331	<i>Stenotrophomonas maltophilia</i>	3.15E+06
ATCC 27028	<i>Citrobacter koseri</i>	5.05E+06
ATCC 49809	<i>Providencia stuartii</i>	3.01E+06
ATCC 49037	<i>Peptostreptococcus anaerobius</i>	5.00E+05
CCUG 29780 / ATCC 12401	<i>Streptococcus agalactiae</i>	5.21E+06
ATCC 15703	<i>Bifidobacterium adolescentis</i>	1.10E+08
ATCC 51697	<i>Enterobacter aerogenes</i>	3.19E+06
ATCC 43071	<i>Proteus mirabilis</i>	1.78E+06
CCUG 34787	<i>Acinetobacter spp.</i>	2.40E+06
CCUG 418	<i>Citrobacter freundii</i>	2.95E+06
CCUG 33629	<i>Corynebacterium diphtheriae</i>	4.48E+06
CCUG 17874	<i>Helicobacter pylori</i>	1.61E+06
CCUG 33548	<i>Listeria monocytogenes</i>	4.77E+06
CCUG 6325	<i>Providencia alcalifaciens</i>	4.91E+06
CCUG 43594 / ATCC 33560	<i>Campylobacter jejuni</i>	3.27E+06
MRVP/ZeptoMetrix	Adenovirus B Type 7A/NY	1.40E+05 TCID <sub>50</sub> /mL
MRVP/ZeptoMetrix	Enterovirus Type 71/NY	4.40E+05 TCID <sub>50</sub> /mL
Clinical Sample – Cepheid	Norovirus GII	2.5 x 10 <sup>7</sup> RNA copies/mL

**Table 5-8. Cell Line Representing Human Genomic DNA**

Organism Name	Source
Bladder Cell Carcinoma (hgDNA)	ATCC HTB-4

### Competitive Interference

A competitive interference study was performed to test whether a high titer of one or more carbapenemase-producing organisms would interfere with the detection of a second target carbapenemase-producing organism that was present at a low titer. High titered samples were formulated at concentrations of  $5 \times 10^6$  CFU/swab and low titered targets were formulated at approximately 2x LoD for the respective strain in rectal swab matrix. One carbapenemase-producing bacterial strain for each gene analyte, i.e., the genes encoding

KPC, NDM, VIM, OXA-48, and IMP, was used in this study. Each carbapenemase-producing bacterial strain type was tested at low titers in conjunction with a high titer of each of the other one or two carbapenemase-producing bacterial strain types (Table 5-9). Samples were tested in replicates of eight. An inhibitory effect was observed for three of the five targets (IMP, VIM, and OXA-48) when a low concentration of each target was present in combination with a high concentration of one or two other targets for samples tested in rectal swab matrix. The three targets (IMP, VIM, and OXA-48) were tested at a higher concentration (4x LoD) in combination with a high concentration of one of two other targets for samples in rectal swab matrix. No inhibitory effect was observed for the three targets (IMP, VIM and OXA-48) at 4x LoD in the presence of clinically relevant co-infections for the Xpert Carba-R Assay. The competitive inhibitory effect on the Carba-R targets (IMP, VIM and OXA-48) is addressed in Section 15, Limitations in the package insert.

**Table 5-9. Combinations of Carbapenemase-producing Bacteria Tested with the Xpert Carba-R Assay**

<b>Combination</b>
High KPC/High NDM/Low VIM
High KPC/High NDM/Low OXA
High KPC/High NDM/Low IMP
High VIM/High OXA/Low KPC
High VIM/High OXA/Low NDM
High VIM/High OXA/Low IMP
High IMP/Low KPC
High IMP/Low NDM
High IMP/Low VIM
High IMP/Low OXA
High OXA/Low VIM
High VIM/Low OXA
High KPC/Low NDM
Negative

### **Potentially Interfering Substances**

The performance of the Xpert Carba-R Assay was evaluated with 24 potentially interfering substances that may be present in rectal swab specimens. Potentially interfering substances (IS) solutions were prepared and tested at concentrations specified in Table 5-10. Positive and negative samples were included in this study. Positive samples consisted of a mix of five carbapenemase-producing organisms harboring KPC, NDM, VIM, IMP-1 and OXA-48 gene sequences seeded into pooled negative rectal swab matrix at approximately 3x LoD. Eight replicate positive samples were tested per substance. Negative samples consisted of pooled negative rectal swab matrix not seeded with carbapenemase-producing organisms. Eight replicate negative samples were tested per substance to determine the effect on the

performance of the sample processing control (SPC). Controls consisted of positive and negative samples with no interfering substances added. The effect of each potentially interfering substance on positive and negative replicates was evaluated by comparing target cycle threshold (Ct) values generated in the presence of the substance to Ct values from controls lacking the substance. The positive and negative replicate samples for 22 potentially interfering substances were correctly identified using the Xpert Carba-R Assay. Interference with the Xpert Carba-R Assay may be observed with barium sulfate at > 0.1% w/v and Pepto-Bismol at > 0.01% w/v in tests with rectal swab matrix samples. See Limitations in the package insert. Rectal swab matrix samples, positive for a mix of five carbapenemase-producing organisms harboring KPC, NDM, VIM, IMP-1 and OXA-48 gene sequences that were tested with fecal fat at 0.25% w/v, did not yield any false negative results, however, delayed cycle threshold values were observed for the VIM target. This potential interference from the presence of 0.25% w/v fecal fat is provided in the Limitations section of the package insert.

**Table 5-10. Potentially Interfering Substances Tested**

Substance/Class	Active Ingredient	Concentration Tested
Non-steroidal anti-inflammatory medication	Naproxen	0.25% w/v
Imaging compound	Barium sulfate	0.25% and 0.1% w/v
Antibiotic (oral)	Cephalexin	0.25% w/v
Antibiotic (oral)	Ciprofloxacin	0.25% w/v
Condom with spermicidal lubricant	Nonoxynol-9	1 condom <sup>a</sup>
Creams/ointment/suppositories	Hydrocortisone	0.25% w/v
Laxative	Sennosides	0.25% w/v
Lipids	Stearic acid/Palmitic acid/Cholesterol (fecal fat)	0.25% w/v
Anti-diarrheal medication	Loperamide hydrochloride/bismuth subsalicylate (Imodium)	0.25% w/v
Anti-diarrheal medication	Loperamide hydrochloride/bismuth subsalicylate (Kaopectate)	0.25% w/v
Topical cream	K-Y Jelly	0.25% w/v
Antacids	Calcium carbonate/aluminum hydroxide/magnesium hydroxide/simethicone (Milk of Magnesia)	0.25% w/v
Enemas	Mineral oil	0.25% w/v
Antibiotic (topical)	Polymixin B/ Neomycin/ Bacitracin (Neosporin)	0.25% w/v
Anti-fungal/anti-itch Vaginal	Nystatin	0.25% w/v
Antacid	Famotidine (Pepcid)	0.25% w/v

Substance/Class	Active Ingredient	Concentration Tested
Anti-diarrheal medication	Loperamide hydrochloride/bismuth subsalicylate (Pepto-Bismol)	0.25%, 0.1%, 0.05%, 0.025%, 0.01% w/v
Topical cream	Petroleum jelly	0.25% w/v
Anti-hemorrhoid creams/ointments	Phenylephrine (Preparation H)	0.25% w/v
Acid reducer; antacid	Omeprazole (Prilosec)	0.25% w/v
Enemas	Saline-enema	0.25% w/v
Antacid	Cimetidine (Tagamet)	0.25% w/v
Anti-fungal/anti-itch Vaginal	Benzocaine, resorcinol (Vagisil)	0.25% w/v
Moist towelettes	Benzalkonium chloride, ethanol (Wet Ones)	1 piece <sup>b</sup>

a. One condom added to 40 mL swab matrix.

b. One piece (5 inch x 7-1/2 inch) added to 40 mL swab matrix.

### Carry-Over Contamination

A study was conducted to demonstrate that single-use, self-contained GeneXpert cartridges prevent carry-over contamination in negative samples run following very high positive samples. The study consisted of a negative sample processed in the same GeneXpert module immediately following a very high positive sample. The high positive sample is composed of inactivated *E. coli* cells containing a plasmid with an insert consisting of a synthetic oligonucleotide of the amplicon sequences from the five Xpert Carba-R target analyte genes (KPC, NDM, VIM, IMP and OXA-48 targets). Positive cells were diluted in pooled negative rectal swab matrix to a concentration of  $1 \times 10^6$  CFU/mL. The testing scheme was repeated 25 times on two GeneXpert modules for a total of 102 tests (25 high positive samples per module and 26 negative samples per module) for the rectal swab matrix. All 50 positive samples correctly reported all Xpert Carba-R targets as **DETECTED**. All 52 negative samples correctly reported all Xpert Carba-R targets as **NOT DETECTED**.

### Clinical Studies

#### Clinical Performance

Performance characteristics of the Xpert Carba-R Assay with rectal swab specimens were determined in a multi-site investigational study. The positive percent agreement (PPA) and negative percent agreement (NPA) of the Xpert Carba-R Assay was evaluated relative to a reference method of culture (MacConkey enrichment broth) and PCR/bi-directional DNA sequence analysis.

Five geographically diverse sites (three across the United States and two in Europe) prospectively collected paired rectal swab specimens from subjects who were hospitalized or in a long-term care facility. Highly soiled rectal swab specimens, were excluded from the study. Due to low prevalence of each of the Xpert Carba-R Assay target genes in the absence of an outbreak, contrived specimens were also included in the study.

One rectal swab of the pair was used for Xpert Carba-R Assay testing. The second rectal swab was inoculated into MacConkey enrichment broth and used for reference method testing. A reference culture laboratory determined the presence of carbapenem non-susceptible organisms by culturing the MacConkey enrichment broth from each of the specimens. The MacConkey enrichment broth was screened for the presence of carbapenem-non-susceptible organisms initially by plating the broth on MacConkey agar plates with a meropenem disk. For specimens that exhibited growth of gram-negative bacteria around the meropenem disk, confirmation of carbapenem non-susceptibility was determined on isolated colonies by using the disk diffusion method (per CLSI document M02) as well as CLSI document M100. DNA extracted from the carbapenem non-susceptible isolates was purified, quantified, and amplified using primers specific to all 5 target genes; amplified regions included more bases than the regions amplified by the Xpert Carba-R Assay. The production of the appropriate size amplification product was confirmed on Agilent 2100 Bioanalyzer (Agilent Technologies, Santa Clara, CA).

If bands shown on the Bioanalyzer corresponded to the expected size of the amplicon from any of the five target genes detected by the Xpert Carba-R Assay, the amplicon for the isolate was sent to an independent laboratory for reference bi-directional sequencing analysis, which was validated for detection of the five targets in the Xpert Carba-R Assay. If no bands were shown on the Bioanalyzer for any of the five target genes, the isolate was not sent for sequence analysis and the reference method result was considered negative for the five target genes.

### **Prospective Specimen Results Obtained with the Xpert Carba-R Assay in Comparison to the Reference Method**

A total of 802 prospective rectal swab specimens were initially enrolled in this clinical study, of which 785 were eligible for inclusion. From the 785 eligible specimens, 755 specimens were included in the final dataset after exclusions based on protocol deviations (including 16 *Stenotrophomonas maltophilia* organisms that were excluded due to their intrinsic resistance to the carbapenems tested).

When tested with prospective rectal swab specimens, the Xpert Carba-R Assay demonstrated a PPA range from 60.0% to 100% for the four assay targets (*bla*<sub>KPC</sub>, *bla*<sub>NDM</sub>, *bla*<sub>VIM</sub>, and *bla*<sub>OXA-48</sub>) relative to the reference method (Table 5-11). The NPA for the *bla*<sub>KPC</sub>, *bla*<sub>NDM</sub>, *bla*<sub>VIM</sub>, *bla*<sub>OXA-48</sub>, and *bla*<sub>IMP</sub> gene sequences ranged from 98.6%-99.9% relative to the reference method (Table 5-11).

For specimens with discordant results (the Xpert Carba-R Assay was positive for a target gene but a carbapenem-non-susceptible organism was not isolated by reference culture), discordant analysis was performed using bi-directional sequencing on DNA extracted directly from the MacConkey enrichment broth. Discrepant testing results are footnoted in Table 5-11.

**Table 5-11. Xpert Carba-R Performance vs. Reference Culture + Sequencing – Prospective Specimens**

Specimen Type	Target	N	TP	FP	TN	FN	PPA% (95 CI)	NPA% (95 CI)
Prospective <sup>g</sup>	IMP	755	0	1 <sup>a</sup>	754	0	N/A	99.9% (99.3-100.0)
	VIM	755	6	8 <sup>b</sup>	737	4	60.0% (31.3-83.2)	98.9% (97.9-99.5)
	NDM	755	7	3 <sup>c</sup>	745	0	100.0% (64.6-100.0)	99.6% (98.8-99.9)
	KPC	755	29	6 <sup>d,e</sup>	720	0	100.0% (88.3-100.0)	99.2% (98.2-99.6)
	OXA-48	755	29	10 <sup>f</sup>	715	1	96.7% (83.3-99.4)	98.6% (97.5-99.2)

N = Number, TP = True Positive, FP = False Positive, TN = True Negative, FN = False Negative

a. 1 discordant specimen was confirmed as FP after discordant analysis.

b. 2 of the 8 FPs were determined to be TPs after discordant analysis.

c. 1 of the 3 FPs was determined to be TP after discordant analysis.

d. 1 of the 6 FPs was determined to be TP after discordant analysis.

e. Site reported that subject was on ertapenem during time of specimen collection.

f. 3 of the 10 FPs were determined to be TPs after discordant analysis.

g. Of the 755 prospective rectal swab specimens evaluated in the study, 636 specimens did not yield a culture isolate by Reference Culture. From the remaining 119 specimens, 112 carbapenem-non-susceptible organisms were recovered by Reference Culture in addition to 7 carbapenem susceptible organisms [*Pseudomonas aeruginosa* (5); *Escherichia coli* (1), and *Enterobacter cloacae* (1)].

Performance of the Xpert Carba-R Assay on the prospective specimens is shown in Table 5-12 by species. Only organisms for which at least one positive specimen was collected are included in Table 5-12.

**Table 5-12. Xpert Carba-R Performance vs. Reference Culture + Sequencing by Organism type – Prospective Specimens**

Species <sup>a</sup>	Target	N	TP	FP	TN	FN	PPA (95% CI)	NPA (95% CI)
<i>Enterobacter cloacae</i>	IMP	4	0	0	4	0	NA	100% (51.0-100.0)
	VIM	4	1	0	3	0	100% (20.7-100.0)	100% (43.9-100.0)
	NDM	4	0	0	4	0	NA	100% (51.0-100.0)
	KPC	4	0	0	4	0	NA	100% (51.0-100.0)
	OXA-48	4	1	0	3	0	100% (20.7-100.0)	100% (43.9-100.0)
<i>E. coli</i>	IMP	10	0	0	10	0	NA	100% (72.3-100.0)
	VIM	10	0	0	10	0	NA	100% (72.3-100.0)
	NDM	10	3	0	7	0	100% (43.9-100.0)	100% (64.6-100.0)
	KPC	10	2	0	8	0	100% (34.2-100.0)	100% (64.6-100.0)
	OXA-48	10	3	0	7	0	100% (43.9-100.0)	100% (64.6-100.0)

Species <sup>a</sup>	Target	N	TP	FP	TN	FN	PPA (95% CI)	NPA (95% CI)
<i>Klebsiella oxytoca</i>	IMP	1	0	0	1	0	NA	100% (20.7-100.0)
	VIM	1	0	0	1	0	NA	100% (20.7-100.0)
	NDM	1	0	0	1	0	NA	100% (20.7-100.0)
	KPC	1	0	0	1	0	NA	100% (20.7-100.0)
	OXA-48	1	1	0	0	0	100% (20.7-100.0)	NA
<i>Klebsiella pneumoniae</i>	IMP	60	0	1	59	0	NA	98.3% (94.0-100.0)
	VIM	60	0	1	59	0	NA	98.3% (94.0-100.0)
	NDM	60	4	1	55	0	100% (51.0-100.0)	98.2% (90.6-99.7)
	KPC	60	27	1	32	0	100% (87.5-100.0)	97.0% (89.6-100.0)
	OXA-48	60	24	3	32	1	96.0% (80.5-99.3)	91.4% (77.6-97.0)
<i>Pseudomonas aeruginosa</i>	IMP	30	0	0	30	0	NA	100% (88.7-100.0)
	VIM	30	5	0	21	4	55.6% (26.7-81.1)	100% (84.5-100.0)
	NDM	30	0	1	29	0	NA	96.7% (83.3-99.4)
	KPC	30	0	1	29	0	NA	96.7% (83.3-99.4)
	OXA-48	30	0	0	30	0	NA	100% (88.7-100.0)

a. *Acinetobacter baumannii* (13) and *Enterobacter amnigenus* (1) were recovered but did not contain target sequences by the Reference Method.

Multiple targets were detected by the Xpert Carba-R Assay in eight prospective specimens. The details are provided in Table 5-13, along with the discrepant sequencing result.

**Table 5-13. Prospective Specimens with Multiple Targets Detected**

Specimen	Targets Detected by Xpert Carba-R Assay	Targets Detected by Reference Sequencing	Discrepant Testing Results - Targets Detected by Reference Sequencing
1	KPC, OXA-48	NEG	NEG
2	VIM, KPC	NEG <sup>a</sup>	NEG <sup>a</sup>
3	VIM, OXA-48	OXA-48	OXA-48
4	KPC, OXA-48	KPC	KPC, OXA-48
5	NDM, OXA-48	NDM	NDM, OXA-48
6	VIM, NDM	NEG <sup>a</sup>	NEG
7	NDM, KPC	KPC	NDM, KPC
8	VIM, KPC	VIM	VIM, KPC

a. An organism was not isolated from reference culture, therefore, reference sequencing was not performed.

## Contrived Specimen Results Obtained with the Xpert Carba-R Assay in Comparison to the Reference Method

A total of 432 contrived specimens prepared in rectal swab matrix were also tested as part of the clinical study.

In addition to *Enterobacteriaceae*, *Pseudomonas aeruginosa*, and *Acinetobacter baumannii* groups tested in the contrived study, 5 other non-*Enterobacteriaceae* strains were also evaluated: *Pseudomonas stutzeri* (1), *Pseudomonas oryzihabitans* (1), *Pseudomonas putida* (2), and *Empedobacter brevis* (1).

When tested with contrived specimens, the Xpert Carba-R Assay demonstrated a PPA range from 95% to 100% across the assay targets (*bla*<sub>KPC</sub>, *bla*<sub>NDM</sub>, *bla*<sub>VIM</sub>, *bla*<sub>OXA-48</sub>, and *bla*<sub>IMP</sub>). The NPA for the *bla*<sub>KPC</sub>, *bla*<sub>NDM</sub>, *bla*<sub>VIM</sub>, *bla*<sub>OXA-48</sub>, and *bla*<sub>IMP</sub> gene sequences was 100% relative to the reference method (Table 5-14).

**Table 5-14. Xpert Carba-R Performance vs. Reference Method–  
Contrived Specimens**

Target	N	TP	FP	TN	FN	PPA% (95 CI)	NPA% (95 CI)
IMP	432	76	0	352	4	95.0% (87.8-98.0)	100.0% (98.9-100.0)
VIM	432	81	0	350	1	98.8% (93.4-99.8)	100.0% (98.9-100.0)
NDM	432	80	0	352	0	100.0% (95.4-100.0)	100.0% (98.9-100.0)
KPC	432	80	0	352	0	100.0% (95.4-100.0)	100.0% (98.9-100.0)
OXA-48	432	79	0	352	1	98.8% (93.3-99.8)	100.0% (98.9-100.0)

## Reproducibility Study

Reproducibility of the Xpert Carba-R Assay was evaluated using two panels of 11 samples, prepared in pooled negative rectal swab matrix. Two operators at each of the three study sites tested one panel of 11 samples in replicates of four per day over six testing days (11 samples x 2 replicates x 2 times/day x 6 days x 2 operators x 3 sites). Three lots of Xpert Carba-R Assay cartridges were used at each of the 3 testing sites. The Xpert Carba-R Assay was performed according to the Xpert Carba-R Assay procedure. Results are summarized in Table 5-15.

**Table 5-15. Summary of Reproducibility Results - % Agreement, Rectal Swab Matrix**

Sample	Matrix <sup>a</sup>	Site 1			Site 2			Site 3			% Total Agreement by Sample
		Op 1	Op 2	Site	Op 1	Op 2	Site	Op 1	Op 2	Site	
Neg	R	100% (24/24)	100% (24/24)	100% (48/48)	100% (24/24)	100% (24/24)	100% (48/48)	100% (24/24)	100% (24/24)	100% (48/48)	100% (144/144)
IMP Mod Pos	R	100% (24/24)	100% (24/24)	100% (48/48)	100% (24/24)	100% (24/24)	100% (48/48)	100% (24/24)	100% (24/24)	100% (48/48)	100% (144/144)
IMP Low Pos	R	91.7% (22/24)	87.5% (21/24)	89.5% (43/48)	83.3% (20/24)	87.5% (21/24)	85.4% (41/48)	87.5% (21/24)	79.2% (19/24)	83.3% (40/48)	86.1% (124/144)
VIM Mod Pos	R	100% (24/24)	100% (24/24)	100% (48/48)	100% (24/24)	100% (24/24)	100% (48/48)	100% (24/24)	100% (24/24)	100% (48/48)	100% (144/144)
VIM Low Pos	R	100% (24/24)	100% (24/24)	100% (48/48)	100% (24/24)	100% (24/24)	100% (48/48)	100% (24/24)	100% (24/24)	100% (48/48)	100% (144/144)
NDM Mod Pos	R	100% (24/24)	100% (24/24)	100% (48/48)	100% (24/24)	100% (24/24)	100% (48/48)	100% (24/24)	100% (24/24)	100% (48/48)	100% (144/144)
NDM Low Pos	R	91.7% (22/24)	95.8% (23/24)	93.8% (45/48)	95.8% (23/24)	95.8% (23/24)	95.8% (46/48)	100% (24/24)	91.7% (22/24)	95.8% (46/48)	95.1% (137/144)
KPC Mod Pos	R	100% (24/24)	100% (24/24)	100% (48/48)	100% (24/24)	100% (24/24)	100% (48/48)	100% (24/24)	100% (24/24)	100% (48/48)	100% (144/144)
KPC Low Pos	R	95.8% (23/24)	100% (24/24)	97.9% (47/48)	100% (24/24)	91.7% (22/24)	95.8% (46/48)	95.8% (23/24)	95.8% (23/24)	95.8% (46/48)	96.5% (139/144)
OXA-48 Mod Pos	R	100% (24/24)	100% (24/24)	100% (48/48)	100% (24/24)	100% (24/24)	100% (48/48)	100% (24/24)	100% (24/24)	100% (48/48)	100% (144/144)
OXA-48 Low Pos	R	95.8% (23/24)	100% (24/24)	97.9% (47/48)	95.8% (23/24)	100% (24/24)	97.9% (47/48)	91.7% (22/24)	100% (24/24)	95.8% (46/48)	97.2% (140/144)

a. R=rectal

The reproducibility of the Xpert Carba-R Assay was also evaluated in terms of the fluorescence signal expressed in Ct values for each target detected. The mean, standard deviation (SD), and coefficient of variation (CV) between-sites, between-lots, between-days, between-operators, and within-assays for each panel member are presented in Table 5-16.

**Table 5-16. Summary of Reproducibility Data, Rectal Swab Matrix**

Sample	Matrix <sup>a</sup>	Assay Channel (Analyte)	N <sup>b</sup>	Mean Ct	Between-Site		Between-Lot		Between-Day		Between-Operator		Within-Assay		Total	
					SD	CV (%)	SD	CV (%)	SD	CV (%)	SD	CV (%)	SD	CV (%)	SD	CV (%)
Neg	R	SPC	144	32.9	0.2	0.5	0.2	0.7	0.0	0.1	0.0	0	0.6	1.8	0.7	2.0
IMP Mod Pos	R	IMP	144	34.5	0.0	0.0	0.2	0.5	0	0.0	0.1	0.2	0.7	2.0	0.7	2.1
IMP Low Pos	R	IMP	140	36.4	0.0	0.0	0.0	0.0	0.2	0.5	0.0	0	1.2	3.3	1.2	3.4
VIM Mod Pos	R	VIM	144	31.0	0.0	0.0	0.3	0.9	0	0.0	0.2	0.5	0.5	1.6	0.6	1.9
VIM Low Pos	R	VIM	144	33.8	0.0	0.0	0.6	1.8	0.3	0.9	0.3	1.0	1.4	4.0	1.6	4.6
NDM Mod Pos	R	NDM	144	33.7	0.0	0.0	0.0	0.0	0.0	0.1	0.0	0.0	0.6	1.7	0.6	1.7
NDM Low Pos	R	NDM	143	36.2	0.2	0.7	0.0	0.0	0.3	0.7	0.0	0.0	0.8	2.3	0.9	2.5
KPC Mod Pos	R	KPC	144	34.2	0.0	0.0	0.3	0.8	0.2	0.6	0.0	0.0	0.4	1.2	0.6	1.6
KPC Low Pos	R	KPC	141	35.8	0.0	0.0	0.5	1.5	0.0	0.0	0.3	0.9	0.7	1.9	0.9	2.6
OXA-48 Mod Pos	R	OXA-48	144	34.3	0.0	0.0	0.2	0.5	0.2	0.5	0.1	0.3	0.5	1.6	0.6	1.7
OXA-48 Low Pos	R	OXA-48	143	36.1	0.0	0.0	0.0	0.0	0.2	0.6	0.0	0.0	0.8	2.3	0.9	2.4

a. R=rectal

b. Results with non-zero Ct values out of 144.

**Conclusions**

The results of the nonclinical analytical and clinical performance studies summarized above demonstrate that the Xpert Carba-R Assay is safe and effective for its intended use and is substantially equivalent to the predicate device.