

Public Health Service

Food and Drug Administration 10903 New Hampshire Avenue Document Control Center - WO66-G609 Silver Spring, MD 20993-0002

October 5, 2016

Beijing Adss Development Co., Ltd
% Ray Wang
General Manager
Beijing Believe Technology Service Co., Ltd.
5-1206, Build 332, Dafangju, No.25 Banbidian Rd.,
Liyuan Town, Tongzhou District, Beijing, 101121 CN

Re: K161925

Trade/Device Name: CO2 Laser Therapy Machine
Regulation Number: 21 CFR 878.4810
Regulation Name: Laser Surgical Instrument For Use In General And Plastic Surgery And In Dermatology
Regulatory Class: Class II
Product Code: GEX
Dated: July 11, 2016
Received: July 13, 2016

Dear Mr. Wang:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR

Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical devicerelated adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to

<u>http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm</u> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm.

Sincerely,

Jennifer R. Stevenson -A

For Binita S. Ashar, M.D., M.B.A., F.A.C.S. Director Division of Surgical Devices Office of Device Evaluation Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number *(if known)* k161925

Device Name CO2 Laser Therapy Machine

Indications for Use (Describe)

The CO2 Laser Therapy Machine is used for human tissue vaporization, coagulation in dermatology and plastic surgery, general surgery, gynecology, podiatry, dental and otorhinolaryngology.

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff *PRAStaff@fda.hhs.gov*

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

Tab #7 510(k) Summary

This 510(k) Summary of 510(k) safety and effectiveness information is being submitted in accordance with requirements of SMDA 1990 and Title 21, CFR Section 807.92.

The assigned 510(k) Number: _K161925_

1. Date of Preparation

10/03/2016

2. Sponsor

Beijing ADSS Development Co., Ltd

F6, Xiandao Bldg., Jinyuan Rd. 36, Daxing Economic Development zone, Beijing, China, 102628

Establishment Registration Number: Not yet registered or the Number Contact Person: Gao Yurong Position: sales manager Tel: 86-10-83625120 Fax:86-10-83625121 Email: gyr@adss.com.cn

3. Submission Correspondent

Mr. Ray Wang

Beijing Believe Technology Service Co., Ltd. 5-1206, Build 332, DaFangJu, No.25 BanBiDian Rd., LiYuan Town, TongZhou District, Beijing, 101121, China

Tel: +86-18910677558 Fax: +86-10-52214696 Email: ray.wang@believe-med.com

4. Identification of Proposed Device

Trade Name: CO2 Laser Therapy Machine Common Name: Powered Laser Surgical Instrument Model(s): FG 900/FG 900-B/FG 900-C

Regulatory Information: Classification Name: Powered Laser Surgical Instrument Classification: II; Product Code: GEX; Regulation Number: 21 CFR 878.4810; Review Panel: General& Plastic Surgery;

Intended Use:

The CO2 Laser Therapy Machine is used for human tissue vaporization, coagulation in dermatology and plastic surgery, general surgery, gynecology, podiatry, dental and otorhinolaryngology.

5. Device Description

The CO2 Laser Therapy Machine is a carbon dioxide laser used in medical and aesthetic industry for treatment of such skin conditions as fine and coarse wrinkles, scars of various origin, uneven pigmentation and dilated pores. Due to the CO_2 laser's high absorption of water, its high-energy beam of laser light interacts with the skin's surface causing the upper layer to peel off and use photothermolysis to stimulate deep cell regeneration and then achieve the target of skin improvement.

The proposed device is mainly used for human tissue vaporization, carbonization, coagulation and exposure to achieve the purpose of treatment.

The CO2 Laser Therapy Machine includes three models in this submission, FG 900, FG 900-B and FG 900-C, all three models have same principle, software, operation etc., only differences are appearance.

The proposed device includes the following components:

Components	Function Description	Applied Model(s)		
Surgery tip	Deliver the laser to area to be treated	FG 900/FG 900-B/FG		
		900-С		
Surgery tip Arm	Articulated arm for holding of Surgery tip	FG 900/FG 900-B/FG		
		900-С		
Touchscreen	The user interface and for controlling of the system	FG 900/FG 900-B/FG		
		900-С		
Emergency	Stop the system in case of emergency situation	FG 900/FG 900-B/FG		

Table 1 Main Components of Proposed Device

Switch		900-C
Key Switch	Start the system	FG 900/FG 900-B/FG
		900-C
Goggles for	Protect the eyes of patient	FG 900/FG 900-B/FG
Patient		900-C
Goggles for	Protect the eyes of Operator	FG 900/FG 900-B/FG
Doctor		900-C
Foot Switch	Activate the laser emission	FG 900/FG 900-B/FG
		900-C

6. Identification of Predicate Device

510(k) Number: K110434 Product Name: TRIXEL CO2 LASER Manufacturer: BEIJING SYNTECH LASER CO., LTD.

7. Non-Clinical Test Conclusion

Non clinical tests were conducted to verify that the proposed device met all design specifications as was Substantially Equivalent (SE) to the predicate device. The test results demonstrated that the proposed device complies with the following standards:

- IEC 60601-1:2012 Medical Electrical Equipment Part 1: General Requirements For Basic Safety And Essential Performance;
- IEC 60601-2-22:2007, Medical Electrical Equipment Part 2-22: Particular Requirements For Basic Safety And Essential Performance Of Surgical, Cosmetic, Therapeutic And Diagnostic Laser Equipment;
- ➢ IEC 60825-1: 2007, Safety of laser products Part 1: Equipment classification and requirements.
- IEC 60601-1-2:2007, Medical electrical equipment- Part 1-2: General requirements for basic safety and essential performance- Collateral standard: Electromagnetic compatibility-Requirements and tests.
- ISO 10993-5:2009, Biological Evaluation of Medical Device, Part 5-Tests for Vitro cytotoxicity
- ISO 10993-10:2002/Amd. 1: 2006, Biological Evaluation of Medical Device, Part 10-Test for irritation and delay-type hypersensitivity AMENDMENT 1
- > Performance Testing for Spot Size Accuracy and Energy Output Accuracy.
- Software Validation & Verification Test
- 8. Clinical Test Conclusion

No clinical study is included in this submission.

9. Substantially Equivalent (SE) Comparison

Proposed Device	Predicate Device	Remark
GEX	GEX	SE
21 CFR 878.4810	21 CFR 878.4810	SE
2	2	SE
hospital	hospital	SE
The equipment is used for human tissue	The equipment is used for human tissue	SE
vaporization, coagulation in dermatology and plastic surgery, general surgery, gynecology, podiatry, dental and otorhinolaryngology.	vaporization, coagulation in dermatology and plastic surgery, general surgery, gynecology, podiatry, dental and otorhinolaryngology.	
	Proposed Device GEX 21 CFR 878.4810 2 hospital The equipment is used for human tissue vaporization, coagulation in dermatology and plastic surgery, general surgery, gynecology, podiatry, dental and otorhinolaryngology.	Proposed DevicePredicate DeviceGEXGEX21 CFR 878.481021 CFR 878.4810221Predicate Device2111

Table 2 General Comparison

ITEM	Proposed Device			Predicate Device				Remark	
Maximum Power	30W		CFL-10: 12W (+20%)				SE		
			UFL-60: 30W (±20%)						
work mode	Surgery	Surgery (Single Pulse, Continuous,		Scanner (half, fast and random)				SE	
		Muti-	Pulse)	Surgery (CW, repeat and pulse)				
Wavelength		10.6	5 um		10.6 um				SE
Mode Structure		TE	M00		TEM00				SE
Beam delivery	7 knuc	klearmke	y join	ts light arm	7 knucklearmkey joints light arm				SE
Light arm		1.3	2 m		CFL-10: 0.97m			Analysis	
				UFL-60: 1.17m					
Aiming Beam	650nn	n red diod	e lase	r (0.5 mW)	650nm red diode laser(<1mW)			SE	
Spot size		0.5	mm		0.5mm (±10%)			SE	
Pulse Setting	Single Pu	ılse	10-1	1000ms	Pulse		1-999	Oms	SE
	Muti-	Time (On	10-1000ms	Repeat	Tim	e On	1-999 ms	SE
	Pulse	Time (Off	10-1000ms		Time	e Off	1-999 ms	SE
	Contin	nuous		0-30W	CW		CFL-	10: 0.1-12W	SE
					UFL-60:		JFL-60:		
					0.1-30W				
Power calibration	Period of 1 year		Period of 1 year			SE			
Control System	Touch screen, footswitch		Touch screen, footswitch			SE			
Laser operation		Footswitch		Footswitch				SE	
Laser	CO2		CO2			SE			
medium/energy									
source									
Cooling System	Air cooling		Air cooling			SE			
Clean Method	,	70% medical alcohol		70% medical alcohol			SE		
Patient		SI	cin		Skin			SE	
Contacted Part			1						
Dimension	FG 900		56*	46*112 cm	CFL-10 Trixel CO2 Laser: 210 x 600x		: 210 x 600x	SE	
	FG 900-E	3	60*	54*32cm	330 (without light arm)				
	FG 900-C	2	46*	42*125cm	UFL-60 Trixel II CO2 Laser: 1300 x				
***	FC 000			40.1	550x 420(without light arm)			0E	
Weight	FG 900			49 kg	CFL-10 Trixel CO2 Laser: 20kg		SE		
	FG 900-E	3		28kg	UFL-60 Trixel II CO2 Laser: 40kg				
Demos transf	FG 900-C			43Kg			CT.		
Power input	AC 110V	/30HZ-60	HZ;		CFL-10 Trixel CO2 Laser: 120 V		2F		
					AC/00HZ	ivol II	CO^{2}	Laser, 120	
			UTL-00 INXELII CO2 Laser: 120 VAC/60Hz						
			VAC/00HZ						

Table 3 Performance Comparison

Item	Proposed Device Predicate Device		Remark		
EMC, Electrical and Laser Safety					
Electrical Safety	Comply with IEC 60601-1, IEC Comply with IEC 60601-1, IEC		SE		
	60601-2-22	60601-2-22			
EMC	Comply with IEC 60601-1-2	Comply with IEC 60601-1-2	SE		
Laser Safety	Comply with IEC 60601-2-22, IEC 60825	Comply with IEC 60601-2-22, IEC 60825	SE		

Table 4 Safety Comparison

Analysis

The proposed device is substantially equivalent to the predicate device. Based on the nonclinical tests performed, the subject device is as safe, as effective, and performs as well as the legally marketed predicate device.

10. Substantially Equivalent (SE) Conclusion

Based on the comparison and analysis above, the proposed device is determined to be Substantially Equivalent (SE) to the predicate device.