May 10, 2002

K020486

510(K) SUMMARY: CARESIDE® GGT SAFETY AND **EFFECTIVENESS**

T. **Applicant Information**

Applicant Name CARESIDE, Inc. Applicant/Manufacturer Address 6100 Bristol Parkway Culver City, CA 90230

310-338-6767 C. Telephone Number

D. Contact Person Renate A. MacLaren, Ph.D.

E. FAX Number 310-670-6986

e-Mail Address rmaclaren@CARESIDE.com F.

Date 510(k) Summary prepared May 10, 2002

II. **Device Information**

CARESIDE® GGT A. Device Name (Trade) В. Device Name (Classification) GGT test system C. Device Classification

Clinical chemistry panel

GGT test system

Regulation Number: 21 CFR 862.1360 Regulatory Class 1 (non-exempt for point of

care use)

Classification Number: 75JPZ

D. Special controls and performance standards None applicable

III. Substantial Equivalence Claim

A. General equivalency claim

The ability to monitor analyte-specific biochemical reactions in dry film and other formats is widely recognized and has gained widespread acceptance for use in chemistry assays.

GGT in vitro diagnostic products, in both dry film and other formats, are currently marketed in the U.S., including GGT products which utilize GGT catalyzed generation of p-nitroaniline from the artificial substrate L-y-glutamyl-p-nitroanilide.

В. Specific equivalency claim

The CARESIDE® GGT test is substantially equivalent in principle, intended use, and clinical performance to the currently marketed Vitros slides for the quantitative measurement of GGT on the Vitros DTSC 60 II.

Name of Predicate Device: Johnson and Johnson's (formerly Eastman Kodak, Inc.)

Vitros GGT Slides for Johnson and Johnson's Vitros

DTSC 60 (formerly Eastman Kodak's DTSC 60 II).

Predicate Device 510K number: K912844/A

Product Code: 75JPZ

IV. Device Description

CARESIDE® GGT cartridges are used with the CARESIDE Analyzer® to measure GGT activity in anti-coagulated whole blood, plasma or serum specimens. The CARESIDE GGT cartridge, a single use disposable *in vitro* diagnostic test cartridge, aids in specimen separation and delivers a measured volume of plasma or serum to a dry film to initiate the measurement of GGT activity. The patented film cartridge contains all reagents necessary to measure GGT activity.

A. Explanation of Device Function

Each CARESIDE[®] GGT cartridge consists of a GGT-specific multi-layer reagent film mounted in a plastic base with a hinged lid. The user introduces the anti-coagulated whole blood, serum, or plasma specimen into the cartridge sample well, closes the lid and inserts the cartridge into the CARESIDE Analyzer.

Once loaded, the CARESIDE Analyzer scans the cartridge barcode, brings the cartridge and the contained specimen to 37°C, and spins the cartridge to move the sample from the sample deposition well into the cartridge channels and chambers. As the cartridge continues to spin, the blood cells are separated from the plasma/serum and the cells accumulate in the separation well. Approximately 8.5 microliters of plasma (or serum, as applicable) remain in the metering passage. Any excess sample flows into an overflow well.

The plasma (or serum, as applicable) is automatically dispensed onto the multi-layer reagent film. The spreading and substrate layer distributes the GGT containing specimen uniformly. The GGT in the specimen reacts with the substrate L- γ -glutamyl-p-nitroanilide to release p-nitroaniline resulting in a change in film color. The rate of change of color intensity, as measured by the amount of reflected light at 425 nanometers, directly relates to the specimen GGT activity.

Test Reaction Sequence:

L-γ-glutamyl-p-nitroanilide + glycylglycine

GGT ↓

L-glutamylglycylglycine + p-nitroaniline

As the cartridges spin, photodiodes measure reflectance of light emitted by wavelength-specific light emitting diodes (LEDs) over a fixed time period. The instrument uses the reflectance measurements and the lot-specific standard curve to calculate GGT activity.

B. Test Summary

γ-Glutamyltransferase measurements are used in the diagnosis and treatment of liver diseases including alcoholic cirrhosis and primary and secondary liver tumors. GGT in the serum originates primarily from the hepatobiliary system even though renal tissue contains the highest levels of GGT. GGT is elevated in all forms of liver disease. GGT elevations are observed earlier and are more pronounced than those of other liver enzymes in cases of obstructive jaundice and metastatic neoplasms. In intra- or posthepatic biliary obstruction, GGT levels may reach levels 5 to 30 times higher than in intra- or post-hepatic biliary obstruction

V. Intended Use

A. Intended Use

The CARESIDE $^{\otimes}$ GGT cartridge is intended for *in vitro* diagnostic use in conjunction with the CARESIDE $^{\otimes}$ Analyzer to quantitatively measure GGT activity in anti-coagulated whole blood, plasma or serum.

B. Indications for Use

For in vitro diagnostic use with the CARESIDE *Analyzer* to quantitatively measure GGT from anti-coagulated whole blood, plasma, or serum specimens to aid in the diagnosis and treatment of liver diseases such as alcoholic cirrhosis and primary and secondary liver tumors.

VI. Technological Characteristics

A. <u>Similarities</u>

	CARESIDE® GGT	Vitros GGT DT Slides
Intended Use	For in vitro Diagnostic use	Same
Indications	Primarily to aid in the	Same
	diagnosis and treatment of	
	patients with liver diseases	
	such as alcoholic cirrhosis and	
	primary and secondary liver	
	tumors.	
Measurement	Quantitative	Same
Method	Dry film based quantitation of	Same
Principle	enzymatic activity by	
	reflectance photometry using	
	L-g-glutamyl transferase and	
	glycylglycine substrate.	
Specimen	Not required	Same
Dilution		
Materials	L-γ-glutamyl-p-nitroanilide	L-γ-glutamyl-p-nitroanilide
	and glycylglycine	and glycylglycine
Detector	Reflectance (425 nm)	Reflectance (400 nm)
Test Time	Approx. 4 min. warm-up (on-	15 minutes slide warm-up (off-
	board) plus 5 minute test time.	line) plus 5 minutes test time.
Sample Type	Anti-coagulated whole blood,	Serum, plasma
	serum or plasma [whole blood	
	applied sample, plasma test	
	sample]	
Specimen	8.5 µl test volume	10 μl
Volume	$(90 \pm 10 \mu l \text{ applied volume})$	
Calibration	Calibration information bar-	Run Vitros DTSC II
	coded on each cartridge.	calibrators whenever a new
	Calibration information may	slide lot is used or when
	change with each lot.	necessary.
Quality Control	2 levels	Same
Reporting Units	U/L	Same
Reaction Temp.	37 °C	Same

B. <u>Differences</u>

	CARESIDE® GGT	Vitros GGT DT Slides
Direct Blood Specimen	Yes, whole blood	No, requires separation of whole blood prior to sample application
Reportable range	20 to 1000 U/L	5 to 1400 U/L
Accurate Pipetting	Not required	Required
Reagent Pre- warming	Not required	Required

C. <u>Comparative Performance Characteristics</u>

	CARESIDE® GGT	Vitros GGT DT Slides	
Detection limit	20 U/L	5 U/L	
Reportable	20 to 1000 U/L	5 to 1400 U/L	
Range			
Accuracy	Mean recovery 101%	Not available	
Precision	Total CV, 344 U/L, 2.6%	Total CV, 166 U/L, 2.4%	
Method	CARESIDE = 1.00 (BM/Hitachi 902 GGT) + 0.91 U/L,		
Comparison	r= 1.00		
Linearity	Linearity by dilution	Not available	
	yielded slope and		
	correlation coefficient		
	within acceptable limits.		
Interference	No significant interference	None stated	
	observed at tested concentration of		
	interferent:		
	Ascorbic Acid 10 mg/dL		
	Bilirubin 10 mg/dL		
	Triglycerides 3000 mg/dL		
Specimen Types	No clinically significant	No clinically significant	
&	difference between sodium	difference between serum,	
Anticoagulants	heparinized whole blood,	heparin plasma, or EDTA	
1 Mary Company of the	serum, and sodium heparin	plasma. Whole blood is	
	plasma.	unsuitable.	

D. <u>Conclusion</u>

The nonclinical and clinical data provided demonstrate that the CARESIDE® GGT product is as safe, effective, and performs as well as or better than the legally marketed predicate device



Food and Drug Administration 2098 Gaither Road Rockville MD 20850

Renate A. MacLauren, Ph.D. Clinical Affairs Manager Careside, Inc. 6100 Bristol Parkway Culver City, CA 90230

JUL 1 2002

Re: k020486

Trade/Device Name: Careside® GGT Regulation Number: 21 CFR 862.1360

Regulation Name: Gamma-glutamyl transpeptidase and isoenzymes test system

Regulatory Class: Class I, reserved

Product Code: JPZ Dated: May10, 2002 Received: May13, 2002

Dear Dr. MacLauren:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4588. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "http://www.fda.gov/cdrh/dsma/dsmamain.html".

Sincerely yours,

Steven I. Gutman, M.D., M.B.A.

Director

Division of Clinical Laboratory-Devices

Steven Butman

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

INDICATIONS FOR USE

K020486

510(k) Number:

Device Name:

- CARESIDE® GGT

Indications for use: For in vitro diagnostic use with the CARESIDE *Analyzer* to quantitatively measure GGT from anti-coagulated whole blood, plasma, or serum specimens to aid in the diagnosis and treatment of liver diseases such as alcoholic cirrhosis and primary and secondary liver tumors.

(Division Sign-Off)

Division of Clinical Laboratory

5100k) Number

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)