

NOV 19 2002

K022687

510(k) SUMMARY

In response to the requirements addressed by the Safe Medical Devices Act (SMDA) of 1990, a summary follows with the safety and effectiveness information upon which the substantial equivalence determination is based.

**510(k) SUMMARY
FOR
AOSEPT Clear Care Cleaning and Disinfecting Solution**

1. **Submitter Information**
CIBA Vision Corporation
11460 Johns Creek Parkway
Duluth, Georgia 30097
Contact Person: Steven Dowdley
Telephone No: 678-415-3897
2. **Device Name**
Classification Name: Soft (hydrophilic) Contact Lens Solution
Proprietary Name: AOSEPT Clear Care Cleaning and Disinfecting Solution
3. **Predicate Devices**
Allergan Ultracare Disinfecting Solution/Neutralizer
4. **Description of the Devices**
The AOSept Clear Care Cleaning and Disinfecting Solution is an aqueous solution contains hydrogen peroxide 3% (stabilized with phosphonic acid), sodium chloride, a phosphate buffer system and a non-ionic surfactant.
5. **Indications for Use**
AOSeptClear Care Cleaning and Disinfecting Solution is indicated for use in simultaneous cleaning, daily protein removal, disinfecting, and storing of soft (hydrophilic) contact lenses as recommended by your eye care practitioner.
6. **Description of Safety and Substantial Equivalence**
A series of preclinical and clinical studies were completed on this product and were previously submitted under submission K003345 and K013512.
7. **Substantial Equivalence**
AOSept Clear Care Cleaning and Disinfection Solution is substantially equivalent to currently marketed AOSept Clear Care Cleaning and Disinfection Solution in cleaning, disinfecting, daily protein removal and storing of soft (hydrophilic) contact lenses as recommended by your eye care practitioner.



Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

NOV 19 2002

CIBA Vision Corp.
c/o Steven Dowdley, RAC
11460 Johns Creek Parkway
Duluth, GA 30097-1556

Re: K022687
Trade/Device Name: AOSEPT Clear Care Cleaning and Disinfecting Solution
Regulation Number: 21 CFR 886.5928
Regulation Name: Soft (hydrophilic) Contact Lens Care Products
Regulatory Class: Class II
Product Code: LYL
Dated: October 4, 2002
Received: October 7, 2002

Dear Mr. Dowdley:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 21 CFR Part 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4613. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address.
<http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,



A. Ralph Rosenthal, M.D.
Director
Division of Ophthalmic and Ear,
Nose and Throat Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

PART III. INDICATIONS FOR USE STATEMENT

510(k) Number: *This is a new 510 (k) Notification. (Number to be assigned)*

K022687

Device Name: AOSEPT Clear Care Cleaning and Disinfecting Solution

Indications for Use:

AOSEPT Clear Care Cleaning and Disinfecting Solution is indicated for use in simultaneous cleaning, daily protein removal, disinfecting, and storing of soft (hydrophilic) contact lenses as recommended by your eye care practitioner.

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use: or over-the-counter:

JS

Myra Smith

Division Sign-Off
Division of Ophthalmic Ear,
Nose and Throat Devices

510(k) Number *K022687*