



**Wiener lab.**

Especialidades para Laboratorios Clínicos

WIENER LABORATORIOS S.A.I.C. - Riobamba 2944 - 2000 Rosario - Argentina  
Phone +54 (341) 432-9191/6 - Fax +54 (341) 432-5454/5555  
Internet: <http://www.wiener-lab.com.ar>

## Section 6 – Summary

DEC 09 2002

### 510(k) Summary

**“This summary of 510(k) safety and effectiveness information is being submitted in accordance with the requirements of SMDA 1990 and 21CFR 807.92”**

**“The assigned 510(k) number is: K023551”**

#### **Introduction**

According to the requirements of 21 CFR 862.1100, the following information provides sufficient details to understand the basis of a determination of substantial equivalence.

---

#### **6-1 Submitter** **Name, Address,** **Contact**

Wiener Laboratorios S.A.I.C.  
Riobamba 2944  
2000 – Rosario – Argentina  
Tel: 54 341 4329191  
Fax: 54 341 4851986  
Contact person: Viviana Cétola  
Date Prepared: September 02, 2002

---

#### **6-2 Device Name**

Proprietary name: Wiener lab. GOT (AST) UV AA Líquida.  
Common name: Aspartate aminotransferase (AST/SGOT) test system.  
Classification name: NADH Oxidation / NAD Reduction, AST / SGOT.  
Device Class II

---



**6-3 Predicate Device**

We claim substantial equivalence to the currently marketed Wiener lab. GOT (AST) UV (Cat. N°1751302).

**6-4 Device Description**

Kinetic Method.

The principle is based on the following reaction system:



The rate of disappearance of NADH and the resulting decrease in absorbance at 340 nm is directly proportional to the activity of GOT (AST).

AST or GOT: Aspartate Aminotransferase.

MDH: Malate Dehydrogenase.

NADH: Nicotinamide-Adenine Dinucleotide (Reduced)

NAD<sup>+</sup>: Nicotinamide-Adenine Dinucleotide.

H<sup>+</sup>: Proton.

**6-5 Intended Use**

The WIENER LAB. GOT (AST) UV AA Líquida test system is a quantitative in vitro diagnostic device intended to be used in the quantitative determination of aspartate amino transferase (AST or GOT) in human serum and plasma on both manual and automated systems. Aspartate amino transferase measurements are used in the diagnosis and treatment of certain types of liver and heart diseases.

**6-6 Equivalencies and Differences**

The WIENER LAB. GOT (AST) UV AA Líquida test system is substantially equivalent to other products in commercial distribution intended for similar use. Most notably it is substantially equivalent to the currently WIENER LAB. GOT (AST) UV test system.

The following table illustrates the similarities and differences between the WIENER LAB. GOT (AST) UV AA Líquida test system and the currently marketed WIENER LAB. GOT (AST) UV test system.

	GOT (AST) UV	GOT (AST) UV AA Líquida
Intended Use	Quantitative determination of Aspartate Aminotransferase in human serum and plasma.	
Test Principle	<p>Kinetic Method.</p> <p>The principle is based on the following reaction system:</p> <p style="text-align: center;">L-aspartate + 2-oxoglutarate</p> <p style="text-align: center;">GOT (AST)</p> <p style="text-align: center;">Oxaloacetate + L-glutamate</p> <p style="text-align: center;">Oxalacetate + NADH + H<sup>+</sup></p> <p style="text-align: center;">MDH</p> <p style="text-align: center;">L-Malate + NAD<sup>+</sup></p> <p>The rate of disappearance of NADH and the resulting decrease in absorbance at 340 nm is directly proportional to the activity of GOT (AST).</p> <p>AST or GOT: Aspartate Aminotransferase, MDH: Malate Dehydrogenase, NADH: Nicotinamide-Adenine Dinucleotide (Reduced), NAD<sup>+</sup>: Nicotinamide-Adenine Dinucleotide, H<sup>+</sup>: Proton.</p>	
Reagents	Buffer: L-aspartate – TRIS (buffer) Substrate: NADH – MDH – LDH – 2-oxoglutarate	
Preparation of Working Reagent	Dissolution of Reagent with stated volume of Buffer.	Reagents may be used separately or as <b>Monoreagent</b> , mixing 4 parts Buffer and 1 part Substrate.
Wavelength of Reading	334 – 340 – 366 nm	340 nm
Linearity	470 U/l	700 U/l
Expected values	Male: up to 38 U/l (37°C) Female: up to 32 U/l (37°C)	
<i>Continued on next page</i>		

	GOT (AST) UV	GOT (AST) UV AA Líquida
Within-run precision	Normal Level Serum: CV = 4.4% High Level Serum: CV = 1.3%	Normal Level Serum: CV = 2.41% High Level Serum: CV = 1.22%
Run-to-run precision	Normal Level Serum: CV = 4.9% High Level Serum: CV = 1.6%	Normal Level Serum: CV = 2.26% High Level Serum: CV = 2.16%

---

**6-7 Conclusion** Above mentioned data show substantial equivalency to the predicate device.

---



DEPARTMENT OF HEALTH & HUMAN SERVICES

Food and Drug Administration  
2098 Gaither Road  
Rockville MD 20850

Dr. Viviana Cetola  
QC/QA Manager  
Wiener Laboratorios S.A.I.C.  
Riobamba 2944  
2000 – Rosario - Argentina

DEC 09 2002

Re: k023551  
Trade/Device Name: Wiener Lab. GOT (AST) UV AA Liquida  
Regulation Number: 21 CFR § 862.1100  
Regulation Name: NADH Oxidation/NAD Reduction, AST/SGOT  
Regulatory Class: II  
Product Code: CIT  
Dated: September 13, 2002  
Received: October 22, 2002

Dear Dr. Cetola:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in Title 21, Code of Federal Regulations (CFR), Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); and good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820). This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

Page 2

If you desire specific information about the application of labeling requirements to your device, or questions on the promotion and advertising of your device, please contact the Office of In Vitro Diagnostic Device Evaluation and Safety at (301) 594-3084. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>.

Sincerely yours,

A handwritten signature in cursive script that reads "Steven Gutman".

Steven I. Gutman, M.D., M.B.A.

Director

Office of In Vitro Diagnostic Device Evaluation and Safety  
Center for Devices and Radiological Health

Enclosure

K023551

Page \_\_\_ of \_\_\_

510(k) Number (if known): K023551

Device Name: Wiener lab.

GOT (AST) UV AA LIQUIDA

**Indications For Use:**

The "Wiener lab. GOT (AST) UV AA Líquida" test system is a quantitative in vitro diagnostic device intended to be used in the quantitative determination of aspartate amino transferase (AST or GOT) in human serum and plasma on both manual and automated systems. Aspartate amino transferase measurements are used in the diagnosis and treatment of certain types of liver and heart diseases.

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use   
(Per 21 CFR 801.109)

OR

Over-The-Counter Use

(Optional Format 1-2-96)

[Signature]  
(Division Sign-Off)  
Division of Clinical Laboratory Devices  
510(k) Number K023551

SK33  
CH  
11