

1023839

JAN 14 2003

Section 3
HemosIL Factor V Deficient Plasma - 510(k) Summary
(Summary of Safety and Effectiveness)

Submitted by:

Instrumentation Laboratory Company
113 Hartwell Avenue
Lexington, MA 02421
Phone: 781-861-4467
Fax: 781-861-4207

Contact Person:

Carol Marble, Regulatory Affairs Manager
Phone: 781-861-4467 / Fax: 781-861-4207

Summary Prepared:

November 15, 2002

Name of the Device:

HemosIL Factor V Deficient Plasma

Classification Name(s):

864.7290	Factor Deficiency Tests	Class II
81GJT	Plasma, Coagulation Factor Deficient	

Identification of Predicate Device(s):

K893533 Hemoliance Factor V Deficient Plasma on ELECTRA Series Analyzers
K002400 IL Test Factor V Deficient Plasma* on ACL Family of Analyzers

*NOTE: Reagent was 510(k) cleared as part of multiple analyzer systems, most recently the ACL Advance.

Description of the Device/Intended use(s):

HemosIL Factor V Deficient Plasma is human plasma immunodepleted of factor V and intended for the *in vitro* diagnostic quantitative determination of factor V activity in citrated plasma, based on the prothrombin time (PT) assay, on IL Coagulation and ELECTRA Systems.

Abnormalities of the extrinsic pathway factors are determined by performing a modified prothrombin time (PT) test. Patient plasma is diluted and added to a plasma deficient in factor V. Correction of the clotting time of the deficient plasma is proportional to the concentration (% activity) of the factor V in the patient plasma, interpolated from a calibration curve.

Statement of Technological Characteristics of the Device Compared to Predicate Device:

HemosIL Factor V Deficient Plasma is substantially equivalent to Hemoliance Factor V Deficient Plasma (on ELECTRA Series Analyzers) and IL Test Factor V Deficient Plasma (on ACL Family of Analyzers) in performance, intended use and safety and effectiveness.

Section 3
HemosIL Factor V Deficient Plasma - 510(k) Summary
(Summary of Safety and Effectiveness)

Summary of Performance Data:

Method Comparison

In method comparison studies evaluating 60 citrated plasma samples (30 normal/30 abnormal), the slopes and correlation coefficients (r) for HemosIL Factor V Deficient Plasma versus the predicate devices are shown below:

NOTE: HemosIL RecombiPlasTin (K012768) was used as the PT reagent in all testing.

**HemosIL Factor V Deficient Plasma vs.
Predicate Hemoliance Factor V Deficient Plasma on ELECTRA**

IL System	Slope	r
E1400C	1.0161	0.9877

**HemosIL Factor V Deficient Plasma vs.
Predicate IL Test Factor V Deficient Plasma on ACL Family**

IL System	Slope	r
ACL 300	1.0005	0.9932
ACL 6000	0.9771	0.9940
ACL 9000	0.9746	0.9945
ACL Futura	1.1331	0.9814

Within Run Precision

Within run and total precision assessed over multiple runs (n=80) using two levels of control gave the following results:

Instrument	Control	Mean % Factor V	Within run CV%	Total CV%
ACL 300	Normal Control	110.6	1.0	3.1
	Low Abnormal Control	30.1	1.8	3.5
ACL 6000	Normal Control	115.1	1.5	3.0
	Low Abnormal Control	31.5	1.9	2.8
ACL 9000	Normal Control	113.7	1.2	1.7
	Low Abnormal Control	29.0	4.0	5.6
ACL Advance	Normal Control	132.4	6.0	7.6
	Low Abnormal Control	27.1	4.1	5.7
ELECTRA 1400C	Normal Control	101.0	2.1	3.5
	Low Abnormal Control	24.0	1.5	2.6



DEPARTMENT OF HEALTH & HUMAN SERVICES

Food and Drug Administration
2098 Gaither Road
Rockville MD 20850

JAN 14 2003

Ms. Carol Marble
IL Regulatory Affairs Director
Instrumentation Laboratory
113 Hartwell Avenue
Lexington, Massachusetts 02421

Re: k023839
Trade/Device Name: HemosIL Factor V Deficient Plasma
Regulation Number: 21 CFR § 864.7290
Regulation Name: Factor deficiency tests
Regulatory Class: II
Product Code: GJT
Dated: November 15, 2002
Received: November 18, 2002

Dear Ms. Marble:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in Title 21, Code of Federal Regulations (CFR), Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Page 2

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); and good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820). This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific information about the application of labeling requirements to your device, or questions on the promotion and advertising of your device, please contact the Office of In Vitro Diagnostic Device Evaluation and Safety at (301) 594-3084. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>.

Sincerely yours,


Steven I. Gutman, M.D., M.B.A.

Director
Office of In Vitro Diagnostic Device
Evaluation and Safety
Center for Devices and Radiological Health

Enclosure

Indications for Use Statement

510(k) Number (if known): K023839


Device Name: HemosIL Factor V Deficient Plasma

Indications for Use:

HemosIL Factor V Deficient Plasma is human plasma immunodepleted of factor V and intended for the *in vitro* diagnostic quantitative determination of factor V activity in citrated plasma, based on the prothrombin time (PT) assay, on IL Coagulation and ELECTRA Systems.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)



(Division Sign-Off)
Division of Clinical Laboratory Devices
510(k) Number K023839

Prescription Use

OR

Over-The-Counter Use