

AUG 27 2003



**Bio-Rad  
Laboratories**

Diagnostics Group  
4000 Alfred Nobel Drive  
Hercules, California 94547  
Telephone: 510 724-7000  
Fax: 510 741-5824

---

## 510(k) SUMMARY

This summary of 510(k) safety and effectiveness information is being submitted in accordance with the requirements of SMDA 1990 and 21 CFR 807.92.

The assigned 510(k) number is: K031043.

**Submitter:** Bio-Rad Laboratories, Inc.  
Clinical Diagnostics Group  
4000 Alfred Nobel Drive,  
Hercules, California 94547  
Phone: (510) 741-5309  
FAX: (510) 741-6471

**Contact Person:** Jackie Buckley  
Regulatory Affairs Representative

**Date of Summary Preparation:** March 31, 2003

**Device Name:** D-10™ Hemoglobin A<sub>1c</sub> Program

**Classification Name:** Assay, Glycosylated Hemoglobin, 81LCP

**Predicate Device:** VARIANT™ II Hemoglobin A<sub>1c</sub> Program  
K984268  
Bio-Rad Laboratories, Inc.

**Statement of Intended Use:** The Bio-Rad D-10 Hemoglobin A<sub>1c</sub> Program is intended for the percent determination of hemoglobin A<sub>1c</sub> in human whole blood using ion-exchange high performance liquid chromatography (HPLC). The D-10 Hemoglobin A<sub>1c</sub> is intended for use only with the Bio-Rad D-10 Hemoglobin Testing System.

For In Vitro Diagnostic Use.

**Description of Device:**

The D-10 Hemoglobin Testing System uses the principles of high performance liquid chromatography (HPLC). The D-10 Hemoglobin A<sub>1c</sub> Program is based on chromatographic separation of Hemoglobin A<sub>1c</sub> on a cation exchange cartridge.

**Technical Characteristics Compared to Predicate:**

D-10 and VARIANT II systems have the same technical characteristics that are summarized in the table below:

<b>Characteristics</b>	<b>D-10 Hemoglobin A<sub>1c</sub></b>	<b>VARIANT II Hemoglobin A<sub>1c</sub></b>
Analyte Measured: Reported	%A <sub>1c</sub>	%A <sub>1c</sub>
Intended Use	The Bio-Rad D10 Hemoglobin A <sub>1c</sub> Program is intended for the percent determination of hemoglobin A <sub>1c</sub> in human whole blood using ion-exchange high performance liquid chromatography (HPLC).	The Bio-Rad VARIANT II Hemoglobin A <sub>1c</sub> Program is intended for the percent determination of hemoglobin A <sub>1c</sub> in human whole blood using ion-exchange high performance liquid chromatography (HPLC).
Assay Principle	Cation exchange high performance liquid chromatography	Cation exchange high performance liquid chromatography
Sample Type	Human anticoagulated whole blood (EDTA)	Human anticoagulated whole blood (EDTA)
Visible Detection	415 nm	415 nm
Standardization	Traceable to the Diabetes Control and Complications Trial (DCCT) reference method and IFCC. Certified via the National Glycohemoglobin Standardization Program (NGSP).	Traceable to the Diabetes Control and Complications Trial (DCCT) reference method and IFCC. Certified via the National Glycohemoglobin Standardization Program (NGSP).

**Testing To Establish Substantial Equivalence:**

Method correlation between D-10 Hemoglobin A<sub>1c</sub> Program and VARIANT II Hemoglobin A<sub>1c</sub> Program was evaluated with 40 anticoagulated whole blood samples ranging from 4.98% to 12.15% HbA<sub>1c</sub>. The results are presented in the following table:

<b>Regression Method</b>	<b>n</b>	<b>r<sup>2</sup></b>	<b>Slope</b>	<b>Intercept</b>
Least Squares	40	0.9945	0.9743	0.3078

**Conclusion:**

When considering the similarities of the intended use, general characteristics of the two assays, the use of the same technology and excellent correlation between the two methods, it can be concluded that the D-10 Hemoglobin A<sub>1c</sub> Program and VARIANT II Hemoglobin A<sub>1c</sub> Program are substantially equivalent.



Food and Drug Administration  
2098 Gaither Road  
Rockville MD 20850

**AUG 27 2003**

Ms. Jackie Buckley  
Regulatory Affairs Representative  
Bio-Rad Laboratories, Inc.  
Diagnostics Group  
4000 Alfred Nobel Drive  
Hercules, CA 94547

Re: k031043  
Trade/Device Name: D-10 Hemoglobin A1c  
Regulation Number: 21 CFR 864.7470  
Regulation Name: Glycosylated Hemoglobin Assay  
Regulatory Class: Class II  
Product Code: LCP  
Dated: July 25, 2003  
Received: July 29, 2003

Dear Ms. Buckley:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in Title 21, Code of Federal Regulations (CFR), Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); and good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820).

Page 2 –

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific information about the application of labeling requirements to your device, or questions on the promotion and advertising of your device, please contact the Office of In Vitro Diagnostic Device Evaluation and Safety at (301) 594-3084. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>.

Sincerely yours,

A handwritten signature in black ink that reads "Steven Gutman". The signature is written in a cursive style.

Steven I. Gutman, M.D., M.B.A.  
Director  
Office of *In Vitro* Diagnostic Device  
Evaluation and Safety  
Center for Devices and  
Radiological Health

Enclosure

# Statement of Indications for Use

510(k) Number: K031043

Device Name: Bio-Rad D-10™ Hemoglobin A<sub>1c</sub> Program

Indications for Use: The Bio-Rad D-10 Hemoglobin A<sub>1c</sub> Program is intended for the percent determination of hemoglobin A<sub>1c</sub> in human whole blood using ion-exchange high performance liquid chromatography (HPLC). The D-10 Hemoglobin A<sub>1c</sub> Program is intended for use only with the Bio-Rad D-10 Hemoglobin Testing System.

For In Vitro Diagnostic Use.

Jean Cooper  
(Division Sign-Off)  
Division of Clinical Laboratory Devices  
510(k) Number K031043

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF

NEEDED

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescriptive Use ✓  
(Per 21 CFR 801.109)

OR Over-The-counter Use \_\_\_\_\_