

AUG 21 2003

K032377

510(k) Summary - Roche/Hitachi Bicarbonate liquid

Introduction According to the requirements of 21 CFR 807.92, the following information provides sufficient detail to understand the basis for a determination of substantial equivalence

Submitter name, address, contact Roche Diagnostics Corporation
9115 Hague Rd
Indianapolis IN 46250
(317) 521-3831

Contact person: Sherri L. Coenen

Date prepared: July 30, 2003

Device Name Proprietary name: Roche Diagnostics Hitachi Bicarbonate liquid

Common name: Bicarbonate Assay

Classification name: Enzymatic bicarbonate/carbon dioxide test system

Device description The Roche/Hitachi Bicarbonate liquid is a ready-to-use liquid enzymatic assay with phosphoenolpyruvate carboxylase and malate dehydrogenase. A decrease in absorbance at 415 nm is proportional to the concentration of bicarbonate in the sample.

Intended use In vitro test for the quantitative determination of bicarbonate (HCO_3^-) in human serum and plasma on Roche automated clinical chemistry analyzers.

Predicate Device We claim substantial equivalence to the currently marketed COBAS Integra Bicarbonate liquid Assay. (K031879).

510(k) Summary - Roche/Hitachi Bicarbonate liquid, continued

Reagent Summary The following table describes the similarities and differences between the Roche/Hitachi Bicarbonate liquid and the predicate device.

Topic	COBAS Integra Bicarbonate liquid (K031879)	Roche/Hitachi Bicarbonate liquid (Modified Device)
Intended Use	The cassette COBAS Integra Bicarbonate liquid (CO2-L) contains an in vitro diagnostic reagent system intended for use on COBAS Integra systems for the quantitative determination of the bicarbonate (HCO_3^-) concentration in human serum and plasma.	In vitro test for the quantitative determination of bicarbonate (HCO_3^-) in human serum and plasma on Roche clinical chemistry analyzers.
Method	Enzymatic, colorimetric test	Same
Sample type	Human Serum and Plasma	Same
Measuring range	0.46 - 50 mmol/L	1.5 - 50 mmol/L
Expected values	22 - 29 mmol/L	Same



AUG 21 2003

Food and Drug Administration
2098 Gaither Road
Rockville MD 20850

Ms. Sherri L. Coenen MT(ASCP)
Regulatory Affairs Consultant
Regulatory Submissions, Centralized Diagnostics
Roche Diagnostics Corporation
9115 Hague Road
P.O. Box 50457
Indianapolis, IN 46250

Re: k032377
Trade/Device Name: Roche/Hitachi Bicarbonate liquid
Regulation Number: 21 CFR 862.1160
Regulation Name: Bicarbonate/carbon dioxide test system
Regulatory Class: Class II
Product Code: KHS
Dated: July 30, 2003
Received: August 1, 2003

Dear Ms. Coenen:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in Title 21, Code of Federal Regulations (CFR), Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); and good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820).

Page 2 –

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific information about the application of labeling requirements to your device, or questions on the promotion and advertising of your device, please contact the Office of In Vitro Diagnostic Device Evaluation and Safety at (301) 594-3084. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>.

Sincerely yours,

A handwritten signature in black ink that reads "Steven Gutman". The signature is written in a cursive style.

Steven I. Gutman, M.D., M.B.A.
Director
Office of *In Vitro* Diagnostic Device
Evaluation and Safety
Center for Devices and
Radiological Health

Enclosure

Indications for Use Statement

510(k) Number (if known): N/A K032377

Device Name: Roche/Hitachi Bicarbonate liquid

Indications For Use:

In vitro test for the quantitative determination of bicarbonate (HCO₃⁻) in human serum and plasma on Roche automated clinical chemistry analyzers. Bicarbonate/carbon dioxide measurements are used in the diagnosis and treatment of numerous potentially serious disorders associated with changes in body acid-base balance.

Carol C Benson for Jean Cooper, DVM
Division Sign-Off

**Office of In Vitro Diagnostic Device
Evaluation and Safety**

510(k) K032377

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF
NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)
Prescription Use X OR Over-The-Counter Use _____
(Per 21 CFR 801.109)

(Optional Format 1-2-96)