

510(k) Summary - Tina-quant® RF II Test System

Introduction According to the requirements of 21 CFR 807.92, the following information provides sufficient detail to understand the basis for a determination of substantial equivalence

Submitter name, address, contact Roche Diagnostics Corporation
9115 Hague Rd
Indianapolis IN 46250
(317) 521-3831

Contact person: Sherri L. Coenen

Date prepared: August 8, 2003

Device Name Proprietary name: Tina-quant® Rheumatoid Factors II Test System

Common name: Tina-quant® RF II

Classification name: Rheumatoid Factor Test System

Device description The Tina-quant® Rheumatoid Factors II Test System is a particle enhanced immunoturbidimetric assay. Latex-bound heat inactivated IgG (antigen) reacts with the RF-antibodies in the sample to form antigen/antibody complexes which, following agglutination, are measured turbidimetrically.

Intended use Immunoturbidimetric assay for the in vitro quantitative determination of rheumatoid factors in human serum and plasma on automated clinical chemistry analyzers. Measurements may be used as an aid in the diagnosis of rheumatoid arthritis.

Predicate Device We claim substantial equivalence to the currently marketed Tina-quant® Rheumatoid Factors II Test System. (K002609).

510(k) Summary - COBAS Integra Creatinine plus ver.2,
continued

**Reagent
Summary**

The following table describes the similarities and differences between the Tina-quant® Rheumatoid Factors II Test System and the predicate device.

Topic	Tina-quant® RF II (K002609)	Tina-quant® RF II (Modified Device)
Intended Use	Immunoturbidimetric assay for the quantitative in vitro determination of rheumatoid factors in human serum on automated clinical chemistry analyzers. Measurements may be used as an aid in the diagnosis of rheumatoid arthritis.	Immunoturbidimetric assay for the quantitative in vitro determination of rheumatoid factors in human serum and plasma on automated clinical chemistry analyzers. Measurements may be used as an aid in the diagnosis of rheumatoid arthritis.
Method	Particle-enhanced immunoturbidimetric assay	Same
Sample type	Serum	Serum Li/Na Heparin, Na ₂ /K ₂ /K ₃ EDTA-Plasma
Measuring range	3 - 120.0 IU/ml	7 - 103 IU/ml
Expected values	< 14 IU/ml	Same



DEPARTMENT OF HEALTH & HUMAN SERVICES

Ms. Sherri L. Coenen
Regulatory Submissions, Centralized Diagnostics
Roche Diagnostics Corporation
9115 Hague Road
P.O. Box 50457
Indianapolis, Indiana 46250-0457

Food and Drug Administration
2098 Gaither Road
Rockville MD 20850

SEP 10 2003

Re: k032535
Trade/Device Name: Tina-quant® RF II Test System
Regulation Number: 21 CFR § 866.5775
Regulation Name: Rheumatoid factor immunological test system
Regulatory Class: II
Product Code: DHR
Dated: August 8, 2003
Received: August 18, 2003

Dear Ms. Coenen:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in Title 21, Code of Federal Regulations (CFR), Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); and good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820). This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

Page 2 –

If you desire specific information about the application of labeling requirements to your device, or questions on the promotion and advertising of your device, please contact the Office of In Vitro Diagnostic Device Evaluation and Safety at (301) 594-3084. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>.

Sincerely yours,

A handwritten signature in black ink that reads "Steven Gutman". The signature is written in a cursive style with a large initial 'S' and 'G'.

Steven I. Gutman, M.D., M.B.A.
Director
Office of In Vitro Diagnostic Device Evaluation and Safety
Center for Devices and Radiological Health

Enclosure

Indications for Use Statement

510(k) Number (if known): N/A K032535

Device Name: Tina-quant® RF II Test System

Indications For Use:

Immunturbidimetric assay for the in vitro quantitative determination of rheumatoid factors in human serum and plasma on automated clinical chemistry analyzers. Measurements may be used as an aid in the diagnosis of rheumatoid arthritis.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF
NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)
Prescription Use OR Over-The-Counter Use
(Per 21 CFR 801.109)

(Optional Format 1-2-96)

Tina Joy
Division Sign-Off

Office of In Vitro Diagnostic Device
Evaluation and Safety

510(k) K032535