

## 510(k) Summary

K030329

**FastPack<sup>®</sup> Free T<sub>4</sub> Immunoassay on the FastPack<sup>®</sup> Analyzer System**

This summary of 510(k) safety and effectiveness is being submitted in accordance with the requirements of SMDA 1990 and 21 CFR 807.92.

- |  |  |
|--|--|
| <b>1. Submitter name, address, contact</b> | Qualigen, Incorporated<br>2042 Corte del Nogal<br>Carlsbad, CA 92009 |
|  | Telephone: (760) 918-9165<br>Fax: (760) 918-9127                     |
|  | Contact Person: Dorothy Deinzer                                      |
|  | Date Prepared: January 28, 2003                                      |
- 
- |                       |  |
|-----------------------|--|
| <b>2. Device name</b> | Proprietary name: FastPack <sup>®</sup> Free T <sub>4</sub> Immunoassay on the FastPack <sup>®</sup> Analyzer System |
|                       | Common name: Chemiluminescence assay for the determination of Free T <sub>4</sub>                                    |
|                       | Classification Name: Quantitative Determination of Free T <sub>4</sub> in Human Serum                                |
- 
- |                            |   |
|----------------------------|---|
| <b>3. Predicate device</b> | Abbott Laboratories IMx Free T <sub>4</sub> (K902834) |
|----------------------------|---|
- 
- |                              |   |
|------------------------------|---|
| <b>4. Device description</b> | <i>FastPack<sup>®</sup> Free T<sub>4</sub> Immunoassay Reagents</i>   |
|                              | The FastPack <sup>®</sup> Free T <sub>4</sub> Immunoassay is a competitive chemiluminescence assay.   |
|                              | <ul style="list-style-type: none"> <li>• Primary incubation: Sample, calibrator, or control (100 µL) is added to the antibody solution (100 µL) to start the sequence. The reaction time is 10 seconds at 37° C.</li> </ul> |

- Secondary incubation: The initial reaction mixture (200  $\mu$ L) is transferred to the magnetic particles with bound testosterone and competition between the sample testosterone and the beads continues for an additional 2.5 minutes at 37° C.
- Removal of unbound materials: The paramagnetic particles are washed three times with wash buffer (0.2 mL/wash cycle) to remove unbound materials.
- Substrate addition and detection: Chemiluminogenic substrate [140  $\mu$ L] is added to the solid phase bound complex to form a chemiluminescent glow, which is measured by the FastPack® Analyzer System at 37° C.

5. **Intended use** The FastPack® Free T<sub>4</sub> Immunoassay is a paramagnetic particle immunoassay intended for the *in vitro* quantitative determination of Free T<sub>4</sub> in human serum. The FastPack® Free T<sub>4</sub> Immunoassay is designed for use with the FastPack® Analyzer System. It is intended strictly for *in vitro* diagnostic use as an aid in the diagnosis and management of thyroid dysfunction.

6. **Comparison to Predicate Device** The following tables compare the FastPack® Immunoassay System for Free T<sub>4</sub> with the Abbott Laboratories Free T<sub>4</sub> method.

Feature	FastPack® System	Abbott IMx® System
Intended Use	For the quantitative measurement of free thyroxine in human serum.	For the quantitative measurement of free thyroxine (FT <sub>4</sub> ) in human serum and heparinized plasma.
Assay Methodology: Storage Condition:	Competitive immunoassay 2-8 °C	Competitive immunoassay 2-8 °C
Data Analysis	Internal data reduction via microcomputer	Internal data reduction via microcomputer
Temperature Control	Required	Required
Test Processing	Automated	Automated
Sample Type:	Serum	Serum, Heparinized Plasma
Detector:	Photomultiplier Tube (PMT)	Photomultiplier Tube (PMT)
Label	Alkaline Phosphatase	Alkaline Phosphatase

Sample Volume	100 $\mu$ L	150 $\mu$ L
Assay Range	0.4 to 6 ng/dL	0.4 to 6 ng/dL
Instrument Required	FastPack <sup>®</sup> Analyzer System	Abbott IMx <sup>®</sup> System
Control Levels	2(not supplied)	3
Antibody	Monoclonal	Polyclonal
Competing component	T <sub>3</sub> covalently coupled to paramagnetic particles	T <sub>3</sub> -Alkaline Phosphatase Conjugate
Solid Phase	Paramagnetic particles	Latex Microparticles
Substrate	ImmuGlow <sup>™</sup> (Indoxyl -3-phosphate and lucigenin)	4-Methylumbelliferyl Phosphate
Detection	Chemiluminescence	Fluorescence
Calibration	Factory generated master curve with single-level calibration adjustment every 14 days.	Full calibration curve (six standards) with change in reagents.
Throughput	Single Sample	Batch
Time to Result	7 minutes	45 minutes to first result
Reagents Supplied as	Box of 50 disposable self contained reagent packs	Free T <sub>4</sub> Reagent Pack for 100 test, 6 vials of calibrator

**Performance Characteristics:**

Feature	FastPack <sup>®</sup> Free T <sub>4</sub>			Abbott IMx <sup>®</sup> Free T <sub>4</sub>		
<i>Precision</i>	Mean	%CV		Mean	%CV	
	ng/dL			ng/dL		
	<i>Between Run</i>			<i>Run to Run</i>		
	1	1.43	7.4	1	1.18	4.47
	2	3.51	6.2	2	3.23	5.59
	<i>Between Analyzer</i>			<i>Between Run</i>		
	1	1.43	2.8	1	1.18	3.80
	2	3.51	1.8	2	3.23	4.79
	<i>Between Reagent Lot</i>					
	1	1.43	11.1			
2	3.51	4.8				
<i>Analytical Sensitivity</i>	0.4 ng T <sub>4</sub> /dL			0.4 ng T <sub>4</sub> /dL		
<i>Method Comparison</i>	versus Abbott IMx Free T <sub>4</sub> :  n = 131 Range of values (IMx): 0.61 to 5.52 ng T <sub>4</sub> /dL Range of values (FastPack): 0.52 to 5.93 ng T <sub>4</sub> /dL y = 1.07x - 0.17 r = 0.95					

<b><i>Interfering Substances</i></b>	No interference up to:	No interference up to:
Bilirubin	40 mg/dL	20 mg/dL
Hemoglobin	1000 mg/dL	1000 mg/dL
Triglycerides	1000 mg/dL	1200 mg/dL
<b><i>Specificity</i></b>	L-T <sub>3</sub> ≤ 0.9%	L-T <sub>3</sub> ≤ 0.5%
	D-T <sub>4</sub> 100%	D-T <sub>4</sub> 100%



DEPARTMENT OF HEALTH & HUMAN SERVICES

Food and Drug Administration  
2098 Gaither Road  
Rockville MD 20850

Ms. Dorothy Deinzer  
Vice President  
Quality and Regulatory Affairs  
Qualigen, Inc.  
2042 Corte Del Nogal  
Carlsbad, CA 92009

JUL 03 2003

Re: k030329  
Trade/Device Name: FastPack<sup>®</sup> Free T<sub>4</sub> Immunoassay  
FastPack<sup>®</sup> Testo/Free T<sub>4</sub> Calibrator  
Regulation Number: 21 CFR 862.1695  
Regulation Name: Free thyroxine test system  
Regulatory Class: Class II  
Product Code: CEC; JIX  
Dated: May 23, 2003  
Received: May 28, 2003

Dear Ms. Deinzer:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in Title 21, Code of Federal Regulations (CFR), Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements; including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); and good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820).

Page 2 --

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific information about the application of labeling requirements to your device, or questions on the promotion and advertising of your device, please contact the Office of In Vitro Diagnostic Device Evaluation and Safety at (301) 594-3084. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>.

Sincerely yours,

A handwritten signature in black ink that reads "Steven Gutman". The signature is written in a cursive style.

Steven I. Gutman, M.D., M.B.A.  
Director  
Office of *In Vitro* Diagnostic Device  
Evaluation and Safety  
Center for Devices and  
Radiological Health

Enclosure

Attachment 4

Indications for Use Statement

510(k) Number K030329

Device Name FastPack® Free T4 Immunoassay

Indications for Use The FastPack® Free T4 Immunoassay is a paramagnetic particle, chemiluminescence immunoassay for the in vitro quantitative determination of free thyroxine (FT4) in human serum. The FastPack® Free T4 Immunoassay is designed for use with the FastPack® Analyzer System. It is intended strictly for in vitro diagnostic use as an aid in the diagnosis and management of thyroid dysfunction.

Jean Coogan
Division Sign-Off

Office of In Vitro Diagnostic Device Evaluation and Safety

510(k) K030329

PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use [checked] (Per 21 CFR 801.109)

OR

Over-The-Counter Use

### Indications for Use Statement

510(k) Number

K030329

Device Name

FastPack<sup>®</sup> Testo/Free T<sub>4</sub> Calibrator

Indications for Use

The FastPack<sup>®</sup> Testo/Free T<sub>4</sub> Calibrator is intended to calibrate the FastPack<sup>®</sup> Analyzer system when used for the quantitative determination of Total Testosterone or Free T<sub>4</sub> in human serum

*Jean Cooper*  
Division Sign-Off

Office of In Vitro Diagnostic Device  
Evaluation and Safety

510(k) K030329

PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use X  
(Per 21 CFR 801.109)

OR

Over-The-Counter Use \_\_\_\_\_