

K031786

AUG - 5 2003

510(k) Premarket Notification: Traditional
Trans-Nasal Esophagoscope with EndoSheath® System

510(k) Summary

Trade Name: Vision-Sciences Trans-Nasal Esophagoscope with EndoSheath® System

Sponsor: Vision-Sciences, Inc.
9 Strathmore Road
Natick, MA 01760
Registration #1223490

Device Common Name: Esophagoscope with sheath

Classification: According to Section 513 of the Federal Food, Drug, and Cosmetic Act, the device classification is Class II.

Predicate Devices: K990354 – Slide-On EndoSheath® for Flexible ENT Scopes
K012543 – EndoSheath® System for Flexible ENT Scopes
K021344 – EndoSheath® System for Flexible Fiberoptic Bronchoscope
Manufactured by:
Vision-Sciences, Inc.
9 Strathmore Road
Natick, MA 01760

Product Description: The device system described in this 510(k) consists of a flexible, fiberoptic esophagoscope and sterile, single use protective sheath.

Indications for Use:

The Vision-Sciences Trans-Nasal Esophagoscope with EndoSheath® System is indicated for use in endoscopic access and examination of the larynx, esophagus and gastro-esophageal junction. The System may also be used to assist in intubation.

Safety and Performance:

Substantial equivalence for the new device was based on design characteristics, comparison to legally marketed predicate devices, and performance testing. Performance testing included sheath burst/leak testing, sheath tensile/elongation testing, sheathed scope articulation testing, sheathed scope image quality evaluation and scope cycle testing

Conclusion:

Based on the indications for use, technological characteristics, performance testing and comparison to predicate devices, the proposed VSI Trans-Nasal Esophagoscope with EndoSheath® System has been shown to be safe and effective for its intended use.

000076

**VSI Trans-Nasal Esophagoscope with EndoSheath® System
Substantial Equivalence Comparison**

Characteristic	Proposed VSI TNE-2000 Trans-Nasal Esophagoscope with EndoSheath® System (Current Submission)	Currently Marketed VSI ENT-2000 with EndoSheath® System (K990354, K012534, K024095)	Currently Marketed VSI Bronchoscope with EndoSheath® System (K021344)	Pentax EE-1540 Trans-Nasal Esophagoscope (510(k) # unknown)	Olympus PEF-V Trans-Nasal Esophagoscope (510(k) # unknown)
Sheath material	Same as VSI predicate devices	Thermoplastic elastomer	Thermoplastic elastomer	N/A - No sheath	N/A - No sheath
Window material	Same as VSI predicate devices	Thermoplastic polymer	Thermoplastic polymer	N/A - No sheath	N/A - No sheath
Luer connector material	Same as VSI predicate devices	N/A - no luer connector	Thermoplastic polymer	N/A - No sheath	N/A - No sheath
Proximal connector tubing material	Same as VSI predicate devices	Thermoplastic polymer	Thermoplastic polymer	N/A - No sheath	N/A - No sheath
Working channel ID materials	N/A - no working channel	N/A - no working channel	2.1 mm	2.0 mm	2.0 mm
Adhesives	N/A - no working channel	N/A - no working channel	Thermoplastic polymer	Unknown	Unknown
Microbial barrier claim	Same as VSI predicate devices	UV curable	UV curable	Unknown	Unknown
Sheath installation method	Yes Slides on and off (no vacuum/pressure source required)	Yes Slides on and off (no vacuum/pressure source required)	Yes Slides on and off (no vacuum/pressure source required)	N/A - No sheath N/A - No sheath	N/A - No sheath N/A - No sheath
Sheath Working Length	27"	12"	24"	N/A - No sheath	N/A - No sheath
Minimum sheath wall thickness	.002"	.002"	.002"	N/A - No sheath	N/A - No sheath
Sheath Packaging	Tyvek/Mylar pouch	Tyvek/Mylar pouch	Tyvek/Mylar pouch	N/A - No sheath	N/A - No sheath
Scope Working Length	685 mm	300 mm	550 mm	600 mm	650 mm
Scope Insertion Tube OD	3.6 mm (w/out sheath) 4.8 mm (with sheath)	3.6 mm	6.0 mm	5.1 mm	5.3 mm
Articulation (Up/Down)	180°/90° (sheathed scope)	135°/135° (sheathed scope)	170°/120° (sheathed scope)	210°/120°	Unknown
Angle of View	90°	75°	90°	140°	Unknown
Depth of Field	3 - 50 mm	3 - 50 mm	3 - 50 mm	Unknown	Unknown

<p>Characteristic</p>	<p>Proposed VSI TNE-2000 Trans-Nasal Esophagoscope with EndoSheath® System (Current Submission)</p> <p>For use in endoscopic access and examination of the larynx, esophagus and gastro-esophageal junction. The System may also be used to assist in intubation.</p>	<p>Currently Marketed VSI ENT-2000 with EndoSheath® System (K990354, K012534, K024095)</p> <p>For use in flexible, endoscopic examination of the upper airway, vocal chords and/or nasal passages.</p>	<p>Currently Marketed VSI Bronchoscope with EndoSheath® System (K021344)</p> <p>Is used during flexible endoscopic examination of the trachea and other major passages of the lungs, to gather specimens, and/or to find and endoscopically remove foreign objects from the lungs.</p>	<p>Pentax EE-1540 Trans-Nasal Esophagoscope (510(k) # unknown)</p> <p>For use in endoscopic examination of the larynx, esophagus and gastro-esophageal junction.</p>	<p>Olympus PEF-V Trans-Nasal Esophagoscope (510(k) # unknown)</p> <p>For use in endoscopic examination of the larynx, esophagus and gastro-esophageal junction.</p>
<p>Indications for Use</p>					



AUG - 5 2003

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

Vision-Sciences, Inc.
c/o Pamela Papineau, RAC
Delphi Medical Device Consulting, Inc.
5 Whitcomb Avenue
Ayer, MA 01432

Re: K031786
Trade/Device Name: Trans-Nasal Esophagoscope with EndoSheath® System
Regulation Number: 21 CFR 874.4710
Regulation Name: Esophagoscope (flexible or rigid) and accessories
Regulatory Class: Class II
Product Code: EOX
Dated: June 6, 2003
Received: June 10, 2003

Dear Ms. Papineau:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (301) 594-4613. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,



A. Ralph Rosenthal, M.D.
Director
Division of Ophthalmic and Ear,
Nose and Throat Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

510(k) Number (if known): K031786

Device Name: Trans-Nasal Esophagoscope with EndoSheath® System

Indications for Use:

The Vision-Sciences Trans-Nasal Esophagoscope with EndoSheath® System is indicated for use in endoscopic access and examination of the larynx, esophagus and gastro-esophageal junction. The System may also be used to assist intubation.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use ✓
(Per 21 CFR 801.109)

OR

Over-the -Counter Use _____



(Division Sign-Off)
Division of Ophthalmic Ear,
Nose and Throat Devices

510(k) Number K031786

000007