

JUL - 7 2004

K041608

510(k) Summary

Applicant Information

Date prepared: June 1, 2004
 Name: Unilens Corp., USA
 Address: 10431 72nd Street, North
 Largo, FL 33777
 Contact person: Josepha Bruno, Director Quality Assurance
 Phone number: (727) 544-2531
 Fax number: (727) 545-1883

Device Information

Device classification: Class II
 Classification number: LPL
 Classification name: Lenses, Soft Contact, Daily Wear
 Trade name: C-Vue⁵⁵ Toric (methafilcon A) Soft (hydrophilic) Multifocal Contact Lens

Equivalent device

The C-Vue⁵⁵ Toric (methafilcon A) Soft (hydrophilic) Multifocal Contact Lens for Daily Wear is substantially equivalent to the predicate device identified below in terms of intended use and design.

Predicate device:

LifeStyle MV2TM Toric (polymacon) Soft (hydrophilic) Multifocal Contact Lens

Device description

The C-Vue⁵⁵ Toric (methafilcon A) Soft (hydrophilic) Multifocal Contact Lens is a front surface asphere consisting of multiple aspheric zones with an eccentric lenticulation for ballast and axis stabilization. The base curve has a toroidal posterior optic zone and a flattened peripheral curve which approximates the curvature of the sclera. The most plus power is in the center of the lens, with the power progressively becoming more minus towards the periphery. The lens material, (methafilcon A), is a copolymer of 2-hydroxyethyl methacrylate and methacrylic acid, crosslinked with ethyleneglycol dimethacrylate.

The C-Vue⁵⁵ Toric Multifocal Contact Lens is a hemispherical shell of the following dimensions:

Chord diameters	14.0 to 15.0mm
Center thickness	0.10 to 0.60mm; varies with power
Base curves	8.4 to 9.3mm
Powers	-20.00 to +20.00 diopters
Add powers	Up to +3.00 diopters
Cylinder	Up to 4.00 diopters
Axis	1° to 180° in 1° steps
Optical zone diameters	5.0 to 10.0mm

The physical/optical properties of the lens are:

specific gravity	1.09	
refractive index (wet)	1.415	
light transmittance (clear and tint)		greater than 90% in the visible light region
water content	55%	
oxygen transmissibility	18.83×10^{-11}	(cm^2/sec)($\text{ml O}_2/\text{ml} \times \text{mm Hg}$)*

* as measured by Schema Versatae Model 920 connected to a polarographic cell

Intended use

The C-Vue⁵⁵ Toric Multifocal Contact Lens is indicated for daily wear for the correction of refractive ametropia (myopia and hyperopia), presbyopia and astigmatism up to 4.00 diopters in aphakic and/or not-aphakic persons with non-diseased eyes.

Eyecare practitioners may prescribe the lens for daily wear in a Planned Replacement Program. The lens may be disinfected using a chemical disinfection system.

Substantial equivalence

The new device will be manufactured according to specified process controls and a quality management system currently in place. The device will undergo manufacturing, packaging and sterilization procedures similar to devices currently manufactured, marketed and distributed by Unilens Corp., USA. The established safety profile (pre-clinical toxicology and manufacturing/chemistry data) of the device is equivalent to the BENZ-MF (methafilcon A), 510(k) K003861. Being similar with respect to indications for use, materials, physical construction and safety and effectiveness to the predicate device, this meets the requirements per section 510(k) of the act regarding substantial equivalence and *does not raise* different questions of safety and effectiveness than the predicate device identified above.



JUL - 7 2004

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

Unilens Corp., USA
c/o Ms. Josepha Bruno
Director, Quality Assurance
10431 72nd Street North
Largo, FL 33777

Re: K041608
Trade/Device Name: C-Vue⁵⁵ Toric Multifocal (methafilcon A) Soft (hydrophilic)
Contact Lens for Daily Wear
Regulation Number: 21 CFR 886.5925
Regulation Name: Soft (hydrophilic) Contact Lens
Regulatory Class: Class II
Product Code: LPL
Dated: June 9, 2004
Received: June 15, 2004

Dear Ms. Bruno:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Page 2 - Ms. Josepha Bruno

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (301) 594-4613. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,



A. Ralph Rosenthal, M.D.
Director
Division of Ophthalmic and Ear,
Nose and Throat Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

Indications for Use

510(k) Number (if known): K041608

Device Name: C-Vue⁵⁵ Toric Multifocal (methafilcon A) Soft (Hydrophilic) Contact Lens

Indications for Use:

The C-Vue⁵⁵ Toric Multifocal Contact Lens is indicated for daily wear for the correction of refractive ametropia (myopia and hyperopia), presbyopia and astigmatism up to 4.00 diopters in aphakic and/or not-aphakic persons with non-diseased eyes.

Eyecare practitioners may prescribe the lens for daily wear in a Planned Replacement Program. The lens may be disinfected using a chemical disinfection system.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

T. ISO
(Division Sign-Off)
Division of Ophthalmic Ear,
Nose and Throat Devices

510(k) Number K041608

[Signature]
Prescription Use
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter-Use _____
(Part 21 CFR 807 Subpart C)