

510(k) Summary

Applicant/Sponsor: Biomet Manufacturing Corp.
P.O. Box 587
Warsaw, Indiana 46581-0587

Contact Person: Lonnie Witham
Telephone: (574) 267-6639
Fax: (574) 372-1683

Proprietary Name: Disposable Cement Spacer Molds for Temporary Knee Prosthesis

Common Name: Bone Cement Spacer Mold

Classification Name: Knee joint, patellofemorotibial, polymer/metal/polymer, semi-constrained, cemented prosthesis (21 CFR 888.3560)

Legally Marketed Devices To Which Substantial Equivalence Is Claimed:
Spacer-K Temporary Knee Prosthesis cleared by the FDA in K032522 and marketed by Exactech, Inc., Gainesville, FL.

Device Description: The femoral cement spacer molds are offered in four sizes (60mm, 65mm, 70mm & 75mm). The tibial cement spacer molds are offered in four sizes (65mm, 70mm, 75mm & 80mm). The disposable cement spacer molds are not implanted. They are filled with polymethylmethacrylate /gentamicin bone cement, or equivalent, either by injecting with a dispenser/gun, or by pouring the prepared cement into the mold. After the cement hardens, the temporary knee prosthesis components are removed from the molds and placed into the joint space. The temporary knee prosthesis remains in place (180 days or less) until the second stage of the two-stage procedure is performed to implant a conventional knee joint prosthesis.

Intended Use: The intended use of the Biomet disposable cement spacer molds is to provide the surgeon with a means to mold a temporary knee prosthesis at the point of care that is substantially equivalent to the Exactech Spacer-k temporary knee prosthesis cleared in K032522. The temporary knee prosthesis made with the Biomet disposable cement spacer molds has the same indication for use as the Exactech Spacer-k.

Indication for Use: Disposable cement spacer molds are indicated for use to mold a temporary total knee replacement (TKR) for skeletally mature patients undergoing a two-stage procedure due to a septic process. The molded temporary knee prosthesis is indicated for an implantation period of 180 days or less. Because of inherent mechanical limitations of the device material (polymethylmethacrylate/gentamicin), the molded temporary prosthesis is only indicated for patients who will consistently use traditional mobility assist devices (e.g. crutches, walkers) throughout the implant period.

Contraindications:

The temporary knee prosthesis made with the disposable cement spacer molds is contraindicated for the following situations:

- The patient's condition is such that a two-stage arthroplasty procedure is contraindicated due to decreased immune response or other relevant systemic clinical conditions.
- Bone loss precluding adequate support of the prosthesis.
- Lack of adequate competence (anatomical and functional) of peripheral ligamentous apparatus and extensor mechanism.
- The procedure is unjustified due to deficiencies in the patient's muscular, nervous or vascular systems.
- Poor bone quality (as in osteoporosis) could cause the prosthesis to migrate or to fracture host bone.
- Infection of the TKR cannot be confirmed.
- The infected TKR devices cannot be removed.
- The infecting pathogens are resistant to gentamicin.
- The patient is sensitive (allergic) to gentamicin, aminoglycosides or PMMA bone cement.
- A systemic or secondary remote infection is expected or confirmed.
- The patient does not have a TKR and the infection is secondary to trauma, septic arthritis or other surgical procedures.
- The patient does not have sufficient bone stock to allow insertion and fixation of the prosthesis.
- The patient has neuromuscular disorders that do not allow control of the knee joint.
- The patient's age, weight, or activity level would cause the surgeon to expect early failure of the system.

Summary of Technologies: The bone cement molds (femoral and tibial) are sterile disposables made of medical grade silicone. The disposable cement spacer molds produce temporary knee prosthesis components that are composed of similar bone cement in similar sizes as the predicate.

Non-Clinical Testing: Comparative testing was performed utilizing a knee joint simulator on both the temporary knee prosthesis made with the Biomet disposable cement spacer molds and the Exactech Spacer-K device. The temporary knee prostheses were found to be equivalent in strength and wear characteristics. Elution testing demonstrated equivalent gentamicin release.

Clinical Testing: No clinical testing was performed.

Page 2 - Mr. Lonnie Witham

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0120. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,



Miriam C. Provost, Ph.D
Acting Director
Division of General, Restorative
and Neurological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

Indications for Use

510(k) Number (if known):

Device Name: Disposable Cement Spacer Molds for Temporary Knee Prosthesis

Indications For Use:

Disposable cement spacer molds are indicated for use to mold a temporary total knee replacement (TKR) for skeletally mature patients undergoing a two-stage procedure due to a septic process. The molded temporary knee prosthesis is indicated for an implantation period of 180 days or less. Because of inherent mechanical limitations of the device material (gentamicin/polymethylmethacrylate), the molded temporary prosthesis is indicated for patients who will consistently use traditional mobility assist devices (e.g. crutches, walkers) throughout the implant period.

Contraindications:

The temporary knee prosthesis made with the disposable cement spacer molds is contraindicated for the following situations:

- The patient's condition is such that a two-stage arthroplasty procedure is contraindicated due to decreased immune response or other relevant systemic clinical conditions.
- Bone loss precluding adequate support of the prosthesis.
- Lack of adequate competence (anatomical and functional) of peripheral ligamentous apparatus and extensor mechanism.
- The procedure is unjustified due to deficiencies in the patient's muscular, nervous or vascular systems.
- Poor bone quality (as in osteoporosis) could cause the prosthesis to migrate or to fracture host bone.
- Infection of the TKR cannot be confirmed.
- The infected TKR devices cannot be removed.
- The infecting pathogens are resistant to gentamicin.
- The patient is sensitive (allergic) to gentamicin, aminoglycosides or PMMA bone cement.
- A systemic or secondary remote infection is expected or confirmed.
- The patient does not have a TKR and the infection is secondary to trauma, septic arthritis or other surgical procedures.
- The patient does not have sufficient bone stock to allow insertion and fixation of the prosthesis.
- The patient has neuromuscular disorders that do not allow control of the knee joint.
- The patient's age, weight, or activity level would cause the surgeon to expect early failure of the system.

Prescription Use YES
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-USE ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Miriam C. Provost
(Division Sign-Off)
Division of General, Restorative
and Neurological Devices

Page 1 of 1

510(k) Number K050210