



APR 29 2005

510(K) Submission for

## UCP Rapid™ Drug Screening Tests (k050540)

### 10. 510(K) SUMMARY

This summary of safety and effectiveness information is being submitted in accordance with the requirements of SMDA 1990 and 21CFR 807.92.

The Assigned 510(k) number is k050540

**Submitter:**

UCP Biosciences, Inc  
1445 Koll Circle, Ste 111  
San Jose, CA 95014  
Tel: 408-392-0064  
Fax: 408-392-0163

**Date:**

February 28, 2005

**Contact Person:**

Dr. Nancy Chen

**Trade Name:**

UCP Rapid™ Drug Screening Test Strips  
UCP Rapid™ Drug Screening Test Devices

**Common Name:**

Amphetamine Test System  
Barbiturate Test System  
Benzodiazepine Test System  
Cocaine and Cocaine Metabolite Test System  
Methamphetamine Test System (Methamphetamine)  
Methamphetamine Test System (MDMA)  
Opiate Test System (Opiates)  
Methadone Test System  
Opiate Test System (Oxycodone)  
Phencyclidine Test System  
Cannabinoid Test System

**Product Code:**

DKZ (Amphetamine Test System)  
DIS (Barbiturate Test System)  
JXM (Benzodiazepine Test System)  
DIO (Cocaine Test System)  
DJC (Methamphetamine Test System)  
DJC (Methamphetamine Test System--MDMA)



510(K) Submission for  
**UCP Rapid™ Drug Screening Tests (k050540)**

DMB (Methadone Test System)  
DJG (Opiates Test System)  
DJG (Opiates Test System--Oxycodone)  
LCM (Amphetamine Test System—Phencyclidine)  
LDJ (Cannabinoid Test System)

**Regulation Section:**

21 CFR 862§ 3100  
21 CFR 862 § 3150  
21 CFR 862 § 3170  
21 CFR 862 § 3250  
21 CFR 862 § 3610  
21 CFR 862 § 3610 (MDMA)  
21 CFR 862 § 3620  
21 CFR 862 § 3250  
21 CFR 862 § 3250 (Oxycodone)  
21CFR 862 § 3100 (Phencyclidine)  
21 CFR 862 § 3870

**Panel:** Toxicology (91)

**Device Classification: II**

**Substantially Equivalent Devices:**

Product: ACON One Step Drugs of Abuse Tests  
Manufactured by ACON Laboratories

**Product Description:**

UCP Rapid™ Drug Screening Tests are competitive binding, lateral flow immunochromatographic assays for qualitatively the detection of drugs and their metabolites at the cut-off levels as indicated. The tests can be performed without the use of an instrument.

**Intended Use:**

UCP Rapid™ Drug Screening Tests are rapid, qualitative, competitive binding immunoassays and intended for qualitatively the detection of Amphetamine, Barbiturate, Cocaine, Methamphetamine, MDMA, Opiates, Methadone, Oxycodone, Phencyclidine, Marijuana (THC) and their metabolites in human urine at the following cut-off concentrations:

Amphetamine	1000 ng/mL
Barbiturates	300 ng/mL
Benzodiazepine	300 ng/mL



510(K) Submission for

## UCP Rapid™ Drug Screening Tests (k050540)

Cocaine	300 ng/mL
Methamphetamine	1000 ng/mL
Methadone	300 ng/mL
MDMA	500 ng/mL
Opiates300	300 ng/mL
Opiates 2000	2000 ng/mL
Oxycodone	100 ng/mL
Phencyclidine	25 ng/mL
Cannabinoids	50 ng/mL

The tests provide only preliminary test results, which should be confirmed by other methods such as gas chromatography/mass spectrometry (GC/MS). Clinical considerations and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are indicated. The tests are not intended to be used in monitoring the drugs levels.

### Comparison to Predicate Devices:

When compared to the predicates, UCP Rapid™ Drug Screening Tests provide the qualitative determination of the same twelve drugs in the same matrix, and utilizes the same cutoff concentrations. Both tests are immunochromatographic, lateral flow assays for the qualitative detection of drugs with visual, qualitative end results. Both tests are intended to provide preliminary analytical test results.

### Safety and Effectiveness Data:

#### Accuracy

A clinical comparison study was conducted using 128 clinical specimens per each drug type including approximately 10% of the specimens containing one type drug at concentrations between -50% cutoff to cutoff ranges, another 10% of the specimens containing one type drug at concentrations between cutoffs to +50% cutoff ranges. The study was compared the test results between UCP Rapid™ Drug Screening Tests with ACON One Step Drugs of Abuse Tests. Total 64 positive clinical urine specimens and 64 negative clinical urine specimens were tested against each drug. All test results were confirmed with GC/MS analysis. UCP Rapid™ Drug Screening Tests demonstrated performance of > 96% for all drugs when performance was compared to a legally marketed device and GC/MS.

### Other Information about Performance Characteristics:

The performance characteristics of UCP Rapid™ Drug Screening Tests was evaluated by precision study, sensitivity study, specificity and cross reactivity study, interference study and stability study. The study results indicate that UCP



510(K) Submission for  
**UCP Rapid™ Drug Screening Tests (k050540)**

Rapid™ Drug Screening Tests perform satisfactorily when used according to the package inserts.

**Conclusion:**

UCP Rapid™ Drug Screening Tests is substantially equivalent to ACON Laboratories One Step Drug of Abuse Tests.



DEPARTMENT OF HEALTH & HUMAN SERVICES

APR 29 2005

Food and Drug Administration  
2098 Gaither Road  
Rockville MD 20850

Dr. Nancy Chen  
UCP Biosciences, Inc.  
1445 Koll Circle, Suite 111  
San Jose, CA 95112

Re: k050540  
Trade/Device Name: UCP Rapid™ Drug Screening Test Strips  
UCP Rapid™ Drug Screening Test Devices  
Regulation Number: 21 CFR 862.3100  
Regulation Name: Amphetamine test system  
Regulatory Class: Class II  
Product Code: DKZ, DIS, JXM, LDJ, DIO, DJC, DJR, DJG, LCM  
Dated: February 28, 2005  
Received: March 2, 2005

Dear Ms. Chen:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in Title 21, Code of Federal Regulations (CFR), Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); and good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820).

Page 2 –

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific information about the application of labeling requirements to your device, or questions on the promotion and advertising of your device, please contact the Office of *In Vitro* Diagnostic Device Evaluation and Safety at (240)276-0484. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>

Sincerely yours,



Jean M. Cooper, MS, D.V.M.

Director

Division of Chemistry and Toxicology

Office of *In Vitro* Diagnostic Device

Evaluation and Safety

Center for Devices and

Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known): k050540

Device Name: UCP Rapid™ Drug Screening Test Strips  
UCP Rapid™ Drug Screening Test Devices

### Indications For Use:

The UCP Rapid™ Drug Screening Test Strips and UCP Rapid™ Drug Screening Test Devices are rapid, qualitative, competitive binding immunoassays for the detection of Amphetamine, Barbiturate, Benzodiazepine, Cocaine, Methamphetamine, MDMA, Opiates, Methadone, Oxycodone, Phencyclidine, Marijuana and their metabolites in human urine at the following cutoff levels:

<u>Test</u>	<u>Calibrator</u>	<u>Cut-off</u>
Amphetamine	D-Amphetamine	1000 ng/mL
Barbiturates	Secobarbital	300 ng/mL
Benzodiazepine	Oxazepam	300 ng/mL
Cocaine	Benzoyllecgonine	300 ng/mL
Methamphetamine	D-Methamphetamine	1000 ng/mL
Methadone	Methadone	300 ng/mL
MDMA	D, L-MDMA	500 ng/mL
Opiates 300	Morphine	300 ng/mL
Opiates 2000	Morphine	2000 ng/mL
Oxycodone	Oxycodone	100 ng/mL
Phencyclidine	Phencyclidine	25 ng/mL
Cannabinoids	11-nor- $\Delta^9$ -THC-9 COOH	50 ng/mL

The tests provide only preliminary data, which should be confirmed by other methods such as gas chromatography/mass spectrometry (GC/MS). Clinical considerations and professional judgment should be applied to any drug of abuse test results, particularly when preliminary positive results are indicated. The tests are not intended to be used in monitoring drug levels.

Prescription U   x   AND/OR  
(Part 21 CFR 801 Subpart D)

Over-The-Counter Use       
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

  
Concurrence of CDRH, Office of Device Evaluation (ODE)  
Division Sign-Off

Page 75 of 125

Office of In Vitro Diagnostic Device  
Evaluation and Safety

510(k)   k050540