

**Section 3****HemosIL Factor II Deficient Plasma - 510(k) Summary  
(Summary of Safety and Effectiveness)****Submitted by:**

Instrumentation Laboratory Co.  
113 Hartwell Avenue  
Lexington, MA 02421  
Phone: 781-861-4467  
Fax: 781-861-4207

**Contact Person:**

Carol Marble, Regulatory Affairs Director  
Phone: 781-861-4467 / Fax: 781-861-4207

**Summary Prepared:**

March 14, 2005

**Name of the Device:**

HemosIL Factor II Deficient Plasma

**Regulatory Information:**

Regulation Section: Factor Deficiency Test (864.7290)  
Classification: Class II  
Product Code: GJT  
Panel: Hematology

**Identification of Predicate Device(s):**

K900133 Hemoliance Factor II Deficient Plasma on ELECTRA Series Analyzers  
K002400 HemosIL Factor II Deficient Plasma\* on ACL Family of Analyzers  
\*NOTE: FDA cleared as part of each ACL instrument 510(k): for example, the ACL Advance (K002400)

**Description of the Device/Intended Use(s):**

HemosIL Factor II Deficient Plasma is human plasma immunodepleted of factor II and intended for the *in vitro* diagnostic quantitative determination of factor II activity in citrated plasma, based on the prothrombin time (PT) assay, on IL Coagulation and ELECTRA Systems.

Abnormalities of the extrinsic pathway factors are determined by performing a modified prothrombin time (PT) test. Patient plasma is diluted and added to a plasma deficient in factor II. Correction of the clotting time of the deficient plasma is proportional to the concentration (% activity) of the factor II in the patient plasma, interpolated from a calibration curve.

**Statement of Technological Characteristics of the Device Compared to Predicate Device:**

The new HemosIL Factor II Deficient Plasma is substantially equivalent to Hemoliance Factor II Deficient Plasma (on ELECTRA Series Analyzers) and HemosIL Factor II Deficient Plasma (on ACL Family of Analyzers) in performance, intended use and safety and effectiveness.

**Section 3 (Cont.)**  
**HemosIL Factor II Deficient Plasma - 510(k) Summary**  
**(Summary of Safety and Effectiveness)**

**Summary of Performance Data:**

**Method Comparison**

In an in-house method comparison study evaluating approximately 60 citrated plasma samples, the correlation statistics for HemosIL Factor II Deficient Plasma versus the predicate devices are shown below:

NOTE: HemosIL RecombiPlasTin was used as the PT reagent in testing.

<b>System</b>	<b>slope</b>	<b>intercept</b>	<b>r</b>	<b>Reference method</b>
ACL 300	1.0357	1.0196	0.9921	HemosIL Factor II Deficient Plasma
ACL 6000	1.0464	-2.1396	0.9954	HemosIL Factor II Deficient Plasma
ACL 9000	1.0458	-1.6123	0.9953	HemosIL Factor II Deficient Plasma
ACL TOP	1.0582	-4.6831	0.9855	HemosIL Factor II Deficient Plasma
E1600C	1.0603	-0.2821	0.9917	Hemoliance Factor II Deficient Plasma

In a separate clinical study (n=61), the following correlation statistics were obtained on an ACL Futura using a specific lot of PT reagent (RecombiPlasTin):

<b>System</b>	<b>slope</b>	<b>intercept</b>	<b>r</b>	<b>Reference method</b>
ACL Futura	1.0602	-3.6113	0.9946	HemosIL Factor II Deficient Plasma

**Within Run Precision**

Within run and between run precision was assessed over multiple runs (n=80) on different instruments using a specific lot of PT reagent (RecombiPlasTin) and both normal and abnormal controls.

<b>Instrument</b>	<b>Control</b>	<b>Mean % Factor II</b>	<b>Within run %CV</b>	<b>Between Run %CV</b>
<b>ACL 9000</b>	Normal Control	99.2	2.5	2.0
	Special Test Control Level 2	33.8	2.5	4.5
<b>ACL Futura</b>	Normal Control	86.8	3.9	4.7
	Special Test Control Level 2	25.9	6.2	3.1
<b>ACL TOP</b>	Normal Control	130.6	2.8	4.3
	Special Test Control Level 2	28.6	3.3	4.8
<b>ELECTRA 1600C</b>	Normal Control	99.4	2.3	5.1
	Special Test Control Level 2	33.3	2.8	6.4



DEPARTMENT OF HEALTH & HUMAN SERVICES

Food and Drug Administration  
2098 Gaither Road  
Rockville MD 20850

MAY - 2 2005

Ms. Carol Marble  
Regulatory Affairs Director  
Instrumentation Laboratory Co.  
113 Hartwell Avenue  
Lexington, MA 02421

Re: k050661  
Trade/Device Name: HemosIL Factor II Deficient Plasma  
Regulation Number: 21 CFR § 864.7290  
Regulation Name: Factor Deficiency Test  
Regulatory Class: II  
Product Code: GJT  
Dated: March 14, 2005  
Received: March 15, 2005

Dear Ms. Marble:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in Title 21, Code of Federal Regulations (CFR), Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); and good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820). This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

Page 2 –

If you desire specific information about the application of labeling requirements to your device, or questions on the promotion and advertising of your device, please contact the Office of In Vitro Diagnostic Device Evaluation and Safety at (240) 276-0484. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>

Sincerely yours,



Robert L. Becker, Jr., MD, Ph.D  
Director

Division of Immunology and Hematology  
Office of *In Vitro* Diagnostic Device  
Evaluation and Safety  
Center for Devices and  
Radiological Health

Enclosure

## Indications for Use Statement

510(k) Number (if known): K050661

Device Name: HemosIL Factor II Deficient Plasma

### Indications for Use:

HemosIL Factor II Deficient Plasma is human plasma immunodepleted of factor II and intended for the *in vitro* diagnostic quantitative determination of factor II activity in citrated plasma, based on the prothrombin time (PT) assay, on IL Coagulation and ELECTRA Systems.

Prescription Use ✓  
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use \_\_\_\_\_  
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

---

Concurrence of CDRH, Office of In Vitro Diagnostic Devices (OIVD)

Josephine Bantata  
Division Sign/Off

Office of In Vitro Diagnostic Device  
Evaluation and Safety

510(k) K050661