

JAN 12 2006

K 053174

Respironics Novamatrix LLC  
LoFlo C5 CO<sub>2</sub> Sensor  
Special 510(k) – Device Modification

## Section iv - 510 (k) Summary

[ Refer to 21 C.F.R § 807.92 ]

Submitted by: Respironics Novamatrix, LLC  
5 Technology Drive  
Wallingford, CT 06484

Contact Person: Kevin Mader  
Q.A. and Regulatory Manager  
Phone: 203-697-6466

Date Prepared: 1/6/2006

Proprietary Name: LoFlo C5 CO<sub>2</sub> sensor

Common Name: CO<sub>2</sub> sensor

Classification Name: Class II, 21 C.F.R 868.1400

Predicate Device: Zoll M Series EtCO<sub>2</sub> LoFlo Option ,[510(k) K042417]

Description of Device: The LoFlo C5 CO<sub>2</sub> sensor is designed for continuous, non-invasive sidestream monitoring of carbon dioxide. Carbon dioxide is measured on-airway using an infrared absorption (IR) technique. The airway adapters and associated nasal cannulas are already legally marketed as accessories to the predicate device. The LoFlo C5 CO<sub>2</sub> sensor is an integrated microprocessor based data acquisition system consisting of CO<sub>2</sub> measurement, control circuitry and a high speed serial interface. The LoFlo C5 CO<sub>2</sub> sensor uses SRAM for data storage and an EEPROM to store system parameters. The firmware resides in a PROM. The operations performed by the LoFlo C5 CO<sub>2</sub> sensor include data acquisition, parameter calculation, zeroing, heater control and corrections to the CO<sub>2</sub> signal for N<sub>2</sub>O, O<sub>2</sub> and barometric pressure.

Intended Use of the Device: This sensor has the same intended use as the predicate device. For reference, the intended use of the LoFlo C5 CO<sub>2</sub> sensor is to provide carbon dioxide monitoring to a host monitoring system during anesthesia / recovery, in the intensive care unit (ICU), and in Emergency Medicine/Transport or Respiratory care.

Technological Characteristics: LoFlo C5 CO<sub>2</sub> sensor is a sidestream CO<sub>2</sub> sensor that interfaces to the patient either via an airway adapter or a nasal cannula. An airway adapter is attached to a mouthpiece or mask, or to the breathing circuit between the endotracheal tube and ventilator circuit wye, if the patient is intubated. A nasal cannula is used for non-intubated patients. The LoFlo C5 CO<sub>2</sub> sensor uses an infrared absorption (IR) technique for monitoring CO<sub>2</sub>.

The LoFlo C5 CO<sub>2</sub> sensor consists of an integrated microprocessor-based data acquisition and measurement system that measures CO<sub>2</sub> and provides a serial



JAN 12 2006

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

Respiroics Novamatrix, LLC  
Mr. Kevin Mader  
Quality Assurance & Regulatory Affairs Manager  
Hospital Division  
Respiroics Novamatrix, LLC  
5 Technology Drive  
Wallingford, Connecticut 06492-1950

Re: K053174  
Trade/Device Name: LoFlo C5 CO<sub>2</sub> Sensor  
Regulation Number: 868.1400  
Regulation Name: Carbon Dioxide Gas Analyzer  
Regulatory Class: II  
Product Code: CCK  
Dated: December 20, 2005  
Received: December 21, 2005

Dear Mr. Mader:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0120. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,



Chiu Lin, Ph.D.  
Director  
Division of Anesthesiology, General Hospital,  
Infection Control and Dental Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

## Section ii Indications for Use

510(k) Number (if known): \_\_\_\_\_

Device Name: LoFlo C5 CO<sub>2</sub> Sensor

Indications for Use:

The intended use of the LoFlo C5 CO<sub>2</sub> sensor is to provide carbon dioxide monitoring to a host monitoring system during anesthesia / recovery, in the intensive care unit (ICU), and in Emergency Medicine/Transport or Respiratory care.

Prescription Use   x    
(Per 21 CFR 801 Subpart D)

AND/OR

Over-The -Counter Use \_\_\_\_\_  
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF  
NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)

*Janette Y. Mitchell M.D.*

A 053174

Page 1 of \_\_\_\_