

MAR 10 2006

**510(k) SUMMARY**

**510(k) number:** K053335

**Submitters name and address:** Meridian Bioscience, Inc.  
3471 River Hills Drive  
Cincinnati, OH 45244  
513-271-3700  
Contact: Susan Rolih

**Date summary prepared:** January 30, 2006

**Name of the device:** Premier Platinum HpSA PLUS  
Enzyme Immunoassay for the detection of *H. pylori* antigen in human stool specimens

**Classification:** LYR, CFR section 866.3110

**Predicate device to which this device is being compared:** Premier Platinum HpSA (Meridian Bioscience, Inc., Cincinnati, OH) (K980076, K983255)

**Device description:** Premier Platinum HpSA PLUS is an in vitro diagnostic, microwell-based, enzyme-linked immunoassay for the detection of *Helicobacter pylori* antigen in human stool. The assay is intended for use in clinical laboratories to test for bacterial colonization to aid diagnosis, or monitor a patient's response during therapy to eradicate infection. The assay consists of Microwells coated with specific antibodies (solid phase/capture antibodies), Enzyme Conjugate (detector antibodies), Sample Diluent, Premier 20X Wash Buffer I, Premier Substrate Solution I, Premier Stop Solution I and Positive Control. Sample Diluent also functions as the Negative Control reagent.

No calibrators are needed to use this device.

**Intended use:** The Premier Platinum HpSA PLUS enzyme immunoassay (EIA) is an in vitro qualitative procedure for the detection of *Helicobacter pylori* antigens in human stool. Test results are intended to aid in the diagnosis of *H. pylori* infection and to monitor response during and post-therapy in patients. Accepted medical practice recommends that testing by any current method, to confirm eradication, be done at least four weeks following completion of therapy.

There is no change to the intended use of this device from its predicate.

**Comparison charts (Premier Platinum HpSA PLUS vs Predicate Device):**

Characteristics	<i>Premier Platinum HpSA PLUS</i>	<i>Premier Platinum HSA (predicate)</i>
<b>Device Type</b>		
In vitro diagnostic device	Yes	Yes
Control	Includes external control reagent	Includes external control reagent
Calibrator	No	No
<b>Intended Use</b>		
Detection of <i>H. pylori</i> antigen	Yes	Yes
Screening test	Yes	Yes
Diagnostic test	No	No
Monitoring therapy	Yes	Yes
<b>Acceptable Sample</b>		
Stool	Yes	Yes

Laboratory Equivalence with (Predicate Device) Combined Totals	Premier Platinum HpSA PLUS	Predicate
Agreement, positive tests	100%	N/A
Agreement, negative tests	94.8%	N/A
Agreement, overall	96.5%	N/A
Performance characteristics		
Precision/Reproducibility (intra-assay)	97%	100%
Linearity/reportable range	N/A	N/A
Analytical limit of detection/sensitivity	≥ 4.67 ng in stool	≥ 184 ng in stool
Assay cutoff	0.100 at OD <sub>450/630</sub>	0.120 at OD <sub>450/630</sub>
Indeterminant range	None	0.100 – 0.120 at OD <sub>450/630</sub>

### Interpretation of test results

The results of bench tests were read using a standard laboratory dual wave length spectrophotometer. Results were interpreted according to the following scale:

Spectrophotometric dual wavelength (450/630 nm)

Negative < 0.100

Positive ≥ 0.100

Results occurring in the 0.100 to 0.120 OD range were tracked to determine if a significant number of results were obtained such that would justify the inclusion of an equivocal range. The absence of equivocal results in the studies showed that this criteria was not necessary.

### Analytical sensitivity – limit of detection

**Study design:** Serial dilutions of purified *H. pylori* flagellar antigen and a *H. pylori* bacterial strain (ATCC 43504) were prepared in stool or Sample Diluent and used to determine the lowest concentration of antigen that would still yield a definitive positive result ( $A_{450/630} \geq 0.100$  on Premier Platinum HpSA PLUS). Final concentrations were calculated from the data points using linear regression analysis. **Conclusions to the study:** The analytical limit for *H. pylori* flagellar antigen is 4.67 ng/mL in stool and 0.69 ng/mL in sample diluent. The limit for *H. pylori* bacterial strain is  $1.0 \times 10^6$  organisms/mL in stool and  $4.4 \times 10^4$  organisms/mL in Sample Diluent. The limits are lower than those reported for the Predicate Device.

### Linearity

Linearity does not apply to the endpoint produced by this device.

### Interfering substance testing

#### Drugs, Nonmicrobial Substances

**Study design:** Selected drugs and other nonmicrobial substances that might be present in stool specimens were added to five known *H. pylori*-positive and five known negative samples. The final concentrations of the additives per 500 uL of human stool are as follows: TUMS - 10 mg, Mylanta - 0.84 mg, Pepto Bismol - 0.35 mg, Tagamet - 1 mg, Prilosec OTC - 1 mg, barium sulfate - 10 mg, whole blood - 100 uL, mucin - 6.7 mg, human hemoglobin (to make dark stool) – 15 mg, steric + palmitic acids (to make fatty stool) – 7.9 mg. The spiked samples were tested in triplicate and in parallel with an unspiked control. Acceptance criteria required that the values within replicates be similar to each other and to the value obtained with the unspiked specimen. None of the potentially interfering substances had a significant effect on positive or negative test results. Values correlated closely with unspiked samples. **Conclusions to the study:** Drugs or nonmicrobial substances that might be present as co-contaminants of stool samples do not adversely affect results obtained with Premier Platinum HpSA PLUS. These data correlate with data published for the predicate device.

#### Microbial/Viral organisms (potentially cross-reactive or interfering species)

**Study design:** The bacteria, yeast and viruses selected were those that might be expected to be present in human stool samples either as part of normal flora or from a disease state. The final concentration of bacteria

or yeast in each sample was  $\geq$  #4 McFarland Standard ( $1.2 \times 10^9$  organisms/mL). The final concentration of viruses in each sample was not determined. Unspiked samples were tested in parallel to provide a reference against which the reactions with spiked samples could be compared. Samples were tested in triplicate. Acceptance criteria required that the values within replicates be similar to each other and to the value obtained with the unspiked specimen. See data in Tables below. None of the potential co-contaminants adversely affected the final positive or negative test results. **Conclusions to the study:** Microbial and viral organisms that might be present as co-contaminants of stool samples do not adversely affect results obtained with Premier Platinum HpSA PLUS. These data correlate with data published for the predicate device.

Effects of various microbial organisms on positive test results.

Sample	42				32			
	Run 1	Run 2	Run 3	Avg	Run 1	Run 2	Run 3	Avg
No Spike	0.465	0.409	0.415	0.430	1.472	1.316	1.099	1.296
No Spike	0.536	0.553	0.494	0.528	1.861	1.979	1.900	1.913
No Spike	0.417	0.468	0.403	0.429	1.464	1.550	1.706	1.573
No Spike	0.445	0.460	0.478	0.472	1.308	1.345	1.562	1.405
No Spike	NT	NT	NT	NT	1.412	1.468	1.330	1.403

NT = not tested

Sample	42				32			
	Run 1	Run 2	Run 3	Avg	Run 1	Run 2	Run 3	Avg
Adenovirus	0.460	0.438	0.412	0.437	1.517	1.511	1.340	1.456
<i>Aeromonas hydrophila</i>	0.508	0.557	0.469	0.511	1.837	1.766	1.860	1.821
<i>Borrelia burgdorferi</i>	0.480	0.306	0.0.293	0.360	1.239	1.307	1.230	1.259
<i>Campylobacter lari</i>	0.421	0.522	0.410	0.451	1.721	1.733	1.848	1.767
<i>Campylobacter fetus</i>	NT	NT	NT	NT	1.505	1.633	1.623	1.587
<i>Campylobacter jejuni</i>	NT	NT	NT	NT	1.544	1.378	1.489	1.470
<i>Campylobacter jejuni 2</i>	NT	NT	NT	NT	1.520	1.526	1.396	1.481
<i>Campylobacter jejuni</i> solution	NT	NT	NT	NT	1.402	1.661	1.612	1.558
<i>Campylobacter lari</i>	NT	NT	NT	NT	1.594	1.512	1.565	1.557
<i>Candida albicans</i>	0.452	0.436	0.423	0.437	1.682	1.642	1.839	1.721
<i>Citrobacter freundii</i>	0.563	0.557	0.527	0.549	2.264	2.264	2.181	2.236
<i>Clostridium difficile</i>	0.543	0.600	0.516	0.553	1.878	1.996	2.047	1.974
<i>Clostridium perfringens</i>	0.514	0.496	0.516	0.509	1.799	1.833	1.922	1.851
<i>Enterobacter cloacae</i>	0.428	0.572	0.553	0.518	2.373	2.374	2.469	2.405
<i>Enterococcus faecalis</i>	0.506	0.517	0.555	0.526	2.033	1.963	2.134	2.043
<i>Escherichia coli</i> 0157:H7	0.569	0.534	0.527	0.543	2.177	2.228	2.243	2.216
<i>Escherichia coli</i> 8739	0.500	0.263	0.499	0.421	1.882	1.890	1.952	1.908
<i>Escherichia coli</i> 9637	0.493	0.536	0.498	0.509	1.850	1.782	1.864	1.832
<i>Escherichia fergusonii</i>	0.533	0.528	0.517	0.526	2.237	2.167	2.124	2.176
<i>Escherichia hermannii</i>	0.437	0.440	0.378	0.418	1.469	1.643	1.714	1.609
<i>Escherichia hermannii</i> EMDI-64	0.488	0.436	0.443	0.456	1.838	1.530	1.748	1.705
<i>Helicobacter pylori</i>	OUT*	OUT	OUT	OUT	OUT	OUT	OUT	OUT
<i>Klebsiella pneumoniae</i>	0.533	0.534	0.519	0.529	2.051	2.269	2.153	2.158
<i>Lactobacillus lactis</i>	0.474	0.456	0.420	0.450	1.486	1.642	1.684	1.604
<i>Listeria monocytogenes</i>	0.441	0.427	0.391	0.420	1.672	1.623	1.582	1.626
<i>Peptostreptococcus anaerobius</i>	0.422	0.341	0.381	0.381	1.559	1.481	1.516	1.519
<i>Proteus vulgaris</i>	0.509	0.487	0.476	0.491	1.880	1.950	1.962	1.931

Sample	42				32			
	Run 1	Run 2	Run 3	Avg	Run 1	Run 2	Run 3	Avg
<i>Pseudomonas aeruginosa</i>	0.217	0.560	0.536	0.438	2.208	2.272	2.158	2.213
<i>Pseudomonas fluorescens</i>	0.448	0.427	0.460	0.445	1.620	1.785	1.749	1.718
Rotavirus	0.384	0.407	0.445	0.412	1.356	1.361	1.387	1.368
<i>Salmonella enterica serovar Hilversum</i>	0.441	0.433	0.429	0.434	1.290	1.661	1.617	1.523
<i>Salmonella enterica subsp. Enterica serovar Hilversum</i>	0.519	0.537	0.534	0.530	2.232	2.281	2.336	2.283

Sample	42				32			
	Run 1	Run 2	Run 3	Avg	Run 1	Run 2	Run 3	Avg
<i>Salmonella enterica subsp. Enterica serovar Minnesota</i>	0.565	0.505	0.600	0.557	2.512	2.529	2.296	2.446
Salmonella Group B	0.413	0.431	0.388	0.411	1.632	1.597	1.784	1.671
<i>Salmonella typhimurium</i>	0.549	0.588	0.569	0.569	2.336	2.180	2.215	2.244
<i>Serratia liquefaciens</i>	0.445	0.429	0.457	0.444	1.508	1.610	1.653	1.590
<i>Serratia liquefaciens</i>	0.482	0.541	0.503	0.509	2.180	1.990	2.125	2.098
<i>Serratia marcescens</i>	0.393	0.470	0.329	0.397	1.019	1.411	1.617	1.349
<i>Shigella boydii</i>	0.540	0.552	0.520	0.537	2.171	2.091	2.235	2.166
<i>Shigella flexneri</i>	0.365	0.549	0.562	0.492	2.321	2.317	2.337	2.325
<i>Shigella dysenteriae</i>	0.536	0.534	0.535	0.535	2.294	2.241	2.275	2.270
<i>Shigella sonnei</i>	0.420	0.465	0.325	0.403	1.598	1.606	1.564	1.589
<i>Staphylococcus aureus</i>	0.618	0.526	0.527	0.557	2.243	2.255	2.189	2.229
<i>Staphylococcus aureus</i> (Cowans 1)	0.505	0.461	0.477	0.481	1.919	1.805	1.808	1.844
<i>Staphylococcus epidermidis</i>	0.587	0.518	0.570	0.558	2.294	2.200	2.121	2.205
<i>Streptococcus faecalis</i>	0.501	0.513	0.469	0.494	1.963	1.965	1.954	1.961
<i>Yersinia enterocolitica</i>	0.500	0.543	0.504	0.516	1.985	1.881	1.490	1.785
<i>Yersinia enterocolitica</i>	0.503	0.490	0.491	0.495	1.974	1.762	1.992	1.909

\*Note: "> 3.000" signifies the signal exceeded the high limit of the plate reader. The limit of the plate reader used in the study was A<sub>450/630</sub> 2.999.

Sample	2-128			
	Run 1	Run 2	Run 3	Avg
Unspiked Sample	0.386	0.382	0.323	0.364
<i>Campylobacter lari</i>	0.399	0.398	0.390	0.396
<i>Campylobacter jejuni</i>	0.468	0.442	0.456	0.455
<i>Campylobacter jejuni</i> solution	0.432	0.427	0.409	0.423
<i>Campylobacter jejuni</i> 2	0.414	0.410	0.393	0.406
<i>Campylobacter fetus</i>	0.398	0.427	0.383	0.403

Sample	2-42			
	Run 1	Run 2	Run 3	Avg
No Spike	0.011	0.009	0.014	0.011
No Spike	0.006	0.005	0.008	0.006
No Spike	0.006	0.007	0.014	0.009
No Spike	0.005	0.007	0.007	0.006
No Spike	0.006	0.008	0.006	0.007

Sample	2-42			
	Run 1	Run 2	Run 3	Avg
Adenovirus	0.007	0.007	0.006	0.007
<i>Aeromonas hydrophila</i>	0.004	0.006	0.006	0.005
<i>Borrelia burgdorferi</i>	0.012	0.006	0.007	0.008
<i>Campylobacter lari</i>	0.008	0.007	0.004	0.006
<i>Campylobacter fetus</i>	0.007	0.008	0.004	0.006
<i>Campylobacter jejuni</i>	0.016	0.006	0.008	0.010
<i>Campylobacter jejuni</i> 2	0.000	0.006	0.007	0.004
<i>Campylobacter jejuni</i> solution	0.004	0.009	0.007	0.007
<i>Campylobacter lari</i>	0.005	0.007	0.006	0.006
<i>Candida albicans</i>	0.006	0.008	0.008	0.007
<i>Citrobacter freundii</i>	0.005	0.016	0.006	0.009
<i>Clostridium difficile</i>	0.010	0.011	0.009	0.010
<i>Clostridium perfringens</i>	0.010	0.008	0.008	0.009
<i>Enterobacter cloacae</i>	0.003	0.003	0.004	0.003
<i>Enterococcus faecalis</i>	0.002	0.002	0.007	0.004
<i>Escherichia coli</i> 0157:H7	0.005	0.034	0.064	0.034
<i>Escherichia coli</i> 8739	0.005	0.005	0.007	0.006
<i>Escherichia coli</i> 9637	0.003	0.007	0.006	0.005
<i>Escherichia fergusonii</i>	0.004	0.002	0.007	0.004
<i>Escherichia hermannii</i>	0.008	0.007	0.008	0.008
<i>Escherichia hermannii</i> EMDi-64	0.006	0.008	0.007	0.007
<i>Helicobacter pylori</i>	> 3.000	> 3.000	> 3.000	> 3.000
<i>Klebsiella pneumoniae</i>	0.004	0.004	0.007	0.005
<i>Lactobacillus lactis</i>	0.007	0.007	0.017	0.010
<i>Listeria monocytogenes</i>	0.008	0.005	0.009	0.007
<i>Peptostreptococcus anaerobius</i>	0.008	0.008	0.009	0.008
<i>Proteus vulgaris</i>	0.004	0.000	0.007	0.004
<i>Pseudomonas aeruginosa</i>	0.011	0.026	0.021	0.019
<i>Pseudomonas fluorescens</i>	0.006	0.007	0.007	0.007
Rotavirus	0.006	0.005	0.005	0.005
<i>Salmonella enterica</i> serovar Hilversum	0.008	0.007	0.009	0.008
<i>Salmonella enterica</i> subsp. <i>Enterica</i> serovar Hilversum	0.002	0.010	0.045	0.019
<i>Salmonella enterica</i> subsp. <i>Enterica</i> serovar Minnesota	0.005	0.008	0.005	0.006
Salmonella Group B	0.007	0.006	0.007	0.007

Sample	2-42			
	Run 1	Run 2	Run 3	Avg
<i>Salmonella typhimurium</i>	0.003	0.007	0.005	0.005
<i>Serratia liquefaciens</i>	0.007	0.007	0.007	0.007
<i>Serratia liquefaciens</i>	0.000	0.003	0.004	0.002
<i>Serratia marcescens</i>	0.007	0.008	0.008	0.008
<i>Shigella boydii</i>	0.006	0.000	0.003	0.003
<i>Shigella flexneri</i>	0.004	0.027	0.004	0.012
<i>Shigella dysenteriae</i>	0.006	0.005	0.050	0.020
<i>Shigella sonnei</i>	0.007	0.002	0.009	0.006
<i>Staphylococcus aureus</i>	0.008	0.004	0.008	0.007
<i>Staphylococcus aureus</i> (Cowans 1)	0.007	0.009	0.007	0.008
<i>Staphylococcus epidermidis</i>	0.008	0.003	0.007	0.006
<i>Streptococcus faecalis</i>	0.006	0.005	0.007	0.006
<i>Yersinia enterocolitica</i>	0.005	0.004	0.008	0.006
<i>Yersinia enterocolitica</i>	0.004	0.003	0.006	0.004

### Performance Evaluation Data Summarized

**Comparison of Premier Platinum HpSA PLUS to Premier Platinum HpSA:** Tests with 291 samples from symptomatic patients collected either prior to or following treatment were used to demonstrate that Premier Platinum HpSA PLUS performed similarly to Premier Platinum HpSA. Thirty three of these samples were originally evaluated in an earlier trial to demonstrate the effectiveness of Premier Platinum HpSA. Test performance including 95% confidence intervals is detailed in the following table.

PP HpSA PLUS	PP HpSA		
	Positive	Negative	Indeterminate
Positive	94	10	3
Negative	0	183	1
	95% CI		
Correlation	277/287 (96.5%)	93.7% - 98.3%	
Agreement	Positive test	94/94 = 100%	
	Negative test	183/193 = 94.8%	

Eight of the ten samples that were positive by Premier Platinum HpSA PLUS, but negative by Premier Platinum HpSA, were positive by CLO, histology or UBT testing. The three samples that were positive by Premier Platinum HpSA PLUS but indeterminate by Premier Platinum HpSA were positive by CLO, histology or UBT testing. The one sample that was negative by Premier Platinum HpSA PLUS but indeterminate by Premier Platinum HpSA was negative by CLO, histology or UBT testing.

### Analysis of samples producing discordant results

Samples producing discordant results between Premier Platinum HpSA PLUS and the predicate were evaluated against test data from other conventional tests such as CLO, Histology, or UBT to determine the trueness of the results. The results of that evaluation are provided shown in the Table below.

Sample Number	PP HpSA PLUS Results	CLO/Histology/UBT Results	PP HpSA Results (Predicate)	Interpretation using CLO/Hist/UBT
UC82	Positive	Negative	Negative	FP--PPHpSAPLUS
2	Positive	Positive	Indeterminate	TP--PPHpSAPLUS
3-44	Negative	Negative	Indeterminate	TN--PPHpSAPLUS
U082	Positive	Positive	Negative	TP--PPHpSAPLUS
U004	Positive	Positive	Negative	TP--PPHpSAPLUS
U120	Positive	Positive	Negative	TP--PPHpSAPLUS
U159	Positive	Positive	Negative	TP--PPHpSAPLUS
U056	Positive	Positive	Negative	TP--PPHpSAPLUS
U137	Positive	Positive	Indeterminate	TP--PPHpSAPLUS
U161	Positive	Positive	Negative	TP--PPHpSAPLUS
P026	Positive	Positive	Negative	TP--PPHpSAPLUS
P040	Positive	Negative	Negative	FP--PPHpSAPLUS
P172	Positive	Positive	Indeterminate	TP--PPHpSAPLUS
P173	Positive	Positive	Negative	TP--PPHpSAPLUS

Legend: FP = false positive, TP = true positive, TN = true negative

### Performance Comparison Table

Performance Characteristics (rounded) in Direct Comparison to Clinical Status or Condition	Premier Platinum HpSA PLUS	Premier Platinum HpSA (Predicate)
Estimated Clinical Sensitivity	N/A	96.1%
Estimated Clinical Specificity	N/A	95.7%
Predictive Value of a Positive Test	N/A	96.1%
Predictive Value of a Negative Test	N/A	95.7%
<b>Laboratory Equivalence with (Predicate Device) Combined Totals</b>		
Agreement, positive tests	100%	N/A
Agreement, negative tests	94.8%	N/A
Correlation	96.5%	95.9%
<b>Performance characteristics</b>		
Precision/Reproducibility	100%	100%
Linearity/reportable range	N/A	N/A
Limit of detection	≥ 4.67 ng in stool	≥ 184 ng in stool
Assay cutoff	0.100 at OD <sub>450/630</sub>	0.120 at OD <sub>450/630</sub>

### Therapeutic Monitoring

**Study design:** A panel of frozen, archival specimens from four patients who were monitored during eradication therapy and tested using the predicate Premier Platinum HpSA (K980076 and K983255) were assessed using the Premier Platinum HpSA PLUS assay. One of the panels represented a low positive state with the predicate at the beginning of eradication therapy (See Figure 1.) The remaining three panels represented strongly positive states. (See Figures 2-4.) That data obtained with Premier Platinum HpSA PLUS was compared to that originally obtained with the predicate. In the case of strongly positive samples, the eradication curves for the two tests are substantively the same. The eradication curve for Premier Platinum HpSA PLUS differs from that of the predicate with low positive samples at the beginning of therapy since it produces stronger test results. However, by week four following treatment, the curves are identical. **Conclusions to the study:** Premier Platinum HpSA performs similarly to the predicate when used to monitor the effectiveness of eradication therapy.

Figure 1.




Figure 2.




Figure 3.





Figure 4.



**Reproducibility**

Assay precision, intra-assay variability and inter-assay variability were assessed with a reference panel prepared from high positive samples (n = 2), low low negative samples (n = 2), and low positive and high negative specimens (n = 1 each). The latter were diluted to near the assay limit of sensitivity. Nine replicates each of the low positive and high negative samples were included in the panel to bring the total cohort to 22 reference specimens. Each reference specimen was coded to prevent its identification during testing. Each was evaluated twice per day for three consecutive days by three different laboratories. In accordance with the IFU, values of < 0.100 are interpreted as negative when results are read at A450/630.

High negative samples (OD values just below 0.100) produced weakly positive results (OD values just above 0.100) in 42 out of 162 times tests. It is expected that high negative samples tested at the cut-off will produce weakly positive results 50% of the time. (See EP12-A, User protocol for evaluation of qualitative performance; approved guideline; NCCLS/CLSI, Vol. 22, no.14, 2002.) Low positive, high positive and low negative samples produced the correct results 100% of the time. Reproducibility was 100% with no intra-assay and inter-assay variability for samples prepared above or below the limit of analytical sensitivity.

Results with reproducibility test panel

Sample ID	Sample Qual. Result	Technologist 1						Technologist 2						Technologist 3					
		Day 1		Day 2		Day 3		Day 1		Day 2		Day 3		Day 1		Day 2		Day 3	
		Run 1	Run 2	run 1	run 2	Run 1	Run 2	run 1	run 2	Run 1	Run 2	run 1	run 2	Run 1	Run 2	run 1	run 2	Run 1	Run 2
1 HP #1	1.510	2.162	1.992	1.830	1.924	1.852	1.754	2.092	1.986	1.941	1.842	1.880	1.745	2.380	2.300	1.962	2.361	2.233	2.364
2 HP #2	1.061	1.445	1.436	1.283	1.318	1.244	1.244	1.537	1.406	1.302	1.307	1.293	1.154	1.858	1.743	1.400	1.698	1.699	1.614
3 Cut off LP #1	0.133	0.217	0.206	0.176	0.230	0.188	0.162	0.274	0.234	0.192	0.203	0.186	0.175	0.292	0.290	0.294	0.289	0.290	0.301
4 Cut off LP #2	0.133	0.206	0.197	0.171	0.214	0.166	0.165	0.271	0.237	0.213	0.189	0.173	0.185	0.310	0.294	0.268	0.273	0.299	0.299
5 Cut off LP #3	0.133	0.179	0.195	0.195	0.207	0.176	0.179	0.197	0.234	0.199	0.185	0.158	0.178	0.318	0.311	0.305	0.298	0.299	0.290
6 Cut off LP #4	0.133	0.185	0.215	0.187	0.210	0.182	0.149	0.265	0.219	0.219	0.175	0.163	0.167	0.295	0.312	0.279	0.308	0.295	0.336
7 Cut off LP #5	0.133	0.180	0.211	0.179	0.210	0.176	0.172	0.243	0.185	0.214	0.164	0.149	0.167	0.286	0.275	0.288	0.318	0.296	0.316
8 Cut off LP #6	0.133	0.201	0.219	0.221	0.220	0.188	0.178	0.239	0.215	0.214	0.146	0.137	0.156	0.300	0.292	0.287	0.315	0.305	0.318
9 Cut off LP #7	0.133	0.207	0.187	0.206	0.203	0.189	0.131	0.217	0.220	0.203	0.140	0.133	0.158	0.293	0.291	0.306	0.287	0.308	0.300
10 Cut off LP #8	0.133	0.186	0.187	0.202	0.223	0.179	0.163	0.233	0.201	0.211	0.118	0.097	0.140	0.309	0.307	0.267	0.298	0.290	0.299
11 Cut off LP #9	0.133	0.193	0.192	0.175	0.173	0.159	0.156	0.238	0.200	0.204	0.140	0.124	0.145	0.331	0.317	0.301	0.298	0.315	0.314
12 Cut off HN #1	0.075	0.094	0.104	0.074	0.121	0.076	0.078	0.131	0.097	0.098	0.059	0.060	0.066	0.096	0.094	0.096	0.094	0.095	0.090
13 Cut off HN #2	0.075	0.101	0.103	0.069	0.105	0.074	0.074	0.152	0.112	0.083	0.090	0.097	0.105	0.081	0.093	0.095	0.093	0.085	0.094
14 Cut off HN #3	0.075	0.073	0.104	0.088	0.092	0.082	0.071	0.161	0.118	0.090	0.105	0.088	0.101	0.074	0.086	0.084	0.081	0.087	0.091
15 Cut off HN #4	0.075	0.058	0.109	0.096	0.079	0.080	0.090	0.142	0.115	0.094	0.105	0.083	0.099	0.096	0.081	0.096	0.087	0.089	0.084
16 Cut off HN #5	0.075	0.080	0.109	0.090	0.089	0.085	0.071	0.147	0.112	0.102	0.079	0.073	0.103	0.092	0.093	0.093	0.089	0.090	0.081
17 Cut off HN #6	0.075	0.083	0.116	0.090	0.101	0.098	0.087	0.137	0.111	0.107	0.085	0.073	0.091	0.090	0.098	0.093	0.091	0.093	0.082
18 Cut off HN #7	0.075	0.088	0.113	0.089	0.113	0.088	0.097	0.128	0.101	0.113	0.064	0.068	0.083	0.082	0.081	0.094	0.080	0.089	0.081
19 Cut off HN #8	0.075	0.098	0.118	0.095	0.089	0.088	0.098	0.140	0.105	0.118	0.060	0.058	0.083	0.089	0.082	0.096	0.085	0.093	0.078
20 Cut off HN #9	0.075	0.090	0.103	0.083	0.112	0.092	0.091	0.126	0.104	0.109	0.046	0.053	0.069	0.085	0.086	0.091	0.093	0.092	0.088
21 LN #1	0.005	0.005	0.006	0.006	0.005	0.005	0.007	0.001	0.002	0.004	0.001	0.001	0.002	0.003	0.004	0.003	0.004	0.004	0.003
22 LN #2	0.005	0.004	0.005	0.006	0.006	0.006	0.009	0.001	0.000	0.004	0.001	0.000	0.002	0.005	0.005	0.004	0.005	0.004	0.003
Average high negative value		0.085	0.109	0.140	0.108	0.087	0.088	0.086	0.100	0.102	0.077	0.093	0.088	0.085	0.084	0.073	0.089	0.090	0.085
Average low positive value		0.195	0.201	0.242	0.216	0.304	0.299	0.190	0.210	0.208	0.162	0.288	0.298	0.178	0.162	0.147	0.163	0.300	0.308
Percent Correlation		100%	100%	100%	95%	100%	100%	59%	100%	100%	100%	95%	100%	100%	100%	100%	100%	100%	100%
Correlation of cut off Specimens												97%							

Legend: LP = low positive, HP = High positive, LN = Low negative, HN = High negative.

## **CONCLUSIONS**

Premier Platinum HpSA PLUS:

1. Can be used reliably for the rapid detection of *H. pylori* in human stool specimens
2. Performs similarly to the existing FDA approved Premier Platinum HpSA (K980076 and K983255).



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
2098 Gaither Road  
Rockville MD 20850

MAR 10 2006

Ms. Susan Rolih  
Vice President, Regulatory Affairs and Quality Assurance  
Meridian Bioscience, Inc.  
3471 River Hills Drive  
Cincinnati, OH 45244

Re: k053335  
Trade/Device Name: Premier Platinum HpSA PLUS  
Regulation Number: 21 CFR 866.3110  
Regulation Name: Campylobacter Fetus Serological Reagents  
Regulatory Class: Class I  
Product Code: LYR  
Dated: January 31, 2006  
Received: February 1, 2006

Dear Ms. Rolih:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in Title 21, Code of Federal Regulations (CFR), Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); and good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820).

Page 2 --

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific information about the application of labeling requirements to your device, or questions on the promotion and advertising of your device, please contact the Office of *In Vitro* Diagnostic Device Evaluation and Safety at (240)276-0484. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>

Sincerely yours,



Sally A. Hojvat, M.Sc., Ph.D.  
Director  
Division of Microbiology Devices  
Office of *In Vitro* Diagnostic Device  
Evaluation and Safety  
Center for Devices and  
Radiological Health

Enclosure

**INDICATIONS FOR USE STATEMENT**  
**Premier Platinum HpSA PLUS**

510(K) Number: K053335

The Premier Platinum HpSA PLUS enzyme immunoassay (EIA) is an in vitro qualitative procedure for the detection of *Helicobacter pylori* antigens in human stool. Test results are intended to aid in the diagnosis of *H. pylori* infection and to monitor response during and post-therapy in patients. Accepted medical practice recommends that testing by any current method, to confirm eradication, be done at least four weeks following completion of therapy.

Prescription Use   X    
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use \_\_\_\_\_  
(21 CFR 807 Subpart C)

---

Concurrence of CDRH, Office of In Vitro Diagnostic Devices (OIVD)

  
\_\_\_\_\_  
Division Sign-Off

**Office of In Vitro Diagnostic Device  
Evaluation and Safety**

510(k) K053335