

K061598

**510(k) Summary**  
**HemosIL Homocysteine and HemosIL Homocysteine Controls**

**Submitted by:**

Instrumentation Laboratory Co.  
113 Hartwell Avenue  
Lexington, MA 02421

**Contact Information:**

Carol Marble, Regulatory Affairs Director  
Phone: 781-861-4467  
Fax: 781-861-4207  
E-mail: cmarble@qa.ilww.com

SEP 22 2006

**Summary Prepared:**

June 7, 2006

**Device Trade Name:**

HemosIL Homocysteine  
HemosIL Homocysteine Controls

**Regulatory Information:**

862.1377 Urinary homocysteine (nonquantitative) test system  
LPS Class II  
862.1660 Single (specified) analyte controls (assayed and unassayed)  
JJX Class I

**Identification of Predicate Device:**

K992858 Abbott IMx Homocysteine

**Device Intended Use and Description:**

HemosIL Homocysteine is an automated latex enhanced immunoassay for the quantitative determination of total L-homocysteine in human citrated plasma on IL Coagulation Systems. Homocysteine (Hcy) values can assist in the diagnosis and treatment of patients suspected of having hyperhomocysteinemia or homocystinuria.

HemosIL Homocysteine Controls are assayed quality controls intended to monitor the accuracy and precision of HemosIL Homocysteine on IL Coagulation Systems.

Hcy levels in patient plasma are measured automatically on IL Coagulation Systems in three stages:

1. Reduction of mixed disulfides and protein-bound forms of Hcy present in the plasma samples to free Hcy.
2. Enzymatic conversion of free Hcy to S-adenosyl-L-homocysteine (SAH) by the SAH hydrolase (SAHH) in the presence of an excess of adenosine.
3. Competitive agglutination reaction between anti-SAH and SAH / conjugate.

The degree of agglutination is inversely proportional to the concentration of total Hcy in the sample and is determined by measuring the decrease of transmitted light caused by the aggregates.

## 510(k) Summary

### HemosIL Homocysteine and HemosIL Homocysteine Controls

#### Statement of Technological Characteristics of the Device Compared to Predicate Device:

HemosIL Homocysteine is substantially equivalent to the commercially available predicate device (Abbott IMx Homocysteine) in performance and intended use.

#### Summary of Performance Data:

##### Method Comparison

In a method comparison study evaluating 76 paired sodium citrate and EDTA patient plasma samples with homocysteine levels ranging from 4.2 to 56.7  $\mu\text{mol/L}$ , the correlation statistics for HemosIL Homocysteine (sodium citrate plasma) versus the predicate device (EDTA plasma) are shown below:

<u>IL System</u>	<u>Slope</u>	<u>Intercept</u>	<u>r</u>
ACL Advance	0.8292	0.3503	0.9915

##### Precision

Within run and total precision assessed over multiple runs using two control levels and a plasma sample gave the following results:

<u>ACL TOP</u>	<u>Mean (<math>\mu\text{mol/L}</math>)</u>	<u>CV% (Within run)</u>	<u>CV% (Total)</u>
Hcy Control Level 1	11.4	2.0	4.8
Hcy Control Level 2	22.4	1.5	3.5
Hcy Plasma Sample	8.1	2.9	5.5

<u>ACL ELITE/ELITE PRO/ 8/9/10000</u>	<u>Mean (<math>\mu\text{mol/L}</math>)</u>	<u>CV% (Within run)</u>	<u>CV% (Total)</u>
Hcy Control Level 1	12.3	2.3	5.1
Hcy Control Level 2	22.8	4.3	6.2
Hcy Plasma Sample	8.4	2.6	5.9

<u>ACL Futura/ACL Advance</u>	<u>Mean (<math>\mu\text{mol/L}</math>)</u>	<u>CV% (Within run)</u>	<u>CV% (Total)</u>
Hcy Control Level 1	10.5	3.5	6.0
Hcy Control Level 2	21.1	2.6	3.5
Hcy Plasma Sample	7.9	3.5	5.6



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
2098 Gaither Road  
Rockville MD 20850

Ms. Carol Marble  
Instrumentation Laboratory Company  
113 Hartwell Ave.  
Lexington  
MA 02421

SEP 22 2006

Re: k061598

Trade/Device Name: ~~HemosIL Homocysteine and~~  
HemosIL Homocysteine Controls  
Regulation Number: 21 CFR 862.1377  
Regulation Name: Urinary homocysteine (non-quantitative) test system  
Regulatory Class: Class II  
Product Code: LPS, JJX  
Dated: August 21, 2006  
Received: August 22, 2006

Dear Ms. Marble:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in Title 21, Code of Federal Regulations (CFR), Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); and good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820).

Page 2 –

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific information about the application of labeling requirements to your device, or questions on the promotion and advertising of your device, please contact the Office of In Vitro Diagnostic Device Evaluation and Safety at (240) 276-0484. Also, please note the regulation entitled, "~~Misbranding by reference to premarket notification~~" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,



Alberto Gutierrez, Ph.D.  
Director  
Division of Chemistry and Toxicology  
Office of In Vitro Diagnostic Device  
Evaluation and Safety  
Center for Devices and  
Radiological Health

Enclosure

## Indications for Use Statement

510(k) Number (if known): K061598

Device Name: HemosIL Homocysteine and  
HemosIL Homocysteine Controls

### Indications for Use:

HemosIL Homocysteine is an automated latex enhanced immunoassay for the quantitative determination of total L-homocysteine in human citrated plasma on IL Coagulation Systems. Homocysteine (Hcy) values can assist in the diagnosis and treatment of patients suspected of having hyperhomocysteinemia or homocystinuria.

HemosIL Homocysteine Controls are assayed quality controls intended to monitor the accuracy and precision of HemosIL Homocysteine on IL Coagulation Systems.

For *in vitro* diagnostic use.

Prescription Use      
(Part 21 CFR 801 Subpart D)

OR Over-The-Counter Use      
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

---

Concurrence of CDRH, Office of In Vitro Diagnostic Devices (OIVD)

  
Division Sign-Off

Office of In Vitro Diagnostic Device  
Evaluation and Safety

K061598