KO62581

510(k) Summary

NOV 2 9 2006

This summary of the 510(k) safety and effectiveness information is being submitted in accordance with the requirements of SMDA 1990 and 21 CFR 807.92.

1. Submitter name, address, contact

Olympus America 3131 W Royal Lane Irving, TX 75063

Telephone:

Fax:

(972) 556-9697 (972) 556-0365

Contact Person:

David Davis

Date Prepared:

October 6, 2006

2. Device name

Proprietary Name:

Olympus TSH Reagent

Olympus TSH Calibrator Olympus TSH Control

Olympus AU3000i™ Immunoassay System

Common Name:

TSH Reagent, Calibrator & Control

Immunoassay Analyzer

Classification Name:

Thyroid-Stimulating Hormone Test System,

Calibrator, secondary

Single (specified) analyte control (assayed and

unassayed)

Analyzer, Chemistry (Photometric, Discrete),

For Clinical Use

3. Predicate device

Reagent:

Roche Elecsys® TSH Assay (K961491)

Calibrator:

Roche Elecsys® TSH CalSet (K060754)

Control:

Roche Elecsys® PreciControl TSH (K962573)

Analyzer:

Roche Elecsys® 2010 analyzer (K961481)

4. Device description

The Olympus TSH assay is a two-step paramagnetic particle enzyme immunoassay. It is based on the sandwich principle and used to quantitate TSH in serum/plasma.

The Olympus TSH assay reagent and sample are added to the assay cuvette in the following sequence:

- I. Samples are incubated first with a monoclonal anti-TSH antibody bound to paramagnetic particles.
- 2. After a washing step, a second monoclonal anti-TSH antibody conjugated with alkaline phosphatase is added. The TSH reacts with the paramagnetic particles and the conjugated antibody to form a sandwich complex. Washing steps remove the unbound material.
- 3. The chemiluminescent substrate is added to the assay cuvette and reacts with the bound alkaline phosphatase (ALP). Light generated by the reaction is measured by the luminometer. The light emission is proportional to the quantity of TSH in the sample.
- 4. Results are calculated from a pre-defined calibration curve. The Olympus AU3000i system automatically calculates the TSH concentration of each sample in mIU/L or μIU/mL.

5. Intended use

The Olympus thyroid stimulating hormone (TSH) assay is a paramagnetic particle, chemiluminescent immunoassay for the quantitative determination of TSH levels in human serum and plasma using the Olympus AU3000i™ Immunoassay System.

Measurements of TSH produced by the anterior pituitary are used in the diagnosis of thyroid or pituitary disorders.

The Olympus TSH Calibrator is used for calibrating the quantitative Olympus TSH Assay on the Olympus AU3000i Immunoassay System.

The Olympus TSH Control is used for quality control of the Olympus TSH test system on the Olympus AU3000i Immunoassay System.

The Olympus AU3000i Immunoassay System is a chemiluminescent discrete photometric chemistry analyzer for the quantitative determination of analytes in human serum and plasma.

6. Comparison to predicate device

The following tables compare the Olympus TSH Test System with the Roche Elecsys $\mbox{\it TSH}$ Test System.

Similarities					
Item AU3000i with TSH Roche Elecsys® TSH					
Intended Use	Quantitative determination of thyroid stimulating hormone (TSH) levels in human serum/plasma	Immunoassay for the in vitro quantitative determination of thyrotropin in human serum and plasma. WHO			
Traceability	WHO				
Detection/Operating Principle	Chemiluminescence	Chemiluminescence			
Assay Methodology	Sandwich immunoassay	Sandwich immunoassay			
Solid Phase	Micro-particle	Micro-particle			
Sample Type	Serum and plasma	Serum and plasma			
Antibody	Monoclonal anti-TSH mouse antibody	Monoclonal anti-TSH mouse antibody			
Analyzer Reagents	On-board Storage Bar-coded (I-D & 2-D) reagent, calibrator & control Cap/septum for increased reagent stability and evaporation control	On-board Storage Bar-coded (1-D & 2-D) reagent, calibrator & control Cap/septum for increased reagent stability and evaporation control			
Reagent, Calibrator and Control storage form	Liquid	Liquid			
Calibration and Control Stability	 Unopened At 2-8°C up to the stated expiration date Opened 28 days / 4 weeks at 2-8°C 	 Unopened At 2-8°C up to the stated expiration date Opened 28 days / 4 weeks at 2-8°C 			
Control levels and target	One	One			
concentration range	Euthyroid/Hyperthyroid threshold	Euthyroid/Hyperthyroid threshold			
Control constituents	Single	Single			
Control Matrix	Human serum with added antigen	Human serum with added antigen			
Analyzer System	Fully automated Random access Computer controlled Stat capability	Fully automated Random access Computer controlled Stat capability			
Analyzer/User Interface	Keyboard/touch screen control Colored touch screen monitor Stationary barcode scanners Graphical user interface On line help	Keyboard/touch screen control Colored touch screen monitor Stationary barcode scanners Graphical user interface On line help			
Analyzer Consumables	Disposable reaction vessels Disposable sample tips	Disposable reaction vessels Disposable sample tips			
Analyzer Sample Detection	Liquid Level Detection Clot Detection	Liquid Level Detection Clot Detection			
Analyzer Host Interface	RS232C bidirectional	RS232C bidirectional			

Differences				
Item	AU3000i with TSH	Roche Elecsys® TSH		
Instrument Required	Olympus AU3000i™	Roche Elecsys and Modular		
	Immunoassay System	analytics immunoassay analyzers.		
Traceability	WHO 3 rd IS 81/565	WHO 2 nd IRP 80/558		
Calibrator levels and target	One	Two		
Concentration	Cal: 13 µIU/mL	Cal 1: 0 μIU/mL		
		Cal 2: 1.5 μIU/Ml		
Solid phase binding principle	Direct coating	Biotin and streptavidin		
Assay Range	0.001 – 130 μIU/mL	0.005 – 100 μIU/mL		
Calibrator Matrix	Bovine serum human pituitary	Horse serum with added		
	TSH	recombinant TSH		
Analyzer Sample Volume	10 – 100 μΙ	10 – 50 μl		
Analyzer Reagent Positions	24 reagent positions	15 reagent positions		
	2 diluent/pretreatment positions	2 diluent/pretreatment positions		
Analyzer Throughput	240 results/hour	88 results/hour		

Performance Characteristics					
Item	AU3000i with TSH		Roche	Roche Elecsys® TSH	
Precision	Sample Total CV%		Sample	Sample Total CV%	
	Low	2.9	Low	8.7	
	Med	6.0	Med	3.3	
	High	3.7	High	3.6	
Functional Sensitivity	0.0013 μΙ	U/mL	0.014 μ	0.014 μIU/mL	
Analytical Sensitivity	0.0002 μΙ	0.0002 μIU/mL		0.005 μIU/mL	
Measurable Range	0.001 - 1	0.001 – 130 μIU/mL		0.005 – 100.0 μIU/mL	
Method Comparison (Passing	Intercept-	Intercept-0.0046		Intercept0.01	
Bablok)	Slope	0.935	Slope	1.01	
	R	0.9931	R	0.944	
Applicable Interfering Substances					
	Bilirubin	$\leq 10\% @ 40$	Bilirub	in $\leq 10\% @ 41$	
	mg/dL Hemolysis ≤ 5% @ 5 g/L		mg/dL		
				Hemolysis ≤ 10% @ 1 g/dL	
		Lipemia ≤ 3% @ 10 g/L		Lipemia ≤ 10% @ 1500 mg/dL	
Specificity	LH No significant		LH	No significant	
	interference			interference	
	FSH No significant			FSH No significant	
	interference		interfer hGH	interference	
		hGH Not Tested			
	nCG N	hCG Not Detected		Not Detected	
	1		1		

DEPARTMENT OF HEALTH & HUMAN SERVICES





Food and Drug Administration 2098 Gaither Road Rockville MD 20850

Mr. David Davis Manager, RA/QA Olympus America, Inc. 3131 West Royal Lane Irving, TX 75063

NOV 2 9 2006

Re: k062581

Trade/Device Name: Olympus TSH (thyroid stimulating hormone) Test System and

Olympus AU3000i™ Immunoassay System

Regulation Number: 21 CFR 862.1690

Regulation Name: Thyroid stimulating hormone test system

Regulatory Class: Class II

Product Code: JLW, JIS, JJX, JJE

Dated: August 28, 2006 Received: August 31, 2006

Dear Mr. Davis:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in Title 21, Code of Federal Regulations (CFR), Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); and good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820).

Page 2 –

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific information about the application of labeling requirements to your device, or questions on the promotion and advertising of your device, please contact the Office of In Vitro Diagnostic Device Evaluation and Safety at (240) 276-0484. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address http://www.fda.gov/cdrh/industry/support/index.html.

Sincerely yours,

Alberto Gutierrez, Ph.D.

Director

Division of Chemistry and Toxicology

Office of In Vitro Diagnostic Device

Evaluation and Safety

Center for Devices and

Radiological Health

Enclosure

Indications for Use

510(k) Number (if know	vn): K062581
Device Name:	Olympus TSH (thyroid stimulating hormone) Test System and Olympus AU3000i™ Immunoassay System
Indications for Use:	
chemiluminescent immi	l stimulating hormone (TSH) assay is a paramagnetic particle, unoassay for the quantitative determination of TSH levels in human serum slympus AU3000i™ Immunoassay System.
Measurements of TSH pituitary disorders.	produced by the anterior pituitary are used in the diagnosis of thyroid or
The Olympus TSH Calil Olympus AU3000i Imm	orator is used for calibrating the quantitative Olympus TSH Assay on the unoassay System.
The Olympus TSH Con Olympus AU3000i Imm	trol is used for quality control of the Olympus TSH test system on the unoassay System.
	Immunoassay System is a chemiluminescent discrete photometric he quantitative determination of analytes in human serum and plasma.
Prescription Use(Part 21 CFR 801 Subp	x OR Over-The-Counter Useart D) (21 CFR 801 Subpart C)
(PLEASE DO NOT WR	ITE BELOW THIS LINE – CONTINUE ON ANOTHER PAGE IF NEEDED)
Concurre	ence of CDRH, Office of In Vitro Diagnostics Devices (OIVD)
	· · · · · · · · · · · · · · · · · · ·

Office of In Vitro Diagnostic Device Evaluation and Safety

ma KO62581