

510(K) SUMMARY OF SAFETY AND EFFECTIVENESS
Grace Medical Partial Ossicular Replacement Prostheses

DEC - 6 2006

Trade Name: The family of Grace Medical Nitinol Partial Ossicular Replacement Prostheses consists of:

- Nitinol Locking Incudo-Stapedial Joint
- Nitinol Locking Angular Piston
- Nitinol Locking Piston
- Nitinol Locking Malleus Piston
- Nitinol Locking PORP
- Footplate Shoe

Common Name: Partial Ossicular Replacement Prostheses
 Classification Name: Partial Ossicular Replacement Prostheses (CFR 21 § 874.3450)

Official Contact: Jeff Cobb
 Vice President of Regulatory Affairs & Quality
 Grace Medical, Inc.
 8500 Wolf Lake Drive, Suite 110
 Memphis, TN 38133
 Telephone: (901) 380-7000
 Telefax: (901) 380-7001
 Date Prepared: November 8, 2006

Predicate Devices

- 1) SMart ISJ Prosthesis manufactured by Gyrus ENT LLC and cleared via 510(k) No. K033554.
- 2) SMart Piston manufactured by Gyrus ENT LLC and cleared via 510(k) No. K003214.
- 3) Angular Prosthesis (Plester) manufactured by Heinz Kurz GmbH Medizintechnik and cleared via 510(k) No. K972492.
- 4) CliP Piston MVP manufactured by Heinz Kurz GmbH Medizintechnik and cleared via 510(k) No. K042503.
- 5) CliP Piston aWengen manufactured by Heinz Kurz GmbH Medizintechnik and cleared via 510(k) No. K021479.
- 6) Adjustable and Fixed-Length PORPs manufactured by Grace Medical, Inc. and cleared via 510(k) No. K061853.
- 7) Partial Ossicular Replacement Prosthesis manufactured by Grace Medical, Inc. and cleared via 510(k) K972815.

Intended Use – The Grace Medical Partial Ossicular Replacement Prostheses have the same primary intended use as the predicate devices.

An ossicular replacement prosthesis is a device intended to be implanted for the functional reconstruction of segments of the ossicular chain and facilitates the conduction of sound wave from the tympanic membrane to the inner ear. Ossicular replacement prostheses are indicated for the functional restoration of the ossicular chain when a conductive hearing loss is present. Indications for use include:

- (a) Chronic middle ear disease,
- (b) Otosclerosis,
- (c) Congenital fixation of the stapes,
- (d) Secondary surgical intervention to correct for a significant and persistent conductive hearing loss from prior otologic surgery, and
- (e) Surgically correctible injury to the middle ear from trauma.

Material(s) – The Grace Medical Partial Ossicular Replacement Prostheses are manufactured from the same materials as the predicate devices:

- (i) Nitinol (ASTM F2063-05)
- (ii) Unalloyed titanium (ASTM F67)
- (iii) Titanium alloy (ASTM F136)
- (iv) Hydroxylapatite (ASTM F1185)
- (v) Fluoroplastic (USP Class VI for plastics, ASTM D1710 and ASTM F754)
- (vi) Medical grade silicone

Comparison Charts

NITINOL LOCKING INCUDO-STAPEDIAL JOINT

| | Nitinol Locking Incudo-Stapedial Joint | SMart ISJ (Gyrus) | SMart Piston (Gyrus) | Angular Prosthesis (Kurz) | Similarities or Differences |
|----------------------|---|-------------------|--------------------------|---------------------------|---|
| Intended Use | An ossicular replacement prosthesis is a device intended to be implanted for the functional reconstruction of segments of the ossicular chain and facilitates the conduction of sound wave from the tympanic membrane to the inner ear. | Same | Same | Same | Same |
| Material(s) | Nitinol | Nitinol | Nitinol Fluoroplastic | Titanium | Substantially Equivalent |
| Method of Attachment | Manual or Heat Crimp-Assist | Same | Same | Manual | Same as SMart ISJ and SMart Piston. Different than but SE to Kurz Angular Prosthesis. |
| How Supplied | Sterile | Sterile | Sterile | Sterile | Same |

NITINOL LOCKING ANGULAR PISTON

| | Nitinol Locking Angular Piston | SMart ISJ (Gyrus) | SMart Piston (Gyrus) | Angular Piston (Kurz) | Similarities or Differences |
|--------------|---|-------------------|--------------------------|-----------------------|-----------------------------|
| Intended Use | An ossicular replacement prosthesis is a device intended to be implanted for the functional reconstruction of segments of the ossicular chain and facilitates the conduction of sound wave from the tympanic membrane to the inner ear. | Same | Same | Same | Same |
| Material(s) | Nitinol Titanium | Nitinol | Nitinol Fluoroplastic | Titanium | Substantially Equivalent |
| Lengths | 3.0mm to 9.0mm | n/a | 3.0mm to 6.0mm | 4.25mm to 4.75mm | Substantially Equivalent |
| How Supplied | Sterile | Sterile | Sterile | Sterile | Same |

NITINOL LOCKING PISTON

| | Nitinol Locking Piston | SMart ISJ (Gyrus) | SMart Piston (Gyrus) | Angular Piston (Kurz) | Similarities or Differences |
|--------------|---|-------------------|--------------------------|-----------------------|-----------------------------|
| Intended Use | An ossicular replacement prosthesis is a device intended to be implanted for the functional reconstruction of segments of the ossicular chain and facilitates the conduction of sound wave from the tympanic membrane to the inner ear. | Same | Same | Same | Same |
| Material(s) | Nitinol Titanium Fluoroplastic | Nitinol | Nitinol Fluoroplastic | Titanium | Substantially Equivalent |
| Lengths | 3.0mm to 9.0mm | n/a | 3.0mm to 6.0mm | 4.25mm to 4.75mm | Substantially Equivalent |
| How Supplied | Sterile | Sterile | Sterile | Sterile | Same |

NITINOL LOCKING MALLEUS PISTON

| | Nitinol Locking Malleus Piston | Clip Piston MVP (Kurz) | Clip Piston aWengen (Kurz) | Similarities or Differences |
|--------------|---|------------------------|----------------------------|-----------------------------|
| Intended Use | An ossicular replacement prosthesis is a device intended to be implanted for the functional reconstruction of segments of the ossicular chain and facilitates the conduction of sound wave from the tympanic membrane to the inner ear. | Same | Same | Same |
| Material(s) | Nitinol Titanium Fluoroplastic | Titanium | Titanium | Substantially Equivalent |
| Lengths | 3.0mm – 9.0mm | 5.00mm – 6.25mm | 4.0mm – 5.5mm | Substantially Equivalent |
| How Supplied | Sterile | Sterile | Sterile | Same |

NITINOL LOCKING PORP

| | Nitinol Locking PORP | SMart ISJ (Gyrus) | SMart Piston (Gyrus) | Grace Adjustable and Fixed-length PORPs | Similarities or Differences |
|-------------------|---|-------------------|--------------------------|--|-----------------------------|
| Intended Use | An ossicular replacement prosthesis is a device intended to be implanted for the functional reconstruction of segments of the ossicular chain and facilitates the conduction of sound wave from the tympanic membrane to the inner ear. | Same | Same | Same | Same |
| Shaft Material(s) | Nitinol Titanium Silicone (sleeve on adjustable models) | Nitinol | Nitinol Fluoroplastic | Titanium Silicone (sleeve on adjustable models) | Substantially Equivalent |
| Head Material(s) | Titanium Hydroxylapatite | No head | No head | Titanium Hydroxylapatite | Same |
| Lengths | Fixed Sizes: 0.5mm to 9.0mm Adjustable model – Trimmed intraoperatively to length | n/a | 3.0mm to 6.0mm | Fixed Sizes: 0.5mm to 9.0mm Adjustable model – Trimmed intraoperatively to length | Substantially Equivalent |
| How Supplied | Sterile | Sterile | Sterile | Sterile | Same |

FOOTPLATE SHOE

| | Footplate Shoe | Grace HA Footplate Shoe K972815 | Similarities or Differences |
|---|---|------------------------------------|--------------------------------|
| Intended Use | An ossicular replacement prosthesis is a device intended to be implanted for the functional reconstruction of segments of the ossicular chain and facilitates the conduction of sound wave from the tympanic membrane to the inner ear. | Same | Same |
| Material(s) | Nitinol Titanium Hydroxylapatite | Hydroxylapatite | Substantially Equivalent |
| Accommodates Shafts with Diameters: | 0.3mm to 1.07mm | 1.07mm | Substantially Equivalent |
| How Supplied | Sterile | Sterile | Same |

Differences between the Grace Medical Modified Partial Ossicular Replacement Prostheses and the predicate devices should not affect the safety or effectiveness.



Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

Grace Medical, Inc.
c/o Jeff Cobb
Regulatory Affairs & Quality
8300 Wolf Lake Drive
Suite 110
Memphis, TN

DEC - 6 2006

Re: K063374

Trade/Device Name: Partial Ossicular Replacement Prostheses, Nitinol Incudo-Stapedial Joint, Nitinol Angular Piston, Nitinol Piston, Nitinol Malleus Piston, Nitinol PORP and Footplate Shoe

Regulation Number: 21 CFR 874.3450

Regulation Name: Partial ossicular replacement prosthesis

Regulatory Class: Class II

Product Code: ETB

Dated: November 8, 2006

Received: November 8, 2006

Dear Mr. Cobb:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0115. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (240) 276-3150 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,



Malvina B. Eydelman, M.D.

Director

Division of Ophthalmic and Ear, Nose
and Throat Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

510(k) Number: K063374

Device Name: Partial Ossicular Replacement Prostheses including: Nitinol Incudo-Stapedial Joint, Nitinol Angular Piston, Nitinol Piston, Nitinol Malleus Piston, Nitinol PORP, and Footplate Shoe

Indications for Use:

An ossicular replacement prosthesis is a device intended to be implanted for the functional reconstruction of segments of the ossicular chain and facilitates the conduction of sound wave from the tympanic membrane to the inner ear. Ossicular replacement prostheses are indicated for the functional restoration of the ossicular chain when a conductive hearing loss is present. Indications for use include:

- (a) Chronic middle ear disease,
- (b) Otosclerosis,
- (c) Congenital fixation of the stapes,
- (d) Secondary surgical intervention to correct for a significant and persistent conductive hearing loss from prior otologic surgery, and
- (e) Surgically correctible injury to the middle ear from trauma.

(PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use X
(Per 21 CFR 801.109)

OR

Over-The-Counter Use _____
(Optional Format 1-2/96)



(Division Sign-Off)
Division of Ophthalmic Ear,
Nose and Throat Devices

510(k) Number K063374