

FEB 12 2007
K-6

510(k) Summary
NovaBone Dental Putty – Bioactive Synthetic Bone Graft

11/21/06

1. Submitter Information:

Name: NovaBone Products, LLC
Address: 13709 Progress Boulevard, #33
Alachua, FL 32615
Telephone: (386) 462-7660
Facsimile: (386) 418-1636
Contact: David M. Gaisser

RECEIVED
2006 NOV 24 A 9:42
FDA/CDRH/ODE/PHD

2. Name of Device:

Trade Name: NovaBone Dental Putty – Bioactive Synthetic Bone Graft
Common Name: Osteoconductive Bone Void Filler
Synthetic Resorbable Bone Graft Material
Classification Name: Endosseous Implant for Bone Filling and/or Augmentation

3. Legally Marketed Predicate Device:

Predicate #1: PerioGlas – Synthetic Bone Graft Particulate
[K053387, K040278, K992416, K962492, K930115]
(Also named as NovaBone per Special 510(k) K000149)
Predicate #2: NovaBone Putty II – Bioactive Synthetic Graft
[K060728]
Predicate #3: Exactech Resorbable Bone Paste - [K020078]
Predicate #4: Grafton DBM - [K051195]

4. Device Description

NovaBone Dental Putty is an osteoconductive, bioactive, bone void filler device. The device is intended for dental intraosseous, oral, and cranio-/maxillofacial bony defects and is supplied sterile. It is composed of a calcium-phosphorus-sodium-silicate (Bioglass) particulate mixed with a synthetic binder that acts as a temporary binding agent for the particulate. The particulate and binder are provided premixed as a pliable cohesive material. On implantation, the binder is absorbed to permit tissue infiltration between the Bioglass particles. The particles are slowly absorbed and replaced by new bone tissue during the healing process. The device is intended for dental intraosseous, oral, and cranio-/maxillofacial bony defects. It is supplied sterile.

5. Intended Use

NovaBone Dental Putty is indicated to be packed into bony voids or gaps to fill and/or augment oral, dental intraosseous, and craniofacial defects. These defects

K063 549

may include: periodontal/infrabony defects; alveolar ridge augmentation (sinusotomy, osteotomy, cystectomy); dental extraction sites (ridge maintenance, implant preparation/placement); sinus lifts; cystic defects; craniofacial augmentation. NovaBone Dental Putty may be used alone in a manner comparable to autogenous bone graft chips or allograft bone particulate (demineralized freeze dried bone), or it may be mixed with either as a bone graft extender.

6. Technological Characteristics

The technological characteristics of the NovaBone Dental Putty device are similar to those of the predicates. The device and the predicates are designed as osteoconductive space-filling devices to be gently packed into defect sites and used as non-structural scaffolds for the body's natural healing and bone regeneration process. The device indications are the same as for the PerioGlas predicate.

The NovaBone Dental Putty device is identical in composition and formulation to the NovaBone Putty II predicate device. The primary component of NovaBone Dental Putty is identical to the bioactive glass (45S5 Bioglass) particulate found in the PerioGlas and NovaBone Putty II predicates. This synthetic material is both biocompatible and osteoconductive. NovaBone Dental Putty also includes a synthetic binder as an inert carrier for ease of handling and delivery, forming a premixed cohesive material. The binder is biocompatible and is absorbed after implantation, opening space between the bioactive glass particles for cell infiltration and bone formation. The bioactive glass particulate remains for a longer post-implantation period, acting as a scaffold for bone ingrowth. This particulate is absorbed and replaced by new bone tissue.

7. Warnings and Precautions

NovaBone Dental Putty does not possess sufficient mechanical strength to support load-bearing defects prior to hard tissue ingrowth. In cases where load support is required, standard internal or external stabilization techniques must be followed to obtain rigid stabilization in all planes.

NovaBone Dental Putty is intended for use by clinicians familiar with bone grafting and internal/external fixation techniques. NovaBone Dental Putty must not be used to gain screw purchase or to stabilize screw placement.

8. Complications

Possible complications are the same as to be expected of autogenous bone grafting procedures. These may include: superficial wound infection, deep wound infection, deep wound infection with osteomyelitis, delayed union, loss of reduction, failure of fusion, loss of bone graft, graft protrusion and / or dislodgement, and general complications that may arise from anesthesia and / or surgery. Complications specific to oral/dental use are those as may be typically observed for similar bone grafting procedures and may include: tooth sensitivity,

16063549

gingival recession, flap sloughing, resorption or ankylosis of the treated root, abscess formation.

9. Conclusion

NovaBone Dental Putty is claimed to be substantially equivalent to the PerioGlas, NovaBone Putty, Exactech Resorbable Bone Paste, and Grafton DBM predicate devices as a non-structural osteoconductive bone void filler for osseous defects. *In vivo* performance data were presented. Additional supporting *in vitro* data were supplied.



Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

Mr. David M. Gaiser
Vice President
NovaBone Products, LLC
One Progress Boulevard, Suite 33
Alachua, Florida 32615

FEB 12 2007

Re: K063549

Trade/Device Name: NovaBone Dental Putty-Bioactive Synthetic Bone Graft
Regulation Number: 872. 3930
Regulation Name: Bone Grafting Material
Regulatory Class: II
Product Code: LYC
Dated: November 21, 2006
Received: December 4, 2006

Dear Mr. Gaiser:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Page 2 – Mr. Gaisser

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0115. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,



Chiu Lin, Ph.D.

Director

Division of Anesthesiology, General Hospital,

Infection Control and Dental Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

K063549

STATEMENT OF INDICATIONS FOR USE

510(k) Number (if known): _____

Device Name: NovaBone Dental Putty – Bioactive Synthetic Bone Graft

Indications For Use:

The intended use of NovaBone Dental Putty is to provide a safe, biocompatible synthetic bone graft material for use in oral, dental intraosseous, and craniofacial defects. It is used alone in a manner comparable to autogenous bone graft chips or allograft bone particulate (Demineralized Freeze Dried Bone) or may be mixed with either (typically 1:1 ratio v/v) as a bone graft extender. Typical uses include:

- Periodontal/Infrabony defects
- Ridge Augmentation (sinusotomy, osteotomy, cystectomy)
- Extraction sites (ridge maintenance/augmentation, implant preparation/ placement)
- Sinus lifts
- Cystic cavities
- Cranio-facial augmentation

For larger defects, a mixture of NovaBone Dental Putty with an equal volume of allograft or autograft bone and bone marrow may improve new bone formation.

Prescription Use XX

OR
(Per 21 CFR 801.109)

Over-The-Counter Use _____

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Susan Purvis
(Division Sign-Off)
Division of Anesthesiology, General Hospital,
Infection Control, Dental Devices

510(k) Number: K063549