

**510(k) SUMMARY as required by 807.92**  
**Summary of Safety & Effectiveness Information**

MAY 29 2007

1. Device Name

Proprietary Name

50 ml Terumo Syringe for Administration of UV Sensitive Medicines

Classification Name

Piston Syringe (80FMF)

21CFR, Section 880.5860

Classification: Class II

2. Reason for Submission

New device

3. Intended Use

The 50 ml Terumo Syringe for Administration of UV Sensitive Medicines are hypodermic syringes for single use, intended for the aspiration of fluids and blood, or for the injection of fluids immediately after filling. Intended for manual use and use with power driven syringe pumps. The opaque color of the barrel allows the administration of UV sensitive medicines.

4. Description

The 50 ml Terumo Syringe for Administration of UV Sensitive Medicines is a hypodermic standard piston syringe for single use, with a 6 % luer lock tip, made of plastic material and a synthetic rubber gasket. The barrel of the syringes has an opaque color to avoid the transmission of UV light. The 6% luer lock tip allows the use with power driven syringe pumps.

5. Performance

The 50 ml Terumo Syringe was tested in accordance with EN ISO 7886-1 (1997) and EN ISO 7886-2 (1997). Other testing included barrel transmission and UV protection of the opaque barrel.

6. Substantial Equivalence

The 50 ml (with scale extension up to 60 ml) Terumo Syringe for administration of UV sensitive medicines is substantially equivalent in intended use, design, technology/principals of operation, materials and performance to the following cleared devices:

1. Terumo Syringe for Administration of UV Sensitive Medicines, manufactured by Terumo Europe N.V. (K053413)
2. 60 mL BD syringe with BD Luer-Lok tip, manufactured by Becton Dickinson & Company (Note: we presume that this syringe is covered by the 510(k) number K980987, but we were unable to verify this based on the published resources)

6. Additional Safety Information

The sterility of the 50 ml Terumo Syringe for Administration of UV Sensitive Medicines is assured by using a validated sterilization method qualified in accordance with EN 550: "Sterilization of Medical Devices: Validation and routine control of ethylene oxide" and ISO 11135: "Medical Devices: Validation and routine control of ethylene oxide sterilization" to a sterility assurance level (SAL) of  $10^{-6}$  as required by EN 556-1: "Sterilization of Medical Devices - Requirements for medical devices to be designated "STERILE" - Part - 1: Requirements for terminally sterilized medical devices".

Ethylene oxide residual levels resulting from EtO sterilization are in compliance with EN ISO 10993-7: " Biological evaluation of medical devices – Part 7: Ethylene oxide sterilization residuals"

The device's blood contacting materials were tested in accordance with the tests recommended in the FDA General Program Memorandum #G95-1 (5/1/95): Use of International Standard ISO-10993, "Biological Evaluation of Medical Devices Part-1: Evaluation and testing". Results of the testing demonstrate that the blood contacting materials are biocompatible.

The expiration dating for the 50 ml Terumo Syringe for Administration of UV Sensitive Medicines has been established at 5 years.

7. Conclusion

In summary, the 50 ml (with scale extension up to 60 ml) Terumo Syringe for Administration of UV Sensitive Medicines are substantially equivalent in intended use, design, technology/principals of operation, materials and performance to the following cleared devices:

1. Terumo Syringe for Administration of UV Sensitive Medicines, manufactured by Terumo Europe N.V. (K053413)
2. 60 mL BD syringe with BD Luer-Lock tip, manufactured by Becton Dickinson & Company (*Note: we presume that this syringe is covered by the 510(k) number K980987, but we were unable to verify this based on the published resources*)

Differences between the devices do not raise any new or different issues of safety or effectiveness.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

MAY 29 2007

Mrs. M.J. Aerts  
Manager Regulatory Affairs  
Terumo Europe N.V.  
Researchpark Zone 2 Haasrode  
Interleuvenlaan 40  
B-3001 Leuven  
BELGIUM

Re: K070264

Trade/Device Name: 50 ml Terumo® Syringe for Administration of UV  
Sensitive Medicines  
Regulation Number: 21 CFR 880.5860  
Regulation Name: Piston Syringe  
Regulatory Class: II  
Product Code: FMF  
Dated: May 21, 2007  
Received: May 23, 2007

Dear Mrs. Aerts:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0115. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (240) 276-3150 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,



Chiu Lin, Ph.D.

Director

Division of Anesthesiology, General Hospital,

Infection Control and Dental Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

K070264

## Indications for Use

510(k) Number (if known):

Device Name: 50 ml Terumo® Syringe for administration of UV sensitive medicines

Indications for Use:


The 50 ml Terumo Syringes for administration of UV sensitive medicines are sterile hypodermic syringes for single use, intended for the aspiration of fluids and blood, or for the injection of fluids immediately after filling. Intended for manual use or for use with power driven syringe pumps. The opaque color of the barrel allows the administration of UV sensitive medicines.

Prescription Use   X   AND/OR Over-The-Counter Use             
(Part 21 CFR 801 Subpart D) (21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE  
OF NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)



(Signature Sign-Off)  
Department of Anesthesiology, General Hospital,  
Injection Control, Dental Devices

510(k) Number: K070264

Page   1   of   1