


K070595

MAY 24 2007

510(k) Summary of Safety and Effectiveness

Submitter:

- SPSmedical Supply Corp.
6789 West Henrietta Road
Rush, NY 14543 U.S.A.
Phone: (585)-359-0130 Fax: (585)-359-0167
- Establishment FDA Registration No.: 1319130
- Date Summary was Prepared May 24th, 2007
- Gary J. Socola
Printed name of person submitting for 510(k)
- 
Signature of person submitting for 510(k)
- Vice President, Scientific Affairs
Title of person submitting for 510(k)

Device Name and Classification

Trade Name: SporView® Steam Self-Contained Biological Indicator

Classification Name: Sterilization Process Biological Indicator

Common Name: Self-Contained Biological Indicator

Device Classification: Class II, Regulation Number 880.2800

Product Code: FRC

Predicate Device: SGM EZTest® - Steam (K963841)

Device Description:

The SPSmedical SporView® biological indicator consists of a self-contained unit that includes bacterial spores of *Geobacillus stearothermophilus* ATCC #7953 inoculated onto a paper filter carrier and a small glass ampoule containing Tryptic Soy Broth with Bromocresol Purple acting as a pH indicator encased in a plastic vial that serves as the culture tube.

Intended Use:

SporView® Steam is a self-contained biological indicator inoculated with viable *Geobacillus stearothermophilus* bacterial spores and is intended for monitoring the efficacy of saturated steam sterilization processes operating at 121°C and 132°C gravity displacement, 132°C flash gravity displacement and 121°C - 134°C prevacuum cycles.

Statement of Similarity to the Legally Marketed Predicate Device:

- Both are intended to monitor steam sterilization cycles from 121°C to 134°C.
- Both utilize the same strain of bacterial spores.
- Both utilize the same carrier material.
- Both use virtually the same media.
- Both are similar in size and shape.
- Both are activated in the same manner.
- Both are incubated at the same temperature

Non-Clinical Testing:

Testing was performed in order to validate the indicators label claims and performance characteristics. Multiple lots of indicators were tested for;

- Resistance
- Spore population
- Media recovery in extended steam sterilization cycles
- Effects of holding time
- Reduced incubation period
- Stability of the color change
- Media Evaporation
- Effects of carrier and package materials

All test results met the defined acceptance criteria.

Conclusion:

Supportive data has demonstrated that the SPSmedical SporView® Steam Self-Contained Biological Indicator is equivalent to the legally marketed predicate device and therefore, as safe and effective as the legally marketed device.



Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

MAY 24 2007

SPSmedical Supply Corporation
C/O Mr. Daniel W. Lehtonen
Responsible Third Party Official
Intertek Testing Services
2307 East Aurora Road
Twinsburg, Ohio 44087

Re: K070595

Trade/Device Name: SporView® Steam Self Contained BI
Regulation Number: 21 CFR 880.2800
Regulation Name: Sterilization Process Indicator
Regulatory Class: II
Product Code: FRC
Dated: May 15, 2007
Received: May 16, 2007

Dear Mr. Lehtonen:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0115. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (240) 276-3150 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,



Chiu Lin, Ph.D.

Director

Division of Anesthesiology, General Hospital,

Infection Control and Dental Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

INDICATIONS for USE STATEMENT

Applicant: SPSmedical Supply Corp.

510(k) Number (if known): K070595

Device Name: SporView® Steam Self Contained BI

Indications For Use:

SporView® Steam is a self-contained biological indicator inoculated with viable *Geobacillus stearothermophilus* bacterial spores and is intended for monitoring the efficacy of saturated steam sterilization processes operating at 121°C and 132°C gravity displacement, 132°C flash gravity displacement and 121°C - 134°C prevacuum cycles.

Prescription Use _____
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use X
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Shirley R. Murphy MD
(Signature)
Department of Anesthesiology, General Hospital,
Infection Control, Dental Devices

510(k) Number Page 1 of 1

K070595