

K07/449

Page 1 of 2



7 September 2007

SEP 14 2007

**510(k) SUMMARY**

This summary of 510(k) safety and effectiveness information is being submitted in accordance with the requirements of 21 CFR §807.92:

(a)(1) Refer to information above and concluding this summary.

**(a)(2) Name of the Device**

Model Number / Name: **BRM3 Brain Monitor**

Classification Name: Electroencephalograph  
Neurology Devices, 21 CFR §882.1400, Class II, GWQ

**(a)(3) Identification of Legally Marketed Devices**

K030489	BRM2 Brain Monitor	Brainz Instruments Ltd
K033010	BRM2 Brain Monitor	Brainz Instruments Ltd
K983229	Lectromed Cerebral Function Monitor	Olympic Medical Corp

**(a)(4) Description of the Device**

The BRM3 Brain Monitor is a three-channel electroencephalograph (EEG) system, as per 21 CFR §882.1400: a device used to measure and record the electrical activity of the patient's brain by placing two or more electrodes on the head.

The BRM3 Brain Monitor system consists of a Neonatal Sensor Set (EEG skin electrodes) or Sensor Adaptor Set (electrode adaptor), Data Acquisition Unit (EEG head stage), Data Cable, Monitor, Power Supply Unit, and Roll-Pole. These components have equivalent configuration and functions to those described in 510(k) K030489 and K033010 for the BRM2 Brain Monitor.

Changes to the BRM3 Brain Monitor consist of a new touch-screen panel PC (Monitor), new EEG head stage (Data Acquisition Unit) and new mobile mounting hardware (Roll-Pole). New software features include display of a further EEG channel, display of an aEEG parameter, and improved signal quality indication and file reviewing capacity. User instructions have been updated accordingly.

**(a)(5) Statement of the Intended Use**

The Brainz Instruments Ltd **BRM3 Brain Monitor** is an electroencephalograph or EEG system, as per 21 CFR §882.1400: a device used to measure and record the electrical activity of the patient's brain obtained by placing two or more electrodes on the head.

- The BRM3 Brain Monitor is intended for use to monitor the state of the brain by acquisition of electroencephalogram (EEG) signals, in clinical environments such as the intensive care unit, and for clinical research.
- The BRM3 Brain Monitor is intended for use only with neonatal patients, defined as from birth to 28 days post-delivery, and corresponding to a post-conceptual age of 24 to 46 weeks.
- The BRM3 Brain Monitor does not provide any diagnostic conclusion about the patient's condition.

<p>www.brainz.com</p>	<p><b>Brainz Instruments Ltd</b> 25 Carbine Rd, Mt Wellington Auckland 1060 New Zealand PO Box 51078, Pakuranga Manukau 2140 New Zealand Tel +64 9 978 8896</p>	<p><b>Brainz Instruments USA Inc</b> Suite 1700 8850 Stanford Blvd Columbia MD 21045 United States Tel +1-410-312-0324 TOLL FREE 1-877-312-0324 Fax +1-410-312-0328</p>	<p><b>Brainz Instruments UK Ltd</b> Fountain Court 2 Victoria Square Victoria Street St Albans Hertfordshire AL1 3TF United Kingdom Tel +44 1727 884 784</p>	<p>see inside</p>
-----------------------	---	---	--	-------------------

**(a)(6) Technological Characteristics Summary**

The technological characteristics of the BRM3 Brain Monitor are equivalent to the predicate devices listed above. The system uses similar components for the touch-screen panel PC display unit, EEG head stage electronics module, and mobile mounting hardware. Technical equivalence has been established to ensure accuracy of the aEEG parameter implementation. Software is updated for compatibility with those changes, and the same skin electrodes are used as previously notified in 510(k) K033010.

**(b)(1) Discussion of the Non-Clinical Tests**

Non-clinical testing of the BRM3 Brain Monitor has been carried out to cover the changes to the system. These include performance testing of the new EEG head stage electronics module, full software verification procedures for updated device functionality, comparative results of the aEEG parameter implementation, and third-party testing of the entire system for compliance to IEC 60601 standards for electromagnetic compatibility, electrical safety, particular requirements for electroencephalographs, and US national requirements. This established correct functionality of the BRM3 Brain Monitor system according to requirements, and is equivalent to the predicate devices.

**(b)(2) Discussion of the Clinical Tests**

Clinical testing was not required to demonstrate the safety and effectiveness of the modifications involved in the BRM3 Brain Monitor. The changes to hardware, electronics, and software have been adequately assessed by establishing technical equivalence and bench testing as above.

**(b)(3) Conclusions Demonstrating Safety, Effectiveness and Performance**

The testing carried out for the BRM3 Brain Monitor indicates that it meets design and performance functional requirements. The software verification demonstrates that the modified device features are effective, and the hardware / electronics configuration functions equivalently to the predicate devices.

This information indicates that the BRM3 Brain Monitor is equivalent to the predicate devices in terms of device safety, effectiveness and performance.

signed:



Chris Mander  
Operations Manager  
**Brainz Instruments Ltd**

date:

7 Sept 2007



SEP 14 2007

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

Brainz Instruments Ltd.  
% Underwriters Laboratories, Inc.  
Morten Simon Christensen  
Staff Engineer & FDA Office Coordinator  
455 East Trimble Road  
San Jose, California 95131-1230

Re: K071449

Trade/Device Name: BRM3 Brain Monitor  
Regulation Number: 21 CFR 882.1400  
Regulation Name: Electroencephalograph  
Regulatory Class: II  
Product Code: GWQ  
Dated: August 13, 2007  
Received: August 14, 2007

Dear Mr. Christensen:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set

forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0120. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (240) 276-3150 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,

  
Mark N. Melkerson  
Director  
Division of General, Restorative,  
and Neurological Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

11 June 2007

**INDICATIONS FOR USE**

510(k) Number:     K071449    

Device Name: **BRM3 Brain Monitor**

Indications for Use:

The Brainz Instruments Ltd **BRM3 Brain Monitor** is an electroencephalograph or EEG system, as per 21 CFR §882.1400: a device used to measure and record the electrical activity of the patient's brain obtained by placing two or more electrodes on the head.


- The BRM3 Brain Monitor is intended for use to monitor the state of the brain by acquisition of electroencephalogram (EEG) signals, in clinical environments such as the intensive care unit, and for clinical research.
- The BRM3 Brain Monitor is intended for use only with neonatal patients, defined as from birth to 28 days post-delivery, and corresponding to a post-conceptual age of 24 to 46 weeks.
- The BRM3 Brain Monitor does not provide any diagnostic conclusion about the patient's condition.

Prescription Use   
(21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use   
(21 CFR 801 Subpart C)

Concurrence of CDRH, Office of Device Evaluation (ODE)



**(Division Sign-Off)**  
**Division of General, Restorative,  
and Neurological Devices**

510(k) Number     K071449