
510(k) Summary

This summary of the 510(k) safety and effectiveness information is being submitted in accordance with the requirements of SMDA 1990 and 21 CFR 807.92.

The assigned 510(k) number is: K073489

- | | |
|-------------------------------------|---|
| 1. Submitter name, address, contact | <p>Olympus America
Inc.
3500 Corporate
Parkway
Center Valley, PA
18034</p> <p>U.S. Telephone: 469-230-0959
U.S. Fax: 972-317-7861</p> <p>Contact Person: Stephanie G. Donnelly</p> <p>Date Prepared: December 10, 2007</p> |
| 2. Device name | <p>Proprietary Name: Olympus IgA Reagent (OSR6X171)</p> <p>Common Name: IgA reagent</p> <p>Classification Name: IgA, Antigen, Antiserum, Control</p> |
| 3. Predicate device | <p>Reagent: Olympus (OSR6X44) IgA Reagent
Submitted K951055</p> |
| 4. Device description | <p>In this Olympus procedure:</p> <ul style="list-style-type: none"> • When a sample is mixed with R1 buffer and R2 antiserum solution, human IgA reacts specifically with anti-human IgA antibodies to yield insoluble aggregates. • Immune complexes formed in solution scatter light in proportion to their size, shape and concentration. • Turbidimeters measure the reduction of incidence light due to reflection, absorption or scatter. • In the Olympus procedure, the decrease in intensity of light transmitted (increase in absorbance) through particles suspended in solution is as a result of complexes formed during the antigen-antibody reaction. |
| 5. Intended use | <p>System reagent for the quantitative determination of IgA immunoglobulins in human serum and plasma on OLYMPUS analyzers.</p> <p>For <i>in vitro</i> diagnostic use.</p> |
-

6.

The following Tables compare the new Olympus IgA (OSR6X171) reagent with the predicate devices outlined in point 3 above.

Similarities		
Item	Olympus IgA (OSR6X171) reagent	Predicate System
Intended Use	System reagent for the quantitative determination of IgA immunoglobulins in human serum and plasma on Olympus analyzers	System reagent for the quantitative determination of IgA immunoglobulins in human serum on Olympus analyzers
Measurement	Quantitative	Same
Instrument Required	Olympus AU400/400 ^e , 600/640/640 ^e and 2700/5400	Same
Reagent handling	Ready for use	Same
Assay Methodology/Operating Principle	Immunoturbidimetric	Same
Reagent storage form	Liquid On-board storage	Same
Calibrator	Olympus Serum Protein Multit-Calibrator (ODR3021)	Same
Calibration Traceability	This method is traceable to the International Reference Preparation CRM 470 (US designation RPPHS lot 91/0619)	Same
Antibody	Goat Anti-IgA antiserum	Same
Expected Values	66-433 mg/dL	Same
Reagent On Board Stability	Opened reagents are stable for 90 days when stored in the refrigerated compartment of the analyzer.	Same

Differences		
Item	Olympus IgA (OSR6X171) reagent	Predicate System
Specimen Type	Serum, EDTA Plasma and Li-Heparinized Plasma	Serum
Calibration Frequency	90 days	7 days

Performance Characteristics				
Item	Olympus IgA (OSR6X171) reagent		Predicate System	
Precision	AU400/400 ^e		AU400/400 ^e	
	Sample	Total CV%	Sample	Total CV%
	1	2.43	1	1.53
	2	2.52	2	1.63
	3	2.95	3	1.26
	AU600/640/640 ^e		AU600	
	Sample	Total CV%	Sample	Total CV%
	1	3.39	1	1.38
	2	3.85	2	1.08
	3	4.01	3	1.81
	AU2700/5400		AU640/640 ^e	
	Sample	Total CV%	Sample	Total CV%
	1	1.50	1	1.8
	2	1.91	2	1.3
	3	1.83	AU2700/5400	
		Sample	Total CV%	
		1	2.64	
		2	2.05	
		3	3.29	

Assay Range	10 to 700 mg/dL	10 to 700 mg/dL
LoQ	10 mg/dL	Not specified
Method Comparison (Linear Regression)	Intercept 15.1 mg/dL Slope 0.923 R ² 0.999 N 111 Range 38-672 mg/dL	Intercept 1.3 mg/dL Slope 0.957 R ² 0.99 N 94 Range 54-660 mg/dL
Interfering Substances	<p>AU400/400^e Bilirubin: Interference less than 2% up to 40 mg/dL Bilirubin Hemolysate: Interference less than 1% up to 500 mg/dL Hemolysate Lipemia: Interference less than 10% up to 1000 mg/dL Intralipid</p> <p>RF: Interferences less than 8% up to 600 IU/mL</p> <p>Not tested</p> <p>AU600/640/640^e Bilirubin: Interference less than 3% up to 40 mg/dL Bilirubin Hemolysate: Interference less than 5% up to 500 mg/dL Hemolysate Lipemia: Interference less than 6% up to 1000 mg/dL Intralipid</p> <p>RF: Interferences less than 8% up to 600 IU/mL</p> <p>Not tested</p> <p>AU2700/5400 Bilirubin: Interference less than 3% up to 40 mg/dL Bilirubin Hemolysate: Interference less than 4% up to 500 mg/dL Hemolysate Lipemia: Interference less than 4% up to 1000 mg/dL Intralipid</p> <p>RF: Interferences less than 4% up to 600 IU/mL</p> <p>Not tested</p>	<p>AU400/400^e Bilirubin: Interference less than 2% up to 40 mg/dL Bilirubin Hemolysate: Interference less than 2% up to 500 mg/dL Hemolysate Lipemia: Interference less than 10% up to 600 mg/dL Intralipid</p> <p>Not Specified</p> <p>Ascorbic Acid: Interferences less than 2% up to 20 mg/dL Ascorbate</p> <p>AU600/640/640^e Bilirubin: Interference less than 5% up to 40 mg/dL Bilirubin Hemolysate: Interference less than 2% up to 500 mg/dL Hemolysate Lipemia: Interference less than 7% up to 1000 mg/dL Intralipid</p> <p>Not Specified</p> <p>Ascorbic Acid: Interferences less than 2% up to 20 mg/dL Ascorbate</p> <p>AU2700/5400 Bilirubin: Interference less than 5% up to 40 mg/dL Bilirubin Hemolysate: Interference less than 3% up to 500 mg/dL Hemolysate Lipemia: Interference less than 10% up to 1000 mg/dL Intralipid</p> <p>Not Specified</p> <p>Ascorbic Acid: Interferences less than 3% up to 20 mg/dL Ascorbate</p>
Prozone Capacity	No high dose effect at IgA concentrations up to 10,000 mg/dL	No high dose effect at IgA concentrations up to 3,200 mg/dL



FEB 11 2008

Food and Drug Administration
2098 Gaither Road
Rockville MD 20850

Olympus America, Inc.
c/o Ms. Stephanie Donnelly
Regulatory Affairs/Quality Assurance Manager
Olympus Life Science Research Europa GmbH
Lismeehan, O, Callaghan's Mills
Co. Claire, Ireland.

Re: k073489

Trade/Device Name: Olympus IgA reagent
Regulation Number: 21 CFR 866.5510
Regulation Name: Immunoglobulins A, G, M, D, E immunological test systems
Regulatory Class: Class II
Product Code: CFN
Dated: December 11, 2007
Received: December 12, 2007

Dear Ms. Donnelly:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820). This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The

Page 2 –

FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific information about the application of labeling requirements to your device, or questions on the promotion and advertising of your device, please contact the Office of In Vitro Diagnostic Device Evaluation and Safety at (240) 276-0450. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding postmarket surveillance, please contact CDRH's Office of Surveillance and Biometric's (OSB's) Division of Postmarket Surveillance at 240-276-3474. For questions regarding the reporting of device adverse events (Medical Device Reporting (MDR)), please contact the Division of Surveillance Systems at 240-276-3464. You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (240) 276-3150 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,



Robert L. Becker, Jr., M.D., Ph.D.

Director

Division of Immunology and Hematology Devices

Office of In Vitro Diagnostic Device Evaluation and Safety

Center for Devices and Radiological Health

Enclosure

Indication for Use

510(k) Number (if known): K073489

Device Name: The Olympus IgA reagent (OSR6X171).

Indication For Use: System reagent for the quantitative determination of IgA immunoglobulins in human serum and plasma on OLYMPUS analyzers.

The spectrum of abnormalities in serum immunoglobulin concentrations is broad. Abnormal concentrations range from a virtual absence of one or more of the three major classes of immunoglobulin (IgG, IgA, and IgM) to polyclonal increases in one or more immunoglobulins. Measurement of these immunoglobulins aids in the diagnosis of abnormal protein metabolism and the body's lack of ability to resist infectious agents.

For *in vitro* diagnostic use.


Prescription Use (21 CFR Part 801 Subpart D)

And/Or

Over the Counter Use (21 CFR Part 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE; CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of In Vitro Diagnostic Device Evaluation and Safety (OIVD)



Division Sign-Off
Office of In Vitro Diagnostic Device
Evaluation and Safety

510(k) K073489