

## GENERAL INFORMATION

JUN 26 2008

**1. Applicant**

Date: May 05<sup>th</sup>, 2008

Name: VIASYS Healthcare GmbH (owned by **Cardinal Health**)

Address: Leibnizstrasse 7  
D-97204 Hoechberg  
Germany

**Contact person in Germany:****(Official Correspondent)**

Address

**Thomas Rust**VIASYS Healthcare GmbH  
Leibnizstrasse 7, 97204 Hoechberg  
Germany

Phone:

+49 931 4972 - 383

FAX:

+49 931 4972 - 62383

E-mail

[Thomas.Rust@viasyshc.com](mailto:Thomas.Rust@viasyshc.com)or at: [Thomas.Rust@cardinalhealth.com](mailto:Thomas.Rust@cardinalhealth.com)**Contact person in USA:****(US Agent)**

Address

**Yvette Lloyd**Cardinal Health  
22745 Savi Ranch Parkway  
Yorba Linda, CA 92887

Phone/Fax:

(714) - 919 - 3247

E-mail

[Yvette.Lloyd@cardinalhealth.com](mailto:Yvette.Lloyd@cardinalhealth.com)**2. Trade Name**FlowScreen  
FlowScreen ECG  
FlowScreen CT**3. Classification Name**Pulmonary calculator (21 CFR 868.1890, Product Code BTY)  
Electrocardiograph (21 CFR 870.2340, Product Code DPS)**4. Establishment Registration Number**

9615102

**5. Facility Address**VIASYS Healthcare GmbH  
Leibnizstrasse 7  
D-97204 Hoechberg  
Germany

## 6. Section 513 Device Classification

### 6.1 Classification

This is a Class II device

### 6.2 Classification Panel

Panel 73, Anesthesiology

Code BTY

Panel 74, Circular System Devices, ECG

Code DPS

## 7. Reason for Premarket Notification

New option for FlowScreen K062011

(Combination of two VIASYS devices, K062011 + K070614)

## 8. Predicate Devices Descriptions

### 8.1 Name

a) FlowScreen

b) CorScreen

### 8.2 Predicate Device Company

Viasys Healthcare GmbH

### 8.3 Predicate Device 510(k)#

a) K062011

b) K070614

## 9. Device Description

**FlowScreen is an active medical device providing following characteristics:**

Mains operation

Colour LCD display for user interface

Alphanumerical keyboard

Colour ink-printer for printout of reports in US-letter and DIN A4 size

Patient information and measurements are stored in an internal database

Data can be stored on an SD memory card

### a) pulmonary functions

- Measurement with ultrasonic handle or pneumotach handle
- Slow spirometry (VCin, VCex, VCmax, ERV, IC, VT, IRV, MV, BF, TI, TE, ...)
- Forced spirometry (FVCin, FVC, FEV1, PEF, FEV1/FVC, FEF 50, FEF 75, PIF, ...)
- Flow-Volume and Volume-Time Loop, pre/post tests
- MVV measurement
- Trending capabilities
- Patient Incentive animations
- Interpretation modules

## b) ECG functions

- Simultaneous acquisition of the 12 standard leads
- Storage of 10 seconds of acquired ECG signal
- Digital filters for base-line drift and mains interference suppression
- Interpretation program Hanover ECG System (HES) providing the following additional information:
  - Representatives templates of each lead including markers on fiducial points
  - Summary of mean measurements
  - Rhythm Analysis statements
  - Signal noise detection and information
  - Specific findings on QRS complex
  - Conduction statements
  - QRS T diagnostic statements
  - Arrhythmia monitoring detection
  - Heart Rate Variability

## 10. Intended Use Statement

The FlowScreen / FlowScreen ECG is a diagnostic system for recording and assessing inspiratory and expiratory pulmonary function (spirometry).

In addition it is intended for measuring a 3/6- or 12-channel surface electrocardiogram (ECG) of a patient. The acquired ECG can be recorded and displayed on a screen or printed on paper. 12-channel ECGs are analysed automatically and suggestions for the interpretation of the 12-channel ECG can be made by the software.

FlowScreen / FlowScreen ECG can be used for non interpretive applications for patients with an age of 4 years and older and a weight of 20 kg or higher. FlowScreen / FlowScreen ECG is intended for use in routine ECG recording by trained physicians in the office or hospital. FlowScreen / FlowScreen ECG is not intended for intracardial use. Automatic interpretation of the ECG is not possible for pediatric patients with an age below 16 years and for pacemaker patients.

FlowScreen CT (Clinical Trial version) incorporates the identical measurements, but individual access rights are defined for different user roles (e.g. Investigator, doctor, study nurse, trainer and service personnel).

The interpretation software is intended to support the physician in evaluation the ECG in terms of morphology and rhythm.

A qualified physician has to reassess all FlowScreen measurements. An interpretation by FlowScreen is only significant if it is considered in connection with other clinical findings. ECG interpretation statements made by the FlowScreen represent partial qualitative and quantitative information on the patient's cardiovascular condition and no therapy or drugs can be administered based solely on the interpretation statements.

The FlowScreen is powered from 100 - 240V / 50 - 60Hz wall outlets. No energy is transferred to the patient.

Federal U.S. law restricts this device to sale by or on the order of a physician.

## 11. Required Components

FlowScreen – monitor  
Ultrasonic handle or pneumotach handle  
ECG-amplifier  
Disposable ECG-electrodes  
User manual

## 12. Summary table of comparison

<b>Pulmonary Function</b>				
#	Parameter	New Device: FlowScreen	Predicate Device: FlowScreen K062011	Result
1	Intended Use	Diagnostic Spirometry (VCin, VCex, etc.)	Diagnostic Spirometry (VCin, VCex, etc.)	identical
2	Patient population	To be used as a screening device to determine whether or not a patient requires further diagnosis for pulmonary function disorders	To be used as a screening device to determine whether or not a patient requires further diagnosis for pulmonary function disorders	identical
3	Performance	Slow Spirometry Forced Spirometry Flow-Volume MVV Pre/Post Tests Trending capabilities Interpretation modules	Slow Spirometry Forced Spirometry Flow-Volume MVV Pre/Post Tests Trending capabilities Interpretation modules	identical
4	Patient user interface	Ultrasonic handle or pneumotach handle	Ultrasonic handle or pneumotach handle	identical
5	Material of patient user interface	<u>Pneumotach handle:</u> Pneumotach (Ultem 1010R, handle (Romira ABS 1001 FRV0), mouthpiece (Polypropylene RG835MO) <u>Ultrasonic handle:</u> handle (Luran S778 TE), mouthpiece (HDPE Eraclene MS 80U)	<u>Pneumotach handle:</u> Pneumotach (Ultem 1010R, handle (Romira ABS 1001 FRV0), mouthpiece (Polypropylene RG835MO) <u>Ultrasonic handle:</u> handle (Luran S778 TE), mouthpiece (HDPE Eraclene MS 80U)	identical
6	Patient contacting accessories	Mouthpiece Nose clip Nose pads	Mouthpiece Nose clip Nose pads	identical
7	Material of patient contacting accessories	Nose clip: Polyacetal Nose pads: Ethylene Vinyl Acetate	Nose clip: Polyacetal Nose pads: Ethylene Vinyl Acetate	identical
8	Dimensional specification	455 x 280 x 380 (W x H x D)	455 x 280 x 380 (W x H x D)	Identical
9	Software	Data acquisition Data calculation Calc. of predicted values Data analysis Data interpretation Data storage Data output Data input	Data acquisition Data calculation Calc. of predicted values Data analysis Data interpretation Data storage Data output Data input	Identical

<b>ECG Function</b>				
<b>#</b>	<b>Parameter</b>	<b>New Device: FlowScreen</b>	<b>Predicate Device: CorScreen K070614</b>	<b>Result</b>
1	Intended Use	3/6- or 12-channel surface ECG recording device	3/6- or 12-channel surface ECG recording device	identical
2	Input dynamic range	+/- 300mV @ DC	+/- 300mV @ DC	identical
3	Frequency response Bandwidth	0,05 – 150 Hz / According to EC11 and IEC 60601-2-51	0,05 – 150 Hz / According to EC11 and IEC 60601-2-51	identical
4	A/D conversion	24 bits	24 bits	identical
5	Leads	12 Standard	12 Standard	identical
6	Paper Speed	25 50 mm/s +/-5% According to EC11	25 50 mm/s +/-5% According to EC11	identical
7	Recorder Sensitivity	5 10 20 mm/s According to EC11	5 10 20 mm/s According to EC11	identical
8	Writing System	Ink-printer US-letter and DIN-A4 size	Ink-printer US-letter and DIN-A4 size	Identical
9	Printed Channels	1/2/6/12	1/2/6/12	Identical
10	Paper	US-Letter and DIN-A4	US-Letter and DIN-A4	Identical
11	Mode of operation	Manual	Manual	Identical
12	Input/output	SD Memory card	SD Memory card	Identical
13	<b>Display</b>			
14	Size	320 x 240 pixels	320 x 240 pixels	Identical
15	No. of displayed channels	1/3/6/12	1/3/6/12	Identical
16	Trace speeds	5 10 25 50 mm/s	5 10 25 50 mm/s	Identical
17	Sensitivity	5 10 20 40 mm/mV	5 10 20 40 mm/mV	Identical

<b>Hardware platform</b>		
<b>FlowScreen</b>	<b>FlowScreen K062011</b>	<b>CorScreen K070614</b>
Power supply	Identical	Identical
Mainboard	Identical	Identical
Connector board	Identical	Identical
Interfaces	Identical	Identical
Color printer	Identical	Identical
Keyboard	Identical	Identical
Color display	Identical	Identical
Enclosure	Identical	Identical
<b>Hardware patient user interfaces</b>		
Patient user interface pneumotach handle	Identical	-----
Patient user interface USS-handle	Identical	-----
ECG Amplifier	-----	Identical
<b>Accessories</b>		
Mouthpiece	Identical	-----
Nose clip	Identical	-----
ECG electrodes	-----	Identical
<b>Software / Firmware / Operating System</b>		
Boot loader (u-boot)	Identical	Identical
Operating system (Linux)	Identical	Identical
<b>Firmware:</b>		
Base module	Identical	Identical
Spirometry module	Identical	-----
ECG module	-----	Identical

**13. Summary of non-clinical performance tests**

The following practices were followed and monitored for development of the FlowScreen / FlowScreen ECG / FlowScreen CT:

The risk analysis method used to assess the impact of the FlowScreen with the new option ECG was a Failure Modes and Effects Analysis (FMEA). The risk analysis and risk control document number is 651003.17.

The design validation tests that were performed as a result of this risk analysis assessment and functional specifications are documented in document "Validation test script" and "validation test log" with document number 671003.32 (Test Script), 671003.33 (Test Log), 671003.40 (Test Script Workflow ECG) and 671003.41 (Test Log Workflow ECG).

The safety test procedures demonstrate satisfaction of all safety requirements and mitigation of all identified hazards. The document numbers for the performed safety tests according IEC 60601-1 are 407.267.2 and 71320123. The document number according IEC60601-2-25 is 407.267.3.

The EMC testing was performed according EN60601-2 with the document numbers 266.182 and 267.281.

Conclusion:

Based on the above, VIASYS HEALTHCARE GMBH concludes, that FlowScreen is determined substantially equivalent to the legally marketed predicate VIASYS devices and is safe and effective for its intended use, and performs at least as well as the predicate devices. The FlowScreen / FlowScreen ECG / FlowScreen CT is a combination of FlowScreen K062011 and CorScreen K070614 on the identical hardware base.



Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

**JUN 26 2008**

Mr. Thomas Rust  
Regulatory Affairs Manager  
VIASYS Healthcare GmbH  
Leibnizstrasse 7  
D-97204 Hoechberg  
GERMANY

Re: K080734  
Trade/Device Name: FlowScreen  
FlowScreen ECG  
FlowScreen CT  
Regulation Number: 21 CFR 868.1890  
Regulation Name: Predictive Pulmonary Function Value Calculator  
Regulatory Class: II  
Product Code: BTY,DPS  
Dated: May 7, 2008  
Received: May 12, 2008

Dear Mr. Rust:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0120. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,



Chiu Lin, Ph.D.

Director

Division of Anesthesiology, General Hospital,

Infection Control and Dental Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

# Indications for Use

510(k) Number (if known): K080734

Device Name: FlowScreen  
FlowScreen ECG  
FlowScreen CT

## Indications for Use:

The FlowScreen / FlowScreen ECG is a diagnostic system for recording and assessing inspiratory and expiratory pulmonary function (spirometry).

In addition it is intended for measuring a 3/6- or 12-channel surface electrocardiogram (ECG) of a patient. The acquired ECG can be recorded and displayed on a screen or printed on paper. 12-channel ECGs are analysed automatically and suggestions for the interpretation of the 12-channel ECG can be made by the software.

FlowScreen / FlowScreen ECG can be used for non interpretive applications for patients with an age of 4 years and older and a weight of 20 kg or higher. FlowScreen / FlowScreen ECG is intended for use in routine ECG recording by trained physicians in the office or hospital. FlowScreen / FlowScreen ECG is not intended for intracardial use. Automatic interpretation of the ECG is not possible for pediatric patients with an age below 16 years and for pacemaker patients.

FlowScreen CT (Clinical Trial version) incorporates the identical measurements, but individual access rights are defined for different user roles (e.g. Investigator, doctor, study nurse, trainer and service personnel).

The interpretation software is intended to support the physician in evaluation the ECG in terms of morphology and rhythm.

A qualified physician has to reassess all FlowScreen measurements. An interpretation by FlowScreen is only significant if it is considered in connection with other clinical findings. ECG interpretation statements made by the FlowScreen represent partial qualitative and quantitative information on the patient's cardiovascular condition and no therapy or drugs can be administered based solely on the interpretation statements.

The FlowScreen is powered from 100 - 240V / 50 - 60Hz wall outlets. No energy is transferred to the patient.

Federal U.S. law restricts this device to sale by or on the order of a physician.


Prescription Use  X   
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use \_\_\_\_\_  
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

  
(Division Sign-Off)

Division of Anesthesiology, General Hospital  
Infection Control, Dental Devices

510(k) Number:  K080734

Page 1 of \_\_\_\_\_