

**Traditional 510(k) for Kimberly-Clark\* STERLING\* Nitrile & STERLING\* Nitrile-Xtra\* Powder-Free Exam Gloves with a Chemotherapy Drug Use Claim****Section 5. 510(k) SUMMARY****SUMMARY OF SAFETY AND EFFECTIVENESS INFORMATION**

Submitter's Name:	Kimberly-Clark Corporation
Submitter's Address:	1400 Holcomb Bridge Road Roswell, GA 30076-2199
Submitter's Phone No:	770-587-8208
Submitter's Fax No:	920-969-5160
Date of Preparation:	July 10, 2008
Name of Device	
Trade Name:	<ul style="list-style-type: none"><li>• Kimberly-Clark* STERLING* Nitrile Powder-Free Exam Glove with a Chemotherapy Drug Use Claim</li><li>• Kimberly-Clark* STERLING* Nitrile-Xtra* Powder-Free Exam Glove with a Chemotherapy Drug Use Claim</li></ul>
Common Name:	Patient examination glove
Classification Name:	Glove, Patient Examination, Specialty – 80 LZC
Legally marketed device to which equivalency is claimed:	<ul style="list-style-type: none"><li>• Kimberly-Clark* STERLING* Nitrile Powder-Free Exam Glove - K051347</li><li>• Perusahaan Getah Asas Sdn. Bhd. Powdered Free Patient Examination Gloves, Blue Colored Non-Sterile (Low Dermatitis Potential and Chemotherapy Drug Protection Labeling Claims) – K042805</li></ul>
Description of the device:	Light gray nitrile, chlorinated, powder-free, textured fingertip, ambidextrous patient examination glove that meets all of the requirements of ASTM D 6319-00a, <i>Standard Specification for Nitrile Examination Gloves for Medical Application</i>
Intended use of device:	The Kimberly-Clark* STERLING* Nitrile and Nitrile-Xtra* Powder-Free Exam Gloves are disposable devices intended for medical purposes that are worn on the examiner's hand to prevent contamination between patient and examiner.

**Traditional 510(k) for Kimberly-Clark\* STERLING\* Nitrile & STERLING\* Nitrile-Xtra\*  
Powder-Free Exam Gloves with a Chemotherapy Drug Use Claim**

Summary of technological characteristics compared to predicate device:	There are no different technological characteristics compared to the predicate devices. They are all powder-free non-sterile nitrile exam gloves, one predicate a gray color and the other a blue color. The Chemotherapy Drug Use Claim is similar to that of Perusahaan Getah Asas Sdn. Bhd. Powdered Free Patient Examination Gloves, Blue Colored Non-Sterile (Low Dermatitis Potential and Chemotherapy Drug Protection Labeling Claims) – K042805		
Brief description of Non-Clinical Tests:	<u>Non-Clinical Tests</u>	<u>Standard</u>	<u>Performance</u>
	Dimensions	ASTM D 6319-00a	Meets
	Physical Properties	ASTM D 6319-00a	Meets
	Freedom from pinholes	ASTM D 6319-00a	Meets
		ASTM D 5151-06	Meets
	Powder Free	ASTM D 6124-06	Meets
		ASTM D 6319-00a	Meets
	ISO Skin Irritation Study	ISO 10993, Part 10	Meets
	Murine Local Lymph Node Assay	ISO 10993, Part 10	Meets
	ISO Systemic Toxicity Study	ISO 10993, Part 11	Meets
	Resistance to Permeation (Protective Clothing)	ASTM F 739-07	See data below
	Resistance to Permeation (Medical Gloves)	ASTM D 6978-05	See data below
	<b>Tested Chemotherapy Drug and Concentration</b>	<b>Average Breakthrough Detection Time (minutes)</b>	
	Cyclophosphamide (20.0 mg/ml)	No breakthrough up to 240 minutes	
Doxorubicin HCl (2.0 mg/ml)	No breakthrough up to 240 minutes		
Etoposide (20.0 mg/ml)	No breakthrough up to 240 minutes		
5-Fluorouracil (50.0 mg/ml)	No breakthrough up to 240 minutes		
Paclitaxel (Taxol) 6.0 mg/ml)	No breakthrough up to 240 minutes		
ThioTEPA (10.0 mg/ml)	Avg. minutes before breakthrough = 54.2		
Cisplatin (1.0 mg/ml)	No breakthrough up to 240 minutes		
Dacarbazine (10.0 mg/ml)	No breakthrough up to 240 minutes		
Ifosfamide (50.0 mg/ml)	No breakthrough up to 240 minutes		
Mitoxantrone (2.0 mg/ml)	No breakthrough up to 240 minutes		
Vincristine sulfate (1.0 mg/ml)	No breakthrough up to 240 minutes		
Brief description of Clinical Tests:	No new clinical tests were required to support this 510(k) application.		
Conclusions drawn from the Non-Clinical and Clinical Tests:	Non-clinical laboratory and animal based biocompatibility test data confirm the Kimberly-Clark* STERLING* Nitrile and Nitrile-Xtra* Powder-Free Exam Gloves Kimberly-Clark* STERLING* meets all applicable performance and biocompatibility requirements.		
Other Information deemed necessary by the FDA:	None		



Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

JUL 14 2008

Mr. Richard V. Wolfe  
Associate Director, Regulatory Affairs  
Kimberly-Clark Corporation  
1400 Holcomb Bridge Road  
Roswell, Georgia 30076

Re: K081089

Trade/Device Name: Kimberly-Clark\* STERLING\* Nitrile and STERLING\* Nitrile-Xtra\* Powder-Free Exam Gloves with a Chemotherapy Drug Use Claim

Regulation Number: 21 CFR 880.6250

Regulation Name: Patient Examination Glove

Regulatory Class: I

Product Code: LZC

Dated: June 27, 2008

Received: July 2, 2008

Dear Mr. Wolfe:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Center for Devices and Radiological Health's (CDRH's) Office of Compliance at (240) 276-*[See Below For Phone Numbers]*. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding postmarket surveillance, please contact CDRH's Office of Surveillance and Biometric's (OSB's) Division of Postmarket Surveillance at 240-276-3474. For questions regarding the reporting of device adverse events (Medical Device Reporting (MDR)), please contact the Division of Surveillance Systems at 240-276-3464. You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (240) 276-3150 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,

*Anthony D. Watson for*  
Chiu Lin, Ph.D.  
Director  
Division of Anesthesiology, General Hospital,  
Infection Control and Dental Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure



INDICATIONS FOR USE

**Applicant:** Kimberly-Clark Corporation

**510(k) Number:** K081089

**Device Name:** Kimberly-Clark\* STERLING\* Nitrile and STERLING\* Nitrile-Xtra\* Powder-Free Exam Gloves with a Chemotherapy Drug Use Claim

**Indications for Use:** Based upon 21CFR§880.6250 "Patient examination glove"


A patient examination glove is a medical device intended for medical purposes that is worn on the examiner's hand to prevent contamination between patient and examiner.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

---

Concurrence of CDRH Office of Device Evaluation (ODE)

Prescription Use \_\_\_\_\_ OR Over-The-Counter   X    
Per 21CFR 801.109 Subpart D Per 21CFR 801.109 Subpart C

  
\_\_\_\_\_  
(Division Sign-Off)  
Division of Anesthesiology, General Hospital  
Infection Control, Dental Devices

510(k) Number:   K 081089