

JAN 16 2009

**510(k) Summary – COBAS INTEGRA 400/800 Bilirubin Total**

---

**Introduction** According to the requirements of 21 CFR 807.92, the following information provides sufficient detail to understand the basis for a determination of substantial equivalence

---

**Submitter name, address, contact** Roche Diagnostics  
9115 Hague Rd  
Indianapolis IN 46250  
(317) 521-4569

Contact person Jennifer Tribbett

Date prepared April 23, 2008

---

**Device Name** Proprietary name Cobas Integra 400/800 Bilirubin Total  
Common name Bilirubin (total or direct) test system  
Classification name Bilirubin (total or direct) test system

---

**Device Description** The Cobas Integra Bilirubin Total reagent is intended for use with the Cobas Integra systems for the quantitative determination of the total bilirubin concentration in human serum

---

**Intended use** In vitro test for the quantitative determination of the total bilirubin concentration in human serum on Cobas Integra Bilirubin is formed in the reticuloendothelial system during the degradation of aged erythrocytes. The heme portion from hemoglobin and from other heme-containing proteins is removed, metabolized to bilirubin and transported as a complex with serum albumin to the liver. In the liver, bilirubin is conjugated with glucuronic acid for solubilization and subsequent transport through the bile duct and elimination via the digestive tract. Diseases or conditions which, through hemolytic processes, produce bilirubin faster than the liver can metabolize it, cause the levels of unconjugated (indirect) bilirubin to increase in the circulation. Liver immaturity and several other diseases in which the bilirubin conjugation mechanism is impaired cause similar elevations of circulating unconjugated bilirubin. Bile duct obstruction or damage to hepatocellular structure causes increases in the levels of both conjugated (direct) and unconjugated (indirect) bilirubin in the circulation.

---

*Continued on next page*

**Predicate Device**

We claim substantial equivalence to the Cobas Integra Bilirubin Total (K951595)

**Substantial equivalency – device comparison**

The table below indicates the similarities and differences between the modified Bilirubin Total and the predicate Bilirubin Total (K951595)

<b>Characteristic</b>	<b>Predicate Cobas Integra Bilirubin Total (K951595)</b>	<b>Modified Cobas Integra Bilirubin Total</b>
<b>Intended Use</b>	The Cobas Integra Cassette Bilirubin Total (BIL-T) contains an in vitro diagnostic reagent system intended for use on Cobas Integra for the quantitative determination of the total bilirubin concentration in serum and plasma (test BIL-T, 0-048)	In vitro test for the quantitative determination of the total bilirubin concentration in human serum (test BIL-T, 0-048) on Cobas Integra
<b>Indications for Use</b>	Bilirubin is formed in the reticuloendothelial system during the degradation of aged erythrocytes. The heme portion from hemoglobin and from other heme-containing proteins is removed, metabolized to bilirubin and transported as a complex with serum albumin to the liver. In the liver, bilirubin is conjugated with glucuronic acid for solubilization and subsequent transport through the bile duct and elimination via the digestive tract. Diseases or conditions which, through hemolytic processes, produce bilirubin faster than the liver can metabolize it, cause the levels of unconjugated (indirect) bilirubin to increase in the circulation. Liver immaturity and several other diseases in which the bilirubin conjugation mechanism is impaired cause similar elevations of circulating unconjugated bilirubin. Bile duct obstruction or damage to hepatocellular structure causes increases in the levels of both conjugated (direct) and unconjugated (indirect) bilirubin in the circulation.	Same

*Continued on next page*

Characteristic	Predicate Cobas Integra Bilirubin Total (K951595)	Modified Cobas Integra Bilirubin Total																																
Assay principle	Diazo method	Same																																
Active Ingredients	<table border="1"> <thead> <tr> <th></th> <th>R1</th> <th>R2</th> <th>Test</th> </tr> </thead> <tbody> <tr> <td>Sulfanilic Acid</td> <td>35</td> <td></td> <td>13.5 mmol/L</td> </tr> <tr> <td>Oxalic Acid</td> <td>40</td> <td></td> <td>15.4 mmol/L</td> </tr> <tr> <td>HEDTA</td> <td>4</td> <td></td> <td>1.5 mmol/L</td> </tr> <tr> <td>Myristyltrimethyl-ammonium bromide</td> <td>2</td> <td></td> <td>0.8 %</td> </tr> <tr> <td>Dodecyltrimethyl-ammonium bromide</td> <td>5</td> <td></td> <td>1.9 %</td> </tr> <tr> <td>Sodium nitrite</td> <td></td> <td>24.5</td> <td>3.1 mmol/L</td> </tr> <tr> <td>pH</td> <td>1.1</td> <td>6.8</td> <td>1.4</td> </tr> </tbody> </table>		R1	R2	Test	Sulfanilic Acid	35		13.5 mmol/L	Oxalic Acid	40		15.4 mmol/L	HEDTA	4		1.5 mmol/L	Myristyltrimethyl-ammonium bromide	2		0.8 %	Dodecyltrimethyl-ammonium bromide	5		1.9 %	Sodium nitrite		24.5	3.1 mmol/L	pH	1.1	6.8	1.4	Same
	R1	R2	Test																															
Sulfanilic Acid	35		13.5 mmol/L																															
Oxalic Acid	40		15.4 mmol/L																															
HEDTA	4		1.5 mmol/L																															
Myristyltrimethyl-ammonium bromide	2		0.8 %																															
Dodecyltrimethyl-ammonium bromide	5		1.9 %																															
Sodium nitrite		24.5	3.1 mmol/L																															
pH	1.1	6.8	1.4																															
Measuring range	<p>Test range 0-340 umol/L (0-20 mg/dL) w/postdilution 0-1020 umol/L (0-60 mg/dL) Postdilution factor 3</p> <p>Sensitivity The sensitivity is defined as the change of analytical response (<math>\Delta A</math>) per unit change in analyte concentration at a pathlength of 1 cm. The sensitivity is <math>2.9 \times 10^{-3} \Delta A</math> per umol/L of total bilirubin (<math>5.0 \times 10^{-2} \Delta A</math> per mg/dL of bilirubin)</p>	<p>1.7 – 340 umol/L (0.099-20 mg/dL) Determine samples having higher concentrations via the rerun function. Dilution of samples via the rerun function is a 1:3 dilution. Results from samples diluted by the rerun function are automatically multiplied by a factor of 3.</p> <p>Lower Detection Limit 1.7 umol/L (0.099 mg/dL) The detection limit represents the lowest measurable analyte level that can be distinguished from zero. It is calculated as the value lying three standard deviations above that of a zero sample (zero sample + 3 SD, within-run precision, n=30)</p>																																

*Continued on next page*

Characteristic	Predicate Cobas Integra Bilirubin Total (K951595)	Modified Cobas Integra Bilirubin Total																											
Precision	<p>Precision was evaluated on Cobas Integra using two human serum pools and following the guidelines of the NCCLS Manual EP5-T2</p> <table border="1" data-bbox="370 463 899 655"> <thead> <tr> <th></th> <th>Level 1</th> <th>Level 2</th> </tr> </thead> <tbody> <tr> <td>Mean</td> <td>24.2 umol/L (1.4 mg/dL)</td> <td>72.5 umol/L (4.2 mg/dL)</td> </tr> <tr> <td>CV w/in run</td> <td>0.46%</td> <td>0.45%</td> </tr> <tr> <td>CV day/day</td> <td>0.48%</td> <td>0.53%</td> </tr> <tr> <td>CV total</td> <td>0.78%</td> <td>0.80%</td> </tr> </tbody> </table>		Level 1	Level 2	Mean	24.2 umol/L (1.4 mg/dL)	72.5 umol/L (4.2 mg/dL)	CV w/in run	0.46%	0.45%	CV day/day	0.48%	0.53%	CV total	0.78%	0.80%	<p>Reproducibility was determined using human samples and controls in an internal protocol (within-run n=20, between-run n=20). The following results were obtained</p> <table border="1" data-bbox="922 497 1442 655"> <thead> <tr> <th></th> <th>Level 1</th> <th>Level 2</th> </tr> </thead> <tbody> <tr> <td>Mean</td> <td>24.2 umol/L (1.4 mg/dL)</td> <td>72.5 umol/L (4.2 mg/dL)</td> </tr> <tr> <td>CV w/in run</td> <td>0.46%</td> <td>0.45%</td> </tr> <tr> <td>CV total</td> <td>0.78%</td> <td>0.80%</td> </tr> </tbody> </table>		Level 1	Level 2	Mean	24.2 umol/L (1.4 mg/dL)	72.5 umol/L (4.2 mg/dL)	CV w/in run	0.46%	0.45%	CV total	0.78%	0.80%
	Level 1	Level 2																											
Mean	24.2 umol/L (1.4 mg/dL)	72.5 umol/L (4.2 mg/dL)																											
CV w/in run	0.46%	0.45%																											
CV day/day	0.48%	0.53%																											
CV total	0.78%	0.80%																											
	Level 1	Level 2																											
Mean	24.2 umol/L (1.4 mg/dL)	72.5 umol/L (4.2 mg/dL)																											
CV w/in run	0.46%	0.45%																											
CV total	0.78%	0.80%																											
Accuracy	<p>Accuracy Total bilirubin values for human sera and plasma samples obtained on Cobas Integra with the cassette bilirubin total were compared to those determined with reagents for total bilirubin on Cobas Mira and the commercially available alternative clinical chemistry system. Samples were measured in duplicate. Sample size (n) represents all replicates. Values ranged from 1.7 to 350.6 umol/L (0.1 to 20.5 mg/dL)</p> <p><u>Cobas Mira</u>  Sample size (n) 208  Corr Coefficient (r) 0.999  (r<sub>s</sub>) 0.985  Lin Regression y=1.00x + 1.2 umol/L  Passing Bablock y=1.00x + 0.6 umol/L</p> <p><u>Alternative System</u>  Sample size (n) 210  Corr Coefficient (r) 0.999  (r<sub>s</sub>) 0.972  Lin Regression y=0.93x + 0.8 umol/L  Passing Bablock y=0.93x + 0.4 umol/L</p>	<p>Method Comparison Total bilirubin values for human serum samples obtained on a Cobas Integra 800 analyzer using the Cobas Integra Bilirubin Total reagent (y) were compared to those determined using Cobas Integra Total Bilirubin Special reagent on a Cobas Integra 800 analyzer (x)</p> <p>Cobas Integra 800 analyzer  Sample size (n) = 49  Passing/Bablok  y=0.956x + 2.113 umol/L  r=0.999</p> <p>Values ranged from 5.49 to 317 umol/L (0.321 to 18.5 mg/dL)</p> <p>Total bilirubin values for human serum samples obtained on a Cobas Integra 400 analyzer using the Cobas Integra Bilirubin Total reagent (y) were compared to those determined with commercially available reagents for total bilirubin on a Roche/Hitachi 911 analyzer (x)</p> <p>Roche/Hitachi 911 analyzer  Sample size (n) = 104  Linear regression  y=0.989x - 0.520 umol/L  r=0.999</p> <p>Passing/Bablok  y=0.991x + 0.219 umol/L</p> <p>Values ranged from 3.29 to 282 umol/L (0.192 to 16.5 mg/dL)</p>																											

Continued on next page

Characteristic	Predicate Cobas Integra Bilirubin Total (K951595)	Modified Cobas Integra Bilirubin Total
Expected Values	<p>Adults and infants &gt;1month 34-17 umol/L (0.2 - 1.0 mg/dL)</p> <p>Newborns (up to 24 h) 34-103 umol/L (2.0 - 6.0 mg/dL)</p> <p>Newborns (up to 48 h) 103-171 umol/L (6.0 - 10.0 mg/dL)</p> <p>Newborns (3 to 5 days) 68 - 137 umol/L (4.0 - 8.0 mg/dL)</p> <p>Note It is recommended that each laboratory establishes and maintains its own reference ranges and that the values given here are used as a guideline only</p>	<p>Same</p> <p>Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range</p>

*Continued on next page*

Characteristic	Predicate Cobas Integra Bilirubin Total (K951595)	Modified Cobas Integra Bilirubin Total
Limitations – interference	<p>Hemolysis Avoid hemolyzed specimens Even slight hemolysis interferes with the test</p> <p>Lipemia Avoid lipemic specimens Even slight lipemia interferes with the test</p> <p>Drugs Of the drugs tested in vitro, propranolol and theophylline cause artificially low total bilirubin values at the tested drug level</p>	<p>Hemolysis No significant interference up to an H index of 10 (approximate hemoglobin concentration 6 umol/L or 10 mg/dL)</p> <p>Lipemia (Intralipid) No significant interference up to an L index of 9 There is a poor correlation between the L index (corresponds to turbidity) and triglycerides concentration</p> <p>Drugs Therapeutic drug interference was tested according to the recommendations of VDPH No interference was found Exception Propranolol and theophylline cause artificially low total bilirubin values at the tested drug level Hydroxocobalamin (Cyanokit) may cause false-high results</p> <p>Other In very rare cases gammopathy, in particular type IgM (Waldenstroms macroglobulinemia) may cause unreliable results</p>



Food and Drug Administration  
2098 Gaither Road  
Rockville MD 20850

JAN 16 2009

Roche Diagnostics  
c/o Jennifer Tribbett  
Regulatory Affairs Principal  
9115 Hague Road  
Indianapolis, IN 46250

Re k081193  
Trade/Device Name Roche COBAS Integra Bilirubin Total (BIL-T)  
Regulation Number 21 CFR 862 1110  
Regulation Name Bilirubin (total or direct) test system  
Regulatory Class Class II  
Product Code CIG  
Dated December 18, 2008  
Received December 19, 2008

Dear Ms Tribbett

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in Title 21, Code of Federal Regulations (CFR), Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to, registration and listing (21 CFR Part 807), labeling (21 CFR Parts 801 and 809), and good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820).

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific information about the application of labeling requirements to your device, or questions on the promotion and advertising of your device, please contact the Office of In Vitro Diagnostic Device Evaluation and Safety at (301) 594-3084. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,



Courtney C. Harper, Ph.D.  
Acting Director  
Division of Chemistry and Toxicology  
Office of *In Vitro* Diagnostic Device  
Evaluation and Safety  
Center for Devices and Radiological Health

# Indication for Use

510(k) Number (if known) K081193

Device Name Roche COBAS Integra Bilirubin Total (BIL-T)

## Indication For Use

In vitro test for the quantitative determination of the total bilirubin concentration in human serum on COBAS Integra systems

Measurement of the level of bilirubin is used in the diagnosis and treatment of liver, hemolytic hematological, and metabolic disorders, including hepatitis and gall bladder block

Prescription Use XXX  
(21 CFR Part 801 Subpart D)

And/Or

Over the Counter Use \_\_\_\_\_  
(21 CFR Part 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE, CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of In Vitro Diagnostic Device Evaluation and Safety (OIVD)

Carol C. Benson  
Division Sign-Off  
Office of In Vitro Diagnostic Device  
Evaluation and Safety

510(k) K081193