

5. 510(k) Summary

Spirometry Option for Innocor

Date of Summary	12/12/2008								
Submitter/Contact Person	<p>H. Carl Jenkins FEB 23 2009</p> <p>The Wood Burditt Group</p> <p>1025 W. Everett Rd., Suite 100</p> <p>Lake Forest, IL 60045</p> <p>(ph) (847) 234-7500 x 205</p> <p>(fax) (847) 574-0728</p> <p>(email) hcjenkins@woodburditt.com</p>								
Applicant	<table border="0"> <tr> <td>Innovision A/S</td> <td>Phone: +45 65 95 91 00</td> </tr> <tr> <td>Lindvedvej 75</td> <td>Fax: +45 65 95 78 00</td> </tr> <tr> <td>DK-5260 Odense S</td> <td>info@innovision.dk</td> </tr> <tr> <td>Denmark</td> <td>www.innovision.dk</td> </tr> </table>	Innovision A/S	Phone: +45 65 95 91 00	Lindvedvej 75	Fax: +45 65 95 78 00	DK-5260 Odense S	info@innovision.dk	Denmark	www.innovision.dk
Innovision A/S	Phone: +45 65 95 91 00								
Lindvedvej 75	Fax: +45 65 95 78 00								
DK-5260 Odense S	info@innovision.dk								
Denmark	www.innovision.dk								
Device Name	Spirometry Option to Innocor								
Common Name	Diagnostic Spirometer								
Classification	<p><i>[Hemodynamic Measurements—Already Cleared K051907]</i></p> <p>Computer, diagnostic, programmable</p> <p>Regulation Number: 21 CFR §870.1425</p> <p>Product Code: DQK</p> <p>Panel Code: Cardiovascular</p> <p>Device Class: IIa</p> <p><i>[Cardiopulmonary Exercise Testing Option – Already Cleared K071911]</i></p> <p>Oxygen uptake computer</p> <p>Regulation Number: 21 CFR §868.1730</p> <p>Product Code: BZL</p> <p>Panel Code: Anesthesiology</p> <p>Device Class: IIa</p>								

5. 510(k) Summary

Spirometry Option for Innocor

	<p><i>[Spirometry Option]</i></p> <p>Diagnostic spirometer</p> <p>Regulation Number: 21 CFR §868.1840</p> <p>Product Code: BZG</p> <p>Panel Code: Anesthesiology</p> <p>Device Class: II</p>
<p>Legally Marketed Predicate Devices</p>	<p>The Spirometry Option for Innocor is substantially equivalent in respect to the intended use, design and method of operation to:</p> <p><i>Predicate Device No. 1</i></p> <p>Name: Innocor</p> <p>510(k) number: K051907</p> <p>Manufacturer: Innovision A/S, Denmark</p> <p><i>Predicate Device No. 2</i></p> <p>Name: Cardiopulmonary Exercise Testing Option</p> <p>510(k) number: K071911</p> <p>Manufacturer: Innovision A/S, Denmark</p> <p><i>Predicate Device No. 3</i></p> <p>Name: Spirobank G</p> <p>510(k) number: K072979</p> <p>Manufacturer: MIR Medical International Research</p>
<p>Device Description</p>	<p>Innocor is a compact point-of-care device intended to be used for non-invasive measurement of a) cardiac output (CO) and other hemodynamic parameters utilizing inert gas rebreathing (IGR) technology, and b) metabolic parameters including oxygen uptake by means of a breath-by-breath gas exchange method.</p> <p>The Cardiopulmonary Exercise Testing Option to Innocor</p>

provides measurements of gas exchange parameters including oxygen uptake (VO_2), carbon dioxide excretion (VCO_2), ventilation (V_E) and end-tidal gas concentrations plus a number of derived parameters. These parameters are determined by simultaneous measurements of the respiratory flow and gas concentrations when breathing ambient air. The respiratory flow is measured by means of a differential pressure type flowmeter (pneumotachometer) placed between the respiratory valve unit and the patient. The gas exchange calculations are carried out online for each breath between the rebreathing tests, providing the opportunity to perform an incremental exercise test on a bicycle ergometer or treadmill and measure the progress of cardiac function, pulmonary function and gas exchange at the same time.

Spirometry is a physiological test that measures how an individual inhales or exhales volumes of air as a function of time. Spirometry is recognized as a valuable screening test of general respiratory health.

The Spirometry Option for Innocor measures the subset of spirometric variables of a patient during a forced expiration testing procedure. These measured variables include FEV_1 (forced expiratory volume in 1 second), FVC (forced vital capacity), $\text{FEV}_1\%$, PEF (peak expiratory flow), MEF 75 (Maximal instantaneous forced expiratory flow where 75% of the FVC remains to be expired), MEF 50, MEF 25, FET (Forced expiratory time) and MVV (Maximum voluntary ventilation). These parameters are determined by measurements of the respiratory flow when breathing ambient air during a spirometry test of a patient (tidal breathing followed by a full inspiration and then finally a maximal forced expiration). The respiratory flow is measured by means of a pneumotachometer.

5. 510(k) Summary

Spirometry Option for Innocor

	<p>The Spirometry Option for Innocor, used in conjunction with the entire Innocor system, provides health care providers with a set of valuable diagnostic tools.</p>																														
<p>Intended Use and Indications</p>	<p>The Spirometry Option for Innocor is intended to be used as a diagnostic spirometer, used in pulmonary function testing, to measure the flow of gas moving in and out of a patient's lungs.</p> <p>In order to produce data regarding the maximum performance with respect to tidal volume and ventilation, the specific parameters measured by the Innocor Spirometry Option are:</p> <table border="1" data-bbox="570 767 1349 1908"> <thead> <tr> <th data-bbox="570 767 764 842"><i>Abbreviation</i></th> <th data-bbox="764 767 1211 842"><i>Name</i></th> <th data-bbox="1211 767 1349 842"><i>Unit</i></th> </tr> </thead> <tbody> <tr> <td data-bbox="570 842 764 957">FEV₁</td> <td data-bbox="764 842 1211 957">Forced expiratory volume in one second</td> <td data-bbox="1211 842 1349 957">L [BTPS]</td> </tr> <tr> <td data-bbox="570 957 764 1071">FVC</td> <td data-bbox="764 957 1211 1071">Forced vital capacity</td> <td data-bbox="1211 957 1349 1071">L [BTPS]</td> </tr> <tr> <td data-bbox="570 1071 764 1146">FEV₁%</td> <td data-bbox="764 1071 1211 1146">FEV₁ / FVC</td> <td data-bbox="1211 1071 1349 1146">%</td> </tr> <tr> <td data-bbox="570 1146 764 1261">PEF</td> <td data-bbox="764 1146 1211 1261">Peak expiratory flow</td> <td data-bbox="1211 1146 1349 1261">l/sec [BTPS]</td> </tr> <tr> <td data-bbox="570 1261 764 1414">MEF 75*</td> <td data-bbox="764 1261 1211 1414">Maximum instantaneous forced expiratory flow where 75% of the FVC remains to be expired</td> <td data-bbox="1211 1261 1349 1414">l/sec [BTPS]</td> </tr> <tr> <td data-bbox="570 1414 764 1567">MEF 50*</td> <td data-bbox="764 1414 1211 1567">Maximum instantaneous forced expiratory flow where 50% of the FVC remains to be expired</td> <td data-bbox="1211 1414 1349 1567">l/sec [BTPS]</td> </tr> <tr> <td data-bbox="570 1567 764 1720">MEF 25*</td> <td data-bbox="764 1567 1211 1720">Maximum instantaneous forced expiratory flow where 25% of the FVC remains to be expired</td> <td data-bbox="1211 1567 1349 1720">l/sec [BTPS]</td> </tr> <tr> <td data-bbox="570 1720 764 1795">FET</td> <td data-bbox="764 1720 1211 1795">Forced expiratory time</td> <td data-bbox="1211 1720 1349 1795">Sec</td> </tr> <tr> <td data-bbox="570 1795 764 1908">MVV</td> <td data-bbox="764 1795 1211 1908">Maximum voluntary ventilation</td> <td data-bbox="1211 1795 1349 1908">L/min [BTPS]</td> </tr> </tbody> </table>	<i>Abbreviation</i>	<i>Name</i>	<i>Unit</i>	FEV ₁	Forced expiratory volume in one second	L [BTPS]	FVC	Forced vital capacity	L [BTPS]	FEV ₁ %	FEV ₁ / FVC	%	PEF	Peak expiratory flow	l/sec [BTPS]	MEF 75*	Maximum instantaneous forced expiratory flow where 75% of the FVC remains to be expired	l/sec [BTPS]	MEF 50*	Maximum instantaneous forced expiratory flow where 50% of the FVC remains to be expired	l/sec [BTPS]	MEF 25*	Maximum instantaneous forced expiratory flow where 25% of the FVC remains to be expired	l/sec [BTPS]	FET	Forced expiratory time	Sec	MVV	Maximum voluntary ventilation	L/min [BTPS]
<i>Abbreviation</i>	<i>Name</i>	<i>Unit</i>																													
FEV ₁	Forced expiratory volume in one second	L [BTPS]																													
FVC	Forced vital capacity	L [BTPS]																													
FEV ₁ %	FEV ₁ / FVC	%																													
PEF	Peak expiratory flow	l/sec [BTPS]																													
MEF 75*	Maximum instantaneous forced expiratory flow where 75% of the FVC remains to be expired	l/sec [BTPS]																													
MEF 50*	Maximum instantaneous forced expiratory flow where 50% of the FVC remains to be expired	l/sec [BTPS]																													
MEF 25*	Maximum instantaneous forced expiratory flow where 25% of the FVC remains to be expired	l/sec [BTPS]																													
FET	Forced expiratory time	Sec																													
MVV	Maximum voluntary ventilation	L/min [BTPS]																													

5. 510(k) Summary

Spirometry Option for Innocor

	<p><i>*MEF 75 is equal to FEF 25 (maximal instantaneous forced expiratory flow where 25% of the FVC has been expired); MEF 50 is equal to FEF 50; MEF 25 is equal to FEF 75.</i></p>
Performance Testing	<p>The Spirometry Option for Innocor has been evaluated against the “Standardisation of Spirometry” document in the series “ATS/ERS Task Force: Standardisation of Lung Function Testing,” issued by The American Thoracic Society (ATS) and the European Respiratory Society (ERS). Performance data demonstrates that the hardware and software meet the ATS/ERS standards, and the Spirometry Option for Innocor is accordingly substantially equivalent to legally marketed predicate diagnostic spirometers.</p>



Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

The Wood Burditt Group LLC.
c/o Mr. H. Carl Jenkins
Regulatory Affairs Counsel
10 E. Scranton Ave., Suite 201
Lake Bluff, IL 60044

FEB 23 2009

Re: K083879
Spirometry Option for Innocor
Regulation Number: 21 CFR 870.1425
Regulation Name: Programmable Diagnostic Computer
Regulatory Class: Class II (two)
Product Code: DQK, BZG and BZL
Dated: December 12, 2008
Received: December 29, 2008

Dear Mr. Jenkins:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Page 2 - Mr. H. Carl Jenkins

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0120. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (240) 276-3150 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,



 Bram D. Zuckerman, M.D.
Director
Division of Cardiovascular Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

4. Indications for Use Statement

Spirometry Option for Innocor

Indications for Use

510(k) Number (if known): K083879

Device Name: Spirometry Option for Innocor

Indications for Use:

The Spirometry Option for Innocor is intended to be used as a diagnostic spirometer, used in pulmonary function testing, to measure the flow of gas moving into and out of a patient's lungs.

In order to produce data regarding the maximum performance with respect to tidal volume and ventilation, the specific parameters measured by the Innocor Spirometry Option are:

Abbreviation	Name	Unit
FEV ₁	Forced expiratory volume in one second	L [BTPS]
FVC	Forced vital capacity	L [BTPS]
FEV ₁ %	FEV ₁ / FVC	%
PEF	Peak expiratory flow	l/sec [BTPS]
MEF 75*	Maximum instantaneous forced expiratory flow where 75% of the FVC remains to be expired	l/sec [BTPS]
MEF 50*	Maximum instantaneous forced expiratory flow where 50% of the FVC remains to be expired	l/sec [BTPS]
MEF 25*	Maximum instantaneous forced expiratory flow where 25% of the FVC remains to be expired	l/sec [BTPS]
FET	Forced expiratory time	Sec
MVV	Maximum voluntary ventilation	L/min [BTPS]

*MEF 75 is equal to FEF 25 (maximal instantaneous forced expiratory flow where 25% of the FVC has been expired); MEF 50 is equal to FEF 50; MEF 25 is equal to FEF 75.

Prescription Use X
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use _____
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Danna R. Cochran
(Division Sign-Off)
Division of Cardiovascular Devices

4. Indications for Use Statement

510(k) Number K083879