

K081231

SEP 26 2008

**ARCHITECT *i*Phenobarbital**

**510(k) Summary (Summary of Safety and Effectiveness)**

This summary of the 510(k) safety and effectiveness information is being submitted in accordance with the requirements of SMDA 1990 and 21 CFR 807.92.

**Applicant Name:**

Carol Jochum  
Senior Regulatory Affairs Specialist  
Abbott Laboratories  
100 Abbott Park Road  
Abbott Park, IL 60064

**Device Name:**

**Reagents:**

Classification Name: Phenobarbital test system  
Trade Name: ARCHITECT *i*Phenobarbital Immunoassay  
Common Name: Phenobarbital test  
Governing Regulation: 862.3660  
Device Classification: Class II  
Classification Panel: Toxicology  
Product Code: DLZ

**Calibrators:**

Classification Name: Calibrator, drug specific  
Trade Name: ARCHITECT *i*Phenobarbital Calibrators (A-F)  
Common Name: Calibrator  
Governing Regulation: 862.3200  
Device Classification: Class II  
Classification Panel: Toxicology  
Product Code: DLJ

**Legally marketed device to which equivalency is claimed:**

AxSYM Phenobarbital (K940596)

**Intended Use of Device:**

The ARCHITECT *i* Phenobarbital assay is an *in vitro* chemiluminescent microparticle immunoassay (CMIA) for the quantitative measurement of phenobarbital, an anticonvulsant and sedative-hypnotic drug, in human serum or plasma on the ARCHITECT *i* System with *STAT* protocol capability. The measurements obtained are used in the diagnosis and treatment of phenobarbital overdose and in monitoring levels of phenobarbital to help ensure appropriate therapy.

**Description of Device:**

The ARCHITECT *i* Phenobarbital assay is a one-step *STAT* immunoassay for the quantitative measurement of phenobarbital in human serum or plasma using CMIA technology, with flexible assay protocols, referred to as Chemiflex. Sample, anti-phenobarbital coated paramagnetic microparticles, and phenobarbital acridinium-labeled conjugate are combined to create a reaction mixture. The anti-phenobarbital coated microparticles bind to phenobarbital present in the sample and to the phenobarbital acridinium-labeled conjugate. After washing, pre-trigger and trigger solutions are added to the reaction mixture. The resulting chemiluminescent reaction is measured as relative light units (RLUs). An indirect relationship exists between the amount of phenobarbital in the sample and the RLUs detected by the ARCHITECT *i* System optics.

**Comparison of Technological Characteristics:**

The ARCHITECT *i* Phenobarbital assay is a chemiluminescent microparticle immunoassay (CMIA) method for the quantitative measurement of phenobarbital, an anticonvulsant and sedative-hypnotic drug, in human serum and plasma. The AxSYM Phenobarbital assay utilizes fluorescence polarization immunoassay (FPIA) technology for the measurement of phenobarbital, an anticonvulsant and sedative-hypnotic drug, in serum or plasma.

**Summary of Non-Clinical Performance:**

The ARCHITECT *i*Phenobarbital assay is substantially equivalent to the AxSYM Phenobarbital assay in terms of precision, linearity and interferences as demonstrated in non-clinical performance data in this 510(k) submission.

**Summary of Clinical Performance:**

The ARCHITECT *i*Phenobarbital demonstrated substantially equivalent performance to the AxSYM Phenobarbital with a correlation coefficient of 1.0.



SEP 26 2008

Abbott Laboratories  
c/o Ms. Carol Jochum  
AP 6C-2, Dept. 049C  
100 Abbott Park Road  
Abbott Park, IL 60064

Re: k081231

Trade/Device Name: ARCHITECT *i*Phenobarbital Assay  
ARCHITECT *i*Phenobarbital Calibrators (A-F)  
Regulation Number: 21 CFR 862.3660  
Regulation Name: Phenobarbital test system  
Regulatory Class: Class II  
Product Code: DLZ, DLJ  
Dated: September 9, 2008  
Received: September 10, 2008

Dear Ms. Jochum:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in Title 21, Code of Federal Regulations (CFR), Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); and good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820).

Page 2 –

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific information about the application of labeling requirements to your device, or questions on the promotion and advertising of your device, please contact the Office of In Vitro Diagnostic Device Evaluation and Safety at (240) 276-0490. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (240) 276-3150 or at its Internet address at <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,

*Jean M. Cooper, M.S., D.V.M.*

Jean M. Cooper, M.S., D.V.M.

Director

Division of Chemistry and Toxicology

Office of *In Vitro* Diagnostic Device

Evaluation and Safety

Center for Devices and

Radiological Health

Enclosure

**Indication for Use**

510(k) Number (if known): K081231

Device Name: ARCHITECT *i*Phenobarbital

**Indication for Use:**

**Reagents**

The ARCHITECT *i*Phenobarbital assay is an *in vitro* chemiluminescent microparticle immunoassay (CMIA) for the quantitative measurement of phenobarbital, an anticonvulsant and sedative-hypnotic drug, in human serum or plasma on the ARCHITECT *i* System with *STAT* protocol capability. The measurements obtained are used in the diagnosis and treatment of phenobarbital overdose and in monitoring levels of phenobarbital to help ensure appropriate therapy.

**Calibrators**


The ARCHITECT *i*Phenobarbital Calibrators are for the calibration of the ARCHITECT *i* System with *STAT* protocol capability when used for the quantitative determination of phenobarbital in human serum or plasma.

Prescription Use   X   And/Or Over-The-Counter Use \_\_\_\_\_  
(Part 21 CFR 801 Subpart D) (21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE; CONTINUE ON ANOTHER PAGE IF NEEDED)

---

Concurrence of CDRH, Office of In Vitro Diagnostic Device Evaluation and Safety  
(OIVD)

  
\_\_\_\_\_  
Division Sign-Off  
Office of In Vitro Diagnostic Device  
Evaluation and Safety  
510(k) K081231