

**5. 510(k) SUMMARY**

JAN 28 2010

This summary of 510(k) safety and effectiveness information is being submitted in accordance with the requirements of SMDA 1990 and 21 CFR §807.92.

The Assigned 510(k) number is K 093143.

**Submitter's Identification:**

ACON Laboratories, Inc.  
10125 Mesa Rim Road  
San Diego, California 92121

Tel.: 858-875-8019  
Fax: 858-875-8099

Date Prepared: October 2009

**Contact Person:**

Richard Lenart  
Regulatory Affairs Manager

**Proprietary Name of the Device:**

Mission<sup>®</sup> Breath Alcohol Detector

**Common Name:**

Breath-alcohol test system

**Regulation Section and Classification:**

21 CFR § 862.3050      Class I, Breath-alcohol test system

**Product Code:**

DJZ Devices, Breath trapping, alcohol

**Medical Specialty:**

Toxicology

**Predicate Device:**

BreathScan Alcohol Detector  
Akers Biosciences Inc., Thorofare, New Jersey, USA 08086.  
510(k) Number: K060761

**Description:**

The Breath Alcohol Detector is a visual semi-quantitative test for the detection of alcohol in the exhaled breath. The Breath Alcohol Detector consists of a plastic tube, two plastic plugs, a blow bag (optional), a tube label and a glass vial encased with reaction crystals. The crystals employ a solid-phase chemistry system based on chemically chromogenic reaction. Alcohol, if present in the exhaled breath, reacts with the chemically coated crystals and produces a color change. This color change is proportional to the concentration of alcohol in the breath, which is an approximation of relative Blood Alcohol Concentration (BAC). The Breath Alcohol Detector is available with or without blow bags in six cut-off levels: 0.02%, 0.04%, 0.05%, 0.06%, 0.08% or 0.10% cut-off levels.

**Intended Use:**

The Breath Alcohol Detector is for the semi-quantitative rapid detection of the presence of alcohol in the exhaled breath. The Breath Alcohol Detector indicates relative Blood Alcohol Concentration (BAC) at 0.02%, 0.04%, 0.05%, 0.06%, 0.08% or 0.10% cut-off levels.

**Technological Characteristics:**

| Feature             | Specifications                             |
|---------------------|--|
| Methodology         | Chromogenic reaction                       |
| Specimen            | Exhaled breath                             |
| Measurement Range   | 0.02%, 0.04%, 0.05%, 0.06, 0.08% and 0.10% |
| Measuring Units     | BAC %                                      |
| Reading Time        | 2 minutes                                  |
| Reading Stability   | 5 minutes                                  |
| Storage Temperature | 2-30°C (36-86°F)                           |
| Shelf Life          | 3 years                                    |
| Dimensions          | 0.94 x 8.0 cm (0.37 x 3.15 inches)         |
| Weight              | 2.0 g (0.07 oz)                            |

**Comparison to Predicate Devices:**

The Mission<sup>®</sup> Breath Alcohol Detector is substantially equivalent to the Akers Biosciences BreathScan Alcohol Detector.

| Feature             | Mission <sup>®</sup> Breath Alcohol Detector  | Akers Bioscience BreathScan Alcohol Detector   |
|---------------------|---|--|
| <b>Similarities</b> |   |  |
| Intended Use        | Detect presence of alcohol in exhaled breath  | Same   |
| Target Population   | Over the counter  | Same   |
| Calibration         | None required   | Same   |
| Methodology         | Chromogenic reaction  | Same   |
| Anatomical Site     | Mouth   | Same   |
| Test Sample         | Exhaled breath  | Same   |
| Result              | Semi-quantitative   | Same   |
| Interpretation      | Visual color change   | Same   |
| Measuring Units     | BAC %   | Same   |
| Mouthpiece          | None required   | Same   |
| Blowing Time        | 12 seconds  | Same   |
| Dimensions          | 0.94 x 8.0 cm (0.37 x 3.15 inches)  | Same   |
| Weight              | 2.0 g (0.07 oz)   | Same   |
| <b>Differences</b>  |   |  |
| Measurement Range   | Separate devices available at different cut-off levels: 0.02%, 0.04%, 0.05%, 0.06%, 0.08% and 0.10% | Separate devices available at different cut-off levels: 0.02%, 0.04%, 0.05%, and 0.08% |
| Blow Bag            | May be used with or without blow bag  | None   |

**Discussion of Non-Clinical Tests Performed for Determination of Substantial Equivalence are as follows:**

Guidance documents included the “NHTSA/DOT Highway Safety Programs; Model Specifications for Screening Devices to Measure Alcohol in Bodily Fluids, (Federal Register/Vol.59, No.147, August 2, 1994/Notices/39382),” and NHTSA/DOT Highway Safety Programs; Model Specifications for Screening Devices to Measure Alcohol in Bodily Fluids, (Federal Register/Vol.73, No.62, March 3, 2008/Notices)”

**Laboratory Testing:**

The performance characteristics of the Mission<sup>®</sup> Breath Alcohol Detector were evaluated by the following studies: precision, analytical specificity-blank reading, analytical specificity-cigarette smoke, analytical specificity-volatile substances, temperature flexibility, vibration effect, lighting effect and device comparison study.

**Discussion of Clinical Tests Performed:**

Clinical studies were conducted with laypersons and trained laboratory technicians using the Mission<sup>®</sup> Breath Alcohol Detector. The study data were presented evaluating the accuracy of the Mission<sup>®</sup> Breath Alcohol Detector compared to an evidentiary breath test, Alco-Sensor IV, manufactured by Intoximeters Inc. which is a DOT/NHTSA approved device (Conforming Product List of Evidential Breath Alcohol Measurement Devices – FR/Vol 72, No 241/December 2007), per the ACON Clinical Study Protocol for the Breath Alcohol Detector. Study results indicate that non-professional, inexperienced laypersons were able to obtain comparable readings when using the Mission<sup>®</sup> Breath Alcohol Detector as compared to the results obtained by the trained technicians. In addition, the participating laypersons were questioned and responded as satisfied with the ease of operation by following the Instructions for Use and the overall performance of the Mission<sup>®</sup> Breath Alcohol Detector.

**Conclusion:**

The laboratory testing and clinical study results demonstrate that the Mission<sup>®</sup> Breath Alcohol Detector is safe, accurate and easy-to-use. It also demonstrates that the Mission<sup>®</sup> Breath Alcohol Detector is substantially equivalent to the Akers Biosciences BreathScan Alcohol Detector, currently sold on the U.S. market.



DEPARTMENT OF HEALTH & HUMAN SERVICES

---

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Mail Center – WO66-0609  
Silver Spring, MD 20993-0002

Acon Laboratories, Inc.  
c/o Mr. Richard Lenart  
Regulatory Affairs Manager  
10125 Mesa Rim Road  
San Diego, CA 92121

Re: k093143  
Trade Name: Mission® Breath Alcohol Detector  
Regulation Number: 21 CFR §862.3050  
Regulation Name: Breath-Alcohol Test System  
Regulatory Class: Class I, reserved  
Product Codes: DJZ  
Dated: December 11, 2009  
Received: December 15, 2009

**JAN 28 2010**

Dear Mr Lenart:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in Title 21, Code of Federal Regulations (CFR), Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); and good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820).

Page 2

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of *In Vitro* Diagnostic Device Evaluation and Safety at (301) 796-5450. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding postmarket surveillance, please contact CDRH's Office of Surveillance and Biometric's (OSB's) Division of Postmarket Surveillance at (301) 796-5760. For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or ( 301 ) 796-5680 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

A handwritten signature in black ink, appearing to read 'CCH', with a long horizontal flourish extending to the right.

Courtney C. Harper, Ph.D.  
Director  
Division of Chemistry and Toxicology  
Office of *In Vitro* Diagnostic Device  
Evaluation and Safety  
Center for Devices and Radiological Health

Enclosure

4. INDICATIONS FOR USE STATEMENT

510(k) Number (if known): K093143

Device Name: Mission<sup>®</sup> Breath Alcohol Detector

Indications for Use:


The Breath Alcohol Detector is for the semi-quantitative rapid detection of the presence of alcohol in the exhaled breath. The Breath Alcohol Detector indicates relative Blood Alcohol Concentration (BAC) at 0.02%, 0.04%, 0.05%, 0.06%, 0.08% or 0.10% cut-off levels. Measurements obtained by this device are used in the diagnosis of alcohol intoxication.

Prescription Use \_\_\_\_\_ AND/OR Over-The-Counter Use X  
(Part 21 CFR 801 Subpart D) (21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

---

Concurrence of CDRH, Office of In Vitro Diagnostic Devices (OIVD)

  
Division Sign-Off

Page 1 of 1

Office of In Vitro Diagnostic  
Device Evaluation and Safety

E-1 of E-1

K093143