

**510(k) Summary of Safety and Effectiveness****MAR 25 2010****General Provisions****Trade Name:** AngioDynamics, Inc. Morpheus® SMART PICC CT and Procedure Kit**Classification Name:** Percutaneous, Implanted, Long-Term, Intravascular Catheters 80 LJS**Name of Predicate Devices**

The following predicate devices have been identified for the 6F Triple Lumen Morpheus® SMART PICC CT and Procedure Kits:

Device Name	510(k) Number	Concurrence Date
AngioDynamics Morpheus® CT PICC and Procedure Kits	K070615	May 04, 2007
	K060887	April 24, 2006
	K041420	July 26, 2004
	K040446	March 05, 2004
	K031626	June 19, 2003
	K030415	April 30, 2003
Bard Access Systems, Inc. 6 Fr TL PowerPICC® Catheter	K053501	January 13, 2006

**Contact Name**

Jodi Lynn Frasier  
 Senior Regulatory Affairs Professional  
 Access Business Unit  
 AngioDynamics, Inc  
 603 Queensbury, NY 12804  
 (518) 798-1215 ext 1676

**Date Summary Prepared**

March 01, 2010

**Classification**

Class II

**Performance Standards**

No performance standards have been established under Section 514 of the Food, Drug and Cosmetic Act for Percutaneous, Implanted, Long-Term, Intravascular Catheters.

*Continued on next page*

## 510(k) Summary of Safety and Effectiveness, Continued

### Intended Use and Device Description

The AngioDynamics 6F Triple Lumen Morpheus® SMART PICC CT and Procedure Kit is intended for short or long term peripheral access to the central venous system for intravenous therapy, power injections of contrast media, and allows for central venous pressure monitoring. For blood sampling, infusion, or therapy use a 4 French or larger catheter. For central venous pressure monitoring, it is recommended that a catheter lumen of 20 gauge or larger be used. The device is available as a procedural kit with either a Stylet or a Nitinol Wire.


### Biocompatibility

The 6F Triple Lumen Morpheus® SMART PICC CT and Procedure Kit have been tested for biocompatibility per ISO 10993. All data demonstrate this device is biocompatible for its intended use.

### Summary Comparing Technological Modifications

The AngioDynamics 6F Triple Lumen Morpheus® SMART PICC CT and Procedure Kits device description is as follows:

- 6 French outside diameter, reversed tapered shaft design
- Catheter usable length is 55 cm
- Shaft inner lumen is a triple lumen design
- Catheter shaft tubing is marked with depth indicators
- The catheter has one power injectable lumen
- The product labeling warns against power injection procedures through the two small lumens, which are clearly identified as non CT.
- Three extension legs to facilitate injection through each lumen of the catheter shaft.



*Continued on next page*

**510(k) Summary of Safety and Effectiveness, Continued**

---

**Summary  
Comparing  
Technological  
Modifications  
(continued)**

*The only modifications that were made are as follows:*

- Expand the existing product line to provide a 6F triple lumen catheter shaft. The existing Morpheus® product line currently provides a 6F and 7F Dual lumen catheter design, while
- The catheter shaft is a single durometer material vs. the dual durometer design of the existing product portfolio. The material is the same base material as the distal end of the currently marketed Morpheus® PICC CT and Procedure Kits, with the exception of the barium loading. The 6F triple device will have a 20% Barium loading which is the same as the proximal end of the currently marketed Morpheus® PICC CT and Procedure Kits.
- An additional extension leg has been added to facilitate injection through the third lumen of the catheter shaft. The extension leg materials are identical to those of the currently marketed dual lumen Morpheus® PICC CT and Procedure Kits.
- A natural colored luer will be utilized for the third lumen of the catheter. This luer has identical specifications as those used on the existing Morpheus® PICC CT and Procedure kits and is of the same material. The only difference is the omission of a colorant (red or blue) so that each lumen has a distinctly colored luer for identification purposes.
- Use of a silicone processing aid, for which leave trace amounts may remain on the finished device.

---

*Continued on next page*

## 510(k) Summary of Safety and Effectiveness, Continued

### Summary of Verification Activities

The table immediately following outlines the verification/validation activities completed on the proposed device and compares that to the predicate device.

Test	Proposed Device Acceptance Criteria	Proposed 6F Triple Catheter (K093406) T=0	Proposed 6F Triple Catheter (K093406) T=1	Predicate Device Acceptance Criteria	AngioDynamics Morpheus® CT PICC and Procedure Kits (K070615, K060887, K041420, K040446, K031626, K030415) T=0	AngioDynamics Morpheus® CT PICC and Procedure Kits (K070615, K060887, K041420, K040446, K031626, K030415) T=1
Tip Whip	Tip "whipping" is defined as the tip moving from left to right during an injection. Tip whipping must be less than or equal to predicate devices	Pass	Pass	Tip "whipping" is defined as the tip moving from left to right during an injection. Tip whipping must be less than or equal to predicate devices	Pass	Pass
Tip Displacement	Tip displacement is defined as the tip backing-up during an injection. Tip displacement must be less than or equal to predicate devices	Pass	Pass	Tip displacement is defined as the tip backing-up during an injection. Tip displacement must be less than or equal to predicate devices	Pass	Pass

*Continued on next page*

## Summary of Verification Activities (continued)

Test	Proposed Device Acceptance Criteria	Proposed 6F Triple Catheter (K093406) T=0	Proposed 6F Triple Catheter (K093406) T=1	Predicate Device Acceptance Criteria	AngioDynamics Morpheus® CT PICC and Procedure Kits (K070615, K060887, K041420, K040446, K031626, K030415) T=0	AngioDynamics Morpheus® CT PICC and Procedure Kits (K070615, K060887, K041420, K040446, K031626, K030415) T=1
CT Injection – Catheter Flow Rate	Flow Rate = 5mL/sec $\pm$ 0.5mL/sec with 95% confidence that 95% of the population meets specification	Pass	Pass	Flow Rate = 3 to 8mL/sec $\pm$ 0.5mL/sec (depending on catheter configuration) with 95% confidence that 95% of the population meets specification	Pass	Pass
CT Injection – CT Lumen Integrity	CT Lumen must withstand a minimum of 10 injections at a minimum flow rate of 5mL/sec with a 95% confidence that 80% of the population meets specification	Pass	Pass	CT Lumen must withstand a minimum of 10 injections at a minimum flow rate of 3 to 8mL/sec (depending on catheter configuration) with a 95% confidence that 80% of the population meets specification	Pass	Pass

Continued on next page

## Summary of Verification Activities (continued)

Test	Proposed Device Acceptance Criteria	Proposed 6F Triple Catheter (K093406) T=0	Proposed 6F Triple Catheter (K093406) T=1	Predicate Device Acceptance Criteria	AngioDynamics Morpheus® CT PICC and Procedure Kits (K070615, K060887, K041420, K040446, K031626, K030415) T=0	AngioDynamics Morpheus® CT PICC and Procedure Kits (K070615, K060887, K041420, K040446, K031626, K030415) T=1
Aspiration	Minimum aspiration rate to be 3cc/min using a 10cc syringe without total collapse (all three catheter lumens) with 95% confidence that 95% of the population meets specification	Pass	Pass	Minimum aspiration rate to be 3cc/min using a 10cc syringe without total collapse (all three catheter lumens) with 95% confidence that 95% of the population meets specification	Pass	Pass
Gravity Flow Rate	Minimum gravity flow rate to be 750 ml/hr for CT lumen and 182 ml/hr for non-CT lumens with 95% confidence that 95% of the population meets specification	Pass	Pass	Minimum gravity flow rate to be 750 ml/hr for CT lumen and 182 ml/hr for non-CT lumens with 95% confidence that 95% of the population meets specification	Pass	Pass

Continued on next page

## Summary of Verification Activities (continued)

Test	Proposed Device Acceptance Criteria	Proposed 6F Triple Catheter (K093406) T=0	Proposed 6F Triple Catheter (K093406) T=1	Predicate Device Acceptance Criteria	AngioDynamics Morpheus® CT PICC and Procedure Kits (K070615, K060887, K041420, K040446, K031626, K030415) T=0	AngioDynamics Morpheus® CT PICC and Procedure Kits (K070615, K060887, K041420, K040446, K031626, K030415) T=1
Stylet Withdrawal Testing	Stylet withdrawal force must be less than 2 lbs with 95% confidence that 95% of the population meets specification	Pass	Pass	Stylet withdrawal force must be less than 2 lbs with 95% confidence that 95% of the population meets specification	Pass	Pass
Wire Withdrawal Testing	Guidewire withdrawal force must be less than 2 lbs with 95% confidence that 95% of the population meets specification	Pass	Pass	Guidewire withdrawal force must be less than 2 lbs with 95% confidence that 95% of the population meets specification	Pass	Pass
Static Burst	Catheter burst pressure must be greater than 150 psi with 95% confidence that 95% of the population meets the specification	Pass	Pass	Catheter burst pressure must be greater than 150 psi with 95% confidence that 95% of the population meets the specification	Pass	Pass

Continued on next page

## Summary of Verification Activities (continued)

Test	Proposed Device Acceptance Criteria	Proposed 6F Triple Catheter (K093406) T=0	Proposed 6F Triple Catheter (K093406) T=1	Predicate Device Acceptance Criteria	AngioDynamics Morpheus® CT PICC and Procedure Kits (K070615, K060887, K041420, K040446, K031626, K030415) T=0	AngioDynamics Morpheus® CT PICC and Procedure Kits (K070615, K060887, K041420, K040446, K031626, K030415) T=1
Dynamic Burst	Catheter dynamic burst pressure must be greater than 300 psi with 95% confidence that 95% of the population meets specification	Pass	Pass	Catheter dynamic burst pressure must be greater than 300 psi with 95% confidence that 95% of the population meets specification	Pass	Pass
Tensile Testing						
Extension Leg to Natural Color Hub	Tensile strength must be greater than 5 lbs with 95% confidence that 99% of the population meets the specification	Pass	Pass	Tensile strength must be greater than 5 lbs with 95% confidence that 99% of the population meets the specification	Pass	Pass
Non- CT Extension Leg to Trifurcate	Tensile strength must be greater than 5 lbs with 95% confidence that 99% of the population meets the specification	Pass	Pass	Tensile strength must be greater than 5 lbs with 95% confidence that 99% of the population meets the specification	Pass	Pass

Continued on next page



## Summary of Verification Activities (continued)

Test	Proposed Device Acceptance Criteria	Proposed 6F Triple Catheter (K093406) T=0	Proposed 6F Triple Catheter (K093406) T=1	Predicate Device Acceptance Criteria	AngioDynamics Morpheus® CT PICC and Procedure Kits (K070615, K060887, K041420, K040446, K031626, K030415) T=0	AngioDynamics Morpheus® CT PICC and Procedure Kits (K070615, K060887, K041420, K040446, K031626, K030415) T=1
CT Extension Leg to Trifurcate	Tensile strength must be greater than 5 lbs with 95% confidence that 99% of the population meets the specification	Pass	Pass	Tensile strength must be greater than 5 lbs with 95% confidence that 99% of the population meets the specification	Pass	Pass
Shaft to Trifurcate	Tensile strength must be greater than 5 lbs with 95% confidence that 99% of the population meets the specification	Pass	Pass	Tensile strength must be greater than 5 lbs with 95% confidence that 99% of the population meets the specification	Pass	Pass
Catheter Shaft (at nominal OD)	Tensile strength must be greater than 5 lbs with 95% confidence that 99% of the population meets the specification	Pass	Pass	Tensile strength must be greater than 5 lbs with 95% confidence that 99% of the population meets the specification	Pass	Pass

Continued on next page

**Summary of Verification Activities (continued)**

Test	Proposed Device Acceptance Criteria	Proposed 6F Triple Catheter (K093406) T=0	Proposed 6F Triple Catheter (K093406) T=1	Predicate Device Acceptance Criteria	AngioDynamics Morpheus® CT PICC and Procedure Kits (K070615, K060887, K041420, K040446, K031626, K030415) T=0	AngioDynamics Morpheus® CT PICC and Procedure Kits (K070615, K060887, K041420, K040446, K031626, K030415) T=1
<b>CVP Pressure Monitoring</b>						
Pressure Monitoring	No difference observed between SVS and catheter in mean pressure, systolic pressure and diastolic pressure as compared to predicate devices	Pass	Pass	No difference observed between SVS and catheter in mean pressure, systolic pressure and diastolic pressure as compared to predicate devices	Pass	Pass
Natural Frequency	The natural frequency must be equal to or greater than the predicate device	Pass	Pass	The natural frequency must be equal to or greater than the predicate device	Pass	Pass

**Summary of Substantial Equivalence**

The 6F Triple Lumen Morpheus® SMART PICC CT and Procedure Kit have been tested and compared to the predicate device. All data gathered demonstrate this device is substantially equivalent. No new issues of safety or efficacy have been raised.

## Additional Information

---

**Material  
Change**

In addition to the previously submitted information, AngioDynamics would like to identify an additional change for the device submitted for the pending 510(k) K093406. A silicone material is used as a processing aid which may leave trace amounts on the device. As a result, this MDX silicone should be identified as a material present and included in the 510(k). All testing previously conducted were on units built with this processing aid and as a result representative of the finished device for which we seek marketing clearance. This was an omission in the original submission.

---

**Biocompatibility  
Data**

MDX Silicone, manufactured by Dow is a commonly used medical grade material. Extensive biocompatibility testing has previously been conducted. The following page provides written authorization to allow the Food & Drug Administration to access these data for purpose of this submission review.

---



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Room - WO66-G609  
Silver Spring, MD 20993-0002

MAR 25 2010

AngioDynamics, Incorporated  
Ms. Jodi Lynn Frasier  
Senior Regulatory Affairs Professional  
Access Business Unit  
603 Queensbury Avenue  
Queensbury, New York 12804

Re: K093406  
Trade/Device Name: Morpheus<sup>®</sup> SMART PICC CT and Procedure Kit  
Regulation Number: 21 CFR 880.5970  
Regulation Name: Percutaneous, Implanted, Long-Term Intravascular Catheter  
Regulatory Class: II  
Product Code: LJS  
Dated: March 1, 2010  
Received: March 2, 2010

Dear Ms. Frasier:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you; however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

A handwritten signature in black ink, appearing to read "Anthony D. Watson" followed by a small flourish.

Anthony D. Watson, B.S., M.S., M.B.A.  
Director  
Division of Anesthesiology, General Hospital,  
Infection Control and Dental Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

## Section 6

### Statement of Indications For Use

---

#### INDICATIONS FOR USE

**510(k) Application:** Special 510K Application

**Device Name:** AngioDynamics, Inc. Morpheus<sup>®</sup> SMART PICC CT and Procedure Kit

**Indications for Use:**

The AngioDynamics, Inc. Morpheus<sup>®</sup> SMART PICC CT and Procedure Kit is indicated for short or long term peripheral access to the central venous system for intravenous therapy, power injections of contrast media, and allows for central venous pressure monitoring. For blood sampling, infusion, or therapy use a 4 French or larger catheter. For central venous pressure monitoring, it is recommended that catheter lumen of 20 gauge or larger be used.

Prescription Use   X    
(Per 21 CFR 801.109)


OR Over-the-Counter Use \_\_\_\_\_

---

Please do not write below this line - continue on another page if needed

---

Concurrence of CDRH, Office of Device Evaluation (ODE)

  
Division Sign-Off

Division of Anesthesiology, General Hospital  
Infection Control, Dental Devices

Page 20 of 48

510(k) Number: K093406