510(k) SUBSTANTIAL EQUIVALENCE DETERMINATION DECISION SUMMARY DEVICE ONLY TEMPLATE

A. 510(k) Number:

k050182

B. Purpose For Submission:

Premarket Notification 510(k) of intention to manufacture and market the Biomedix, Inc. Q. STEPS Biometer G/C Dual Monitoring System

C. Analyte:

Whole Blood Glucose Whole Blood Cholesterol

D. Type of Test:

Quantitative, utilizing Glucose Oxidase technology. Quantitative, utilizing Cholesterol Oxidase technology.

E. Applicant:

Biomedix, Inc.

F. Proprietary and Established Names:

Q. STEPSTM Biometer G/C Dual Monitoring System.

G. Regulatory Information:

1. Regulation section:

Regulation	Standard Product Nomenclature	Panel	Product	Class
Number			Code	
862.1345	System, Test, Blood Glucose, Over	Chemistry	NBW	II
	The Counter	(75)		
862.1345	Glucose Oxidase, Glucose	Chemistry	CGA	II
		(75)		
862.1175	Enzymatic Esterase-Oxidase,	Chemistry	CHH	Ι
	Cholesterol	(75)		
862.1660	Single (Specified) Analyte Controls	Chemistry	JJX	I
	(Assayed and Unassayed)	(75)		

H. Intended Use:

1. Intended use(s):

See Indications for use below.

2. <u>Indication(s) for use:</u>

The Q. STEPS Biometer G/C Dual Monitoring System is intended for use with Q.STEPS Glucose and Cholesterol Test Strips with Q. STEPS Biometer G/C by healthcare professionals and home users. Q. STEPS Biometer G/C System provides a quantitative measurement of Glucose and Cholesterol in whole blood from the fingertips. The Glucose measurements are used in helping the management of carbohydrate metabolism disorders including diabetes mellitus, idiopathic hypoglycemia and pancreatis islet cell tumors. Cholesterol measurements are used in the management of disorders involving excess cholesterol in the blood, lipid and lipoprotein metabolism disorders.

3. Special condition for use statement(s):

Provides plasma equivalent results.

4. Special instrument Requirements:

Q. STEPS Biometer G/C Dual Monitoring System.

I. Device Description:

The Q. STEPSTM Biometer G/C Dual Monitoring System uses enzymatic electrochemical biosensor technology to measure whole blood glucose and cholesterol levels. When finger blood is applied to the test spot of the biosensor (test strip) a reduction-oxidation reaction occurs. The oxidase of D-Glucose or Cholesterol which is catalyzed by Glucose Oxidase or by Cholesterol Oxidase respectively, causes an electron transfer at the electrode (silver) surface; and therefore, the magnitude of the current produced is proportional to the glucose or cholesterol concentration in the blood. The Biometer G/C uses that current to quantify the glucose and the cholesterol levels in the blood, and then display on the readout of the monitor.

J. Substantial Equivalence Information:

1. Predicate device name(s):

One Touch Basic Blood Glucose Monitoring System PTS Panels Lipid Panel Test Strips

2. Predicate K number(s):

k031472 k023558

3. Comparison with Predicate:

Substantial Equivalence Comparison Glucose

Similarities

T .	Similariues	D 11 / 1024 (T2
Item	Devices k050182	Predicate k031472
Intended Use	Q. STEPS TM Biometer G/C Dual	One Touch Basis System
	Monitoring System (Glucose	Intended to be used with One
	Side Sensor)	Touch ® basic and One
	Intended to be used with	Touch® Test Strips for the
	Q.STEPS TM Biometer G/C and	quantitative measurement of
	Q. STEPS TM Glucose Test Strips	glucose in fresh capillary
	for the quantitative measurement	whole blood from fingertip.
	of glucose in fresh capillary	
	whole blood from the fingertip.	For professionals and Diabetes patients use.
	For professionals and Diabetes	
	patients use.	
Test Principle	Based on Glucose Oxidase	Based on Glucose Oxidase
	oxidation reduction to convert	oxidation reduction to convert
	glucose into gluconic acid.	glucose into gluconic acid.
Labeling Instruction	Test should be run with liquid	Test should be run with liquid
Regarding Response to	quality control material, Q.	quality control material, One
Unusual Results	STEPS [™] Control Solution	Touch® Control Solution
	whenever a new vial of test strip	whenever a new vial of test
	is opened or unusual blood test	strip is opened or unusual
	result is obtained.	blood test result is obtained.
Meter Functional Test	A check strip is provided to	A check strip is provided to
	ensure the system is working	ensure the system is working
	properly	properly
Test Strips Storage	Must be stored in the original	Must be stored in the original
Condition	vial with the cap tightly closed.	vial with the cap tightly
	Stored in a cool dry place not	closed. Stored in a cool dry
	above 86°F (36°C) and away	place not above 86°F (36°C)
	from heat and direct sunlight, not	and away from heat and direct
	refrigerated.	sunlight, not refrigerated.
Strip Shelf Life	After opening the vial, 4 months	After opening the vial, 4 months
		monuis

Differences

Item	Devices k050182	Predicate k031472
Methodology	Amperometric	Photometric
Test Recall Memory	99 test results	75 test results memory
		capacity.
Blood Sample Volume	Minimum is 5µ1/30 seconds	Minimum is 10 μl / 45
Reaction Time		seconds
Physical	The test spot is on the side of the	The test spot is on the center
Characteristics	test strip and is a shape of half-a-	of the test strip and is in the
	circle.	shape of a circle.

Substantial Equivalence Comparison Cholesterol

Similarities

Item	Devices k050182	Predicate k023558
Intended Use	The Q. STEPS TM Test Strip is	The Lipid Panel Test Strips
	intended to measure cholesterol	are intended to measure
	in whole blood on the G/C Dual	cholesterol, HDL and
	Monitoring System.	triglycerides in whole blood
		on a BIOScanner Plus
	For Professionals and Diabetes	(CardioChek Brand) analyzer.
	patients use.	
		For Professionals and home
		user.
Matrix	Finger Whole Blood	Finger Whole Blood
Result Display	Directly displays results without	Directly displays results
	requiring calculation.	without requiring calculation.
Enzymatic reaction	Cholesterol Oxidase and Esterase	Cholesterol Oxidase and
	Reaction	Esterase Reaction
Physical	Test strip with circular spot	Test strip with circular spot
Characteristics		
Calibration Chip	It contains a lot specific	It contains a lot specific
	electronically erasable,	electronically erasable,
	programmable read-only memory	programmable read-only
	(EEPROM) chip in the same	memory (EEPROM) chip in
	package with the strips. The	the same package with the
	EEPROM chip has the curve	strips. The EEPROM chip has
	information programmed into it	the curve information
	and based on a multipoint curve	programmed into it and based
	and is established for each lot.	on a multipoint curve and is
	The user inserts this chip into the	established for each lot. The
	meter with each new lot of test	user inserts this chip into the
	strips.	meter with each new lot of test
		strips.
Device Storage	2-30°C	2-30°C

Differences

Item	Devices k050182	Predicate k023558
Methodology	Amperometric	Photometric
Measuring Range	150-350 mg/dL	100-400 mg/dL
Test Read Time	30 seconds	45 seconds
Sample Volume	Approximately 15 μL or 1 drop	Approximately 135µL or 3
Application	is added to the test spot	drops are added to the test spot
Test	Cholesterol only	Lipid Panel

K. Standard/Guidance Document Referenced (if applicable):

- 1) General Principles of Software Validation; Final Guidance for Industry and FDA Staff, FDA, January 11, 1997.
- 2) IEEE Standard 1012-1986, IEEE Standard for Software Verification and Validation Plans, The institute of Electrical and Electronic Engineers, Inc. 1997.
- 3) Guidance for the Content of Premarket Submissions for Software Contained in Medical Devices, Department of Health and Human Services, Food and Drug Administration, July 2000.
- 4) Quality Systems Model for quality assurance in design, development, production, installation and servicing, ISO 9001:2000, International Organization for Standardization.

L. Test Principle:

The Test Principle used by this device is enzyme electrochemical sensor technology. Biomedix's biosensor uses a separate disposable dry reagent strip a for glucose and cholesterol determinations. When a drop of blood from the fingertip is applied to the half-circle test spot on the glucose test strip or the circular test spot on the cholesterol test strip, the reduction and oxidation reaction causes electron transfer at the electrode surfaces. Current is generated and detected by the Q.STEPS Biometer G/C Dual Monitoring System. The magnitude of the current generated is proportional to the analyte concentration in the blood.

M. Performance Characteristics (if/when applicable):

- 1. Analytical performance:
 - a. Precision/Reproducibility:

Glucose

The sponsor indicated glucose spiked venous whole blood was used to perform the Within-Run Precision Study. Commercially purchased venous blood was pooled together and spun down to separate the erythrocytes from plasma. The venous blood was then adjusted to $45 \pm 3\%$ hematocrit with the plasma. Five different concentrations of glucose spiked whole blood solutions were made by adding the concentrated glucose solution into the whole blood. These glucose spiked whole blood samples were used for the within-run precision study. Each sample was measured twice per day for 20 days. As shown in the table below, the values of C.V were between 2.7% and 8.0%

Within-run	precision	results with	venous whole	blood or	n Biometer G/C
VVICILIII I WII	PI CCIBIOII	i courts with	vellous whole	Dioou of	i Diometer of

Lot #	Glucose	Total Number		Within Run	
	Concentration	of Samples (N)	Mean by	Standard	Coefficient
	(mg/dL) measured		Biometer	Deviation	of Variation
	by YSI		(mg/dL)	(mg/dL)	%
1	50	20	60	2.5	4.3
	80	20	87	4.5	5.2
	120	20	110	8.5	7.8
	200	20	224	11.0	4.9
	400	20	353	16.5	4.7
2	50	20	49	3.1	6.3
	80	20	91	7.3	8.0
	120	20	124	9.6	7.8
	200	20	217	13.9	6.4
	400	20	372	10.3	2.8
3	50	20	65	2.8	4.9
	80	20	84	3.5	4.1
	120	20	125	3.1	2.6
	200	20	240	6.6	2.7
	400	20	432	16.8	4.2

The average within-run total variation coefficient with whole blood was between 3.75 and 5.86%.

Summary results of within-run precision with whole blood.

	Glucose	Average Mean	Average	Average
	Concentration	by Biometer	Standard	Coefficient of
	(mg/dL)	G/C	Deviation by	Variation
	measured by	(mg/dL)	Biometer G/C	[<8%]
	YSI		(mg/dL)	%
All Three Lots	50	58	2.8	4.82
All Three Lots	80	87	5.1	5.86
All Three Lots	120	120	7.0	5.83
All Three Lots	200	227	10.5	4.63
All Three Lots	400	386	14.5	3.75

Cholesterol

The specimens used in the precision study were EDTA preserved whole blood samples that were obtained commercially. First, the whole blood samples were pooled together, spun down and plasma removed. The serum-base cholesterol stock solutions with the desired cholesterol level were then mixed with red blood cells to become a whole blood sample. The cholesterol-spiked blood was prepared at four target concentrations of 190, 200, 240, and 260 mg/dL. The venous blood was then adjusted to a final hematocrit of

45%. 15 μl of these blood samples were tested with the Q.STEPS Biometer G/C System, 20 times for each concentration. The values of within-run (average) total variation coefficient with whole blood were between 1.0 and 4.0%.

Whole blood precision

whole blood precision					
Lot#	Cholesterol	Total		Within-run	
	Concentration	Number of	Mean by	Standard	Coefficient
	by Cobas	Samples (n)	Biometer	Deviation	of Variation
	(mg/dL)		(mg/dL)	(mg/dL)	%CV
Lot 1	190	20	185	6.11	3.30
	200	20	210	7.32	3.49
	240	20	238	7.94	3.33
	260	20	258	4.87	1.89
Lot 2	190	20	196	4.08	2.08
	200	20	203	6.43	3.17
	240	20	230	5.33	2.32
	260	20	263	5.30	2.01
Lot 3	190	20	181	7.35	4.05
	200	20	206	8.35	4.06
	240	20	226	8.71	3.85
	260	20	270	6.59	2.44

The sponsor determined the within run and day-to-day precision on the Q.STEPS System using Standard Cholesterol Solutions at 200mg/dL and 240 mg/dL. 25µl of the standard solution was applied to two test strips per concentration. Readings were taken twice per day for 10 consecutive days. Three different lots of test strips were tested. The within-run C.V. varied form 1.95% to 2.47%

Within-run precision results of Q.STEPS G/C with serum specimen

Lot#	Cholesterol	Total	Mean	With	in-run	Between	Total p	recision
	Concentration	Number	mg/dL	SD	% CV	run	SD	%CV
Lot 1	200	10	202	4.97	2.47	0	6.27	3.11
	240	10	243	5.34	2.19	0	6.00	2.47
Lot 2	200	10	202	449	2.22	0	5.91	2.93
	240	10	242	4.88	2.02	0	5.46	2.26
Lot 3	200	10	203	4.45	2.19	0	5.59	2.75
	240	10	243	4.71	1.95	0	5.13	2.12

Clinical Sites Precision Studies

Venous whole blood samples were collected from patients at three different clinical trial sites. 15µl of the whole blood was applied to each cholesterol test strip. Each blood sample was applied 20 times.

Precision	data	abtained	from	aliniaal	gitog
Precision	aata	optained	irom	ciinical	l sites

Location	Total number	Within-run			
	of samples (n)	Mean	SD	CV	
		(mg/dL)	(mg/dL)	%	
Site 1	20	192	5.91	3.08	
	20	239	6.38	2.67	
Site 2	20	209	6.02	2.89	
	20	239	5.15	2.15	
Site 3	20	206	6.35	3.09	
	20	241	6.86	2.85	

b. Linearity/assay reportable range:

Glucose

The sponsor pooled $980\mu l$ of venous blood and adjusted the hematocrit to $45\% \pm 3$. The pooled sample was then spiked with a Glucose stock solution resulting in blood glucose concentrations of 25 mg/dL, 50 mg/dL, 100 mg/dL, 200 mg/dL 300 mg/dL, 400 mg/dL, 500 mg/dL and 600 mg/dL. Each sample was then measured by the YSI STAT PLUS analyzer, as a reference, and then $10~\mu l$ of the sample was placed on the Q.STEPS Side Sensor Test Strip and Biometer G/C. Each concentration level of whole blood glucose spiked solution was applied to four test strips. The average reading was taken for the measurements.

Results

For each lot of test strips, the glucose concentration determined by the Q.STEPS Biometer G/C was plotted against the glucose concentration determined by the YSI 2300

Linearity studies of whole blood

Lot #	Tested Range	\mathbb{R}^2	Slope	Intercept
Lot 1	50-400 mg/dL	0.99	1.025	3.9
Lot 2	25-600 mg/dL	0.98	1.192	-13
Lot 3	25-600 mg/dL	0.99	1.033	-2.6

Cholesterol

The sponsor spiked 15 µl of venous blood with cholesterol stock solutions resulting in blood cholesterol concentrations of 150 mg/dL, 200 mg/dL, 240 mg/dL, 300 mg/dL and 350 mg/dL. The cholesterol concentration of each sample was first measured by the Cobas Mira PLUS, as a reference, and then 15µl of the sample was placed on the Q.STEPS Cholesterol Test Strip and the Biometer G/C. Each level of whole blood cholesterol spiked solution was applied to four test strips. The average readings were taken for the measurement of linear regression.

Results

For each lot of test strips, the Cholesterol concentration determined by the Q.STEPS Biometer was plotted against the Cholesterol concentration determined by the Cobas Mira PLUS.

Linearity studies of whole blood

Lot #	Tested Range	R^2	Slope	Intercept
Lot 1	150-350 mg/dL	0.96	1.03	-6.47
Lot 2	150-350 mg/dL	0.99	1.01	-0.09
Lot 3	150-350 mg/dL	0.98	1.05	-15.76

c. Traceability (controls, calibrators, or method):

The traceability of the Glucose and Cholesterol calibrators and controls are verified against commercially available Standard Reference Material from the National Institute of Standards and Technology (formerly NBS).

- 1. Commercially available (serum based) Standard Cholesterol Solutions
- 2. "Current Status of Blood Cholesterol Measurement in Clinical Laboratories in the United States: A report from the Laboratory Standardization Panel of the National Cholesterol Education Program". Clin. Chem 34/1, 1998, 193-201.
- 3. "Recommendations for Improving Cholesterol Measurement: A report from the Laboratory Standardization Panel of the National Cholesterol Education Program". NIH Publication No. 90-2964, February 1990.
 - d. Detection limit:

Glucose 50 - 400 mg/dL

Cholesterol 150 - 350 mg/dL

e. Analytical specificity:

Glucose

According to the sponsor twenty-three commonly tested interferent substances were examined. No interference was observed in bilirubin, creatinine and citrate at physiological levels. Based on the tested concentrations, 4-acetamidophenol, ascorbic acid, dopamine, L-Dopa, Ibuprofen, methyldopa, and uric acid interfered with some glucose measurements.

Tested interference substances, concentrations, and their effects.

Substance	Physiological Therapeutic Levels (mg/dL)	Toxic (mg/dL)	Substance Tested Conc. (mg/dL)	Interference At Glucose 80 mg/dL	Interference At Glucose 200 mg/dL)	Effect of low gluc conc. (80mg/dL)	Effect of high gluc conc. (200 mg/dL)
4Acetamindophenol	1.0-2.0	15	2.0-6.0	No inter up to 2.0 mg/dL	No inter up to 4.0 mg/dL	A	A
Ascorbic Acid	0.8- 1.2	-	1.0-4.0	No inter up to 2.0 mg/dL	Inter at 1.0 mg/dL	▼	▼
Bilirubin	0.1-1.2	-	6.8-20	No inter up to 20 mg/dL	No inter up to 20 mg/dL	_	_
Cholesterol	<200	-	300	No inter	No inter	_	_
Citrate sodium salt	1.7-3.0	-	8.9- 26.6	No inter	No inter	_	_
Creatinine	0.6-1.2	-	21.9- 65.7	No inter	No inter	_	_
Dextrin	-	-	0.1-0.2	No inter	No inter	_	_
L-Dopa	0.02-0.03	-	0.23- 6.9	No inter up to 2.3 mg/dL	No inter up to 2.3 mg/dL	A	A
Dopamine	0.4-1.6	-	0.02-2.6	No inter up to 0.7 mg/dL	No inter up to 0.52 mg/dL	A	A
EDTA	61	240	100- 400	No inter	No inter	_	_
D-Galactose		-	22.1 – 66.2		No inter	_	_
Ibuprofen	0.5-4.2	-	27-81	No inter	No inter up to 54 mg/dL	_	•
K3Fe(CN)6	0.07	2.86	0.2-0.5	No inter	No inter	_	_
Maltose	-	-	18-54	No inter	No inter		
D-Mannose	5.0-7.5 mg/day	-	10-30	No inter	No inter	_	_
Mega8	-	-	0.2-0.5	No inter	No inter	-	—
Methyldopa	0.1-0.5	>1.0	1.2-3.5	No inter up to 2.0 mg/dL	No inter up to 1.0 mg/dL	▼	▼
Salicylic acid	15-30	>40	8.7-	No inter	No inter	▼	▼

Substance	Physiological	Toxic	Substance	Interference	Interference	Effect of	Effect
	Therapeutic	(mg/dL)	Tested	At Glucose	At Glucose	low gluc	of high
	Levels		Conc.	80 mg/dL	200 mg/dL)	conc.	gluc
	(mg/dL)		(mg/dL)			(80mg/dL)	conc.
							(200
			261	2.5	- 70		mg/dL)
			26.1	up to 25	up to 50		
				mg/dL	mg/dL		
Tetracycline	0.4	-	1.8-5.4	No inter	No inter	_	-
Tolazamide	2.0-2.5	-	1.6-5.4	No inter	No inter	_	-
Triglyceride	<190	-	100-	No inter	No inter	_	_
			300				
Uric Acid	M: 2.1-7.8	-	2.7-8.1	No inter	No inter	A	A
	F: 2.0- 6.4			up to	up to		
				2.7	5.4		
				mg/dL	mg/dL		
D-Xylose	-	-	15.3-	No inter	No inter	_	_
			45.9				

Cholesterol

The sponsor tested ten common endogenous or exogenous substances for cholesterol meters. Hemoglobin was also tested due to the occurrence of hemolyzed specimens. The sponsor utilized a Student t test for statistical analysis. Paired differences were analyzed to determine statistically significant differences between 3 lots of cholesterol test strips. A p value less than .05 was considered to be statistically significant. The therapeutic and toxic concentrations of the interferents are listed in the below table.

Tested interference substances, concentrations, and their effects.

Substance	Physiological Therapeutic Levels (mg/dL)	Toxic (mg/dL)	Substance Tested Conc. (mg/dL)	Interference At Chol 200 mg/dL	Interference At Chol 240 mg/dL)	Effect of low chol conc. (200 mg/dL)	Effect of high chol conc. (240 mg/dL)
4Acetamindophenol	1.0-2.0	>15	1.5-6.0	No inter up to 3.0 mg/dL	No inter up to 3.0 mg/dL	A	A
Ascorbic Acid	0.8- 1.2	-	1.0-4.0	No inter up to 2.0 mg/dL	Inter at 1.0 mg/dL	▼	▼
Bilirubin	0.1-1.2	-	6.8-20	No inter up to 20 mg/dL	No inter up to 20 mg/dL	_	_
Dopamine	0.4-1.6	-	0.02- 2.6	No inter up to 8.0 mg/dL	No inter up to 4.0	•	▼
EDTA	61	240	100-	No inter	No inter	_	_

Substance	Physiological Therapeutic Levels (mg/dL)	Toxic (mg/dL)	Substance Tested Conc. (mg/dL)	Interference At Chol 200 mg/dL	Interference At Chol 240 mg/dL)	Effect of low chol conc. (200 mg/dL)	Effect of high chol conc. (240 mg/dL)
			400	up to 400 mg/dL	up to 400 mg/dL		
Hemoglobin	10 mg/dL	-	22.1 – 66.2	No inter up to 200 mg/dL	No inter up to 200 mg/dL	_	_
Ibuprofen	0.5-4.2	-	27-81	No inter up to 10 mg/dL	No inter up to 20 mg/dL	A	A
Methyldopa	0.1-0.5	>1.0	1.2-3.5	No inter up to 2.0 mg/dL	No inter up to 1.0 mg/dL	•	▼
Salicylic acid	15-30	>40	8.7- 26.1	No inter up to 25 mg/dL	No inter up to 50 mg/dL	•	▼
Triglyceride	<190	-	475 – 1900	No inter up to 1900 mg/dL	No inter up to 1900 mg/dL	_	_
Uric Acid	M: 2.1-7.8 F: 2.0- 6.4	-	5.0- 20.0	No inter up to 10.0 mg/dL	No inter up to 5.0 mg/dL	▼	▼

Hematocrit Study

Glucose

The sponsor obtained venous whole blood samples that were pooled, and then spun down to separate the red cells from the plasma. The plasma was adjusted to the desired target hematocrit concentration levels of approximately 20%, 25%, 30%, 40%, 45%, 50%, and 60%. Each hematocrit level had 5 target glucose concentrations (50 mg/dL, 80 mg/dL, 120 mg/dL, 200 mg/dL and 400 mg/dL), which were prepared by spiking with appropriate volumes of glucose stock solutions.

The sponsor's acceptance criteria for this study was determined to be hematocrit levels that exhibit glucose concentrations within \pm 20% of that glucose reading of the same specimen at 45% hematocrit level.

Glucose		Hematocrit Percent Levels						
Conc.								
(mg/dL)								
	20%	25%	30%	35%	40%	50%	55%	60%
50	16.1	19.7	12.7	8.1	-0.7	-6.4	-19.7	-22.1
80	25.3	14.5	19.6	2.1	-2.2	-9.8	-7.3	-7.4
120	21.7	11.4	10.5	5.6	-3.7	-9.7	-10.2	-18.5
200	22.9	23.0	14.2	-1.3	-2.7	-1.4	-7.0	-17.7
400	10.0	-4.7	-3.0	4.7	-0.3	-7.2	-11.0	-18.2
			Не	matocrit P	ercent Lev	rels		
50	20.9	27.4	8.2	6.9	3.6	-2.5	-5.9	-10.0
80	8.5	3.8	3.1	2.2	-5.6	-8.0	1.3	-8.3
120	13.5	11.6	8.2	6.7	3.4	-2.0	-10.1	-15.9
200	11.3	15.3	11.1	7.2	1.5	-5.3	-14.4	-14.5
400	1.8	-11.7	4.0	7.8	6.3	-4.0	-5.7	-11.1

Cholesterol

The sponsor obtained venous whole blood samples that were pooled, and then spun down to separate the red cells from the plasma. The plasma was adjusted to the desired target hematocrit concentration levels of approximately 20%, 25%, 30%, 40%, 45%, 50%, and 60%. Each hematocrit level had 2 target cholesterol concentrations (200 and 240 mg/dL), which were prepared by spiking with appropriate volumes of cholesterol stock solutions. 15 μ l of the adjusted whole blood sample was applied to the cholesterol test strips. The data is presented below.

Percent difference in hematocrit

Cholesterol		Hematocrit Percent Levels						
Conc.								
(mg/dL)								
	20%	25%	30%	35%	40%	50%	55%	60%
200	13.4	13.4	11.5	8.7	6.1	5.1	2.3	-12.4
240	22.1	22.1 22.1 13.7 7.7 5.6 -2.6 3.9 -11.6					-11.6	
		Hematocrit Percent Levels						
200	20.5	10.5	15.0	5.0	1.5	-5.0	3.0	-10.3
240	10.5	10.5	10.0	5.0	2.9	-2.9	-6.7	-10.0

f. Assay cut-off: Not Applicable

2. Comparison studies:

a. Method comparison with predicate device:

Glucose

Method comparison studies were performed at three different clinical sites. The correlation studies were made between the Biometer (k050182) new device, One Touch (k031472) predicate device, and the YSI (reference method) utilizing finger stick whole blood samples. The sponsor indicated that each home user performed their own fingerstick and performed the test on the Q.STEPS Biometer G/C. The measurement was also read by the professionals.

The professionals then performed another fingerstick on the same home user and performed a glucose test on the same device with the same lot of test strips. All of the results were masked from each other. The professional then tested the same user with the One Touch and YSI methods. The comparison results are presented the table below. The results are expressed as the mean absolute bias and regression analysis.

Comparison of fingertip whole blood results from three different clinical sites.

Site	Test	Results	Total	Mean of	Regression	Regression	Regression
	Strip Lot	Comparisons	Patients	Absolute	Analysis	Analysis	Analysis
	#	Compunsons	Tested	Bias (%)	Slope	Intercept	Coefficient
	"		Tested	Dias (70)	Stope	тистесрі	Variation
							(r)
Physician Office	Lot 1	Home User vs. YSI	171	8.0	1.01	3.12	0.96
(Site 1)		Home User vs. Professional	171	6.6	1.00	-0.98	0.96
		Professional vs. Predicate Device (One Touch)	175	9.3	0.91	15.03	0.95
		Professional vs. YSI	163	8.4	0.97	8.15	0.96
Physician Office	Lot 2	Home User vs. YSI	126	9.5	1.00	-1.59	0.93
(Site 2)		Home User vs. Professional	126	7.8	1.00	-1.01	0.95
		Professional vs. Predicate Device (One Touch)	130	8.0	0.98	5.8	0.92
		Professional vs. YSI	126	7.7	0.97	2.1	0.95
Physician Office	Lot 3	Home User vs. YSI	40	7.7	0.96	4.07	0.96
(Site 3)		Home User vs. Professional	40	5.4	1.00	-1.44	0.96
		Professional vs. Predicate Device (One Touch)	40	8.7	0.96	10.12	0.94
		Professional vs. YSI	40	6.4	0.96	6.39	0.95

Cholesterol

Method comparison studies were compared with the Abel-Kendall reference method performed in a CDC- certified Cholesterol Reference Method Network Laboratory (CRMLN). The Abell-Kendall method is performed with serum only. External studies were also done at three different clinical sites by Lay-users. A total of 456 patients from the three different sites participated in the clinical trial.

The lay-users performed their fingerstick and performed the cholesterol test on the Q.STEPS Biometer G/C System and the results were recorded by the lay-user and professional. The professional also performed a finger stick cholesterol test on the same lay-user using the same Biometer G/C System with the same lot of test strips. The professional from each site also drew a tube of venous blood from 14-16 lay-users in order to send the serum to CRMLN for comparison. Comparison results are presented in the table below.

Comparison Results from Three Different Clinical Sites

Site	Test	Results	Total	Regre	Regressio	Regression
	Strip	Comparisons	Patients	ssion	n	Analysis
	Lot	1	Tested	Analy	Analysis	Coefficient
	#		105000	sis	Intercept	Variation
	"			Slope	тистеері	(r)
Physician Office	Lot 1	Home User Performed / Read vs. Professional Performed/ Read	195	0.94	13.57	0.93
(Site 1)		Home User Performed/ Read vs. Home User Performed/ Read	195	1.00	0.00	1.00
		Home User Performed/ Read vs. CRMLN	36	0.93	15.10	0.92
		Professional Performed/ Read vs. CRMLN	36	0.87	26.15	0.91
Physician Office	Lot 2	Home User Performed / Read vs. Professional Performed/ Read	146	1.00	-1.84	0.95
(Site 2)		Home User Performed/ Read vs. Home User Performed/ Read	146	1.00	0.00	1.00
		Home User Performed/ Read vs. CRMLN	23	0.96	11.15	0.98
		Professional Performed/ Read vs. CRMLN	22	0.93	17.51	0.97
Physician Office	Lot 3	Home User Performed / Read vs. Professional Performed/ Read	115	0.91	16.59	0.91
(Site 3)		Home User Performed/ Read vs. Home User Performed/ Read	115	1.00	0.00	1.00
		Home User Performed/ Read vs. CRMLN	28	0.94	16.55	0.98
		Professional Performed/ Read vs. CRMLN	28	0.95	19.09	0.98

The bias between the CRMLN and Q.STEPS Biometer G/C Dual Monitoring System at the medical decision points of 200 mg/dL and 240 mg/dL, about 4.65% were positively misclassified and 2.33% were negatively classified.

b. Matrix comparison:

Not Applicable

3. Clinical studies:

a. Clinical sensitivity:

See above

b. Clinical specificity:

See above

c. Other clinical supportive data (when a and b are not applicable): See Comparison Studies referenced above.

4. Clinical cut-off:

Not Applicable

5. Expected values/Reference range:

Glucose

Patient glucose ranges for non-diabetic, non-pregnant adults:

	1 -6
Glucose fasting	70 mg/dL - 110 mg/dL
	3.9 mmol/L – 6.1 mmol/L
1 hour after meal ¹	< 160 mg/dL (<8.9 mmol/L)

1. Krall, LP and Deaser, RS: Joslin Diabetes Manual. Lea and Febiger. (1989) 138.

Cholesterol

According to NCEP classification, patient cholesterol ranges for non-diabetic, non-pregnant adults¹:

Desirable blood cholesterol	<200 mg/dL (<5.17 mmol/L)
Borderline-high blood cholesterol	200 mg/dL -239 mg/dL
	(5.17 mmol/L – 6.18 mmol/L)
High blood cholesterol	>240 mg/dL (>6.21 mmol/L)

1. National Cholesterol Education Program. Summary of the Second report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (Adults Treatment Panel II). JAMA 1993; 269:3015-23 Davidsohn & Henry, Clinical Diagnosis by Laboratory Methods. Todd-Sandford.

N. Proposed Labeling:

The labeling is sufficient and it satisfies the requirements of 21 CFR Part 809.10

O. Conclusion:

The submitted information in this premarket notification is complete and supports a substantial equivalence decision.