510(k) SUBSTANTIAL EQUIVALENCE DETERMINATION DECISION SUMMARY ASSAY ONLY TEMPLATE

A. 510(k) Number:

k060452

B. Purpose for Submission:

New Device

C. Measurand:

17-alpha-Hydroxyprogesterone

D. Type of Test:

Quantitative, Enzyme Immunoassay

E. Applicant:

Neo-Genesis

F. Proprietary and Established Names:

Accuwell 17-Alpha-Hyderoxyprogesterone Enzyme Immunoassay with Models 6015XX-ECAH

G. Regulatory Information:

Product Code	Classification	Regulation Section	Panel
JLX- Radioimmunoassay,	Class I, meets	<u>21 CFR 862.1395,</u>	75 Clinical
<u>17-Hydroxyprogesterone</u>	the limitation of	<u>17-Hydroxy-</u>	<u>chemistry</u>
	exemptions,	progesterone test	<u>(CH)</u>
	21 CFR 862.9	<u>system.</u>	
	(c), (2)	21 CFR 862.9 (c), (2)	

H. Intended Use:

1. <u>Intended use(s)</u>:

See indications for use below.

2. <u>Indication(s) for use:</u>

The Accuwell 17α -Hydroxyprogesterone EIA Kit is designed for the quantitative measurement of 17α -Hydroxyprogesterone (17-OHP) concentrations in neonatal blood samples that have been collected onto Whatman 903® specimen collection paper. The results are used to screen newborns for classical congenital adrenal hyperplasia (CAH).

- 3. <u>Special conditions for use statement(s):</u> For prescription use only
- 4. Special instrument requirements:

Microplate reader capable of reading a wavelength of 650 nm

I. Device Description:

The Accuwell 17a-Hydroxyprogesterone EIA Kit contains the following:

- Microplates coated with rabbit anti-17-OHP.
- Enzyme conjugate which contains 17-OHP derivative conjugated to horseradish peroxidase.
- Conjugate diluent which contains tris-buffer, bovine serum and preservatives.
- Wash buffer concentrate which is a concentrated solution of phosphate buffered saline containing a surfactant.
- Color developer which contains 3,3',5,5'-Tetramethylbenzidine in a diluted organic solvent with citrate buffer and hydrogen peroxide.
- Stopping reagent which is a dilute solution of sodium fluoride and a red dye.
- Multi-Analyte standards (0, 10, 25, 50, 100 and 250 ng/mL) serum equivalents. Each card of standards is made from human whole blood, adjusted to a hematocrit of 55% and each level is then spotted onto Whatman 903 specimen collection paper.
- Controls (15, 40 and 90 ng/mL) serum equivalent. Each control card is made from human whole blood, adjusted to a hematocrit of 55% and each level is then spotted onto Whatman specimen collection paper.

Human source material was tested and found negative for HIV 1 and 2, HBV and HCV using FDA approved methods.

J. Substantial Equivalence Information:

1. <u>Predicate device name(s):</u>

Wallac Oy AutoDelfia Neonatal 17 α -Hydroxyprogesterone Time-resolved fluoroimmunoassay

2. <u>Predicate 510(k) number(s):</u>

k042425

3. <u>Comparison with predicate:</u>

	Similarities										
Item	Device	Predicate									
Indications for use	Screening for increased	Screening for increased									
	levels of 17-OHP in	levels of 17-OHP in									
	newborns	newborns									
Sample Requirements	Neonatal blood collected	Neonatal blood collected									
	on Whatman 903 filter	on Whatman 903 filter									
	paper	paper									
Specimen	1/8 inch spot punched	1/8 inch spot punched									
	from a standard	from a standard									
	collection card of	collection card of									
	Whatman 903 paper	Whatman 903 paper									

	Differences			
Item	Device	Predicate		
Assay Type	Enzyme-immunoassay	Time-resolved		
		fluoroimmunoassay		
Detection Method	Peroxidase which	Europium ions		
	remains bound to micro-	dissociated from the		
	wells reacts with	labeled antiserum form		
	peroxide and TMB	highly fluorescent		
	subsequently converting	chelates with		
	the TMB from colorless	components of an		
	to a blue color which is	enhancement solution.		
	measured	Fluorescence in each		
		well is then measured		
Standard Range	3.8 to250 ng/mL	Up to 190 ng/mL		
	approximate serum	approximate serum		
	equivalent	equivalent		
Controls	Two levels	Three levels		
	concentration 14 and 40	concentration 15, 40 and		
	ng/mL	90 ng/mL		

K. Standard/Guidance Document Referenced (if applicable):

Area of Study	Reference Procedure	Reference Title
Stability Testing	CEN	Stability of In Vitro Diagnostic reagents (13640)
Interference	CLSI/NCCLS EP7-A	Interference Testing in Clinical Chemistry
Method Comparison	CLSI/NCCLS EP9-A	Method Comparison and Bias Estimation using Patient Samples
Precision	CLSI/NCCLS EP5-A	Evaluation of Precision Performance of Quantitative Measurement Methods

L. Test Principle:

The samples can be processed with two different methods; (1) a 1/8 inch disc from each blood spot is placed directly into designated wells of the coated microplate with an overnight incubation, (2) An elution is done using an 1/8 inch disc from each blood spot which is placed in the designated well of an uncoated microplate and saline is added to elute the blood from the dried blood spot. A portion of the eluate is then transferred into a designated well of the coated microplate with an overnight incubation or for 3 hours.

The Accuwell 17 α -hydroxyprogesterone immunoassay (EIA) is based on the competition principle. An unknown amount of antigen present in the sample and a fixed amount of enzyme labeled antigen compete for the binding sites of the antibodies coated onto the wells. After incubation the wells are washed to stop the competition reaction and the color developer is added. The color development is terminated by the addition of the stopping reagent and the absorbance measured at 650 nm and is inversely proportional to the amount of the antigen in the sample. Results of samples can be determined directly using the standard curve.

M. Performance Characteristics (if/when applicable):

- 1. Analytical performance:
 - a. Precision/Reproducibility:

Precision was assessed using three neonatal dried bloodspots prepared by the Centers for Disease Control and Prevention Neonatal Screening Quality Control Program. Samples were prepared using three different assay procedures: Over night Eluate (ON Eluate), Three hour Eluate (3 Hr Eluate) and Over Night Direct (ON Direct). Each of these samples was assayed in duplicate once a day for 20 days. Results are presented in the table below:

	ON Eluate			3 Hr Eluate			ON Direct		
Sample ID	CDC	CDC	CDC	CDC	CDC	CDC	CDC	CDC	CDC
	451	452	453	451	452	453	451	452	453

	Number	40	40	40	40	40	40	40	40	40
	Mean	28.2	56.1	108.1	26.2	53.0	104.6	29.5	57.9	108.6
	SD	3.0	4.9	12.3	3.4	4.1	11.9	4.2	5.8	13.2
	%CV	10.7	8.8	11.4	12.9	7.7	11.3	14.1	10.0	12.2
Within-	Sr	1.0	2.1	2.9	1.9	2.2	5.4	1.5	4.4	11.7
Run SD										
Between-	Sdd	2.9	4.5	12.1	2.9	3.4	10.7	3.9	3.9	6.2
Day SD										
Within	ST	3.1	5.0	12.5	3.4	4.1	12.0	4.2	5.8	13.3
Device										
SD										
Daily	В	3.0	4.8	12.3	3.1	3.8	11.4	4.1	5.0	10.4
Mean										
Standard										
Error										

b. Linearity/assay reportable range:

The linearity of the Accuwell 17-OHP measurement was demonstrated by spiking four serum samples with purified analyte to concentrations of 3.8, 10, 25 and 300 ng/mL. Samples were inter-diluted to make a total of twelve levels with a range between 3.8 - 300 ng/mL. Each sample was combined with lysed, commercially obtained, human, red blood cells to obtain an equivalent hematocrit of 55%. The samples were the spotted onto filter paper and dried. Each sample was tested in quadruplicate for each assay procedure and the average recovery was calculated and plotted against the targeted recoveries. Linear regression of comparison data yielded the following relationship:

ON Eluate - y = 1.0083x - 1.3865 r = 0.9998

3 Hr Eluate - y = 1.1015x - 4.7189 r = 0.997

ON Direct - y = 1.1146x - 4.9283 r = 0.997

The reportable range for the Accuwell 17-OHP measurements is 3.8-250 ng/mL

c. Traceability, Stability, Expected values (controls, calibrators, or methods):

The standards and controls are prepared using commercially available human serum stripped of the analyte then purified 17-OHP analyte is added gravimetrically to various concentrations. Each serum concentration is then mixed with lysed red blood cells to a hematocrit concentration equivalent to 55%. The blood mixtures are then spotted onto filter paper in $\frac{1}{2}$ inch circles. The spotted cards are randomly selected throughout the lot and the concentrations tested. Each new lot of cards is tested against a standard curve prepared from the previously approved lots.

The stability of the standards and controls were determined by studies that support the expiration date of 12 months.

d. Detection limit:

The limit of the blank was determined by calculating 2 SD from the mean of n=26 sample aliquots (zero standard) in one run on a single day using each of the three assay procedures. The analytical sensitivity was estimated to be 2.2 ng/mL for the Over night Eluate, 1.5 ng/mL 3 Hour Eluate and 2.3 ng/mL Over Night Direct.

e. Analytical specificity:

Studies were performed to assess common or known substances that could interfere with the method. Hemoglobin up to 1200 mg/dL, bilirubin, conjugated and unconjugated up to 300 mg/dL, and lipids up to 1875 mg/dL caused no detectable interference.

Cross-reactivity was assessed by testing compounds whose structures could potentially cause interference. The lists of compounds tested are listed below:

Substance	Cross	Substance	Cross
	Reactivity		Reactivity
21-Desoxycortisol	2.4 %	Cholesterol	<0.01 %
16α-	1.2 %	Corticosterone	<0.01 %
Hydroxyprogestrone			
11-Desoxycortisol	0.6 %	Cortisol Glucuronide	<0.01 %
Progesterone	0.5 %	Cortisone	<0.01 %
5-Pregnen-3β,17-Diol-	0.2 %	Cortisone 21-Sulphate	<0.01 %
20-one 3 sulfate		Sodium Salt	
Cortisol	0.06 %	Dehydroisoandrosterone	<0.01 %
Desoxycorticosterone	0.06 %	Dehydroepiandrosterone	<0.01 %
		Sodium Sulphate	
1-Dehydrotestosterone	<0.01 %	Dexamethasone	<0.01 %
5β-Dihydrocortisol	<0.01 %	Estriol	<0.01 %
5β-Dihydrocortisone	<0.01 %	Estrone	<0.01 %
S6β-Hydroxycortisol	<0.01 %	Prednisolone	<0.01 %
t 11-	<0.01 %	Prednisone	<0.01 %
Dehydrocorticosterone			
d16α-	<0.01 %	Pregnenolone	<0.01 %
Hydroxypregnenolone			
e17α-Estradiol	<0.01 %	Pregnenolone Sulphate,	<0.01 %
S		Sodium Salt	
17α-	<0.01 %	Spironolactone	<0.01 %
wHydroxypregnenolone			
e17β-Estradiol	<0.01 %	Testosterone	<0.01 %
r20α-	<0.01 %	Tetrahydrocortisol	<0.01 %
eHydroxyprogesterone			
Aldosterone	<0.01 %	Tetrahydrocorisone	<0.01 %

p e

f. Assay cut-off:

See expected values below.

- 2. Comparison studies:
 - a. Method comparison with predicate device:

A retrospective study was conducted to compare the results obtained with the Accuwell 17-OHP EIA to those obtained by a currently marketed neonatal 17-OHP screening device. Test samples were submitted to the study as blinded neonatal dried blood spots collected in sequence under routine screening conditions from a U.S. department of public health laboratory. Original screening results for each sample using the predicate test kit were also obtained from the submitting laboratory for method comparison.

A summary of results of the method comparisons are provided in Tables 9 – 20, below. (Samples interpreted below as "Follow-Up" were not confirmed as positive for CAH.)

Table 9:Predicate vs. Accuwell O/N DirectValues at the 99th, 97.5th, and 95th Percentiles for aPopulation of Babies \geq 2500 gm and 0-1 Days of Age at Sample Collection

		99t	h Percen	tile		97.5th Percentile					95th Percentile			
	Accuwell	Normal	Follow- up			Accuwell	Normal	Follow- up			Accuwell	Normal	Follow- up	
Predicate		< 27.5	≥27.5	Totals:	Predicate		< 25.8	≥25.8	Totals:	Predicate		< 21.8	≥21.8	Totals:
Normal	<43.8 N=130	114	16	130	Normal	<43.8 N=130	112	18	130	Normal	<43.8 N=130	100	30	130
Follow- Up	\geq 43.8 N = 3	0	3	3	Follow- Up	\geq 43.8 N = 3	0	3	3	Follow- Up	\geq 43.8 N = 3	0	3	3
	Totals:	114	19	133		Totals:	112	21	133		Totals:	100	33	133

The Accuwell 17-OHP EIA (ON Direct) results for this population were: mean =16.5 ng/mL, with a range of 3.7 to 85.6 ng/mL; while the predicate kit yielded: mean =22.6 ng/mL, with a range of 10.1 to 86.0 ng/mL

The correlation was found to be: y (Predicate) = 0.970 (Accuwell ON Direct) + 6.569, R = 0.9570

		99t	h Percen	tile		97.5th Percentile						95t	h Percen	tile
	Accuwell	Normal	Follow- up			Accuwell	Normal	Follow- up			Accuwell	Normal	Follow- up	
Predicate		< 23.9	≥23.9	Totals:	Predicate		< 22.3	≥ 22.3	Totals:	Predicate		< 20.2	\geq 20.2	Totals:
Normal	<43.8 N=130	106	24	130	Normal	<43.8 N=130	103	27	130	Normal	<43.8 N=130	101	29	130
Follow- Up	\geq 43.8 N = 3	0	3	3	Follow- Up	\geq 43.8 N = 3	0	3	3	Follow- Up	\geq 43.8 N = 3	0	3	3
	Totals:	106	27	133		Totals:	103	30	133		Totals:	101	32	133

Table 10: Predicate vs. <u>Accuwell O/N Eluate</u>

Values at the 99th, 97.5th, and 95th, Percentiles for a

Population of Babies >2500 gm and 0-1 Days of Age at Sample Collection

The Accuwell 17-OHP EIA (ON Eluate) results for this population were: mean =16.7 ng/mL, with a range of 5.3 to 65.6 ng/mL; while the predicate kit yielded: mean =22.6 ng/mL, with a range of 10.1 to 86.0 ng/mL.

The correlation was found to be: y (Predicate) = 1.049 (Accuwell ON Direct) + 5.066, R =

	Population	on of Ba	abies <u>></u> 2	500 gn	n and 0-1	Days of A	Age at S	Sample (Collecti	on				
		99t	h Percen	tile		97.5th Percentile					95th Perce			
	Accuwell	Normal	Follow-			Accuwell	Normal	Follow-			Accuwell	Normal	Follov	
		1.0111101	up					up				1.0111101	up	
Predicate		< 23.0	\geq 23.0	Totals:	Predicate		< 21.7	\geq 21.7	Totals:	Predicate		< 19.8	≥ 19.	
Normal	<43.8 N=130	104	26	130	Normal	<43.8 N=130	103	27	130	Normal	<43.8 N=130	99	31	
Follow- Up		0	3	3	Follow- Up	\geq 43.8 N = 3	0	3	3	Follow- Up		0	3	
	Totals:	104	29	133		Totals:	103	30	133		Totals:	99	34	

Table 11: Predicate vs. Accuwell 3 Hour EluateValues at the 99th, 97.5th, and 95th Percentiles for a

The Accuwell 17-OHP EIA (3Hr Eluate) results for this population were: mean =17.3 ng/mL, with a range of 5.1 to 71.1 ng/mL; while the predicate kit yielded: mean =22.6 ng/mL, with a range of 10.1 to 86.0 ng/mL.

The correlation was found to be: y (Predicate) = 0.998 (Accuwell ON Direct) + 5.410, R = 0.9291

Population of I	Valu Babies <u>></u> 25	es at the 9 00 gm and	5 th Perce I 2-3 Da j	entile for ys of Ag	a e at San	nple Collection
		tile				
		Accuwell				
	Predicate		< 18.5	≥18.5	Totals:	
	Normal	<29.2 N=118	113	5	118	
	Follow- Up	≥ 29.2 N = 15	0	15	15	
		Totals:	113	20	133	

Table 12: Predicate vs. Accuwell O/N Direct

The Accuwell 17-OHP EIA (ON Direct) results for this population were: mean =11.3 ng/mL, with a range of 3.4 to 81.1 ng/mL; while the predicate kit yielded: mean =16.2 ng/mL, with a range of 6.9 to 81.0 ng/mL.

The correlation was found to be: y (Predicate) = 1.045 (Accuwell ON Direct) + 4.392, R = 0.9785

	95th Percentile										
	Accuwell	Normal	Follow- up								
Predicate		≥18.7	Totals:								
Normal	<29.2 N=118	114	4	118							
Follow- Up	≥ 29.2 N = 15	0	15	15							
	Totals:	114	19	133							

Table 13: Predicate vs. Accuwell O/N Eluate Values at the 95th Percentile for a Population of Babies ≥2500 gm and 2-3 Days of Age at Sample Collection

The Accuwell 17-OHP EIA (ON Eluate) results for this population were: mean =12.2 ng/mL, with a range of 3.4 to 67.5 ng/mL; while the predicate kit yielded: mean =16.2 ng/mL, with a range of 6.9 to 81.0 ng/mL.

The correlation was found to be: y (Predicate) = 1.156 (Accuwell ON Eluate) + 2.030, R = 0.9729

Table 14: Predicate vs. Accuwell 3 Hour EluateValues at the 95th Percentile for a

Population of Babies **>2500 gm and 2-3 Days** of Age at Sample Collection

	95th Percentile								
	Accuwell Normal		Follow- up						
Predicate		< 17.7	≥17.7	Totals:					
Normal	<29.2 N=118	115	3	118					
Follow- Up	≥ 29.2 N = 15	0	15	15					
	Totals:	115	18	133					

The Accuwell 17-OHP EIA (3Hr Eluate) results for this population were: mean =12.8 ng/mL, with a range of 4.5 to 69.6 ng/mL; while the predicate kit yielded: mean =16.2 ng/mL, with a range of 6.9 to 81.0 ng/mL.

The correlation was found to be: y (Predicate) = 1.150 (Accuwell ON Direct) + 1.495, R = 0.9691

Table 15:Predicate vs. Accuwell O/N Direct

Values at the 97.5th and 95th Percentiles for a Population of **Babies** \geq 2500 gm and \geq 4 Days of Age at Sample Collection

	97.5th Percentile						95t	h Percen	tile
	Accuwell	Normal	Follow- up			Accuwell	Normal	Follow- up	
Predicate		< 14.8	≥14.8	Totals:	Predicate		< 11.5	≥11.5	Totals:
Normal	< 21.9 N=192	181	11	192	Normal	< 21.9 N=192	170	22	192
Follow- Up	≥ 21.9 N = 5	0	5	5	Follow- Up	$ \ge 21.9 $ N = 5	0	5	5
	Totals:	181	16	197		Totals:	170	27	197

The Accuwell 17-OHP EIA (ON Direct) results for this population were: mean =7.3 ng/mL, with a range of 2.2 to 22.6 ng/mL; while the predicate kit yielded: mean =9.5 ng/mL, with a range of 2.7 to 31.0 ng/mL.

The correlation was found to be: y (Predicate) = 1.090 (Accuwell ON Direct) + 1.565, R = 0.9538

Table 16: Predicate vs. <u>Accuwell O/N Eluate</u>	
Values at the 97.5 th and 95^{th} Percentiles for a	
Population of Babies \geq 2500 gm and \geq 4 Days of Age at Sample Collection	tion

	97.5th Percentile						95t	h Percen	tile
	Accuwell	Normal	Follow- up			Accuwell	Normal	Follow- up	
Predicate		< 15.2	≥15.2	Totals:	Predicate		< 13.6	≥13.6	Totals:
Normal	< 21.9 N=192	181	11	192	Normal	< 21.9 N=192	178	14	192
Follow- Up	≥ 21.9 N = 5	1	4	5	Follow- Up	≥ 21.9 N = 5	0	5	5
	Totals:	182	15	197		Totals:	178	19	197

The Accuwell 17-OHP EIA (ON Eluate) results for this population were: mean =7.8 ng/mL, with a range of 1.4 to 26.5 ng/mL; while the predicate kit yielded: mean =9.5 ng/mL, with a range of 2.7 to 31.0 ng/mL.

The correlation was found to be: y (Predicate) = 1.070 (Accuwell ON Eluate) + 1.172, R = 0.9142

97.5th Percentile							95t	h Percen	tile
	Accuwell	Normal	Follow- up			Accuwell	Normal	Follow- up	
Predicate		< 15.6	≥15.6	Totals:	Predicate		< 14.1	≥14.1	Totals:
Normal	< 21.9 N=192	176	16	192	Normal	< 21.9 N=192	174	18	192
Follow- Up	≥ 21.9 N = 5	0	5	5	Follow- Up	≥ 21.9 N = 5	0	5	5
	Totals:	176	21	197		Totals:	174	23	197

Table 17: Predicate vs. Accuwell 3 Hour Eluate
Values at the 97.5th and 95th Percentiles for aPopulation of Babies \geq 2500 gm and \geq 4 Days of Age at Sample Collection

The Accuwell 17-OHP EIA (3Hr Eluate) results for this population were: mean =8.3 ng/mL, with a range of 2.2 to 28.4 ng/mL; while the predicate kit yielded: mean =9.5 ng/mL, with a range of 2.7 to 31.0 ng/mL.

The correlation was found to be: y (Predicate) = 1.024 (Accuwell 3Hr Eluate) + 1.081, R = 0.9350

		99th (1	00th) Per	centile			97.5	th Percer	ntile
	Accuwell	Normal	Follow- up			Accuwell	Normal	Follow- up	
Predicate		< 46.0	\geq 46.0	Totals:	Predicate		< 43.9	\geq 43.9	Totals:
Normal	< 45.1 N=28	28	0	28	Normal	< 40.3 N=26	26	0	26
Follow- Up	\geq 45.1 N = 2	1	1	2	Follow- Up	\geq 40.3 N = 4	3	1	4
	Totals:	29	1	30		Totals:	29	1	30

Table 18: Predicate vs. Accuwell O/N DirectValues at the 99th and 97.5th Percentiles for aPopulation of Babies 1400 - 2500 gm and 1-6 Days of Age at Sample Collection

The Accuwell 17-OHP EIA (ON Direct) results for this population were: mean =19.3 ng/mL, with a range of 6.5 to 47.8 ng/mL; while the predicate kit yielded: mean =25.5 ng/mL, with a range of 9.0 to 65.0 ng/mL.

The correlation was found to be: y (Predicate) = 1.22 (Accuwell ON Direct) + 1.971, R = 0.9556

-	reputation of Buoles 1100 2000 gin and 1 0 Buys of Age a Sumple Concerton								
99th (100th) Percentile							97.5	5th Perce	ntile
	Accuwell	Normal	Follow- up			Accuwell	Normal	Follow- up	
Predicate		< 31.7	≥ 31.7	Totals:	Predicate		< 28.4	≥ 28.4	Totals:
Normal	<45.1 N=28	24	4	28	Normal	< 40.3 N=26	23	3	26
Follow- Up	≥ 45.1 N = 2	1	1	2	Follow- Up	≥ 40.3 N = 4	1	3	4
	Totals:	25	5	30		Totals:	24	6	30

Table 19: Predicate vs. <u>Accuwell O/N Eluate</u>
Values at the 99 th and 97.5 th Percentiles for a
Population of Babies 1400 - 2500 gm and 1-6 Days of Age at Sample Collection

The Accuwell 17-OHP EIA (ON Eluate) results for this population were: mean =19.2 ng/mL, with a range of 6.7 to 40.2 ng/mL; while the predicate kit yielded: mean =25.5 ng/mL, with a range of 9.0 to 65.0 ng/mL.

The correlation was found to be: y (Predicate) = 1.234 (Accuwell ON Direct) + 1.812, R = 0.8770

	1			0		e			
	99th (100th) Percentile					97.5th Percentile			
	Accuwell	Normal	Follow- up			Accuwell	Normal	Follow- up	
Predicate		< 33.3	≥ 33.3	Totals:	Predicate		< 27.7	≥ 27.7	Totals:
Normal	<45.1 N=28	26	2	28	Normal	<40.3 N=26	23	3	26
Follow- Up	≥ 45.1 N = 2	1	1	2	Follow- Up	≥ 40.3 N = 4	2	2	4
	Totals:	27	3	30		Totals:	25	5	30

Table 20: Predicate vs. <u>Accuwell 3Hr Eluate</u>Values at the 99th and 97.5th Percentiles for a

Population of Babies 1400 - 2500 gm and 1-6 Days of Age at Sample Collection

The Accuwell 17-OHP EIA (3Hr Eluate) results for this population were: mean =19.3 ng/mL, with a range of 6.6 to 39.2 ng/mL; while the predicate kit yielded: mean =25.5 ng/mL, with a range of 9.0 to 65.0 ng/mL.

The correlation was found to be: y (Predicate) = 1.174 (Accuwell 3Hr Eluate) + 2.795, R = 0.8421

b. Matrix comparison:

Not applicable.

3. Clinical studies:

a. Clinical Sensitivity:

Not applicable.

b. Clinical specificity:

Not applicable.

c. Other clinical supportive data (when a. and b. are not applicable):

See the method comparison study above and the expected values below.

4. Clinical cut-off:

See expected values below.

5. Expected values/Reference range:

In a separate study of 444 infants whose samples were collected for testing at 1-3days of age and whose birth weights were ≥ 2500 gm, 17-OHP concentrations were measured using each of the Accuwell 17-OHP methods. Samples were collected on S&S filter paper no. 903 and stored for up to 7 months prior to Accuwell testing. Based on the results obtained, the sponsor provides the following guidelines in the labeling for samples taken at 1-3 days of age from babies with birth weights ≥ 2500 gm. The data show the cut-off values at which re-testing or other follow-up should be considered.

Accuwell O/	17-0	HP (ng/	/mL serum)		
		Percentiles			
Age (days)	n	95 th	97.5 th	99 th	
1-3	444	20.8	25.0	27.8	

Accuwell Expected Values at the 95 th , 97.5 th and 99 th
Percentiles for a Population of Babies >2500 gm and 1-3
Days of Age at Sample Collection

Accuwell O/	17-0	HP (ng/	mL serum)	
	Percentiles			
Age (days)	n	95 th	97.5 th	99 th
1-3	19.1	21.7	23.9	

Accuwell 3 hr Eluate		17-OHP (ng/mL serum)		
		Percentiles		
Age (days)	n	95 th	97.5 th	99 th
1-3	444	19.3	21.5	24.3

In the same study, a small sample (n = 91) of infants with birth weights < 2500 gm were tested using the Accuwell 17-OHP assay methods. The study showed the following:

Accuwell Expected Values at the 95th, 97.5th and 99th Percentiles for a Population of Babies <2500 gm and 1-3 Days of Age at Sample Collection

Accuwell O/N Direct		17-OHP (ng/mL serum)			
			Percentiles		
Age (days)	Ν	95 th	97.5 th	99 th	
1-3	91	33.9	43.9	46.4	

Accuwell O/N Eluate		17-OHP (ng/mL serum)		
		Percentiles		
Age (days)	n	95 th	97.5 th	99 th
1-3	91	25.7	31.7	40.9

Accuwell 3 Hr Eluate		17-OHP (ng/mL serum)		
		Percentiles		
Age (days)	n	95 th	97.5 th	99 th
1-3	91	26.3	32.5	40.5

In a study of 391 infants whose samples were collected for testing at \geq 4 days of age and whose birth weights were \geq 2500 gm, 17-OHP concentrations were measured using each of the Accuwell 17-OHP methods. Samples were collected on S&S filter paper no. 903 and stored for up to 7 months prior to Accuwell testing. Based on the results obtained, the following guidelines can be given for samples taken at \geq 4 days of age from babies with birth weights \geq 2500 gm. The data show the cut-off values at which re-testing or other follow-up should be considered.

Accuwell Expected Values at the 95^{th} , 97.5^{th} and 99^{th} Percentiles for a Population of Babies ≥ 2500 gm and 1-3 Days of Age at Sample Collection

Accuwell O/N Direct		17-OHP (ng/mL serum)		
		Percentiles		
Age (days)	n	95 th	97.5 th	99 th
<u>></u> 4	391	11.5	14.8	21.1

Accuwell O/N Eluate		17-OHP (ng/mL serum)				
		Percentiles				
Age (days)	n	95 th	97.5 th	99 th		
<u>></u> 4	391	13.6	15.2	19.5		
Accuwell 3 hr Eluate		17-OHP (ng/mL serum)				
		Percentiles				
Age (days)	n	95 th	97.5 th	99 th		
<u>></u> 4	391	14.1	15.6	19.4		

The values shown here should be used only as a guideline, and were calculated based on a specific population of samples as defined by the manufacturer. The sponsor recommends in the labeling that each laboratory should establish its own cut-off values for 17-OHP in its own population.

N. Proposed Labeling:

The labeling is sufficient and it satisfies the requirements of 21 CFR Part 809.10.

O. Conclusion:

The submitted information in this premarket notification is complete and supports a substantial equivalence decision.