510(k) SUBSTANTIAL EQUIVALENCE DETERMINATION DECISION SUMMARY

A. 510(k) Number: k062035

- **B.** Purpose for Submission: New Device
- C. Measurand: Rheumatoid Factor
- **D. Type of Test:** Quantitative, nephelometry

E. Applicant:

Dade Behring Inc.

F. Proprietary and Established Names:

Dimension VistaTM Rheumatoid Factors Flex ® reagent cartridge Dimension VistaTM Protein 2 Calibrator Dimension VistaTM Protein 2 Control L Dimension VistaTM Protein 2 Control H

G. Regulatory Information:

- 1. <u>Regulation section:</u>
 - 21 CFR § 866.5775 Rheumatoid Factor Immunological test system
 - 21 CFR § 862.1150 Calibrator
 - 21 CFR § 862.1660 Quality Control Material
 - 2. <u>Classification:</u>

Class II

3. <u>Product code:</u>

DHR - System, Test, Rheumatoid Factor

JIX - Calibrator, multi-analyte mixture

JJY - Multi-analyte controls, all kinds (assayed and unassayed)

- 4. <u>Panel:</u>
 - Immunology (82)

H. Intended Use:

1. Intended use(s):

<u>Flex® reagent cartridge</u>: The RF method is an in vitro diagnostic reagent for the quantitative determination of rheumatoid factors (RF) in human serum and lithium heparin plasma on the Dimension VistaTM System. Measurements are used as an aid in the diagnosis of rheumatoid arthritis.

Protein 2 Calibrator: PROT2 CAL is an in vitro diagnostic product for the calibration of the C-reactive protein (CRP), high sensitivity CRP (hsCRP) and Rheumatoid Factors (RF) methods on the Dimension VistaTM System. Protein 2 Controls L and H: PROT2 CON L and H are assayed intra-laboratory quality controls for the assessment of precision and analytical bias in the determination of C-reactive protein (CRP) and Rheumatoid Factor (RF) on the Dimension VistaTM System.

 2. <u>Indication(s) for use:</u> Same as Intended Use.

- 3. <u>Special conditions for use statement(s):</u> Prescription use only
- 4. <u>Special instrument requirements:</u> Dimension VistaTM System

I. Device Description:

<u>Dimension Vista TM Rheumatoid Factors Flex ® regent cartridge (RF)</u>: RF reagent consists of polystyrene particles coated with an immunocomplex consisting of human-IgG/anti-human IgG from sheep and RF supplement reagent consisting of a phosphate buffer and polyethylene glycol.

Dimension Vista TM Protein 2 Calibrator:

PROT2 CAL is a liquid, human serum based product containing C-reactive protein and Rheumatoid Factors.

Dimension Vista TM Protein 2 Control L and H:

PROT2 CON L and H are liquid, multi-analyte, human serum based products containing C-reactive protein and Rheumatoid Factors in two concentration ranges: low and high.

J. Substantial Equivalence Information:

- Predicate device name(s) and 510(k) numbers:
ROCHE Tina-Quant RF II assay
N Rheumatology Standard SL
N/T Rheumatology Control SLk964527
k962373
- 2. <u>Comparison with predicate:</u>

Dimension Vista TM Rheumatoid Factors Flex ® regent cartridge

Similarities				
Item	Device	Predicate		
Intended Use	For the quantitative	For the quantitative in		
	determination of	vitro determination of		
	rheumatoid factors (RF)	rheumatoid factors in		
	in human serum and	human serum and plasma		
	lithium heparin plasma			
Indications for Use	Measurements are used	Same		
	as an aid in the diagnosis			
	of rheumatoid arthritis.			
Measurement type	Quantitative	Same		
Stability: Unopened	90 days 2-8°C	Same		

Differences				
Item	Item Device			
Method	Immunonephelometry	Immunoturbidimetry		
Capture Antibody	Immune complexes of	Human IgG		
	human immunoglobulin			
	and sheep anti-human			
	IgG antibody			
Instrument system	Dimension Vista TM	ROCHE/Hitachi 917		
System		analyzer		

Differences			
Item	Device	Predicate	
Plasma	lithium heparin plasma	lithium heparin plasma and EDTA plasma	
Measuring Range	10-600 IU/mL	7-650 IU/mL	
Expected values	<15 IU/mL	<14 IU/mL	
Stability: Open	21 days 2-8°C	90 days 2-8°C	

Dimension Vista TM Protein 2 Calibrator

Similarities				
Item	Device	Predicate		
Intended Use	For the calibration of the	For establishment of		
	specified methods	reference curves for		
		specific		
		immunonephelometric		
		determinations		
Composition	Pooled Human sera	Same		
Reagent Preparation	Liquid, Ready-for-Use	Same		
Traceability	1 st British Standard	Same		
	64/002			

Differences				
Item	Device	Predicate		
Analytes	C-reactive protein (CRP),	Rheumatoid factors (RF),		
	high sensitivity CRP	anti-Streptolysin O		
	(hsCRP) and Rheumatoid	(ASL) and C-reactive		
	Factors (RF)	protein (CRP)		
Instrument system	Dimension Vista TM	BN Systems		
	System			
Stability: Open	12 days 2-8°C	14 days 2-8°C		
Materials Provided	6 vials, 2 mL per vial	3 vials, 1 mL per vial		

Dimension Vista TM Protein 2 Control L and H

Similarities				
Item	Device	Predicate		
Intended Use	Assayed intra-laboratory quality controls for the assessment of precision and analytical bias in the determination of specific assays	Assayed controls for accuracy and precision in the quantitative determination of specific analytes		
Composition	Pooled human sera	Same		
Reagent Preparation	Liquid, Ready to use	Same		
Concentration Range	Low and High	Same		

Differences				
Item	Device	Predicate		
Analytes	C-reactive protein (CRP) and Rheumatoid Factor (RF)	Rheumatoid factors (RF), anti-streptolysin O (ASL) and C-reactive protein (CRP)		
Instrument system	Dimension Vista [™] System	BN Systems and the TurbiTime System		
Stability Open	12 days 2-8°C	14 days 2-8°C		
Materials Provided	6 vials, 2 mL per vial	3 vials, 1 mL per vial		

K. Standard/Guidance Document Referenced (if applicable):

CLSI Evaluation of Precision Performance of Clinical Devices Approved Guideline: EP5-A2; CLSI Method Comparison and Bias Estimation Using Patient Samples Approved Guideline EP9-A2; CLSI Interference Testing in Clinical Chemistry Approved Guideline EP7-A2.

L. Test Principle:

Polystyrene particles coated with an immunocomplex consisting of human-IgG/antihuman IgG from sheep are aggregated when mixed with samples containing RF. These aggregates scatter a beam of light passed through the sample. The intensity of the scattered light is proportional to the concentration of the respective protein in the sample. The result is evaluated by comparison with a standard of known concentration.

M. Performance Characteristics (if/when applicable):

- 1. Analytical performance:
 - a. Precision/Reproducibility:

Precision testing was performed in accordance with CLSI EP5-A2. Four samples derived from serum pools representing RF concentrations that span the measuring range, and the PROT2 CON L and H were analyzed in duplicate, for 20 days, two times per day. The repeatability and within-lab standard deviations (SD) and percent coefficient of variation (%CV) were calculated. Repeatability imprecision was <6%, and within-lab imprecision was $\leq 6.2\%$.

Motorial	Mean	Repeatability		Within-Lab	
Material	IU/mL	SD	%CV	SD	%CV
PROT2 CON L	73.7	1.5	2.0	1.6	2.2
PROT2 CON H	183.9	3.5	1.9	3.9	2.1
Serum pool* VL	14.5	0.9	5.9	0.9	6.2
Serum pool VL	30.5	0.8	2.6	1.8	5.9
Serum pool L	90.0	2.1	2.3	2.3	2.6
Serum pool H	544.3	21.7	4.0	23.9	4.4

*Tested for 5 days

b. Linearity/assay reportable range:

Linearity across the assay range of 10-600 IU/mL was confirmed by testing a

calibrator with high concentration of RF. The calibrator was serially diluted in 7, two-fold increments (measuring range from 8.36 to 836.86 IU/mL). Each dilution was tested in replicates of three. Data were analyzed using linear regression analysis with the following results:

Slope	Y-Intercept IU/mL	Correlation Coefficient	n
1.040	-2.348	1.0	21

A recovery study using 1^{st} British Standard 64/002 (5 replicates) showed that % recovery ranged from 96.4 to 98.6% with a mean recovery of 97.9%

c. Traceability, Stability, Expected values (controls, calibrators, or methods): Calibration material values were referenced to the 1st British Standard 64/002.

Limit of Detection: 10 IU/mL represents the lower limit of the reportable range of the RF assay.

Analytical sensitivity was determined by assaying analyte-free System Diluent 20 times and determining the mean value plus two SDs: 0.814 IU/mL.

- e. Analytical specificity:
 - i. <u>Interference Studies</u>: Interference testing was performed according to CLSI EP7-A2, to determine the effect of various exogenous substances on the Dimension VistaTM RF assay. The following exogenous interferents were tested by spiking into test samples containing two concentrations of RF (one at the clinical decision point 12-15 IU/mL and one higher), and compared to a control sample without interferent: bilirubin (conjugated and unconjugated, 60 mg/dL), hemoglobin (1000 mg/dL), creatinine (30 mg/dL), albumin (6 g/dL), urea (500 mg/dL), cholesterol (500 mg/dL), uric acid (20 mg/dL). Minimal interference (<8%) was observed. Interference by triglycerides (native and clarified) was tested on 5 test samples containing RF concentrations between 12 and 25 IU/mL, and minimal interference (<7%) was observed. Additionally, 41 potentially interfering drugs were also assayed and shown to exhibit minimal interference (<10%).
 - ii. <u>Cross-reactivity</u> with autoimmune antibodies common to other systemic autoimmune diseases was not tested.
 - iii. <u>Antigen Excess</u>: the effect of antigen excess was evaluated using a serum sample above the assay range. No effect was seen up to 2877 IU/mL. Samples above this level are reported as exceeding assay range (>600 IU/mL).
- f. Assay cut-off:

The assay cut-off of 15 IU/mL was determined by measuring 300 normal blood donors from Europe and taking the mean between the 95^{th} and 97.5^{th} percentile value.

- 2. Comparison studies:
 - a. *Method comparison with predicate device:* The Dimension Vista[™] RF assay was compared to the ROCHE Tina-Quant

d. Detection Limit

RF II assay on the Hitachi 917 Analyzer by assaying 120 serum samples with concentrations ranging form 10-520 IU/mL on both systems. Passing Bablock regression analysis of these results yielded the following statistics:

Comparative	Slope	Y-Intercept (IU/mL)	Correlation	n
Method	(95%CI)	(95%CI)	Coefficient	
ROCHE Tina-Quant RF II/Hitachi 917	0.953 (0.898 to 1.024)	-9.71 (-16.192 to -6.021)	0.95	120

b. Matrix comparison:

Fourteen matched serum and lithium heparin plasma samples, containing concentrations of RF spanning the measuring range (10-515 IU/mL), were assayed and compared. Regression analysis showed no significant bias between the two matrices and yielded the following information:

% recovery in Li-Hep	Slope	Y-Intercept	Correlation	n
plasma vs. serum	(95%CI)	(IU/mL)	Coefficient	
Mean = -2.5% (Range = -9.7% to 2.3%)	0.98 (0.94-1.01)	-0.04	1.00	14

3. <u>Clinical studies</u>:

- *a. Clinical Sensitivity:* Not provided
- *b. Clinical specificity:* Not provided
- *c. Other clinical supportive data (when a. and b. are not applicable):* Not applicable.
- 4. <u>Clinical cut-off</u>: See Assay cut-off

5. Expected values/Reference range:

Testing of serum from 300 European blood donors resulted in a 95th percentile of <10 IU/ml and a 97.5th percentile of 19.5 IU/ml. The expected value of the normal healthy adults is <15 IU/ml.

N. Proposed Labeling:

The labeling is sufficient and it satisfies the requirements of 21 CFR Part 809.10.

O. Conclusion:

The submitted information in this premarket notification is complete and supports a substantial equivalence decision.