510(k) SUBSTANTIAL EQUIVALENCE DETERMINATION DECISION SUMMARY INSTRUMENT ONLY TEMPLATE

A. 510(k) Number:

K063407

B. Purpose for Submission:

This is an original-traditional 510(k) for a 16 parameter automated hematology analyzer and leukocyte differential cell counter with reagents, calibrator, and controls.

C. Manufacturer and Instrument Name:

Shenzhen Mindray Bio-Medical Electronics Co. LTD, BC-3200 Auto Hematology Analyzer

D. Type of Test or Tests Performed:

WBC, Lymph #, Mid#, Gran#, Lymph%, Mid%, Gran%, RBC, HGB, MCV, MCH, MCHC, RDW, HCT, MPV, and WBC Histogram, RBC Histogram, and PLT Histogram

E. System Descriptions:

1. Device Description:

The BC-3200 Auto Hematology Analyzer is a quantitative, automated hematology analyzer and leukocyte differential cell counter for In Vitro Diagnostic use in clinical laboratories. It is only to be used by trained medical professionals to identify the normal patient, with all normal system-generated parameters, and toflag or identify patient results that require additional studies. The analyzer provides analysis results of 16 parameter of human blood and three histograms.

2. Principles of Operation:

WBCs are counted and sized by the impedence method. This method is based on the measurement of changes in electrical resistance produced by a particle, which in this case is a blood cell suspended in a conductive diluent as it passes through an aperture of known dimensions. HGB is determined by the colorimetric method. RBCs and PLTs are counted by the impedence method also. In addition, for RBCs and Plts, volumetric metering is used. An accurate cell count cannot be obtained unless the precise volume of diluted sample that passes through the aperture during the count cycle is known. The analyzer uses a volumetric metering unit to control the count cycle and to ensure that a precise volume of sample is analyzed for the measurement.

3. Modes of Operation:

Closed Vial Whole Blood Mode, Whole Blood Mode for veinous blood, and Predilute Mode for capillary blood.

4. Specimen Identification:

Barcode or manual keyboard entry.

5. Specimen Sampling and Handling:

Samples are manually loaded into a sample compartment one at a time. The BC-3200 utilizes an automatic sampling and mixing device for sample processing. The Mindray calibrator is called SC-CAL PLUS.

6. Calibration:

The device has two calibration programs: manual calibration and auto calibration using commercial calibrators. The Mindray calibrator is called SC-CAL PLUS.

7. Quality Control:

The device has two QC programs: L-J Analysis and X-B Analysis. The Mindray three level control is called BC-3D.

8. Software:

FDA has reviewed applicant's Hazard Analysis and Software Development processes for this line of product types:

Yes <u>X</u> or No Mindray has provided software documentation at a moderate level of concern that conforms to the FDA software guidance document.

F. Regulatory Information:

1. <u>Regulation section:</u>

21 CFR 864.5220, Automated differential cell counter

2. Classification:

Class II

3 <u>Product code:</u>

GKZ

4. <u>Panel:</u>

Hematology (81)

G. Intended Use:

1. Indication(s) for Use:

The BC-3200 Auto Hematology Analyzer is a quantitative, automated hematology analyzer and leukocyte differential cell counter to be used in clinical laboratories for In Vitro Diagnostic Use. The intended use of the BC-3200 Auto Hematology analyzer is to identify the normal patient, with all normal systemgenerated parameters, and to flag or identify patient results that require additional studies.

2. <u>Special Conditions for Use Statement(s):</u>

N/A

H. Substantial Equivalence Information:

1. Predicate Device Name(s) and 510(k) numbers:

COULTER® A^CT diff 2[™] Analyzer, K0990352

2. <u>Comparison with Predicate Device:</u>

	Similarities								
Item	Device	Predicate							
Intended Use	The BC-3200 auto hematology analyzer is a quantitative, automated hematology analyzer and leukocyte differential counter for In Vitro Diagnostic Use in clinical laboratories.	Same							

	Similarities	
Item	Device	Predicate
Sample Types	Whole Blood Mode and Prediluted Mode	Same
Sample Processing	Utilizes an automatic sampling, diluting, and mixing device for sample processing.	Same
Calibration	Provides 2 calibration programs: manual calibration and auto calibration using commercial calibrators.	Same
Aperture Alert	Minimize the possibility of reporting erroneous results caused by a partial or transient aperture clog or by other aperture disturbance.	Same

	Differences	
Item	Device	Predicate
Operating Modes	Closed Vial Whole Blood Mode	Closed Vial Whole Blood Mode and Open Vial Whole Blood Mode
Throughput	1 minute/analysis	60 seconds or les
Quality Control	Provides 2 QC programs: L-J Analysis and X-B Analysis	Provides 1 QC program: L-J Analysis
Recommended Controls	BC-3D: Low, Normal, & High	4C PLUS: Abnormal Low, Normal, Abnormal High
Sample Volume Aspirated	13μL whole blood 20μL prediluted blood	18μL whole blood 20μl prediluted blood
Parameters	WBC, RBC, HGB, HCT, MCV, MCH, MCHC, PLT, Lymph%, Lymph#, Mid%, Mid#, Gran%, Gran#, RDW, MPV	WBC, RBC, HGB, HCT, MCV, MCH, MCHC, PLT, Ly%, Ly#, MO%, MO#, GR%, GR#, RDW, MPV

I. Special Control/Guidance Document Referenced (if applicable):

Class II Special Controls Guidance Document: Premarket Notifications for Automated Differential Cell Counters for Immature or Abnormal Blood Cells; Final Guidance for Industry and FDA. December 4, 2001

J. Performance Characteristics:

- 1. Analytical Performance:
 - a. Accuracy:

Correlation is determined by comparing the results (both CBC and DIFF) obtained by the BC-3200 to those by the Coulter AC \cdot T diff 2TM and by comparing the DIFF results obtained by the BC-3200 to those by manual differential.

Table 1: Correlation to Coulter A^{C} ·T diff 2^{TM}

		Mear	ſ	Difference			Correlation
Parameter	Samples (n)	BC-3200	A°.T. diff 2	ratio (D%)	Slope (a)	Intercept (b)	Correlation coefficient (r)
WBC	103	10.4	10.3	2.4	1.0097	-0.0282	0.9994
Lymph#	98	1.9	2.1	11.8	0.9918	-0.1864	0.9890
Mid#	98	0.7	0.5	40.5	2.1022	-0.3798	0.9187
Gran#	98	6.1	6.0	3.7	0.9886	0.146	0.9978
Lymph%	98	25.8	29.3	11.5	0.7935	2.5772	0.9751
Mid%	98	9.0	6.7	43.0	0.7569	3.8798	0.4644
Gran%	98	65.2	64.0	3.4	0.9046	7.347	0.9707
RBC	103	4.31	4.27	1.7	0.9916	0.0702	0.9971
HGB	103	12.6	12.5	1.2	0.9951	0.0853	0.9982
НСТ	103	37.6	37.2	2.2	1.0041	0.2953	0.9950
MCV	103	87.8	87.5	1.2	0.9549	4.3174	0.9824
MCH	103	29.2	29.5	1.6	0.9426	1.4345	0.9791
MCHC	103	33.3	33.7	1.8	0.7759	7.1720	0.6784
RDW	103	13.1	13.5	4.7	0.4393	7.1667	0.9569
PLT	103	226	230	8.0	0.8882	21.837	0.9961
MPV	102	8.5	8.9	4.7	0.7037	2.2287	0.9334

Table 2: Correlation to Manual Differential

		Ν	lean			Correlation
Parameter	Samples	BC- 3200	Manual differential	Slope	Intercept	Coefficient r
Lymph%	196	26.8	30.4	0.7575	3.7958	0.95
Mid%	196	9.2	9.0	0.3739	5.822	0.57
Gran%	196	64.0	60.6	0.8456	12.721	0.94

b. Precision/Reproducibility:

Reproducibility is stated in terms of both Standard Deviation (SD) and Coefficient of Variation (CV%). Reproducibility was determined by replicate testing (n = 11) with samples of low, normal and high concentrations, three samples for each concentration. For each sample, results of the 2nd to 11th runs were adopted to calculate the SD and CV%.

1	WBC	RBC	HGB	MCV	PLT
•	$ imes 10^3$ / μL	$ imes 10^{6}$ / μ L	(g/dL)	(fl)	$\times 10^3 / \mu L$
mean	4.1	2.88	9.2	64.6	162
SD	0.07	0.04	0.1	0.40	5.06
CV(%)	1.63	1.45	0.8	0.62	3.12
2	WBC	RBC	HGB	MCV	PLT
2	$ imes 10^3$ / μL	$ imes 10^{6}$ / μ L	(g/dL)	(fl)	$\times 10^3 / \mu L$
mean	3.2	3.02	9.3	72.9	155
SD	0.03	0.03	0.1	0.21	7.02
CV(%)	0.99	1.06	1.0	0.28	4.53
3	WBC	RBC	HGB	MCV	PLT
5	$ imes 10^3$ / μL	$\times 10^{6} / \mu L$	(g/dL)	(fl)	$\times 10^3 / \mu L$
mean	3.1	1.91	5.6	61.0	61
SD	0.06	0.03	0.1	0.24	5.11
CV(%)	1.84	1.76	1.1	0.39	8.39

Table 3: Imprecision, Low Concentration Samples

Γ /μL 4
4
5
0
Г
/µL
9
6
5
Г
/µL
1
3
0

Table 4: Imprecision, Normal Concentration Samples

Table 5: Imprecision, Normal Concentration Samples

	WBC	RBC	HGB MCV	PLT
1	$ imes 10^3$ / μL	$ imes 10^{6}$ /µL	(g/dL) (fl)	$\times 10^3 / \mu L$
mean	10.1	4.60	13.1 83.3	244
SD	0.12	0.03	0.09 0.38	8.05
CV(%)	1.18	0.73	0.7 0.45	3.30
2	WBC	RBC	HGB MCV	PLT
2	$\times 10^3$ / μL	$ imes 10^{6}$ / μL	(g/dL) (fl)	$ imes 10^3 / \mu L$
mean	9.8	5.34	15.2 83.1	249
SD	0.10	0.04	0.12 0.27	4.86
CV(%)	0.99	0.78	0.8 0.33	1.95
3	WBC	RBC	HGB MCV	PLT
5	$ imes 10^3$ / μL	$\times 10^{6}/\mu L$	(g/dL) (fl)	$ imes 10^3$ / μL
mean	11.3	5.27	15.0 85.9	231
SD	0.13	0.04	0.06 0.21	8.53
CV(%)	1.11	0.73	0.4 0.25	3.70

Inter-laboratory Precision:

Two laboratories, each having one BC-3200 installed, were selected for the test. Three samples of various concentrations (respectively low, normal and high) were prepared, each with sufficient volume to run twice on both of the BC-3200s. Each BC-3200 was operated by one operator, who conducted the test from beginning to the end. Each sample was divided into two aliquots, and the two aliquots were analyzed respectively by the two selected laboratories within the same day of preparation. Each aliquot was run twice on the BC-3200 and both runs were conducted within a short interval. No outlier was found during the test. Based on the data acquired, repeatability variance (S^{r2}), between laboratory variance (S^{j2}), and reproducibility variance (S $_{R}^{2}$) of the following parameters, WBC, RBC, PLT, HGB, Lymph%, Mid% and Gran%,, were calculated for ea concentration.

Table 6: Within-run Precision and Total Precision							
		Low	Normal	High			
	WBC>	«10 ³ / μL	·				
Mean		2.13	8.10	20.68			
Repeatability variance	S_r^2	0.0025	0.0098	0.0613			
Between Laboratory variance	S_L^2	0.0000	0.0151	0.0000			
	S_R^2	0.0025	0.0249	0.0613			
Reproducibility variance	S_R	0.0500	0.1578	0.2476			
	CV%	2.35%	1.95%	1.20%			
	Gra	an(%)					
Mean		32.53	60.98	81.30			
Repeatability variance	S_r^2	1.1050	0.0221	0.0637			
Between Laboratory variance	S_L^2	1.7588	0.7703	0.0932			
	S_R^2	2.8638	0.7924	0.1569			
Reproducibility variance	S_R	1.6923	0.8902	0.3961			
	CV%	5.20%	1.46%	0.49%			
	Lym	ph (%)					
Mean		12.65	28.83	51.30			
Repeatability variance	S_r^2	0.2073	0.0613	3.0439			
Between Laboratory variance	S_L^2	0.0000	1.3307	6.4781			

Table 6: Within-run Precision and Total Precision

Mid (%) Mid (%) Mean 6.05 10.20 16 Repeatability variance S_r^2 0.0490 0.0098 0.6	220 858 2% .18 655
CV% 3.60% 4.09% 6.0 Mid (%) Mid (%) 6.05 10.20 16 Repeatability variance S_r^2 0.0490 0.0098 0.6 Between Laboratory S_r^2 0.0205 0.0751 1.3	.18
CV% 3.60% 4.09% 6.0 Mid (%) Mid (%) 6.05 10.20 16 Repeatability variance S_r^2 0.0490 0.0098 0.6 Between Laboratory S_r^2 0.0205 0.0751 1.3	.18
Mean 6.05 10.20 16 Repeatability variance S_r^2 0.0490 0.0098 0.6 Between Laboratory S^2 0.0205 0.0751 1.3	
Mean 6.05 10.20 16 Repeatability variance S_r^2 0.0490 0.0098 0.6 Between Laboratory S^2 0.0205 0.0751 1.3	
Repeatability variance S_r^2 0.04900.00980.6Between Laboratory S_r^2 0.02050.07511.3	655
Between Laboratory S ² 0.0205 0.0751 1.3	
	786
S_R^2 0.0695 0.0849 2.0	441
	297
CV% 4.36% 2.86% 8.8	4%
RBC (×10 ⁶ /µL)	
	80
Repeatability variance S_r^2 0.00040.00650.0	085
Between Laboratory	000
S_R^2 0.0011 0.0078 0.0	085
Reproducibility variance S_R 0.03320.08830.0	922
CV% 1.34% 1.81% 1.5	9%
HGB (g/L)	
	.13
Repeatability variance S_r^2 0.00000.00250.0	123
Between Laboratory variance S_L^2 0.00500.06010.0	952
K	075
Reproducibility variance S_R 0.07070.25020.3	279
	1%
MCV (fl)	
Mean 77.28 86.73 96	.33
Repeatability variance S_r^2 0.11030.01230.0	907
Between Laboratory variance S_L^2 2.25621.52522.7	160
S_R^2 2.3665 1.5375 2.8	067
	753
CV% 1.99% 1.43% 1.7	4%

PLT (×10 ³ /μL)								
Mean	94.75	258.25	468.50					
Repeatability variance	S_r^2	13.2453	16.2699	12.5033				
Between Laboratory variance	S_L^2	14.5024	69.9901	65.7484				
	S_R^2	27.7477	86.2600	78.2517				
Reproducibility variance	S_R	5.2676	9.2876	8.8460				
	CV%	5.56%	3.60%	1.89%				

Table 6 Appendix

WBC	Form A				Form B			Form C	
Laboratory	Low	Normal	High	Low	Normal	High	Low	Normal	High
1	2.2	8.1	20.5	2.15	8.2	20.75	0.07	0.14	0.35
	2.1	8.3	21	2.13	0.2	20.75	0.07	0.14	0.55
2	2.1	8	20.6	2.1	Q	20.6	0	0	0
	2.1	8	20.6	2.1	0	20.0	0	0	0

Gran(%)		Form A			Form B			Form C		
Laboratory	Low	Normal	High	n Low Normal		High	Low	Normal	High	
1	32.5	60.2	81.1	31.45	60.35	81.05	1.48	0.21	0.07	
	30.4	60.5	81	51.45	00.55	81.05	1.40	0.21	0.07	
2	33.5	61.6	81.8	33.6	61.6	81.55	0.14	0	0.35	
	33.7	61.6	81.3	33.0	01.0	01.55	0.14	0	0.55	

Lymph (%)	Form A				Form B			Form C		
Laboratory	Low	Normal	High	Low	Normal	High	Low	Normal	High	
1	12.8	29.9	51.7	12.75	29.65	53.3	0.07	0.35	2.26	
	12.7	29.4	54.9	12.75	29.05	55.5	0.07	0.55	2.20	
2	12.1	28	50	12.55	28	49.3	0.64	0	0.99	
	13	28	48.6	12.33	20	47.5	0.04	0	0.99	

Mid (%)	Form A				Form B			Form C		
Laboratory	Low	Normal	High	Low	Normal	High	Low	Normal	High	
1	6.1	9.9	15.8	6.2	10	15.25	0.14	0.14	0.78	
	6.3	10.1	14.7	0.2	10	13.23	0.14	0.14	0.78	
2	6.1	10.4	16.5	5.9	10.4	17.1	0.28	0	0.85	
	5.7	10.4	17.7	5.9	10.4	17.1	0.28	0	0.85	

RBC		Form A			Form B			Form C		
Laboratory	Low	Normal	High Low		Normal	High	Low	Normal	High	
1	2.44	4.78	5.84	2.455	4.845	5.765	0.02	0.09	0.11	
	2.47	4.91	5.69	2.433	4.045	5.705	0.02	0.09	0.11	
2	2.51	4.99	5.89	2.495	4.94	5.84	0.02	0.07	0.07	
	2.48	4.89	5.79	2.493	4.94	5.84	0.02	0.07	0.07	

HGB	Form A			Form B			Form C			
Laboratory	Low	Normal	High	Low	Normal	High	Low	Normal	High	
1	6.3	13.9	18.8	6.3	13.9	18.9	0	0	0.14	
	6.3	13.9	19	0.5	13.9	10.9	0	0	0.14	
2	6.4	14.3	19.4	6.4	14.25	19.35	0	0.07	0.07	
	6.4	14.2	19.3	0.4	14.23	19.55	U	0.07	0.07	

MCV	Form A	Form B	Form C						
Laboratory	Low	Normal	High	Low	Normal	High	Low	Normal	High
1	76.5	85.9	95.2	76.2	85.85	95.15	0.42	0.07	0.07
	75.9	85.8	95.1	70.2	85.85	95.15	0.42	0.07	0.07
2	78.5	87.5	97.8	78.35	87.6	97.5	0.21	0.14	0.42
	78.2	87.7	97.2	/0.55	07.0	97.5	0.21	0.14	0.42

PLT	Form A			Form B			Form C			
Laboratory	Low	Normal	High	Low	Normal	High	Low	Normal	High	
1	88	265	466	91.5	264.5	462.5	4.95	0.71	4.95	
	95	264	459	91.5	204.3	402.3	4.95	0.71	4.95	
2	97	248	474	98	252	474.5	1.41	5.66	0.71	
	99	256	475	70	232	4/4.3	1.41	5.00	0.71	

c. Linearity:

Linearity was determined by running diluted samples. RBC,HGB are diluted by blood plasma of the sample, while WBC and PLT are diluted by specified diluent. Concentrations from 0 to 100% were tested, each concentration twice. The average of the two runs is taken as the result, together with the concentration, to calculate per the linear regression equation.

Table 7: WDC Elifearity										
Dilution%	Test 1	Test 2	Mean	Ideal	Error	Proportional Error				
100	117.1	115.9	116.50	120.01	3.51	2.9				
80	99.8	100.1	99.95	96.01	-3.94	-4.1				
60	73.4	72.1	72.75	72.00	-0.75	-1.0				
40	47.8	48.6	48.20	48.00	-0.20	-0.4				
20	23.1	23.1	23.10	23.99	0.89	3.7				
10	12.1	12.0	12.05	11.99	-0.06	-0.5				
5	6.0	6.2	6.10	6.00	-0.10	-1.7				
2.5	3.0	2.9	2.95	2.99	0.04	1.3				
1.25	1.3	1.3	1.30	1.49	0.19	12.8				
0.625	0.5	0.5	0.50	0.74	0.24	32.4				
0.3125	0.2	0.1	0.15	0.36	0.21	58.3				
0	0	0	0.00	-0.01	-0.01	/				
Slope	1.2002									
Intercept	-0.0129									

Table 7: WBC Linearity

Table 8: RBC Linearity

Dilution%	Test 1	Test 2	Mean	Ideal	Error	Proportional Error				
100	8.46	8.43	8.445	8.519	0.074	0.9				
80	6.91	6.86	6.885	6.819	-0.066	-1.0				
60	5.12	5.17	5.145	5.119	-0.026	-0.5				
40	3.42	3.46	3.440	3.419	-0.021	-0.6				
20	1.71	1.69	1.700	1.719	0.019	1.1				
10	0.89	0.87	0.880	0.869	-0.011	-1.3				
5	0.46	0.46	0.460	0.444	-0.016	-3.6				
2.5	0.21	0.22	0.215	0.232	0.017	7.3				
1.25	0.10	0.13	0.115	0.125	0.010	8.0				
0	0.00	0.00	0.000	0.019	0.019	/				
Slope		0.0850								
Intercept	0.0191									

Tuble 7. HOD Emeanly										
Dilution%	Test 1	Test 2	Mean	Ideal	Error	Proportional Error				
100	25.6	25.6	25.60	25.40	-0.20	-0.8				
80	20.5	20.1	20.30	20.33	0.03	0.1				
60	15.1	14.9	15.00	15.26	0.26	1.7				
40	10.1	10.1	10.10	10.19	0.09	0.9				
20	5.2	5.0	5.10	5.11	0.01	0.2				
10	2.7	2.6	2.65	2.58	-0.07	-2.7				
5	1.4	1.4	1.40	1.31	-0.09	-6.9				
2.5	0.7	0.7	0.70	0.68	-0.02	-2.9				
1.25	0.4	0.4	0.40	0.36	-0.04	-11.1				
0	0.0	0.0	0.00	0.04	0.04	/				
Slope	0.2536									
Intercept	0.0425									

Table 9: HGB Linearity

Table 10: PLT Linearity

Dilution%	Test 1	Test 2	Mean	Ideal	Error	Proportional Error				
100	1014	1008	1011.0	1040.3	29.3	2.8				
80	850	858	854.0	832.5	-21.5	-2.6				
60	631	650	640.5	624.8	-15.7	-2.5				
40	425	419	422.0	417.0	-5.0	-1.2				
20	221	208	214.5	209.3	-5.2	-2.5				
10	109	101	105.0	105.4	0.4	0.4				
5	53	53	53.0	53.5	0.5	0.9				
2.5	23	17	20.0	27.5	7.5	27.3				
1.25	8	5	6.5	14.5	8.0	55.2				
0	0	0	0.0	1.6	1.6	/				
Slope		10.3871								
Intercept	1.5618									

d. Carryover:

Carryover was determined by first running the high concentration sample for three consecutive times (i1, i2, i3) and then the low concentration sample three consecutive times (j1, j2, j3), and finally calculating per the following equation: Carryover (%) = $[(j1 - j3)/(i3-j3)] \times 100\%$. The test was then repeated using the high level control.

	High	Concent	ration	Low (Concent	ration	
Parameter	San	ple (Wl	nole	San	nple (Wl	Carryover	
Falametei	Blood)				Blood)	%	
	i1	i2	i3	j1	j2	j3	
WBC(× 10^3 / μ L)	19.7	20.4	20.0	1.9	1.9	1.9	0%
$RBC(\times 10^6 / \mu L)$	6.34	6.24	6.2	1.87	1.96	1.85	0.46%
HGB(g/dL)	25.4	25.0	24.8	3.3	3.2	3.2	0.46%
PLT(×10 ³ / μ L)	404	390	396	31	34	33	0%

Table 11: Carryover, High Concentration Sample

Table 12: Carryover, High Level Control

	High Concentration			Low Concentration				
Parameter	Sample (High Level			Sample (Specified			Carryover	
Parameter	Control)			Diluent)			%	
	i1	i2	i3	j1	j2	j3		
WBC(×10 ³ / μ L)	21.7	21.3	21.7	0.0	0.0	0.0	0%	
$RBC(\times 10^6 / \mu L)$	5.88	5.79	5.79	0.00	0.00	0.00	0%	
HGB(g/dL)	18.8	18.7	18.9	0.0	0.0	0.0	0%	
$PLT(\times 10^3 / \mu L)$	453	438	429	0	0	0	0%	

e. Interfering Substances:

N/A

2. <u>Other Supportive Instrument Performance Data Not Covered Above:</u>

Ability to Flag Abnormal Histograms:

BC-3200's ability to flag abnormal WBC histograms was determined by comparing 200 sample results obtained by the BC-3200 to those obtained by manual differential.

Table 13: Ability to Flag Abnormal WBC Histograms

ι.	.5. Admity to Flag Admonthan WDC Instograms						
	Manual	BC-3200					
	differential	Positive (39)	Negative (161)				
	Positive (40) Negative (160)	TP 22 FP 17	FN 18 TN 143				
		False Positive Ratio	False Negative				
	Agreement%	(%)	Ratio (%)				
	82.5	10.6	45				

Reference Ranges

A Normal Ranges Study was conducted to assess the Reference Ranges for the BC-3200 analyzer. Whole-blood samples were collected from 121 donors.

				90%Confidence	90%Confidence
Parameter	Units	Sex	Mean	Low Limit	High Limit
WBC	$\times 10^3$ cells /µL	M/F	6.86	3.47	10.25
RBC	$\times 10^6$ cells /µL	M/F	4.56	3.54	5.58
Hgb	g/ dL	M/F	13.40	10.27	16.52
Hct	ratio	M/F	40.12	30.98	49.26
MCV	fL	M/F	88.18	80.82	95.55
MCH	pg	M/F	29.36	26.57	32.15
MCHC	g/ dL	M/F	33.33	32.09	34.56
Plt	×10 ³ cells/µL	M/F	209.92	119.62	300.22
RDW	%	M/F	12.81	11.53	14.10
MPV	fL	M/F	8.47	7.07	9.87
LY	%	M/F	27.33	18.11	36.55
MO	%	M/F	9.45	5.23	13.67
GR	%	M/F	63.26	51.62	74.89

 Table 14: Reference Ranges

K. Proposed Labeling:

The labeling is sufficient and it satisfies the requirements of 21 CFR Part 809.10.

L. Conclusion:

The submitted information in this premarket notification is complete and supports a substantial equivalence decision.