# 510(k) SUBSTANTIAL EQUIVALENCE DETERMINATION DECISION SUMMARY ASSAY ONLY TEMPLATE

# **A.** 510(k) Number:

K070005

## **B.** Purpose for Submission:

To seek clearance for HemosIL RecombiPlasTin 2G

HemosIL RecombiPlasTin 2G is modified from the current PT reagent HemosIL RecombiPlasTin (K043184). The ratio of raw material phospholipids used in the current PT reagent is modified to address the assay's labeled sensitivity to the recently released antibiotic CUBICIN (Daptomycin for injection).

## C. Measurand:

- PT
- Fibrinogen

# **D.** Type of Test:

Absorbance or light-scatter clot detection

# E. Applicant:

Instrumentation Laboratory Co.

# F. Proprietary and Established Names:

HemosIL RecombiPlasTin 2G

# **G.** Regulatory Information:

- 1. Regulation section:
  - CFR 864.7750
  - CFR 864.7340

## 2. Classification:

Π

#### 3. Product code:

- GJS
- GIS

#### 4. Panel:

Hematology

#### H. Intended Use:

## 1. Intended use(s):

HemosIL RecombiPlasTin 2G is a high sensitivity thromboplastin reagent based on recombinant human tissue factor (RTF) for the quantitative determination in human citrated plasma of Prothrombin Time (PT) on IL Coagulation and ELECTRA Systems and Fibrinogen on IL Coagulation Systems only. The product is used for the evaluation of the extrinsic coagulation pathway and the monitoring of Oral Anticoagulant Therapy (OAT).

#### 2. Indication(s) for use:

HemosIL RecombiPlasTin 2G is a high sensitivity thromboplastin reagent based on recombinant human tissue factor (RTF) for the quantitative determination in human citrated plasma of Prothrombin Time (PT) on IL Coagulation and ELECTRA Systems and Fibrinogen on IL Coagulation Systems only. The product is used for the evaluation of the extrinsic coagulation pathway and the monitoring of Oral Anticoagulant Therapy (OAT).

# 3. Special conditions for use statement(s):

## 4. Special instrument requirements:

- PT on IL Coagulation and ELECTRA Systems
- Fibrinogen on IL Coagulation Systems only

## I. Device Description:

The RecombiPlasTin 2G reagent is formulated to be insensitive to therapeutic levels of heparin. In the PT test, the addition of the tissue thromboplastin (RecombiPlasTin 2G reagent) to the patient plasma in the presence of calcium ions initiates the activation of the extrinsic pathway. This results ultimately in the conversion of fibrinogen to fibrin, with formation of a solid gel. For the IL Coagulation Systems only, the Fibrinogen is quantitated (PT-based method) by relating the absorbance or light-scatter during clotting to a calibrator.

# J. Substantial Equivalence Information:

1. <u>Predicate device name(s):</u> HemosIL RecombiPlasTin

2. Predicate 510(k) number(s): K043184

3. Comparison with predicate:

	Similarities					
<u>Item</u>	<u>Device:</u> HemosIL RecombiPlasTin 2G	Predicate: HemosIL RecombiPlasTin K043184				
Indications for Use/Intended Use	Same	For the quantitative determination in human citrated plasma of PT on IL Coagulation and ELECTRA Systems and Fibrinogen on IL System only.				
Test Principle	Same	Reagent is formulated to be insensitive to therapeutic levels of heparin.  - PT: the addition of the tissue thromboplastin to the patient plasma in the presence of calcium initiates the activation of the extrinsic pathway. This results ultimately in the conversion of fibrinogen to fibrin, with the formation of a solid gel.  - Fibrinogen: it is quantitated (PT based method) by relating the absorbance or light-scatter during clotting to a calibrator.				
Reagent: RecombiPlasTin Diluent	Same	Aqueous solution of calcium chloride, polybrene and a preservative (<0.1% sodium azide).				

	Differences				
Item  Reagent: RecombiPlasTin	Device: HemosIL RecombiPlasTin 2G Same except for change in phospholipid raw material ratios to address interference on PT results from CUBICIN®.	Predicate: HemosIL RecombiPlasTin K043184  Lyophilized recombinant human tissue factor, synthetic phospholipids with stabilizers, preservative (<0.1% sodium azide) and buffer.			
Labeling: Limitations/ interfering substances	PT results maybe affected by many commonly administered drugs and further studies should be made to determine the source of unexpected abnormal results.	PT results maybe affected by many commonly administered drugs and further studies should be made to determine the source of unexpected abnormal results. Due to a significant dose response, evaluate the PT/INR status of patients being treated with the antibiotic CUBICIN (Daptomycin for Injection) using an alternate HemosIL PT reagent			

# K. Standard/Guidance Document referenced (if applicable):

# L. Test Principle:

Reagent is formulated to be insensitive to therapeutic levels of heparin. In the PT test, the addition of the tissue thromboplastin to the patient plasma in the presence of calcium initiates the activation of the extrinsic pathway. This results ultimately in the conversion of fibrinogen to fibrin, with the formation of a solid gel.

For the IL Coagulation Systems only, the Fibrinogen is quantitated (PT based method) by relating the absorbance or light-scatter during clotting to a calibrator.

## M. Performance Characteristics (if/when applicable):

The following are instruments, reagents, calibrator, and controls used in the performance testing:

## - Representative instruments:

• ACL Family of Coagulation Analyzers:

o ACL 10000 K060162

• ACL Futura/ACL Advance Coagulation Analyzers:

o ACL Advance K002400

ACL TOP Coagulation Analyzer	K033414
• ELECTRA Series of Coagulation Analyzers:	
o ELECTRA 1600C	K931206
Additional instruments used for instrument equivalence	e testing:
o ACL 100	K881367
o ACL 1000	K912087
o ACL 6000	K961991
- Reagents:	
<ul> <li>HemosIL RecombiPlasTin (Predicate)</li> </ul>	K043184
HemosIL Factor II Deficient Plasma	K050661
HemosIL Factor V Deficient Plasma	K023839
HemosIL Factor VII Deficient Plasma	K024082
HemosIL Factor X Deficient Plasma	K031122
- Calibrator and controls:	
HemosIL Calibration Plasma	K041905
HemosIL Normal Control	K021023
HemosIL Low Abnormal Control	K021022
HemosIL High Abnormal Control	K021024
HemosIL Low Fibrinogen Control	K033414

# 1. Analytical performance:

# a. Precision/Reproducibility:

HemosIL Special Test Control Level 2

PT and Fibrinogen precision studies were performed over 10 days, 2 runs per day, 4 replicates per run, using several different representative IL Coagulation System analyzers.

K040359

Three levels of Quality Control were used for the PT study, and two levels of Quality Control were used for the Fibrinogen study. All within-run, betweenrun, and total run results ranged from 0.5-3.4% CV for PT study and from 0.0-11.8% CV for Fibrinogen study.

PT Seconds						
Analyzer	Control Level	N	Mean	CV% Within Run	CV% Between Run	CV% Total
	Normal	80	11.9	0.8	1.1	2.2
ACL TOP	Low Abnormal	80	22.0	0.8	2.0	3.1
	High Abnormal	80	34.0	0.9	2.9	3.1
	Normal	80	12.5	1.1	0.7	1.9
ACL Advance	Low Abnormal	80	23.8	1.6	0.5	1.9
	High Abnormal	80	36.0	1.8	1.0	2.4
	Normal	80	11.7	0.6	0.5	1.5
ACL 10000	Low Abnormal	80	21.6	1.0	1.0	1.9
	High Abnormal	80	32.9	1.1	1.4	2.6
	Normal	80	10.4	0.7	0.8	2.1
ACL 6000	Low Abnormal	80	20.2	0.9	0.9	2.1
	High Abnormal	80	31.1	1.3	1.7	2.1
	Normal	80	10.8	1.1	0.6	1.5
ACL 3000	Low Abnormal	80	20.5	1.0	1.1	1.5
	High Abnormal	80	31.6	1.3	1.6	2.0
	Normal	80	10.8	1.1	0.8	1.4
ACL 300	Low Abnormal	80	20.9	1.4	0.9	1.8
	High Abnormal	80	32.1	1.5	1.5	2.2
	Normal	80	12.2	1.3	0.9	1.9
E1600C	Low Abnormal	80	21.1	1.2	1.7	2.6
(PT Only)	High Abnormal	80	31.3	1.3	1.8	3.4

Fibrinogen mg/dL						
Analyzer	Control Level	N	Mean	CV% Within Run	CV% Between Run	CV% Total
	Normal	80	296	1.4	1.3	2.4
ACL TOP	Low Fibrinogen	80	135	2.9	1.8	3.6
	Normal	80	229	3.0	0.0	3.1
ACL Advance	Low Fibrinogen	80	157	3.7	0.0	4.5
	Normal	80	319	4.2	2.3	5.0
ACL 10000	Low Fibrinogen	80	149	5.9	0.0	6.9
	Normal	80	292	4.0	2.2	4.9
ACL 6000	Low Fibrinogen	80	134	6.4	3.2	9.1
	Normal	80	297	3.9	3.8	6.2
ACL 3000	Low Fibrinogen	80	128	6.5	9.8	8.3
	Normal	80	327	4.6	2.5	7.4
ACL 300	Low Fibrinogen	80	144	5.7	6.4	11.8

Precision study results for PT and Fibrinogen are within acceptance criteria.

# b. Linearity/assay reportable range:

## • Linearity study

Fibrinogen linearity studies were performed on representative IL Coagulation Systems. Results met the acceptance criteria of r > 0.95 supporting their claims in the device's product insert of:

ACL Family and ACL Futura/Advance: 80-700 mg/dL ACL TOP: 80-700 mg/dL

 $\circ$  Extrinsic factor linearity studies were also performed on the factors II, V, VII, and X. Results met the acceptance criteria of r>0.95.

## Factor recovery

 $\circ$  Extrinsic factor recovery studies were performed on the factors II, V, VII, and X. Results met the acceptance criteria of r >0.95.

c. Traceability, Stability, Expected values (controls, calibrators, or methods):

Stability studies (onboard, reconstituted, and shelf-life) also met the acceptance criteria for validating their claims.

- d. Detection limit:
- e. Analytical specificity:
  - CUBICIN (Daptomycin for injection)

The sponsor provided the justification for the use of plasma spiking studies to test for CUBICIN (Daptomycin for injection) interference; the expected peak dosing for CUBICIN in circulation blood is  $75\mu g/mL$ . The interference study was performed on representative members of the IL Coagulation system families using the following specifications:

- o Normal sample  $\pm 1$  second from unspiked sample
- $\circ$  Coumadin sample  $\pm$  10% recovery of the unspiked sample

The study met the specifications and indicated no clinically significant CUBICIN dose response with the modified HemosIL RecombiPlasTin 2G.

• Heparin, hemoglobin, triglycerides, and bilirubin

Interference study results met the acceptance criteria for validating their claims.

f. Assav cut-off:

## 2. <u>Comparison studies:</u>

- a. Method comparison with predicate device:
  - In-house Study:

A method comparison study was performed to compare the performance of the modified HemosIL RecombiPlasTin 2G vs. the predicate HemosIL RecombiPlasTin on representative members of the different IL Coagulation System Families.

198 citrated frozen plasma samples were used in the study, 130 normal donors and 68 abnormal patient samples. Abnormal patient samples include: OAT, Lupus, liver disease, low fibrinogen, high fibrinogen, unfractionated heparin, low molecular weight heparin, DIC, factor II deficient, factor V deficient, factor VII deficient, factor X deficient.

One sample result was removed from calculations due to instrument's error. A single sample outlier was also removed from calculations due to severe factor VII deficient concentration. The removal of these two results brings their total to N=196.

System	Assay	Slope	Intercept	r
ELECTRA		0.7637	3.0427	0.9789
ACL Family	PT	0.7910	2.7860	0.9885
ACL Futura/	(Seconds)			
ACL Advance		0.8075	2.8899	0.9913
ACL TOP		0.8010	2.7138	0.9916
ACL Family		0.9350	6.1043	0.9866
ACL Futura/	Fibrinogen			
ACL Advance	(mg/dL)	0.9711	10.933	0.9783
ACL TOP		1.0129	-3.6298	0.9969

Additional studies on other instruments, ACL 300 and ACL 3000 using the ACL 6000 as the reference device, were performed to demonstrate instrument equivalency.

The product insert provides a statement indicating that customers can expect to see a discrepancy in slopes for PT seconds when comparing reagents with different ISI assignments.

## • Field Site Study:

The following IL coagulation instruments were grouped as instrument equivalency:

- Linear Analyzer Group: ACL TOP, ACL Advance/Futura, and ELECTRA Series
- o Centrifugal Analyzer Group: ACL 8/9/10000/Elite/ElitePro, ACL 100-7000 Series

#### - First field site study:

PT seconds, INR and Fibrinogen method comparison study was performed on an ACL TOP representing the linear analyzer group. 207 patient samples were used in this study. The patient samples included 119 normal samples and 88 abnormal samples, including 61 patients on OAT and 27 patients with various disease states.

For Fibrinogen, thirteen samples were not included in the calculations due to no available results from either the reference or test instrument, N=194.

System	Assay	Slope	Intercept	r
	PT (Sec.)	0.8137	4.0035	0.9934
ACL TOP	PT INR	1.0838	-0.1071	0.9945
	Fibrinogen (mg/dL)	0.9805	5.3466	0.9946

The sponsor also provided the INR study using only samples from patients on OAT. The resultant statistics indicate no significant difference in the correlation statistics with the results obtained from OAT and normal samples. The slope for the subset of OAT samples alone was 1.0855.

## - Second field site study:

PT seconds, INR and Fibrinogen method comparison study was performed on an ACL 10000 representing the centrifugal analyzer group. 88 patient samples were used in this study. The patient samples included 20 normal samples and 68 patients on OAT.

Assay	Slope	Intercept	r
PT (Sec.)	0.7935	2.1733	0.9887
PT INR	0.9446	0.0367	0.9881
Fibrinogen	0.9431	7.5763	0.9832
	PT (Sec.) PT INR	PT (Sec.) 0.7935 PT INR 0.9446 Fibrinogen	PT (Sec.) 0.7935 2.1733 PT INR 0.9446 0.0367 Fibrinogen

#### b. Matrix comparison:

A fresh versus frozen sample was performed and the data showed no significant difference in performance for PT seconds, INR or fibrinogen. The patient samples included 20 normal samples and 68 patients on OAT.

	Fresh vs. Frozen			
Parameter	Slope	Intercept	$\mathbb{R}^2$	
INR	1.0120	0.0077	0.9844	
PT seconds	1.0119	0.1063	0.9843	
Fibrinogen	0.9527	29.35	0.9819	

## 3. Clinical studies:

- a. Clinical Sensitivity:
- b. Clinical specificity:

c. Other clinical supportive data (when a. and b. are not applicable):

# 4. Clinical cut-off:

# 5. Expected values/Reference range:

Normal range studies were performed for PT and Fibrinogen with 130 donors on representative analyzers. For Fibrinogen, a single normal sample was eliminated from calculations as an outlier, bringing their total to N=129.

PT	N	Range
ELECTRA (PT only)	130	9.8 - 12.2 (seconds)
ACL Family	130	9.1 - 12.1 (seconds)
ACL Futura/ACL Advance	130	9.9 - 12.9 (seconds)
ACL TOP	130	9.4 - 12.5 (seconds)
Fibrinogen	N	Range
ACL Family	129	308 - 613 (mg/dL)
ACL Futura/ACL Advance	129	222 - 340 (mg/dL)
ACL TOP	129	276 - 471 (mg/dL)

The product insert advises customers that "These results were obtained using a specific lot of reagent. Due to many variables which may affect clotting times, each laboratory should establish its own normal range."

# N. Proposed Labeling:

The labeling is sufficient and it satisfies the requirements of 21 CFR Part 809.10.

#### O. Conclusion:

1. The submitted information in this premarket notification is complete and supports a substantial equivalence decision.