

**510(k) SUBSTANTIAL EQUIVALENCE DETERMINATION
DECISION SUMMARY
ASSAY AND INSTRUMENT COMBINATION TEMPLATE**

A. 510(k) Number:

k102740

B. Purpose for Submission:

To obtain a substantial equivalence determination for a 510k for methicillin-resistant *Staphylococcus aureus* (MRSA) from nasal specimens

C. Measurand:

MRSA DNA

D. Type of Test:

The test utilizes NASBA™ (nucleic acid sequence-based amplification) coupled with molecular beacons (sequence-specific fluorescent probes) to detect the presence of MRSA DNA.

E. Applicant:

bioMérieux, Inc

F. Proprietary and Established Names:

NucliSENS EasyQ® MRSA Assay

G. Regulatory Information:

1. Regulation section:

21 CFR 866. 1640 Antimicrobial susceptibility test powder

2. Classification:

Class II

3. Product code:

NQX- nucleic acid amplification test, DNA, methicillin resistant *Staph aureus*,
direct specimen

OOI- Nucleic acid amplification systems, real time

4. Panel:

83 Microbiology

H. Intended Use:

1. Intended use(s):

The NucliSENS EasyQ[®] MRSA assay is a qualitative *in vitro* diagnostic test for the direct detection of methicillin-resistant *Staphylococcus aureus* (MRSA) from nasal swabs in patients at risk for nasal colonization. The NucliSENS EasyQ[®] MRSA assay is performed on the NucliSENS EasyQ[®] platform.

The test utilizes NASBA[™] (nucleic acid sequence-based amplification) coupled with molecular beacons (sequence-specific fluorescent probes) to detect the presence of MRSA DNA.

The NucliSENS EasyQ[®] MRSA assay is used as a screening tool to aid in the prevention and control of MRSA infections in health care institutions.

NucliSENS EasyQ[®] MRSA is not intended to diagnose, guide or monitor treatment for MRSA infections, or provide results of susceptibility to methicillin. A negative result does not preclude MRSA nasal colonization. Concomitant cultures are necessary to recover organisms for epidemiological typing or for further susceptibility testing.

2. Indication(s) for use:

The NucliSENS EasyQ[®] MRSA assay is a qualitative *in vitro* diagnostic test for the direct detection of methicillin-resistant *Staphylococcus aureus* (MRSA) from nasal swabs in patients at risk for nasal colonization. The NucliSENS EasyQ[®] MRSA assay is performed on the NucliSENS EasyQ[®] platform.

The test utilizes NASBA[™] (nucleic acid sequence-based amplification) coupled with molecular beacons (sequence-specific fluorescent probes) to detect the presence of MRSA DNA.

The NucliSENS EasyQ[®] MRSA assay is used as a screening tool to aid in the prevention and control of MRSA infections in health care institutions.

NucliSENS EasyQ[®] MRSA is not intended to diagnose, guide or monitor treatment for MRSA infections, or provide results of susceptibility to methicillin. A negative result does not preclude MRSA nasal colonization. Concomitant cultures are necessary to recover organisms for epidemiological typing or for further susceptibility testing.

3. Special conditions for use statement(s):

For Prescription Use only

4. Special instrument requirements:

NucliSENS EasyQ[®] MRSA assay protocol software (included in the NucliSENS EasyQ[®] MRSA assay kit)

NucliSENS EasyQ[®] Analyzer

NucliSENS EasyQ[®] Incubator II

NucliSENS EasyQ[®] Director Software 2.6

I. Device Description:

The NucliSENS EasyQ[®] MRSA assay kit contains reagents (nucleic acid primers, fluorescent probes, enzymes, buffers and other reagents), controls, software, accessories, and labeling (Instructions for Use) designed for use with the NucliSENS EasyQ System. The NucliSENS EasyQ System consists of NucliSENS EasyQ System instruments and software (NucliSENS EasyQ[®] Analyzer, NucliSENS EasyQ[®] Director software, and NucliSENS EasyQ[®] Incubator).

For the EasyQ MRSA assay kit, the indicated specimen is a nasal swab, collected using Copan nylon flocked swab. After transport to the laboratory in a dry collection tube, the swab is manually processed for lysis of cells and extraction of nucleic acids using lysis tubes provided in the EasyQ MRSA assay kit, in conjunction with a vortex instrument fitted with a specified vortex plate adaptor to accommodate the lysis tubes. The specimen lysates and controls serve as template material for the subsequent amplification reaction in the EasyQ System.

The NucliSENS EasyQ system is indicated for use solely for nucleic acid sequence-based amplification (NASBA) assays. Following specimen processing and nucleic acid extraction, an operator manually prepares amplification reactions, including a key semi-automated incubation step in the EasyQ Incubator.

Automated amplification and detection are performed simultaneously in a dedicated fluorescent reader, the NucliSENS EasyQ Analyzer. The Analyzer is attached to a computer, monitor, and printer. Amplification is detected by means of sequence-specific probes which each fluoresce upon binding of the probe to its target. Result calculation and reporting is performed automatically via dedicated operator software (NucliSENS EasyQ Director Software, in combination with assay-specific “assay protocols” provided on a CD-ROM in the reagent kit.) Results are based on the analysis of fluorescence signal curves measured by the NucliSENS EasyQ[®] Analyzer.

For the EasyQ MRSA assay, results of the Director Software analysis for each test are provided to users as 1 of 3 possible qualitative results for each test specimen: MRSA positive, MRSA negative, or invalid. Control tests yield qualitative output of valid or invalid.

The EasyQ MRSA kit includes 2 external controls: a blank control and a positive control to perform quality control for the amplification and detection portions of the test. The Instructions for Use also recommend use of specimen positive control (MRSA) and specimen negative control (MSSA) strains, which have been validated for use with the EasyQ MRSA test, to provide quality control for the specimen sample preparation processing steps of the test. Operators are instructed to invalidate all tests in a run if any control in the run yields an invalid result.

The NucliSENS EasyQ® MRSA is developed for amplification of two targets simultaneously: the *SCCmec* cassette junction and the *mecA* gene. The NucliSENS EasyQ® MRSA will report a MRSA positive result only if both targets, the *SCCmec* cassette junction and the *mecA* gene, are detected, thus reducing the risk of false positive detection for “*mecA* dropout”- methicillin-sensitive *Staphylococcus aureus* strains that lost the *mecA* gene but still harbor the cassette junction. Primers and beacons in the kit are complementary to the *SCCmec* cassette junction sequences of MREJ types i-v, vii, and xii.

Components of NucliSENS EasyQ® MRSA

Component	Composition	Number of tubes and tube description	Caps Color code
LYS Tub (lysis tube)	Proclin solution with glass beads*	48 tubes 1150 µl per tube	Uncolored cap tubes
REAGENT spheres sachet (amplification/detection reagents)	Lyophilized white sphere: <ul style="list-style-type: none"> primers, molecular probes inhibition control: non-infectious DNA NASBA™ reagents (nucleotides, DTT, KCl, MgCl₂) carbohydrate Lyophilized orange sphere: <ul style="list-style-type: none"> Restriction enzyme Bovine Serum Albumin** carbohydrate 	6 tubes in foil pack with silica gel desiccant 2 spheres per tube (white and orange)	Blue cap tubes
ENZ I sphere sachet (enzyme reagents)	Lyophilized blue sphere : <ul style="list-style-type: none"> AMV-RT, RNase H, T7 RNA polymerase Bovine Serum Albumin** 	6 tubes in foil pack with silica gel desiccant 1 sphere per tube (blue)	Red cap tubes
REAGENTdil	Reagent spheres diluent (Tris/HCl, KCl, DMSO,PVP)	2 tubes. 600 µl per tube	Blue cap tubes
ENZdil	Enzyme sphere diluent (Sorbitol, Bovine Serum Albumin**, Proclin)	2 tubes. 350 µl per tube	Red cap tubes
NW	NASBA™ water with proclin for dissolution of the CONTROL+ sphere and for the NASBA™ blank	2 tubes. 1500 µl per tube	White cap tubes

Component	Composition	Number of tubes and tube description	Caps Color code
CONTROL+ sphere sachet	NASBA™ positive control DNA: lyophilized sphere containing non infectious DNA plasmids	4 tubes in foil pack with silica gel desiccant 1 sphere per tube (white)	White cap tubes
TBSTR CPSTR	8 tubes strip: 6 pieces 8 cap strip: 12 pieces		Not applicable

Instruments/Software

Product	Catalogue number
NucliSENS EasyQ® System It includes the following: NucliSENS EasyQ® Analyzer NucliSENS EasyQ® Incubator V2 NucliSENS EasyQ® Computer Monitor Printer NucliSENS EasyQ® Director Software 2.6 Mini Microcentrifuge Operator Manuals	4700016

Additional materials

Product	Catalogue number
Micro tubes; 1.5 ml (500 x)	200294
0.2 ml 8-tube strips (125 pcs)	285048
0.2 ml 8-tube caps (125 pcs)	285051
EasyQ® Tube tray holder	45685049
EasyQ® Cap Installing Tool (Capping aid)	45685050
Pipetting aid	45685207

Interpretation of the assay results

SCC <i>mec</i> cassette junction result (FAM signal)	pos	pos	neg	neg
<i>mecA</i> result (Cyanine Dye signal)	pos	neg	pos or neg	pos or neg
IC result (ROX signal)	pos or neg	pos or neg	pos	neg
Test validity	Valid	Valid	Valid	Invalid
Test result	MRSA POS	MRSA NEG	MRSA NEG	Invalid

IC : Inhibition control, **pos** :positive, **neg** : negative

Note: IC is co-amplified with the cassette junction. Therefore, when the cassette junction is positive but IC is negative, the test result by design is valid (2 co-amplified targets compete for the same primer(s)).

J. Substantial Equivalence Information:

1. Predicate device name(s):

Assay- BD GeneOhm MRSA (formerly IDI-MRSA)

2. Predicate 510(k) number(s):

k033415

3. Comparison with predicate:

Similarities		
	Device	Predicate
Item	NucliSENS EasyQ® MRSA	BD GeneOhm MRSA (k033415)
Intended Use	Qualitative diagnostic screening for MRSA in nasal swab specimens	Same
Technological Principles	Manual specimen processing with automated amplification, detection, analysis and reporting of results	Same
Probes	Molecular Beacons	Same
Detection	Fluorogenic target-specific hybridization	Same
Analysis and reporting of results	Automated using diagnostic software	Same
DNA Target Sequences	Sequence incorporating the insertion site of Staphylococcal Cassette Chromosome <i>mec</i> (<i>SCCmec</i>) at <i>orfX</i> junction, with 6 variant sets primers and probes to target 7 MREJ types i, ii, iii, iv, v, vii, xii (target sequences for MREJ i and ii are identical)	Same

Differences		
Amplification	NASBA (Nucleic Acid Sequence-Based Amplification)	PCR (Polymerase Chain Reaction)
Instrument	NucleiSENS EasyQ® system	Cepheid SmartCycler®
Additional DNA Target	<i>mecA</i> gene	None
Controls	Internal control in each test, plus 2 external controls (positive and blank) supplied in kit, plus 2 specimen controls (MRSA and MSSA) recommended for each run	One internal reagent control and external positive and negative controls required per run

EasyQ system components comparison

Similarities		
System Component	NucliSENS EasyQ® MRSA	NucliSENS EasyQ Instrumentation and software (k093383)
Software	NucliSENS Director 2.6 software	Same

Similarities

System Component	NucliSENS EasyQ® MRSA	NucliSENS EasyQ Instrumentation and software (k093383)
Incubator	Bath, incubator/water NucliSENS EasyQ incubator II	Same
Computer	NucliSENS EasyQ computer 110V	Same
Monitor	Monitor 110V	Same
Printer	Printer 110V	Same

Differences

Specimen processing	Manual	NucliSENS miniMAG
Assay Reagents	Specific for MRSA	Specific for <i>Enterovirus</i>
Analyzer	Use of three filter pairs	Use of two filter pairs
Software: Assay Protocol	Specific for MRSA	Specific for <i>Enterovirus</i>
Vortex Genie 2, Vortex adapter plate	Required	Not required

K. Standard/Guidance Document Referenced (if applicable):

Not applicable

L. Test Principle:

The NucliSENS EasyQ® MRSA assay includes extraction, amplification and detection of MRSA DNA. Essential processes of the assay are the release of the MRSA-specific nucleic acids, NASBA™ and generation of fluorescence signals during the course of the reaction.

For fluorescence-based detection, specific molecular beacons are added to the NASBA™ reaction. These can instantaneously bind to the complementary single-stranded RNA amplicons that are generated in the amplification process. In DNA NASBA™, the bacterial DNA is first digested with a restriction enzyme and then copies of RNA derived from a target DNA sequence are manifold multiplied by isothermal amplification by the combined activities of 3 enzymes: Avian Myeloblastosis Virus Reverse Transcriptase (AMV-RT), Bacteriophage T7 RNA Polymerase (T7 RNA Pol), and RNase H, along with both ribo- and deoxyribo-nucleotides (rNTPs and dNTPs) and the cofactors dithiothreitol (DTT), MgCl₂, KCl and DMSO. To start the DNA NASBA™ process, the restriction enzyme Hae III cuts the target DNA near the binding site of primer 1 (P1). Afterwards, the reaction mixture is heated to 95°C to denature the linearised DNA and to inactivate the restriction enzyme. Subsequently, the reaction mixture is cooled down to 41°C and the primer P1 anneals to the positive strand of the target DNA. This starts the “linear phase” of the NASBA™ reaction where double-stranded DNA amplicons with a T7 promoter tail are generated. In the cyclic phase of the NASBA™ reaction, the T7 promoter tail on the cDNA target amplicon serves as a binding site for T7 RNA Pol,

which then transcribes from the cDNA amplicon using rNTP's to generate a negative strand RNA copy. A single original positive strand target sequence (with T7 tail applied) will serve as a template for multiple RNA copies as the reaction proceeds.

Molecular beacons are hairpin structure-forming single-stranded oligonucleotides labeled at one end with a fluorescent dye named the fluorophore and at the other end with a quencher. In the absence of complementary target amplicons, the hairpin structure closes, bringing the fluorophore and quencher into close proximity, and as a result no fluorescent signal is emitted.

In the presence of target amplicons, the beacon hybridizes to its complementary target instead of forming an internal hairpin structure. Once the beacon opens the quencher and fluorophore separate, quenching ceases and fluorescence is emitted.

Internal Control DNA is used to monitor nucleic acid amplification. Primer binding sites on the Internal Control DNA and MRSA DNA are identical, whereas differently labeled molecular beacons are used for the detection of the different kinds of amplicons.

Reactions are performed in closed tubes in a NucliSENS EasyQ® Analyzer in which fluorescence is continuously measured. NucliSENS EasyQ® Director Software, in combination with NucliSENS EasyQ® MRSA assay protocols software, provides automated analysis of the resulting fluorescence curves and reporting of assay results.

M. Performance Characteristics (if/when applicable):

1. Analytical performance:

a. Precision/Reproducibility:

The precision and reproducibility studies were performed on the complete assay procedure from the swab specimen to the final result. Two lots of NASBA reagents were used for both the within- and between-laboratory precision studies. A frozen pool of nasal clinical swabs was used to avoid any variability due to the use of individual clinical matrix.

A panel of 3 bacterial inputs was tested on each of 12 days, using the first lot of reagents for the first 6 days and the second lot of reagents for the last 6 days. The study involved 2 operators who each performed one run per day. Two different instruments were used. The table below demonstrated the summary of the precision studies. The data was acceptable.

Within-Laboratory Variation

Sample Input ¹	% Agreement	95% CI	MaxLambda mecA			MaxLambda cassette junction		
			Mean	SD	%CV	Mean	SD	%CV
High Negative (HN)	90.3% (65/72)	81.0-96.0%	2.669	0.0592	2.2	1.034	0.0723	7.0
Low Positive (LP)	98.6% (71/72)	92.5-99.9%	2.641	0.0782	3.0	1.578	0.2178	13.8
Moderate Positive (MP)	100% (72/72)	95.0-100%	2.680	0.0564	2.1	1.996	0.1498	7.5
Total Agreement	96.3% (208/216)	92.8-98.4%						

¹HN = approx. 0.03 x LoD
 LP= approx. 1.00 x LoD
 MP= approx. 4.82 x LoD

The between-laboratory reproducibility study was performed at 3 sites, over six days with two operators, one instrument at each site and two lots of reagent. A panel of 3 bacterial inputs was tested at each laboratory on each of 6 days, 3 days with one reagent lot and 3 days with a second reagent lot. The study involved 2 operators who each performed one run per day. Three different instruments were used (one instrument per site), with the following results:

Between –Laboratory Variation

Sample Input ¹					MaxLambda mecA			MaxLa
	Site 1	Site 2	Site 3	Total Agreement	Mean	SD	%CV	Mean
	% Agreement	% Agreement	% Agreement	% Agreement and 95% CI				
High Negative (HN)	94.4% (34/36)	100% (36/36)	86.1% (31/36)	93.5% (101/108) 87.1 - 97.4%	2.650	0.0662	2.5	1.024
Low Positive (LP)	85.7% (31/36)	75.0% (27/36)	100% (36/36)	87.0% (94/108) 79.2 - 92.7%	2.631	0.0837	3.2	1.347
Low Positive 2 (LP2)	100% (36/36)	97.2% (35/36)	100% (36/36)	99.1% (107/108) 94.9 - 99.9%	2.557	0.0995	3.9	1.735
Moderate Positive (MP)	100% (36/36)	100% (36/36)	100% (36/36)	100% (108/108) 96.6 - 100%	2.670	0.1119	4.2	1.816
Total Agreement	95.1% (137/144)	93.1% (134/144)	96.5% (139/144)	94.9% (410/432) 92.4 - 96.8%				

¹HN = approx. 0.03 x LoD
 LP= approx. 1.00 x LoD
 LP2= approx. 2.00 x LoD
 MP= approx. 4.82 x LoD

Between-Lot Variation (between –Laboratories)

Sample Input ¹	Lot 1	Lot 2	Total Agreement	MaxLambda mecA			MaxLambda cassette junction		
				Mean	SD	%CV	Mean	SD	%CV
	% Agreement	% Agreement	% Agreement and 95% CI						
High Negative (HN)	94.4% (51/54)	92.6% (50/54)	93.5% (101/108) 87.1 - 97.4%	2.650	0.0662	2.5	1.024	0.0558	5.4
Low Positive (LP)	90.7% (49/54)	83.3% (45/54)	87.0% (94/108) 79.2 - 92.7%	2.631	0.0837	3.2	1.347	0.2707	20.1
Low Positive 2 (LP2)	98.1% (53/54)	100% (54/54)	99.1% (107/108) 94.9 – 99.9%	2.557	0.0995	3.9	1.735	0.2405	13.9
Moderate Positive (MP)	100% (54/54)	100% (54/54)	100% (108/108) 96.6 -100%	2.670	0.1119	4.2	1.816	0.2550	14.0
Total Agreement	95.8% (207/216)	94.0% (203/216)	94.9% (410/432) 92.4 – 96.8%						

¹HN = approx. 0.03 x LoD
 LP= approx. 1.00 x LoD
 LP2= approx. 2.00 x LoD
 MP= approx. 4.82 x LoD

Between lot Variation (Within-Laboratory)

Sample Input ¹	Lot 1	Lot 2	Total Agreement	MaxLambda mecA			MaxLambda cassette junction		
				Mean	SD	%CV	Mean	SD	%CV
	% Agreement	% Agreement	% Agreement and 95% CI						
High Negative (HN)	91.7% (33/36)	88.9% (32/36)	90.3% (65/72) 81.0 – 96.0%	2.669	0.0592	2.2	1.034	0.0723	7.0
Low Positive (LP)	100% (36/36)	97.2% (35/36)	98.6% (71/72) 92.5 - 99.9%	2.641	0.0782	3.0	1.578	0.2178	13.8
Moderate Positive (MP)	100% (36/36)	100% (36/36)	100% (72/72) 95.0 -100%	2.680	0.0564	2.1	1.996	0.1498	7.5
Total Agreement	97.2% (105/108)	95.4% (103/108)	96.3% (208/216) 92.8 – 98.4%						

¹HN = approx. 0.03 x LoD
 LP= approx. 1.00 x LoD
 MP= approx. 4.82 x LoD

The data was acceptable.

b. Linearity/assay reportable range:

Not applicable

c. Traceability, Stability, Expected values (controls, calibrators, or methods):

The controls of the EasyQ MRSA assay included:

- An internal Inhibition Control (IC) - An internal control plasmid is incorporated in each reaction tube to detect inhibition to ensure proper reporting of an individual specimen
- Two external controls
 - One NASBA positive control to detect NASBA™ reagents and instrument failures
 - One NASBA blank control to detect contamination with strains of MRSA or DNA and contamination due to carryover of amplicons specimen controls (MRSA and MSSA)
- Two specimen processing controls:
 - One specimen positive control by rotating one of the six specified MRSA strains of different MREJ types (ii, iii, iv, v, vii and xii) into each run to check for integrity of lysis and NASBA™ reagents and absence of instrument failure
 - One specimen negative control (MSSA ATCC 29213) to detect contamination with strains of MRSA or DNA and contamination due to carry-over of amplicons in the lysis and NASBA™ reagents or in the environment

The total valid run rate was 90% (114/127); there were 13 invalid runs; 11 (84.6%) were due to specimen negative control (SNC) failure.

Specimen Controls (SPC and SNC) Stability Study

The stability study was performed using one type of SPC (MRSA MREJ type ii) and the SNC (MSSA strain). The use of one SPC only to verify the stability of the different MREJ types was accepted since the MREJ type represents the sequence of the right part of the *SCCmec* cassette, which is not linked to any phenotypic characteristics that could have an impact on the stability of MRSA. Moreover, MRSA MREJ type ii is the strain most frequently isolated in clinical practice. This was confirmed after sequencing the MRSA strains obtained during the European clinical trials of the NucliSENS EasyQ MRSA test, 87 % of the 171 MRSA strains isolated were MREJ type ii.

Based on qualitative results, SPC (MRSA MREJ type ii) and SNC (MSSA), prepared in mucin, were shown stable up to 6 months, when frozen at -19°C /- 31°C.

Reagent Stability Study

NASBA Reagent

Real-time stability data were provided from a study in real time on 3 NASBA reagent lots (PTB7, PVB1 and PVB2) produced in manufacturing conditions for a period of 19 months when stored at 2-8°C, testing MREJ wild type (WT)

plasmids corresponding to MREJ types ii, iii, iv, v, vii and xii, each plasmid is also associated with *mecA* WT plasmid. NASBA reagents included the NASBA controls: positive control (PC), external blank control (BLK), are stable for at least 19 months when stored at 2/8°C.

In-use and multiple use combined stabilities:

- All reconstituted spheres are still functional within 2 weeks after first use and one (reagent and restriction enzyme spheres, enzyme sphere) or two freeze-thawing cycles (positive control sphere) upon storage at – 19 / - 31°C
- Reconstituted sphere (reagent and restriction enzyme spheres, enzyme sphere and positive control sphere) are still functional within 30 minutes at room temperature.
- Six openings for the diluents (reagent and enzyme) and 16 openings for the NASBA water were performed on 6 non consecutive days separated by 12 days +/- 2 days to cover a total period of 2 months when stored at 2-8°C

MRSA Lysis Tubes

The study was conducted on three batches, stored at 2/8°C after thermal shocks.

Lot PTB4

Time point		T0	T'0	M3	M6	M7	M9	M12	M13	M18	M19	
Experiments BTL number		BTL 061242	BTL 061988	BTL 062959	BTL 064496	BTL 065020	BTL 066135	BTL 067156	BTL 067268	BTL 068418	BTL 068672	
Characteristics		Criteria										
Aspect of reagents	Color	No modification	Passed	Passed								
	Texture											
	Cleanliness											
	Odor											
	Precipitation											
Run controls	Negative control	3/3 valid	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3
	Positive control	3/3 valid	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3
Acceptance criteria	SNC low input	9/9 valid	9/9	9/9	9/9	9/9	9/9	9/9	9/9	9/9	9/9	9/9
	SPC low input	21/24 valid	23/24	24/24	24/24	24/24	24/24	24/24	24/24	24/24	24/24	24/24
	SPC high input	9/9 valid	9/9	9/9	9/9	9/9	9/9	9/9	9/9	9/9	9/9	9/9
	Result	Passed/failed	Passed	Passed								

Lot PVB2

Time point		T0	T'0	M3	M6	M7	M9	M12	M13	M18	M19	
Experiments BTL number		BTL 061896	BTL 065260	BTL 063585	BTL 065166	BTL 065832	BTL 066556	BTL 067393	BTL 067650	BTL 068674	BTL 068880	
Characteristics		Criteria										
Aspect of reagents	Color	No modification	Passed	Passed								
	Texture											
	Cleanliness											
	Odor											
	Precipitation											
Run controls	Negative control	3/3 valid	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	
	Positive control	3/3 valid	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	
Acceptance criteria	SNC low input	9/9 valid	9/9	9/9	9/9	9/9	9/9	9/9	9/9	9/9	9/9	
	SPC low input	21/24 valid	24/24	24/24	24/24	24/24	24/24	24/24	23/24	23/24	23/24	24/24
	SPC high input	9/9 valid	9/9	9/9	9/9	9/9	9/9	8/9 *	9/9	9/9	9/9	9/9
	Result	Passed/failed	Passed	Passed								

* M9 months (run 100106jo05, well F1). The error message was: "candidate *mecA* was not detected in the sample, no amplification for co-amplified candidate and IC".

Lot PVB1.1

Time point		T0	T'0	M3	M6	M7	M9	M12	M13	M18	M19
Experiments BTL number		BTL 061867	BTL 062686	BTL 063586	BTL 065154	BTL 065831	BTL 066360	BTL 067392	BTL 067660	BTL 068673	BTL 068883
Characteristics		Criteria									
Aspect of reagents	Color	No modification	Passed								
	Texture										
	Cleanliness										
	Odor										
Run controls	Precipitation	3/3 valid	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3
	Negative control										
Acceptance criteria	Positive control	3/3 valid	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3	3/3
	SNC low input	9/9 valid	9/9	9/9	9/9	9/9	9/9	9/9	9/9	9/9	9/9
	SPC low input	21/24 valid	23/24	24/24	24/24	23/24	24/24	24/24	24/24	23/24	24/24
	SPC high input	9/9 valid	9/9	9/9	9/9	9/9	9/9	9/9	9/9	8/9 *	9/9
Result		Passed/failed	Passed								

* M9 months (run 100106jo02, well G1). The error message was: "candidate *mecA* was not detected in the sample, no amplification for co-amplified candidate and IC".

The additional real time stability data of the MRSA lysis tubes up to 19 months demonstrated the stability of 19 months for NucliSENS EasyQ[®] MRSA assay. An additional study was conducted to compare the new pre-marked lysis tube to the former unmarked lysis tube. The study demonstrated that the unmarked tube were 85% positive results and the new pre-marked tubes were 90% positive when testing at slightly lower than the LoD (191 CFU/swab).

Specimen, lysates stability

Specimen can be stored:

- Up to 4 days at 2-8 °C before testing (with or without a maximum of 4h transport at room temperature).
- Up to 4 days at room temperature before testing.
- Up to 1 day at room temperature followed by 3 days at 2-8 °C.

Result summary for clinical swabs stored under various conditions and duration

Temp	Storage Conditions		Result	
	Hours		MRSA (+)	MRSA (-)
RT	24		10	0
	48		9*	0
	72		10	0
	96		10	0
2-8°C	24		10	0
	48		10	0
	72		10	0
	96		10	0
4 hrs RT AND 96 hrs at 2-8°C	24		10	0
	48		10	0
	72		10	0
	96		10	0
24 hrs RT AND additional 72 hrs at 2- 8°C	24 (RT)		10	0
	48 (2-8°C)		10	0
	72 (2-8°C)		10	0
	96 (2-8°C)		10	0

* One result was discarded for the analysis because the second swab corresponding to the same patient was characterized MRSA positive

Lysates can be stored:

- Up to 24 hours at 2-8 °C.
- Up to 1 month at -19 °C/-31 °C. A maximum of 5 freeze/thaw cycles are allowed within the month.

Result summary for lysates stored under different conditions and duration

Storage Conditions			Result	
Time	Temperature	freeze/thaw cycles#	MRSA (+)	MRSA (-)
Time zero*	on ice	NA	10	0
4 hrs	2-8°C	NA	9	1
24 hrs	2-8°C	NA	10	0
4 days	-19°C/-31°C	2	10	0
1 month	-19°C/-31°C	2	10	0
1 month	-19°C/-31°C	4	10	0
1 month	-19°C/-31°C	5	10	0

*Time zero: Lysate kept on ice or cooling block while preparing the NASBA

#Freeze/Thaw cycle: freeze at indicated temperature for ≥ 4 days, followed by thaw at RT.

d. *Detection limit:*

Limit of Detection Study

The analytical sensitivity of the NucliSENS EasyQ[®] MRSA assay was determined using various strains of MRSA representing the seven types of MREJ detected by the test (MREJ types i, ii, iii, iv, v, vii and xii) using frozen matrices of clinical swabs. Bacterial suspensions with a bacterial input ranging roughly from 10 to 1000 CFU per swab were tested in 24 replicates of each input using a single lot of reagent. The limit of detection was estimated as the input corresponding to the 95% Hit Rate predicted by statistical analysis using a statistical regression model (Probit). The estimated type-specific LoDs were then confirmed on two reagent lots with 15 replicates each.

The claimed type-specific LoD summarized in the following table corresponds to the higher concentrations claimed from these two studies, i.e. either real input determined by plating and counting bacteria colony forming units in the confirmation study, or predicted concentrations by the probit regression model.

LoD for the seven MRSA types

MREJ type	SCC _{mec} type	LoD (CFU/swab)
i	I	410
ii	II	227
iii	III	493
iv	III	234
v	IV	200
vii	III	482
xii	V	182

Analytical Inclusivity

A total of 193 well characterized MRSA isolates were tested at an input level of approximately 3 times the LoD, in the presence of frozen specimen matrix from nasal swab specimens previously shown to be MRSA negative, using the NucliSENS EasyQ® MRSA assay, in the Inclusivity Study. The strains included pulse field gel electrophoresis (PFGE) types: USA100, 200, 300, 400, 500, 600, 700, 800, 1000, and 1100. Nine were not detected, resulted a detection rate of 95.3% (184/193). Three of the undetected strains were MREJ type iii.

Ten heterogeneous oxacillin resistant strains were tested with acceptable results.

Six "empty cassette" variants (i.e. *mecA* dropouts) were tested with *mecA* gene negative clinical matrices at both input levels of 7.5×10^6 CFU/swab and 1.0×10^4 CFU/swab. They all gave a positive signal for the cassette junction target and a negative signal for the *mecA* gene target, leading to a MRSA negative result. However, the risk of false positive results exists with swabs that possess the *mecA* gene (e.g. carried by methicillin/oxacillin resistant coagulase negative *Staphylococci*) and MSSA *mecA* drop-out strains.

e. **Analytical specificity:**

Cross Reactivity

A total of 100 phylogenetically related *S. aureus* and members of the nasal commensal flora were tested at an input of 7.5×10^6 CFU/swab. There were 26 methicillin-resistant coagulase negative staphylococci (MRCN), 29 methicillin-sensitive coagulase-negative staphylococci (MSCN), 25 methicillin-sensitive coagulase-positive staphylococci (MSCP) of which 20 were MSSA, and 5 were *S. intermedius*; 20 organisms other than *Staphylococci* belonging to the nasal commensal flora (1 Gram-positive rods,

10 Gram-negative rods, 7 non-staphylococci Gram-positive cocci, 1 Gram-negative cocci, 1 yeast).

An additional study was conducted on eight well characterized *K. pneumoniae* subsp. *pneumoniae*, *K. pneumoniae* subsp. *ozaenae* with *mecA* gene negative clinical matrices at both input levels of 7.5×10^6 CFU/swab and 1.0×10^4 CFU/swab. Results demonstrated that at an input level of 7.5×10^6 CFU/swab, the 8 strains showed a positive signal for the cassette junction target and a negative signal for the *mecA* target. When tested at an input level of 1.0×10^4 CFU/swab, the eight strains showed a negative signal for both targets. Therefore, the risk of false positive results exists with swabs that possess the *mecA* gene (e.g. carried by methicillin resistant coagulase negative *Staphylococci* in the clinical sample) with *K. pneumoniae* subsp. *pneumoniae*, *K. pneumoniae* subsp. *ozaenae* due to a positive signal for the cassette junction target at a high input of 7.5×10^6 CFU/swab.

Biological Interference Study

To assess whether the NucliSENS EasyQ[®] MRSA assay could be affected by the presence of other microorganisms that are normal nasal flora (biological interference), an analytical study was carried out to test potential biological interference on both MRSA positive swabs, as well as on MRSA negative swabs using frozen matrices of clinical swabs.

To determine the potential interference of common nasal microorganisms on MRSA positive swabs, dry flocked swabs were spiked with a strain of MRSA at 3 times the LoD and with an excess of a potential interfering microorganism (approximately 10^6 CFU/swab). A control test was performed without interfering microorganisms using dry flocked swabs spiked with MRSA only and had to be MRSA positive validate the experiment.

To determine the potential interference of common nasal microorganisms on MRSA negative swabs, dry flocked swabs were spiked with an excess of the potentially interfering microorganism. A control test was performed without interfering microorganisms using dry flocked swabs spiked with sterile water only and had to be MRSA negative to validate the experiment.

The spiked flocked swabs were then introduced into the lysis tubes containing a thawed matrix of MRSA negative swabs and then submitted to the sample preparation and NASBA amplification
Each type of spiked swab was tested in four replicates.

Two different MRSA strains (MREJ types ii and iv) were used in the study:

- MRSA MREJ type ii was used because it is the most frequently isolated strain in clinical practice (as confirmed by sequencing of the MRSA strains obtained during the European clinical trials for the NucliSENS

EasyQ MRSA test, 87 % of the 171 MRSA strains isolated were MREJ type ii).

- MRSA MREJ type iv was used because the Internal control (IC) is also a type iv and competes with this target for the primers. Therefore, MREJ iv strain is expected to be more susceptible to the potential effect of biological interference than any other MREJ type.

Three strains were used as potentially interfering microorganisms in this study:

- Methicillin Sensitive *Staphylococcus aureus* (MSSA)
- Methicillin Resistant *Staphylococcus epidermidis* (MRSE)
- Methicillin Sensitive *Staphylococcus epidermidis* (MSSE)

The targeted input for the MRSA strains was 3 times the LoD for both type ii and type iv. The real input on the swab was determined by plate counting of the suspension used to spike the swab.

For both the MREJ type ii and the MREJ type iv spiked swabs, the 4 replicates all gave “MRSA positive” results in the presence of the interfering microorganisms (MSSA, MRSE or MSSE). For the negative swabs, the 4 replicates all gave “MRSA negative” results in the presence of the interfering microorganisms (MSSA, MRSE or MSSE).

Chemical Interference Study

To assess whether the NucliSENS EasyQ[®] MRSA assay could be affected by the presence of well identified chemical substances, which can be naturally present in the nasal cavity or can be artificially introduced into the nose, an analytical study was carried out to test potential chemical interference on both MRSA positive swabs, as well as on MRSA negative swabs using frozen matrices of clinical swabs.

To determine the potential chemical interference on MRSA positive swabs, dry flocked swabs were spiked with a strain of MRSA at 3 times the LoD and with an excess of a potential interfering chemical substance (50 ul). A control test was performed without interfering chemical substance using dry flocked swabs spiked with MRSA only and had to be MRSA positive validate the experiment.

The spiked flocked swabs were then introduced into the lysis tubes containing a thawed matrix of MRSA negative swabs and then submitted to the sample preparation and NASBA amplification. Each type of spiked swab was tested in 4 replicates. If a chemical substance yielded 1) invalid results or false negative results for the swabs spiked with MRSA or 2) invalid results for the swabs without MRSA, then 20 additional replicates were performed to confirm the initial results.

Two different MRSA strains (MREJ types ii and iv) were used in the study:

- MRSA MREJ type ii was used because it is the most frequently isolated strain in clinical practice (as confirmed by sequencing of the MRSA strains obtained during the European clinical trials for the NucliSENS EasyQ MRSA test, 87 % of the 171 MRSA strains isolated were MREJ type ii).
- MRSA MREJ type iv was used because the Internal control (IC) is also a type iv and competes with this target for the primers. Therefore, MREJ iv strain is expected to be more susceptible to the potential effect of biological interference than any other MREJ type.

The following naturally or artificially introduced chemical substances were tested in the study:

- Whole blood at different levels
- Mucus (using mucin)
- Astelin (azelastine HCl) antihistamine nasal spray
- Nasonex (mometasone furoate), steroidal nasal spray
- Nasacort (acetonide triamcinolone), glucocorticoid nasal spray
- Vicks Sinex (oxymetazoline HCl), topical decongestants nasal spray
- Neo-Synephrine (phenylephrine HCl 1.0%), decongestant nasal spray
- Sino Fresh (eucalyptus globulus, kalium bichromicum), antiseptic homeopathic nasal spray
- NasalCrom (cromolyn sodium), antiallergic nasal spray
- Zicam cold remedy (zincum gluconicum), homeopathic nasal spray

The targeted input for the MRSA strains was 3 times the LoD for both MREJ type ii and MREJ type iv. The real input on the swab was determined by plate counting of the suspension used to spike the swab.

Five chemical substances were found to be interferents: Astelin, Sino Fresh, Nasal chrom, Zicam and Blood (15 and 50 µl).

Carry-Over Contamination

Carryover study was performed to assess the carry-over risk of the NucliSENS EasyQ MRSA assay during sample preparation and amplification/detection. Five NucliSENS EasyQ MRSA NASBA runs were performed on MRSA high positive (positive for both the cassette junction and *mecA*) and MRSA negative specimen (negative for both the cassette junction and *mecA*) in an alternating pattern, resulting in a total of 105 MRSA high positive samples and 105 MRSA negative samples.

This study was performed using frozen matrices of clinical swabs. In order to obtain MRSA high positive specimen, lysis tubes containing a matrix of MRSA negative clinical swabs were spiked with a high concentration of MRSA MREJ type ii (5×10^6 CFU/lysis tube). In order to obtain MRSA negative specimen, lysis tubes containing

a matrix of MRSA negative clinical swabs, that were proven to be negative for the cassette and *mecA*, were spiked with sterile water.

All 105 replicates of the MRSA high positive sample were tested MRSA positive. All 103 out of 103 valid replicates of the MRSA negative sample were tested MRSA negative (negative for both the cassette and *mecA*). Two MRSA negative replicates were tested “invalid” due to failure of detecting IC. No cross-over contamination was observed amongst the 103 negative samples surrounded by high positive samples.

f. Assay cut-off:

Studies were conducted to determine the following thresholds (cut-offs) for the interpretation of the run controls and the specimen controls. The final thresholds for the controls were further validated in the prospective U.S. clinical trial.

Threshold value	SPC	SNC	PC	BLK
MaxLambda cassette junction	1.2	1.08	1.5	1.08
MaxLambda <i>mecA</i>	2.0	1.1	2.0	1.1
MaxLambda IC	2.0	3.5	2.0	3.5

2. Comparison studies:

a. Method comparison with predicate device:

Not applicable

b. Matrix comparison:

Not applicable

3. Clinical studies:

Prospective specimens were collected at seven geographically diverse U.S. sites, of which two were for specimen collections only, and were sent to the other five sites for testing. The testing sites also included two pediatric sites. The NucliSENS EasyQ MRSA Assay was compared to a reference enriched culture. The reference culture was trypticase soy broth with 6.5% NaCl, incubated for 48 hrs if negative; any suspected *S. aureus* was confirmed by Gram stain, latex agglutination or DNase. The *mecA* mediated oxacillin resistance was tested by disk diffusion using a 30µg cefoxitin disk and report of ≤21 mm (R), ≥22 mm (S).

There were a total of 1355 nasal specimens, with 1238 reported results. There were 42 initial invalid results, resulting in an invalid rate of 3.4% (42/1238) from 114 valid runs. Fourteen remained invalid after repeat testing from frozen lysate.

The overall invalid run rate was 10.2% (13/127).

Clinical Sensitivity/Specificity:

The performance of the NucliSENS EasyQ® MRSA assay (Adult + pediatric) all sites combined is summarized in the following tables. Invalids were not included in the performance calculation of adults, and pediatric+ adults; there were no invalids in pediatric samples:

NucliSENS EasyQ® MRSA Performance (Adult + Pediatric) comparing to Reference Culture

NucliSENS	Reference Culture*		Total
	MRSA+	MRSA-	
MRSA+	137	35	172
MRSA-	6	1046	1052
Total	143	1081	1224
MRSA:		95%CI:	
Sensitivity: (137/143)	95.8%	(91.1 - 98.4%)	
Specificity: (1046/1081)	96.8%	(95.5 - 97.7%)	
Positive Predictive value	79.7%	(72.9 - 85.4%)	
Negative Predictive value	99.4%	(98.8 - 99.8%)	

* *mecA* mediated oxacillin resistance using 30µg cefoxitin disk

NucliSENS EasyQ® MRSA Performance (Adults) comparing to Reference Culture

NucliSENS	Reference Culture*		Total
	MRSA+	MRSA-	
MRSA+	107	26	133
MRSA-	6	724	730
Total	113	750	863
MRSA:		95%CI:	
Sensitivity: (107/113)	94.7%	(88.8 - 98.0%)	
Specificity: (724/750)	96.5%	(95.0 - 97.7%)	
Positive Predictive value	80.5%	(72.7 - 86.8%)	
Negative Predictive value	99.2%	(98.2 - 99.7%)	

* *mecA* mediated oxacillin resistance using 30µg cefoxitin disk

NucliSENS EasyQ® MRSA Performance (Pediatric) comparing to Reference Culture

NucliSENS	Reference Culture*		Total
	MRSA+	MRSA-	
MRSA+	30	9	39
MRSA-	0	322	322
Total	30	331	361
MRSA:		95%CI:	
Sensitivity: (137/143)	100%	(88.4 – 100%)	
Specificity: (1046/1081)	97.3%	(94.9 -98.8%)	
Positive Predictive value	76.9%	(60.7 – 88.9%)	
Negative Predictive value	100%	(98.9 - 100%)	

* *mecA* mediated oxacillin resistance using 30µg cefoxitin disk

c. *Other clinical supportive data (when a. and b. are not applicable):*

Not applicable

4. Clinical cut-off:

Not applicable

5. Expected values/Reference range:

A total of 1355 nasal specimens were collected from 1355 patients at 7 enrolling sites across the United States. There were a total of 1238 final evaluable specimens by reference culture. The study population was grouped into subjects hospitalized for greater than 3 days, hospitalized for 3 days or less, outpatients, Emergency Room (ER) patients, patients from long term care facilities and healthcare workers. Hospitalized patients included, but are not limited to patients previously known to be colonized with MRSA or with MRSA infection, transferred from a long term care institution/nursing home, history of prolonged hospitalization, previous hospital stay, contact with a known MRSA carrier, admitted to intensive care unit, cardiology unit, orthopedic surgery unit or hemodialysis unit, history of intravenous drug use, on contact precautions, and HIV infected. The number and percentage of positive and negative specimens with the reference culture method are presented in the following table:

Group	% Positive	% Negative	Total
Hospitalized >3 days	16.7% (65/390)	83.3% (325/390)	31.5% (390/1238)
Hospitalized ≤ 3 days	10.4% (49/471)	89.6% (422/471)	38.0% (471/1238)
Long Term Care/ Nursing Home	16.7% (3/18)	83.3% (15/18)	1.5% (18/1238)
Outpatients	9.3% (21/227)	90.7% (206/227)	18.3% (227/1238)
ER	4.5% (4/88)	95.5% (84/88)	7.1% (88/1238)
Healthcare workers	4.6% (2/44)	95.4% (42/44)	3.6% (44/1238)
Total	11.6% (144/1238)	88.4% (1094/1238)	100% (1238/1238)

N. Instrument Name:

NucliSENS EasyQ[®] MRSA assay protocol software included in the NucliSENS EasyQ[®] MRSA assay kit

NucliSENS EasyQ[®] Director Software 2.6 (Cleared in k093383)

NucliSENS EasyQ[®] Analyzer (Cleared in k093383)

NucliSENS EasyQ[®] Incubator II (Cleared in k093383)

O. System Descriptions:

1. Modes of Operation:

Batch or Random access through tube strip (Analyzer)

2. Software:

FDA has reviewed applicant's Hazard Analysis and software development processes for this line of product types:

Yes X or No

3. Specimen Identification:

Entered by user (Analyzer)

4. Specimen Sampling and Handling:

Specimens are processed according to assay instructions

5. Calibration:

The analyzer contains an auto-calibration function that reads the fluorescence from the reference materials on the plate carrier. The reading is compared to the value in the system memory. The measured value can differ from the memory value.

6. Quality Control:

The EasyQ MRSA assay includes:

- Inhibition Control (IC) - An internal control plasmid is incorporated in each reaction tube to detect inhibition to ensure proper reporting of an individual specimen
- NASBA positive control to detect NASBA[™] reagents and instrument failures

NASBA blank control to detect contamination with strains of MRSA or DNA and contamination due to carryover of amplicons specimen controls (MRSA and MSSA)

- Specimen positive control to check for integrity of lysis and NASBA™ reagents and absence of instrument failure
- Specimen negative control to detect contamination with strains of MRSA or DNA and contamination due to carry-over of amplicons in the lysis and NASBA™ reagents or in the environment

P. Other Supportive Instrument Performance Characteristics Data Not Covered In The “Performance Characteristics” Section above:

Not applicable

Q. Proposed Labeling:

The labeling is sufficient and it satisfies the requirements of 21 CFR Part 809.10.

R. Conclusion:

The submitted information in this premarket notification is complete and supports a substantial equivalence decision