

**510(k) SUBSTANTIAL EQUIVALENCE DETERMINATION
DECISION SUMMARY
ASSAY ONLY TEMPLATE**

A. 510(k) Number:

k103162

B. Purpose for Submission:

New device

C. Measurand:

Calibration verification and assay range verification for Elecsys Parathyroid Hormone CalCheck 5

D. Type of Test:

Not applicable

E. Applicant:

Roche Diagnostics

F. Proprietary and Established Names:

Elecsys Parathyroid Hormone CalCheck 5

G. Regulatory Information:

1. Regulation section:

21 CFR 862.1660, Quality Control Material (Assayed and Unassayed)

2. Classification:

Class I, reserved

3. Product code:

JJX – Single (Specified) Analyte Controls (Assayed and Unassayed)

4. Panel:

Clinical Chemistry (75)

H. Intended Use:

1. Intended use(s):

See indications for use below

2. Indication(s) for use:

The Elecsys Parathyroid Hormone (PTH) CalCheck 5 is an assayed control for use in calibration verification and for use in the verification of the assay range established by the Elecsys PTH reagent on the indicated Elecsys and cobas e immunoassay analyzers

3. Special conditions for use statement(s):

For prescription use only

4. Special instrument requirements:

Roche Elecsys 2010 / cobas e 411 and MODULAR ANALYTICS E170 / cobas e 601 / cobas e 602 immunoassay analyzers

I. Device Description:

Elecsys PTH CalCheck 5 set is an *in vitro* diagnostic device that contains 5 levels of PTH in a human serum matrix (all human derived products are prepared exclusively from donors tested and found free of HBsAg and antibodies to HCV and HIV). The vials contain lyophilized material which is reconstituted to 1.0 mL each and are packaged in sets of 5. Each CalCheck contains 0.08% methylisothiazolinone and 0.1% oxyprion as preservatives.

CalCheck Level	Approximate Target Range (pg/mL)
1	≤ 1
2	55 - 65
3	2250 - 2750
4	3800 - 4200
5	4800 - 5200

J. Substantial Equivalence Information:

1. Predicate device name(s):

Elecsys Parathyroid Hormone CalCheck

2. Predicate K number(s):

k993642

3. Comparison with predicate:

Similarities		
Item	Device	Predicate
Indications for Use	Same	For use in the verification of the calibration established by the Elecsys PTH reagent on the indicated Elecsys and cobas e immunoassay analyzers.
Analyte	Same	PTH
Format	Same	Lyophilized
Handling	Same	Reconstitute each CalCheck level with exactly 1.0 mL distilled water. Allow to stand closed for 15 minutes, then gently mix by inversion.
Matrix	Same	Human serum matrix
Assay Measuring Range	Same	1.20 – 5000 pg/mL

Differences		
Item	Device	Predicate
Levels	Five	Three
Stability	Unopened: • Store at 2-8°C until expiration date Reconstituted: • 20-25°C: 5 hours	Unopened: • Store at 2-8°C until expiration date Reconstituted: • 20-25°C: 4 hours
Check Target Values (pg/mL)	Check 1: ≤ 1 Check 2: 60 Check 3: 2500 Check 4: 4000 Check 5: 5000	Check 1: < 5 Check 2: 60 Check 3: 3000

K. Standard/Guidance Document Referenced (if applicable):

Guidance for Industry and FDA Staff – Assayed and Unassayed Quality Control Material

L. Test Principle:

Not applicable

M. Performance Characteristics (if/when applicable):

1. Analytical performance:

a. *Precision/Reproducibility:*

Not applicable

b. *Linearity/assay reportable range:*

Not applicable

c. *Traceability, Stability, Expected values (controls, calibrators, or methods):*

Traceability and Value Assignment

A Master calibration curve is generated on the MODULAR ANALYTICS E170 / cobas e 601 / cobas e 602 analyzer and is used to calculate pg/mL PTH in each CalCheck for value assignment. Each master calibrator is standardized using a commercial PTH radioimmunoassay (RIA).

For each Elecsys PTH CalCheck 5 lot manufactured, the CalCheck devices are assayed in duplicate on at least three MODULAR ANALYTICS E170 / cobas e 601 / cobas e 602 analyzers (master platform). The assigned value of each CalCheck is defined as the median value obtained over at least 6 determinations (duplicate runs on at least 3 analyzers) of the respective CalCheck. The assigned range is calculated as $\pm 30\%$ of the assigned value (chosen as the sum of $\pm 10\%$ for between-analyzer platform tolerance, $\pm 10\%$ for stability tolerance, and $+10\%$ CV of the assay). The label states that each laboratory should establish appropriate acceptance criteria when using this product for its intended use.

To ensure the values assigned using the master platform are transferrable and valid for the other instrument platforms, the value assignment process is repeated for each additional analyzer and the assigned values are compared to the values obtained on the master platform. The additional analyzer assigned values must be within 10% of the values assigned using the master platform. After this acceptance criterion is met, the assigned values from the master platform are deemed valid for the Elecsys 2010 / cobas e 411 immunoassay analyzers.

Stability

Open-vial and accelerated stability study protocols and acceptance criteria were described. The data support the package insert claims that open-vial reconstituted Elecsys PTH CalCheck 5 is stable for up to 5 hours at 20-25°C, and that closed vial storage at 2-8°C provides stability for 18 months. The expiration date listed on each CalCheck bottle label reflects this stability time of 18 months. A real-time stability study is ongoing.

d. Detection limit:

Not applicable

e. Analytical specificity:

Not applicable

f. Assay cut-off:

Not applicable

2. Comparison studies:

a. Method comparison with predicate device:

Not applicable

b. Matrix comparison:

Not applicable

3. Clinical studies:

a. Clinical Sensitivity:

Not applicable

b. Clinical specificity:

Not applicable

c. Other clinical supportive data (when a. and b. are not applicable):

Not applicable

4. Clinical cut-off:

Not applicable

5. Expected values/Reference range:

There are 5 levels of the CalCheck 5 kit which span the 1.2-5000 pg/mL measuring range of the PTH assay (see Device Description under section I above). The final assigned values are lot specific and are provided in the labeling. There is a note in the labeling stating that the level 1 CalCheck recovers below the detection limit.

If CalCheck 5 generated a result greater than 5000 pg/mL, the sponsor recommends that the user perform a 1:1 dilution of CalCheck 5 with CalCheck 4 to determine whether it is within the assay measuring range or not. Dilution studies were performed where CalCheck 4 and CalCheck 5 were mixed in a 1:1 ratio (as instructed in the draft package insert for dilutions) and measured in two-fold determination. The PTH CalCheck 5 dilution values were within the assay measuring range after dilution.

N. Proposed Labeling:

The labeling is sufficient and it satisfies the requirements of 21 CFR Part 809.10.

O. Conclusion:

The submitted information in this premarket notification is complete and supports a substantial equivalence decision.