

**510(k) SUBSTANTIAL EQUIVALENCE DETERMINATION  
DECISION MEMORANDUM  
ASSAY ONLY TEMPLATE**

**A. 510(k) Number:**

K141338

**B. Purpose for Submission:**

To obtain a substantial equivalence determination for the Liat™ Strep A Assay performed on the Liat™ Analyzer for the detection of *Streptococcus pyogenes* (Group A *Streptococcus*).

**C. Measurand:**

Group A *Streptococcus* DNA

**D. Type of Test:**

Real-time PCR assay for qualitative detection of Group A *Streptococcus* DNA in throat swab specimens.

**E. Applicant:**

IQuum, Inc.

**F. Proprietary and Established Names:**

Liat™ Strep A Assay

**G. Regulatory Information:**

1. Regulation section:

21 CFR 866.2680

2. Classification:

Class II

3. Product code:

PGX

4. Panel:

83 - Microbiology

**H. Intended Use:**

1. Intended use(s):

The Liat™ Strep A Assay, performed on the Liat™ Analyzer, is a qualitative *in vitro* diagnostic test for the detection of *Streptococcus pyogenes* (Group A β-hemolytic *Streptococcus*, Strep A) in throat swab specimens from patients with signs and symptoms of pharyngitis.

The Liat™ Strep A Assay utilizes nucleic acid purification and polymerase chain reaction (PCR) technology to detect *Streptococcus pyogenes* by targeting a segment of the *Streptococcus pyogenes* genome.

2. Indication(s) for use:

Same as Intended Use

3. Special conditions for use statement(s):

For prescription use only.

4. Special instrument requirements:

IQuum Liat™ Analyzer

**I. Device Description:**

The Liat™ Strep A Assay, performed on the Liat™ Analyzer, is a rapid, automated *in vitro* diagnostic test for the qualitative detection of *Streptococcus pyogenes* (Group A β-hemolytic *Streptococcus*, Strep A) DNA in throat swab specimens in Amies medium.

The Liat™ Strep A Assay targets a conserved region of the *S. pyogenes* genome. An Internal Process Control (IPC) is also included to monitor the adequacy of process steps involved in nucleic acid extraction and amplification/detection, as well as for the presence of inhibitors. In order for a sample to be called negative for *S. pyogenes*, the IPC must be detected. If the IPC is not detected, the result is reported as invalid and the operator is instructed to repeat the test. The time-to-result is approximately 15 minutes.

The assay employs a single-use disposable Liat™ Tube that holds the sample preparation and PCR reagents, and in which the nucleic acid extraction and amplification/detection processes take place. The reagents are housed in unit dose pre-packed tube segments separated by frangible seals that are ruptured sequentially by actuators within the Liat™ Analyzer to effect

sample processing, including bacterial cell lysis, DNA recovery and removal of inhibitors, PCR amplification and detection.

To perform the Liat™ Strep A Assay the operator transfers an aliquot of a throat swab sample in Amies medium into a Liat™ Strep A Assay tube, scans the relevant tube and sample identification barcodes and then inserts the tube into the Liat™ Analyzer for automated processing and result interpretation. An embedded microprocessor coordinates the functions of the Liat™ Analyzer to move the sample from one segment of the tube to another, and to control the reaction volume, temperature and duration of each step.

Positive and Negative External Controls are provided separately from the assay reagents in the Liat™ Strep A Assay Quality Control Kit. The Positive Control comprises dried, inactivated *S. pyogenes* bacteria that are rehydrated by the operator with Amies medium. The Negative Control comprises Amies medium alone.

**J. Substantial Equivalence Information:**

1. Predicate device name(s):

Lyra™ Direct Strep Assay

2. Predicate 510(k) number(s):

K133883

3. Comparison with predicate:

Similarities		
Item Name	Liat™ Strep A	Lyra™ Direct Strep
Intended Use	<p>The Liat™ Strep A Assay, performed on the Liat™ Analyzer, is a qualitative in vitro diagnostic test for the detection of <i>Streptococcus pyogenes</i> (Group A β-hemolytic <i>Streptococcus</i>) in throat swab specimens from patients with signs and symptoms of pharyngitis.</p> <p>The Liat Strep A Assay utilizes nucleic acid purification and polymerase chain reaction (PCR) technology to detect <i>Streptococcus pyogenes</i> by targeting a segment of the <i>Streptococcus pyogenes</i> genome.</p>	<p>The Lyra Direct Strep Assay is a Real-Time PCR <i>in vitro</i> diagnostic test for the qualitative detection and differentiation of Group A β-hemolytic <i>Streptococcus</i> (<i>Streptococcus pyogenes</i>) and pyogenic Group C and G β-hemolytic <i>Streptococcus</i> nucleic acids isolated from throat swab specimens obtained from patients with signs and symptoms of pharyngitis, such as sore throat. The assay does not differentiate between pyogenic Groups C and G β-hemolytic <i>Streptococcus</i>.</p> <p>All negative test results should be confirmed by bacterial culture, because negative results do not preclude Group A, C or G Strep infection and should not</p>

<b>Similarities</b>		
<b>Item Name</b>	<b>Liat™ Strep A</b>	<b>Lyra™ Direct Strep</b>
		be used as the sole basis for treatment. The assay is intended for use in hospital, reference, or state laboratory settings. The device is not intended for point-of-care use.
Regulation	21 CFR 866.2690	Same
Product Code	PGX	Same
Analyte	Group A <i>Streptococcus</i> ( <i>S. pyogenes</i> )	Group A, C/G <i>Streptococcus</i>
Sample Type	Throat swab	Same
Internal Control	Yes	Yes
Strep A Target	Conserved region of Group A <i>Streptococcus</i> genome	Conserved regions within the genomes of Group A and C/G streptococci
Assay Method	PCR for detecting the presence / absence of bacterial DNA in clinical specimens	Same
Detection Technique	Different reporter dyes for target analyte and Internal Control	Same
Assay Result	Qualitative	Same

<b>Differences</b>		
<b>Item Name</b>	<b>Liat™ Strep A</b>	<b>Lyra™ Direct Strep</b>
Sample Processing	Automated bacterial lysis and silica-magnetic bead-based nucleic acid extraction and purification	Manual heat lysis without nucleic acid purification
Bacterial Lysis	Chaotrope and enzymatic digestion	Heat
Equipment Required	Liat™ Analyzer	<ul style="list-style-type: none"> <li>• ABI 7500 Fast Thermocycler</li> <li>• Plate centrifuge for 96 well plate</li> <li>• Heat block</li> <li>• Thermometer</li> <li>• Timer</li> <li>• Micropipette</li> </ul>
Assay Automation	Yes: computer controlled sample processing and PCR amplification/detection	No: manual sample processing and PCR set-up

Differences		
Item Name	Liat™ Strep A	Lyra™ Direct Strep
Reagents / Kit Components	<ul style="list-style-type: none"> <li>• Unitized Liat™ Strep A Assay Tube</li> <li>• Transfer pipette</li> </ul>	<ul style="list-style-type: none"> <li>• Unitized Process Buffer for heat lysis</li> <li>• Bulk PCR Master Mix</li> <li>• Bulk Rehydration Solution for Master Mix</li> </ul>
Reagent Format	<ul style="list-style-type: none"> <li>• Unitized ready for use</li> <li>• Manual reagent manipulation not required</li> </ul>	<ul style="list-style-type: none"> <li>• Bulk reagents</li> <li>• Manual pipetting required</li> </ul>
Result Interpretation	Automated	Manual
Time-to-result	~15 minutes	>70 minutes

**K. Standard/Guidance Document Referenced (if applicable):**

Not applicable

**L. Test Principle:**

To perform the Liat™ Strep A Assay a throat swab sample in Amies medium is first added to the Liat™ Tube using a transfer pipette. The tube is then loaded into the Liat™ Analyzer for automated processing, the first step of which is mixing of the sample with the Internal Process Control (IPC). Lysis of both the IPC and the *S. pyogenes* target organism, if present, are then induced by chaotropic and proteolytic reagents and the nucleic acids are recovered by binding to the surface of silica-coated magnetic beads. The beads and bound nucleic acid are captured using a magnetic field and the lysate is removed, after which a wash step is conducted to remove potential inhibitors and the captured nucleic acids are eluted under low-salt conditions into a small volume of buffer.

Target amplification and detection are performed using TaqMan (hydrolysis) probe-based real-time PCR. Primers and probes for both the *S. pyogenes*-specific target and the IPC are included in the reaction mixture. The fluorescence intensity for each optical channel is monitored and results are interpreted by an automated algorithm based on a combination of cycle threshold (Ct) and endpoint fluorescence values.

The sample preparation and real-time PCR amplification processes are conducted within the Liat™ Analyzer in approximately 15 minutes.

**M. Performance Characteristics (if/when applicable):**

1. Analytical performance:

a. *Precision/Reproducibility:*

The reproducibility of the Liat™ Strep A Assay was evaluated at 3 sites. Two operators at each site tested a 4 member panel in triplicate on 5 different days, for a total of 360 assay runs (4 panel members X 3 replicates X 2 operators X 5 days X 3 sites). Nine Liat™ Analyzers and 3 Liat™ Strep A Assay tube lots were included in the study. Each reproducibility panel comprised a negative, a high negative, a low positive and a medium positive sample prepared in throat swab matrix (**Table 1**).

**Table 1. Reproducibility Panel Members**

	Strain	Multiple of LOD	CFU/Test	CFU/mL	Expected Result for <i>S. pyogenes</i> <sup>1</sup>
Negative	Not applicable	0	0	0	Negative
High Negative	ATCC BAA-946	0.03X	0.03	0.15	Negative
Low Positive		1X	1	5	Positive
Moderate Positive		3X	3	15	Positive

<sup>1</sup> Used in determination of percent agreement

The percent agreement with expected results for each panel member from each site is shown in **Table 2** for the *S. pyogenes* target and in **Table 3** for the Internal Process Control (IPC). The mean and coefficient of variation (%CV) for the Ct and endpoint fluorescence values are also provided for all samples with valid Ct scores (*i.e.*, this analysis was not performed for the *S. pyogenes* target in Negative or High Negative target in Negative or High Negative samples).

All (90/90; 100%) Moderate Positive and 89/90 (98.9%) Low Positive samples produced positive assay results. All (90/90; 100%) Negative and High Negative samples produced negative assay results. Overall there was 99.7% (359/360) agreement with expected results.

At the Low Positive target level, the %CV for the *S. pyogenes* Ct and endpoint fluorescence values ranged from 2.5 to 3.6% and from 29.2 to 43.8%, respectively depending on the site. For the IPC in *S. pyogenes* negative samples, the %CV for Ct ranged from 1.7 to 1.9% and for endpoint fluorescence between 12.2 and 13.3%. Overall, the PCR Ct and endpoint fluorescence values for both the *S. pyogenes* target and IPC demonstrated acceptable precision.

**Table 2. Reproducibility Study Results for *S. pyogenes* Stratified by Site**

	Site A					Site B					Site E					Total				
	Agmt (%)	Ct		Endpoint		Agmt (%)	Ct		Endpoint		Agmt (%)	Ct		Endpoint		Agmt (%)	Ct		Endpoint	
		Mean	%CV	Mean	%CV		Mean	%CV	Mean	%CV		Mean	%CV	Mean	%CV		Mean	%CV	Mean	%CV
<b>Neg</b>	30/30 (100)	NA	NA	NA	NA	30/30 (100)	NA	NA	NA	NA	30/30 (100)	NA	NA	NA	NA	90/90 (100)	NA	NA	NA	NA
<b>HN</b>	30/30 (100)	NA	NA	NA	NA	30/30 (100)	NA	NA	NA	NA	30/30 (100)	NA	NA	NA	NA	90/90 (100)	NA	NA	NA	NA
<b>LP</b>	29/30 (96.7)	29.4	2.5	1.76	29.2	30/30 (100)	29.8	3.6	1.53	43.8	30/30 (100)	29.2	2.9	1.79	31.8	89/90 (98.9)	29.5	3.1	1.69	35.1
<b>MP</b>	30/30 (100)	27.2	2.0	3.15	9.7	30/30 (100)	27.9	2.4	2.79	14.3	30/30 (100)	26.8	2.0	3.19	7.86	90/90 (100)	27.3	2.7	3.04	12.1
<b>Total</b>	119/120 (99.2)	--	--	--	--	120/120 (100)	--	--	--	--	120/120 (100)	--	--	--	--	359/360 (99.7)	--	--	--	--

**Table 3. Reproducibility Study Results for Internal Process Control Stratified by Site**

	Site A					Site B					Site E					Total				
	Agmt (%)	Ct		Endpoint		Agmt (%)	Ct		Endpoint		Agmt (%)	Ct		Endpoint		Agmt (%)	Ct		Endpoint	
		Mean	%CV	Mean	%CV		Mean	%CV	Mean	%CV		Mean	%CV	Mean	%CV		Mean	%CV	Mean	%CV
<b>Neg</b>	30/30 (100)	29.0	1.7	2.93	13.3	30/30 (100)	29.0	1.9	2.89	12.4	30/30 (100)	29.1	1.8	2.75	12.2	90/90 (100)	29.0	1.8	2.86	12.8
<b>HN</b>	30/30 (100)	28.8	2.1	3.01	13.3	30/30 (100)	29.1	2.0	2.86	14.9	30/30 (100)	29.1	2.2	2.72	16.4	90/90 (100)	29.0	2.1	2.86	15.2
<b>LP</b>	30/30 (100)	28.9	1.6	2.95	10.9	30/30 (100)	28.8	1.6	2.93	10.4	30/30 (100)	29.1	1.4	2.62	10.0	90/90 (100)	28.9	1.6	2.84	11.6
<b>MP</b>	30/30 (100)	28.5	1.9	2.69	11.6	30/30 (100)	28.8	1.7	2.66	10.6	30/30 (100)	28.7	1.8	2.23	15.5	90/90 (100)	28.7	1.8	2.52	14.8
<b>Total</b>	120/120 (100)	--	--	--	--	120/120 (100)	--	--	--	--	120/120 (100)	--	--	--	--	360/360 (100)	--	--	--	--

Legend to **Tables 3 and 4:**

NA: Not Applicable; Neg: Negative; HN: High Negative; LP: Low Positive; MP: Moderate Positive; Agmt: Agreement with expected result; %CV: Percent Coefficient of Variation

b. *Linearity/assay reportable range:*

Not applicable

c. *Traceability, Stability, Expected values (controls, calibrators, or methods):*

**Controls**

The Liat™ Strep A Assay incorporates 3 controls:

- (1) Internal Process Control (IPC);
- (2) External Positive Control;
- (3) External Negative Control.

*Internal Process Control*

The IPC is designed to avoid reporting of false negative results due to excessive sample inhibition or system operation outside the normal range. The IPC is comprised of an inactivated bacterium that is pre-packaged in each Liat™ Tube. When conducting an assay, the sample is first mixed with the IPC which then goes through all the test process steps to monitor both the sample processing and PCR performance. The IPC DNA is detected in a separate optical channel from the *S. pyogenes* target nucleic acid. If the IPC Ct and fluorescence endpoint values do not fall within specified limits and *S. pyogenes* target DNA is not detected, the test is reported as invalid and the operator is instructed to repeat the test.

*External Positive and Negative Controls*

Positive and Negative Controls are provided separately from the assay reagents in the Liat™ Strep A Assay Quality Control Kit and must be tested each time a new lot of reagents is used. Should either the Positive or Negative Control fail to produce the expected result during the “Add Lot” process the system software precludes use of the Liat™ tube lot until the cause of the failure is resolved.

The Positive Control and consists of dried, inactivated *S. pyogenes* bacteria. The target level for the Positive Control is designed to be close to the LOD of the assay.

The Negative Control consists of Amies medium without *S. pyogenes* target organisms and is provided in unit dose quantity.

*Evaluation of Control Performance*

The ability of the Liat™ Strep A IPC and Positive Control to monitor the performance of the assay was verified under various simulated failure modes (**Table 3**). All of the replicates tested under each simulated failure mode were reported as invalid thereby demonstrating that the IPC and Positive Control are effective in monitoring substantial failures in reagent and process integrity.

**Table 3. Simulated failure modes used to evaluate assay controls**

	<b>Process Failure</b>	<b>Reagent Failure</b>
Sample Preparation	Failure to capture silica magnetic beads during nucleic acid extraction	Premature rupture of frangible seal between Lysis and Wash Buffers
PCR Amplification and Detection	Deviation in PCR annealing temperature	Premature rupture of frangible seal between Wash and Elution Buffers

The ability to monitor for contamination was evaluated by spiking Negative Controls with low levels of *S. pyogenes* DNA. In each case the system software indicated that the control failed.

During the Clinical Study to validate the performance of the Liat™ Strep A Assay, 15 Positive and 15 Negative Controls were tested across 6 lots of Liat™ Tubes and 11 Liat™ Analyzers. All (100%) of the controls produced the expected results as did 20/20 (100%) Negative Controls and 22/22 (100%) Positive Controls or known positive samples that were used to monitor assay performance during the course of the in-house analytical studies.

### **Sample Stability**

Sample stability was evaluated by testing *S. pyogenes* positive and negative throat swabs after storage at 2-8°C, 25°C or -20°C. *S. pyogenes* positive samples were spiked with organisms at 3X LOD. At each test point, all of the samples produced the expected result, supporting the Package Insert claims for sample stability (48 hours at 2-25°C).

*d. Analytical Sensitivity (Detection limit):*

The Limit of Detection (LOD) of the Liat™ Strep A Assay was determined by limiting dilution studies using titered stocks of 4 strains of *S. pyogenes* that were spiked into throat swab matrix. The LOD was determined as the lowest concentration that was detected  $\geq 95\%$  of the time (*i.e.*,  $\geq 19/20$  replicates tested positive) and ranged from 5-20 CFU/mL (**Table 4**).

**Table 4. Limits of Detection in Throat Swab Matrix**

<b><i>S. pyogenes</i> Strain</b>	<b>LOD (CFU/mL)</b>
ATCC BAA-946	5
ATCC 12370	10
ATCC BAA-1066	10
ATCC 700294	20

The analytical reactivity of the Liat™ Strep A Assay was evaluated with 5 strains of *S. pyogenes* in addition to those tested in the LOD study (**Table 5**). Titered stocks of

each strain were diluted in throat swab matrix and tested in triplicate. All 5 strains produced positive results at or below a concentration of 80 CFU/mL.

**Table 5.** Analytical reactivity of the Liat™ Strep A Assay

<i>S. pyogenes</i> Strain	Reactivity Titer CFU/mL
ATCC 700949	20
ATCC 21548	40
ATCC 10403	80
ATCC 700497	20
ATCC 700499	40

BLAST analysis also showed the target region for the Liat™ Strep A Assay to be conserved across different strains of *S. pyogenes*.

*e. Analytical specificity:*

**Cross-Reactivity Study**

The analytical specificity of the Liat™ Strep A Assay was evaluated by testing a panel of 72 potentially cross-reactive organisms and viruses (**Table 6**). Except where noted, bacteria and yeast were tested at 1.1-4.4 x 10<sup>6</sup> Colony Forming Units/mL and viruses were tested at 1.00-4.45 x 10<sup>5</sup> TCID<sub>50</sub>/mL. Each virus or species of organism was tested in triplicate in *S. pyogenes* negative throat swab matrix. In each case, all three assay replicates produced negative results indicating no evidence of cross-reaction with the Liat™ Strep A Assay.

BLAST analysis of the Liat™ Strep A primer, probe and amplicon sequences also showed no homology with other organisms or viruses that is likely to impact assay performance.

**Table 6. Organisms and Viruses Tested for Potential Cross-reaction or Interference with the Liat™ Strep A Assay**

<i>Arcanobacterium haemolyticum</i>	<i>Legionella micdadei</i>	<i>Streptococcus canis</i>
<i>Acinetobacter baumannii</i>	<i>Legionella pneumophila</i>	<i>Streptococcus constellatus</i>
<i>Bacillus cereus</i>	<i>Listeria monocytogenes</i>	<i>Streptococcus dysgalactiae</i>
<i>Bacteroides oralis</i>	<i>Moraxella catarrhalis</i> (2 strains)	<i>Streptococcus equi</i>
<i>Bordetella bronchiseptica</i>	<i>Moraxella lacunata</i>	<i>Streptococcus gallolyticus</i>
<i>Bordetella parapertussis</i>	<i>Mycoplasma pneumoniae</i> <sup>3</sup>	<i>Streptococcus intermedius</i>
<i>Bordetella pertussis</i>	<i>Neisseria gonorrhoeae</i> <sup>3</sup>	<i>Streptococcus mitis</i>
<i>Campylobacter rectus</i>	<i>Neisseria lactamica</i>	<i>Streptococcus mutans</i>
<i>Burkholderia cepacia</i>	<i>Neisseria meningitidis</i>	<i>Streptococcus oralis</i>
<i>Candida albicans</i>	<i>Neisseria mucosa</i>	<i>Streptococcus pneumoniae</i>
<i>Chlamydia pneumoniae</i> <sup>1</sup>	<i>Neisseria sicca</i>	<i>Streptococcus salivarius</i>

<i>Chlamydia trachomatis</i> <sup>2</sup>	<i>Neisseria subflava</i>	<i>Streptococcus sanguis</i>
<i>Corynebacterium diphtheriae</i>	<i>Proteus mirabilis</i>	<i>Treponema denticola</i> <sup>3</sup>
<i>Corynebacterium pseudodiphtheriticum</i>	<i>Proteus vulgaris</i>	<i>Veillonella parvula</i>
<i>Enterococcus faecalis</i>	<i>Pseudomonas aeruginosa</i>	<i>Yersinia enterocolitica</i>
<i>Enterococcus faecium</i>	<i>Pseudomonas fluorescens</i>	Adenovirus, Type 1
<i>Escherichia coli</i>	<i>Serratia marcescens</i>	Adenovirus, Type 7
<i>Haemophilus influenza</i>	<i>Staphylococcus aureus</i>	Cytomegalovirus
<i>Haemophilus parahaemolyticus</i>	<i>Staphylococcus epidermidis</i>	Epstein-Barr Virus <sup>4</sup>
<i>Haemophilus parainfluenzae</i>	<i>Staphylococcus haemolyticus</i>	Hepatitis B Virus <sup>4</sup>
<i>Klebsiella pneumonia</i>	<i>Stenotrophomonas maltophilia</i>	Herpes Simplex Virus 1
<i>Lactobacillus acidophilus</i>	<i>Streptococcus agalactiae</i>	Human Papilloma Virus, Type 11 <sup>4</sup>
<i>Lactococcus lactis</i>	<i>Streptococcus anginosus</i>	Human Papilloma Virus, Type 6 <sup>4</sup>
<i>Legionella jordanis</i>	<i>Streptococcus bovis</i>	

Concentration tested (per mL):<sup>1</sup> 1.40 x 10<sup>5</sup> TCID<sub>50</sub>; <sup>2</sup> 1.25 x 10<sup>6</sup> Elementary Bodies; <sup>3</sup> 1.25-1.63 x 10<sup>6</sup> genomic copies; <sup>4</sup> 4.45 x 10<sup>4</sup> TCID<sub>50</sub>; <sup>5</sup> 2.15-5.00 x 10<sup>5</sup> genomic copies

### Cross-Contamination Study

In order to evaluate the risk of false-positive results due to contamination between samples (run-to-run) and between analyzers (instrument-to-instrument), 20 high positive (3.13 x 10<sup>6</sup> CFU/mL) and 20 *S. pyogenes* negative samples were tested on each of two Liat™ Analyzers (i.e., 40 high positive and 40 negative samples in total). All 40 high positive and all 40 negative samples produced the expected results, indicating that the risk of contamination between runs and instruments is acceptably low.

#### f. Assay cut-off:

The Liat™ Strep A Assay result algorithm analyses the characteristics of the fluorescence amplification curves for the *S. pyogenes* target and Internal Process Control (IPC) to disposition samples as “Strep A Detected,” “Strep A Not Detected,” “Indeterminate” or “Invalid.”

The algorithm employs cut-offs for both the cycle threshold (Ct) value and endpoint amplitude in addition to other parameters. The cut-offs were determined through analysis of a combination of negative clinical samples and samples that were spiked with different strains at the LOD target level.

In order for a sample to be called negative for *S. pyogenes* (“Strep A Not Detected”), the IPC must be detected. If the IPC is not detected, the result is reported as “Assay Invalid. Repeat Assay.” In cases in which *S. pyogenes* is detected but the curve shape is abnormal, the result is reported as “Indeterminate.”

As indicated in the Package Insert, if the test result is “Indeterminate” or “Invalid”, the assay should be repeated with the same patient specimen or, if possible, a new specimen from the same patient. Specimens that have repeat “Indeterminate” or “Invalid” results should be sent for confirmatory testing by another method.

g. Assay Interference:

**Interfering Substances**

The Liat™ Strep A was evaluated with 28 substances that may be encountered in throat swabs (**Table 7**). Testing was performed in triplicate with medically and/or physiologically relevant concentrations of each potential interfering substance in throat swab matrix and in the presence and absence of *S. pyogenes* ATCC 700294 at 60 CFU/mL (*i.e.*, 3X LOD). In all cases the expected results were obtained, although with bovine mucin at the concentration tested all three assay replicates exhibited delayed Ct values and lower endpoint fluorescence for the *S. pyogenes* target. This is reflected in the Package Insert in a footnote to the list of potential interfering substances that were evaluated.

**Table 7. Interfering Substance Panel**

<b>Substance</b>	<b>Concentration</b>
Acetaminophen (Tylenol)	100 µg/mL
Adult Robitussin Peak Cold, Maximum strength, Cough+Chest	5% v/v
Adult Robitussin Peak Cold, Nighttime, Multi-symptom cold	5% v/v
Amoxicillin	25 µg/mL
Blood (human)	5% v/v
Brompheniramine Maleate	60 ng/mL
Cepacol Sore Throat	5 mg/mL
Cepacol Ultra Sore Throat Spray	5% v/v
Children's Dimetapp Cold & Cough	5% v/v
Children's Robitussin Cough & Cold	5% v/v
Children's Dimetapp Nighttime Cold & Congestion	5% v/v
Chloraseptic Max	5% v/v
Chlorpheniramine Maleate	25 ng/mL
Cool Mint Listerine, antiseptic	5% v/v
Crest Pro-Health	5% v/v
Dextromethorphan HBr	20 ng/mL
Diphenhydramine HCl	350 ng/mL
Doxylamine Succinate	300 ng/mL
Erythromycin	15 µg/mL
Guaifenesin (Guaiacol glyceryl)	5 mg/mL
Halls Mentho-lyptus Cherry	5 mg/mL
Halls Mentho-lyptus Sugar Free	5 mg/mL
Ibuprofen (Advil)	25 µg/mL
Mucin: bovine submaxillary gland, type I-S	25 mg/mL <sup>1</sup>
Penicillin G	1.2 mg/mL
Sucrets Complete	5 mg/mL
Tussin Adult Chest Congestion	5% v/v
Tylenol Cold Sore Throat	5% v/v

<sup>1</sup> Delayed Ct and reduced endpoint fluorescence for *S. pyogenes*

## **Interfering Microorganisms**

The organisms and viruses listed in **Table 6**, above, were evaluated for the potential to interfere with the Liat™ Strep A Assay by testing them at the concentrations shown in the table in the presence of a low level of *S. pyogenes* ATCC 700294 (*i.e.*, 60 CFU/mL = 3X LOD). Each organism or virus was tested in triplicate. All assay replicates produced positive results and there was therefore no evidence of microbial interference with the Liat™ Strep A Assay.

### 2. Comparison studies:

#### *a. Method comparison with predicate device:*

Not applicable

#### *b. Matrix comparison:*

Not applicable

### 3. Clinical studies:

#### *a. Clinical Sensitivity:*

The performance of the Liat™ Strep A Assay was evaluated in a study that was conducted at six clinical sites in the U.S. between December 2013 and April 2014. Specimens were collected from subjects who were 3 years of age or older and who presented with symptoms that are characteristic of pharyngitis (*i.e.*, sore throat plus at least one of the following: pharyngeal redness, pharyngeal or tonsillar exudate, tonsillar swelling, tender cervical lymphadenopathy or fever). Two or three throat swabs were collected from each subject. One swab for use with the Liat™ Strep A Assay was placed in 1mL Amies transport medium. The additional swabs were used for standard of care diagnosis. The order of swab collected was at the discretion of the clinical sites.

The swabs in Amies medium were tested with the Liat™ Strep A Assay at the clinical sites on the day of collection. Residual specimen was then shipped to a central laboratory for reference culture where an aliquot of the transport medium was plated on Trypticase Soy Agar with 5% Sheep Blood and incubated at 35-37°C in an atmosphere of 5-7% CO<sub>2</sub>. Isolated colonies that exhibited β-hemolysis were typed by latex agglutination (Remel Strepex®). Culture plates that did not exhibit β-hemolytic colonies after 48 hours were recorded as negative.

Specimens that produced discordant results between the reference culture and the Liat™ Strep A Assay were subject to additional characterization by an alternative PCR method, followed by bi-directional sequencing.

Throat swabs were collected from a total of 799 subjects during the course of the study. Two hundred and twenty-nine (229) specimens were excluded from the analysis of

performance due to failure to comply with inclusion criteria (2), delayed reference culture (16), absence of a Liat™ Strep A result (1), use of the incorrect swab type (130), incorrect Liat™ Strep assay procedure (20) or protocol deviation (60). A total of 570 specimens was therefore included in the analysis of performance (**Table 8**).

**Table 8. Liat™ Strep A Assay Clinical Performance vs Reference Culture**

		Reference Culture		
		Positive	Negative	Total
<b>Liat™ Strep A</b>	<b>Positive</b>	170	23 <sup>1</sup>	<b>193</b>
	<b>Negative</b>	3 <sup>2</sup>	374	<b>377</b>
	<b>Total</b>	<b>173</b>	<b>397</b>	<b>570</b>
Sensitivity		170/173 = 98.3% (95% CI: 95.0-99.4%)		
Specificity		374/397 = 94.2% (95% CI: 91.5-96.1%)		
Positive Predictive Value		170/193 = 88.1% (95% CI: 82.8-91.9%)		
Negative Predictive Value		374/377 = 99.2% (95% CI: 97.7-99.7%)		

<sup>1</sup> 23/23 positive by alternative PCR and bi-directional sequencing

<sup>2</sup> 3/3 positive by alternative PCR and bi-directional sequencing and upon repeat testing with Liat™ Strep A

Of the 570 unique specimens included in the analysis, 7 initially produced invalid Liat™ Strep A test results (invalid rate = 7/577 = 1.2%). All 7 produced valid results on retest and those results are included in the analysis of performance.

The performance of the LiatStrep A Assay stratified by clinical site is shown in **Table 9**.

**Table 9. Liat™ Strep A Assay performance stratified by site**

Site	Samples (%)	Culture Positive (% Prevalence)	Percent (95% CI)			
			Sensitivity	Specificity	PPV	NPV
<b>A</b>	246 (43.2)	75 (30.5)	96.0 (88.9-98.6)	91.2 (86.0-94.6)	82.8 (73.5-89.3)	98.1 (94.6-99.4)
<b>B</b>	85 (14.9)	32 (37.6)	100 (89.3-100)	94.3 (84.6-98.1)	91.4 (77.6-97.0)	100 (92.9-100)
<b>C</b>	52 (9.1)	11 (21.2)	100 (74.1-100)	100 (91.4-100)	100 (74.1-100)	100 (91.4-100)
<b>E <sup>1</sup></b>	4 (0.7)	0 (0)	NA	100 (51.0-100)	NA	100 (51.0-100)
<b>I</b>	40 (7.0)	11 (27.5)	100 (74.1-100)	93.1 (78.0-98.1)	84.6 (57.8-95.7)	100 (87.5-100)
<b>L</b>	143 (25.1)	44 (30.8)	100 (92.0-100)	97.0 (91.5-99.0)	93.6 (82.8-97.8)	100 (96.2-100)
<b>Total</b>	<b>570</b> <b>(100)</b>	<b>173</b> <b>(30.4)</b>	<b>98.3</b> <b>(95.0-99.4)</b>	<b>94.2</b> <b>(91.5-96.1)</b>	<b>88.1</b> <b>(82.8-91.9)</b>	<b>99.2</b> <b>(97.7-99.7)</b>

<sup>1</sup> Enrollment at Site E was low but the samples were collected according to the study protocol and are therefore included in the analysis.

*b. Clinical specificity:*

Refer to Section 3a, above.

4. Clinical cut-off:

Not applicable

5. Expected values/Reference range:

The clinical study included 570 specimens from six U.S. sites that were collected between December 2013 and April 2014. The overall prevalence of *S. pyogenes* as determined by culture was 30.4% (173/570), and as determined by the Liat™ Strep A Assay the overall prevalence was 33.9% (193/570). The prevalence by age and sex of the subjects is shown in **Table 10**.

**Table 10. Prevalence of *S. pyogenes* stratified by Age and Gender**

Age/Gender	Liat™ Strep A		Total	Prevalence (%)
	Positive	Negative		
≤5 years	59	82	141	41.8
6-21 years	130	271	401	32.4
22-59 years	4	21	25	16.0
≥60 years	0	3	3	0
Male	100	168	268	37.3
Female	93	209	302	30.8

**N. Proposed Labeling:**

The labeling is sufficient and it satisfies the requirements of 21 CFR Part 809.10.

**O. Conclusion:**

The submitted information in this premarket notification is complete and supports a substantial equivalence decision.