

**510(k) SUBSTANTIAL EQUIVALENCE DETERMINATION  
DECISION SUMMARY**

**A. 510(k) Number:**

K141520

**B. Purpose for Submission:**

To obtain a Substantial Equivalence Determination for a new 510(k) application for the Alere™ i Influenza A&B assay on the Alere™ i Instrument.

**C. Measurand:**

Target RNA sequences for the highly conserved regions of the PB2 gene (encodes an RNA polymerase – polymerase basic protein 2) segment of Influenza A virus and the PA gene (encodes another RNA polymerase – polymerase acidic protein) segment of Influenza B virus.

**D. Type of Test:**

Qualitative isothermal nucleic acid amplification assay for the amplification and detection of specific Influenza A and B RNA gene targets.

**E. Applicant:**

Alere Scarborough, Inc.

**F. Proprietary and Established Names:**

Alere™ i Influenza A&B  
Alere™ i Instrument

**G. Regulatory Information:**

1. Regulation section:

21 CFR section 866.3980, Respiratory Viral Panel Multiplex Nucleic Acid Assay

2. Classification:

Class II

3. Product code(s):

OCC - Respiratory virus panel nucleic acid assay system  
OZE - Influenza A and influenza B multiplex nucleic acid assay  
OOI - Real time nucleic acid amplification system

4. Panel:

Microbiology (83)

## H. Intended Use:

### 1. Intended use(s):

The Alere™ i Influenza A & B assay performed on the Alere™ i Instrument is a rapid molecular *in vitro* diagnostic test utilizing an isothermal nucleic acid amplification technology for the qualitative detection and discrimination of influenza A and B viral RNA in nasal swabs from patients with signs and symptoms of respiratory infection. It is intended for use as an aid in the differential diagnosis of influenza A and B viral infections in humans in conjunction with clinical and epidemiological risk factors. The assay is not intended to detect the presence of influenza C virus.

Negative results do not preclude influenza virus infection and should not be used as the sole basis for diagnosis, treatment or other patient management decisions.

Performance characteristics for influenza A were established during the 2012-2013 influenza season when influenza A/H3 and A/H1N1 pandemic were the predominant influenza A viruses in circulation. When other influenza A viruses are emerging, performance characteristics may vary.

If infection with a novel influenza A virus is suspected based on current clinical and epidemiological screening criteria recommended by public health authorities, specimens should be collected with appropriate infection control precautions for novel virulent Influenza viruses and sent to state or local health department for testing. Viral culture should not be attempted in these cases unless a BSL 3+ facility is available to receive and culture specimens.

### 2. Indication(s) for use:

Same as Intended Use

### 3. Special conditions for use statement(s):

For prescription use only

### 4. Special instrument requirements:

Alere™ i Instrument

## I. Device Description:

### Overview

The Alere™ i Influenza A & B assay system utilizes an isothermal nucleic acid amplification technology and is comprised of:

- Sample Receiver – single use, disposable containing the elution buffer
- Test Base – single use, disposable comprising two sealed reaction tubes, each containing a lyophilized pellet
- Transfer Cartridge – single use, disposable for transfer of the eluted sample to the Test Base
- Alere™ i Instrument – repeat use reader, and
- Alere™ i Influenza A & B positive and negative external controls – for quality control purposes

- Nasal Swabs – sterile swabs for use with the Alere™ i Influenza A & B Test.

The first reaction tube (Tube 1) in the Test Base contains the reagents required for amplification of the influenza A target viral nucleic acid, and the second reaction tube (Tube 2) in the Test Base contains reagents required for amplification of the influenza B target viral nucleic acid and an internal control synthetic RNA sequence. Alere™ i Influenza A & B utilizes pairs of templates (similar to primers in a PCR reaction) for the specific isothermal amplification of target RNA sequences from influenza A (highly conserved region of the PB2 gene), influenza B (highly conserved region of the PA gene), and the internal control, and fluorescently labeled molecular beacons designed to specifically identify the amplified products. Alere™ i Influenza A & B is performed within the confinement of the Test Base, and no other part of the Alere™ i Instrument has contact with the sample during the amplification process.

To perform the assay, the Sample Receiver and the Test Base are inserted into the Alere™ i Instrument and the elution buffer is automatically heated by the instrument. The sample is added to the Sample Receiver and transferred via the Transfer Cartridge to the Test Base, resuspending the lyophilized pellet contained within the Test Base and initiating target amplification. Heating, mixing and detection by fluorescence are provided by the instrument, with results automatically reported.

Results (positive, negative, or invalid) are displayed by the Alere™ i Instrument separately for influenza A and influenza B. Results are also stored in an on-board archive and are assigned to a sample ID that has been entered into the Alere™ i Instrument by the operator, and the date/time the test was performed. Data can be retrieved and downloaded by the operator at any time after testing. An external Alere™ Universal Printer can be attached via USB to the Alere™ i Instrument to print test results.

## **Alere™ i Influenza A & B System Components**

### **Sample Receiver**

The Sample Receiver is a foil sealed, single-use, disposable plastic component containing 2.5ml of an elution/lysis buffer. Rapid elution of influenza virus from swabs and virus lysis occur through the action of the elution/lysis buffer. The elution/lysis buffer contains salts at the appropriate concentrations that are required for the test reaction, which upon transfer of the sample using the Transfer Cartridge to the Test Base, is used to reconstitute the lyophilized pellets. The elution/lysis buffer also contains a chelating agent to remove any heavy metals found in the sample, as they are potential inhibitors of the test reaction.

### **Transfer Cartridge**

The Transfer Cartridge is a single-use, disposable plastic component for transfer of the eluted sample to the Test Base. The aspiration of the sample is achieved by placing the Transfer Cartridge over the Sample Receiver and pressing downwards until a click is heard and the orange indicator on the top of the Transfer Device is up. The dispense of the sample is achieved by placing the Transfer Cartridge into the Test Base and pressing on the Transfer Cartridge until the orange indicator on the top of the Transfer Cartridge descends. When the indicator descends, so does the plunger, thus driving the sample out of the two nozzles and into the reaction tubes in the Test Base and mixing with the lyophilized reagent in the reaction tubes.

### **Test Base**

The Test Base is a single-use, disposable plastic component comprising two sealed reaction tubes, the first reaction tube (Tube 1) is for the influenza A reaction and the second tube (Tube 2) is for the influenza B and the Internal Control reactions. Each reaction tube contains a lyophilized pellet which provides all of the necessary reagents to drive the respective reaction(s) except for salts. Tube 2 also contains an internal

control synthetic RNA sequence. The salts that are required for the test reactions are provided at their appropriate concentrations within the elution/lysis buffer contained in the Sample Receiver. Upon transfer of the sample from the Sample Receiver to the Test Base using the Transfer Cartridge, the sample in elution/lysis buffer is used to reconstitute the lyophilized pellets. The Alere™ i Influenza A & B reaction is performed within the Test Base. As soon as the sample is transferred, it comes into contact with the lyophilized pellets, reconstitutes the lyophilized pellets and initiates the assay reaction. Each Test Base contains a barcode with assay specific information (assay ID, expiration date and lot number). The barcode is read automatically from the Test Base upon insertion into the Alere™ i Instrument, providing the information necessary to run the assay.

## **Alere™ i Instrument**

The Alere™ i Instrument is a portable bench-top unit designed to perform in laboratory and point-of care environments. Alere™ i Influenza A & B components are color-coded to match the corresponding holder on the instrument and designed to only fit in the corresponding holder in one direction. The Alere™ i Instrument uses image analysis to automatically detect the presence of the Test Base and to check the test type and expiry date on the barcode present on the Test Base. This triggers the test sequence. Heating, mixing, and detection by fluorescence are provided by the instrument. The Alere™ i Instrument is equipped with two optical fluorescence modules, FAM™ and ROX™. Each optical module consists of a LED excitation source and photodiode receiver configured to the fluorescent channel. Fluorescence signals are optically filtered by a confocal lens and mirror arrangement and processed by software running on an on-board computer. The results are calculated automatically and presented on the instrument graphical user interface and available for export in JSON format or printing.

## **Alere™ i Influenza A & B External Controls**

Alere™ i Influenza A & B External Controls are designed for use with Alere™ i Influenza A & B. The Positive Control swab is coated with inactivated influenza A virus and inactivated influenza B virus dried onto a swab. The Negative Control swab is coated with inactivated *Streptococcus* Group C dried onto a swab.

## **Nasal Swabs**

Sterile swabs for use with the Alere™ i Influenza A & B test are provided as a part of the Alere™ i Influenza A & B System.

## **Alere™ i Influenza A & B Workflow**

An overview of the Alere™ i Influenza A & B workflow is summarized below:

1. User ID and Patient ID are entered via the touch screen prior to initiating the testing sequence. Patient information can be entered by two means:
  - by using the Instrument touch screen
  - by using a barcode reader
2. The Alere™ i Influenza A & B Test Base barcode contains assay information (test type, expiry date, and lot number). Upon proper insertion into the Test Base holder of the instrument, using image analysis the instrument automatically detects the presence of the Test Base and checks the test type, expiry date, and lot number encoded in the barcode. Insertion of an expired Test Base will generate an error message on the instrument screen and the user is not able to proceed with testing. Upon confirmation of the Test Base insertion the user is instructed to insert the Sample Receiver.

3. The Alere™ i Instrument detects the insertion of the Sample Receiver into the Sample Receiver holder and automatically progresses to the warm-up state. The instrument software implements temperature control of both the Test Base and the Sample Receiver holders (heater blocks). Each of the heater blocks is managed independently with high precision temperature monitoring and control.
4. Upon completion of the warm-up state, the Alere™ i Instrument instructs the user to remove the foil seal from the Sample Receiver, insert (direct nasal swab specimen) sample to the elution buffer contained in the Sample Receiver, and mix as described on the screen.
5. The Alere™ i Instrument then instructs the user to transfer the sample in the Sample Receiver to the Test Base using the Transfer Cartridge.
6. The Alere™ i Instrument uses image analysis to confirm that the Transfer Cartridge is present in the Test Base, instructs the user to close the lid, and automatically proceeds to the measurement state upon detection of test bay lid closure. The lid operates using a magnetic switch system that detects lid closure and an electrical signal triggers testing commencement. Heating, rotation and detection by fluorescence is provided by the instrument.
7. Upon completion of the testing sequence, results are displayed on the screen and are presented to the user as positive, negative, or invalid separately for influenza A and influenza B. Upon acknowledgement of the test results, the user is instructed to dispose of the used test components.

The turnaround time for testing a sample is approximately 15 minutes.

## **Quality Control**

### **Internal Control (Procedural Control)**

Alere™ i Influenza A & B contains an internal control (IC) that was designed to control for sample inhibition, amplification, and assay reagent function for both influenza A and influenza B. The IC is an RNA oligonucleotide contains 5' and 3' ends that are complementary to the influenza B target template set's recognition regions, but with a spacer region (the region between the two recognition regions) that differs from the influenza B target's spacer region. The IC is present in the second reaction tube (Tube 2) of the Test Base along with all of the necessary reagents to drive the influenza B and the IC reactions except for salts. Detection of the IC occurs via a molecular beacon that specifically detects the amplified product generated from the IC RNA oligonucleotide. The two molecular beacons present in the second reaction tube of the Test Base are labeled with fluorophores with distinct excitation/emission profiles (Influenza B – FAM and IC – ROX), enabling detection of both beacons simultaneously using a two-color detector.

In positive samples where influenza target amplification is strong, the IC result is ignored. The influenza B target amplification in Tube 2 of the Test Base also serves as the “Control” to confirm that the clinical sample was not inhibitory and that influenza A reagent performance in Tube 1 was robust. At a low frequency, clinical samples can contain inhibitors that may generate invalid results.

### **Procedural Controls**

In addition to the test specific Internal Control described previously, the Alere™ i instrument also employs additional “Procedural Controls” to ensure that a test is performed appropriately. Each Test Base has a barcode which identifies the test type – the instrument automatically detects the Test Base and asks the user to confirm the test type. Each of the additional steps (sample heating, sample transfer, and closing the lid, etc.) are monitored by the instrument and if the user takes too long to complete a certain step, the instrument will time out, generate an error and a result is not reported.

“Procedural Control Valid” displayed on the instrument screen indicates that the assay reagents maintained their functional integrity and the sample did not significantly inhibit assay performance.

### External Positive and Negative Controls

Good laboratory practice suggests the use of positive and negative controls to ensure that test reagents are working and that the test is correctly performed. Aleré™ i Influenza A & B kits contain Positive and Negative Control Swabs. These swabs monitor the entire assay. Test these swabs once with each new shipment received and once for each new operator. Further controls may be tested in order to conform to local, state and/or federal regulations, accrediting groups, or laboratory’s standard Quality Control procedures.

### Results Interpretation

The Aleré™ i Instrument software automatically interprets test results. Aleré™ i Influenza A & B assay operates in the Test Base with two separate reaction tubes, Tube 1 and Tube 2, for influenza A and influenza B plus IC respectively. The assay uses the following three fluorescent channels: Flu A (FAM) in Tube 1, Flu B (FAM) in Tube 2, and IC (ROX) in Tube 2. The ROX channel in Tube 1 is not used. Each channel of each reaction tube is first individually assessed using a defined channel decision algorithm and assigned a channel result as “Asserted”, “Not Asserted” or “Indeterminate”. The individual channel results are then combined to arrive at the influenza A and influenza B results. All nine possible assay result combinations and the respective underlying assay decision logic rules are presented in Table 1 below:

**Table 1: All Possible Assay Result Combinations and Respective Underlying Assay Decision Logic Rules**

Assay Result Reported	Underlying Assay Decision Logic Rules		
	Flu A Channel (FAM) in Tube 1	Flu B Channel (FAM) in Tube 2	IC Channel (ROX) in Tube 2
Flu A Positive/ Flu B Negative	Asserted	Not Asserted	Asserted
Flu A Positive/ Flu B Invalid	Asserted	Not Asserted	Not Asserted
	Asserted	Not Asserted	Indeterminate
	Asserted	Indeterminate	Asserted
	Asserted	Indeterminate	Not Asserted
	Asserted	Indeterminate	Indeterminate
Flu A Positive/ Flu B Positive	Asserted	Asserted	Asserted
	Asserted	Asserted	Not Asserted
	Asserted	Asserted	Indeterminate
Flu A Negative/ Flu B Positive	Not Asserted	Asserted	Asserted
	Not Asserted	Asserted	Not Asserted
	Not Asserted	Asserted	Indeterminate
Flu A Negative/ Flu B Invalid	Not Asserted	Indeterminate	Asserted
Flu A Negative/ Flu B Negative	Not Asserted	Not Asserted	Asserted
Flu A Invalid/Flu B Positive	Indeterminate	Asserted	Asserted
	Indeterminate	Asserted	Not Asserted
	Indeterminate	Asserted	Indeterminate
Flu A Invalid/Flu B Negative	Indeterminate	Not Asserted	Asserted
Flu A Invalid/Flu B Invalid	Not Asserted	Not Asserted	Not Asserted
	Not Asserted	Not Asserted	Indeterminate
	Not Asserted	Indeterminate	Not Asserted
	Not Asserted	Indeterminate	Indeterminate
	Indeterminate	Not Asserted	Not Asserted
	Indeterminate	Not Asserted	Indeterminate
	Indeterminate	Indeterminate	Asserted
	Indeterminate	Indeterminate	Not Asserted
Indeterminate	Indeterminate	Indeterminate	

All possible assay result combinations and the respective interpretations and follow-up actions are presented

in Table 2 below:

**Table 2: All Possible Assay Result Combinations and Respective Interpretations and Follow-up Actions**

Assay Result Reported	Interpretation of Results and Follow-up Actions
Flu A Positive/ Flu B Negative	Flu A Viral RNA Detected; Flu B Viral RNA Not Detected. This result does not rule out co-infections with other pathogens or identify any specific influenza A virus subtype.
Flu A Positive/ Flu B Invalid	Flu A Viral RNA Detected; The presence or absence of Flu B Viral RNA can not be determined. This result does not rule out co-infections with other pathogens or identify any specific influenza A virus subtype.
Flu A Positive/ Flu B Positive	Flu A Viral RNA Detected; Flu B Viral RNA Detected. Dual infections of Flu A and Flu B are rare. Repeat testing using new test components. Contact Technical Support during normal business hours if multiple samples provide this result. This result does not rule out co-infections with other pathogens or identify any specific influenza A virus subtype or specific influenza B virus lineage.
Flu A Negative/ Flu B Positive	Flu A Viral RNA Not Detected; Flu B Viral RNA Detected. This result does not rule out co-infections with other pathogens or identify any specific influenza B virus lineage.
Flu A Negative/ Flu B Invalid	Flu A Viral RNA Not Detected; The presence or absence of Flu B Viral RNA can not be determined. Infection due to Flu B can not be ruled out. Repeat testing of the sample using new test components. If repeated Flu B Invalid results are obtained, results should be confirmed by another method prior to reporting the results.
Flu A Negative/ Flu B Negative	Flu A Viral RNA Not Detected; Flu B Viral RNA Not Detected.
Flu A Invalid/Flu B Positive	The presence or absence of Flu A Viral RNA can not be determined. Flu B Viral RNA Detected; This result does not rule out co-infections with other pathogens or identify any specific influenza B virus lineage.
Flu A Invalid/Flu B Negative	The presence or absence of Flu A Viral RNA can not be determined. Flu B Viral RNA Not Detected; Infection due to Flu A can not be ruled out. Repeat testing of the sample using new test components. If repeated Flu A Invalid results are obtained, results should be confirmed by another method prior to reporting the results.
Flu A Invalid/Flu B Invalid	The presence or absence of Flu A and Flu B Viral RNAs can not be determined. Repeat testing of the sample using new test components. If repeated Flu A and Flu B Invalid results are obtained, results should be confirmed by another method prior to reporting the results.

**J. Substantial Equivalence Information:**

1. Predicate device name(s):

IQuum Liat™ Influenza A/B Assay

2. Predicate 510(k) number:

K111387

3. Comparison with predicates:

## Similarities and Differences

Item	Device	Predicate
	Alere™ i Influenza A & B (K141520)	IQuum Liat™ Influenza A/B Assay (K111387)
<b>Assay Analyte</b>	Influenza A and Influenza B	Same
<b>Intended Use</b>	<p>The Alere™ i Influenza A &amp; B assay performed on the Alere™ i Instrument is a rapid molecular <i>in vitro</i> diagnostic test utilizing an isothermal nucleic acid amplification technology for the qualitative detection and discrimination of influenza A and B viral RNA in nasal swabs from patients with signs and symptoms of respiratory infection. It is intended for use as an aid in the differential diagnosis of influenza A and B viral infections in humans in conjunction with clinical and epidemiological risk factors. The assay is not intended to detect the presence of influenza C virus</p> <p>Negative results do not preclude influenza virus infection and should not be used as the sole basis for diagnosis, treatment or other patient management decisions.</p> <p>Performance characteristics for influenza A were established during the 2012-2013 influenza season when influenza A/H3 and A/H1N1 pandemic were the predominant influenza A viruses in circulation. When other influenza A viruses are emerging, performance characteristics may vary.</p> <p>If infection with a novel influenza A virus is suspected based on current clinical and epidemiological screening criteria recommended by public health authorities, specimens should be collected with appropriate infection control precautions for novel virulent Influenza viruses and sent to state or local health department for testing. Viral culture should not be attempted in these cases unless a BSL 3+ facility is available to receive and culture specimens.</p>	<p>The IQuum Liat™ Influenza A/B Assay performed on the Liat™ Analyzer is an automated multiplex real-time RT-PCR assay for the rapid <i>in vitro</i> qualitative detection and discrimination of influenza A virus and influenza B virus RNA in nasopharyngeal swab specimens from patients with signs and symptoms of respiratory infection in conjunction with clinical and epidemiological risk factors. The test is intended for use as an aid in the differential diagnosis of influenza A and influenza B in humans and is not intended to detect influenza C.</p> <p>Negative results do not preclude influenza virus infection and should not be used as the sole basis for treatment or other patient management decisions. Conversely, positive results do not rule-out bacterial infection or co-infection with other viruses. The agent detected may not be the definite cause of disease.</p> <p>Performance characteristics for influenza A were established when influenza A/H1 and A/H3 were the predominant influenza A viruses in circulation. When other influenza A viruses are emerging, performance characteristics may vary.</p> <p>If infection with a novel influenza A virus is suspected based on current clinical and epidemiological screening criteria recommended by public health authorities, specimens should be collected with appropriate infection control precautions for novel virulent Influenza viruses and sent to state or local health department for testing. Viral culture should not be attempted in these cases unless a BSL 3+ facility is available to receive and culture specimens.</p>
<b>Intended Environment for Use</b>	Professional use, in a medical laboratory or point-of-care	Same
<b>Instrumentation</b>	Alere™ i Instrument	Liat™ Analyzer
<b>Self-contained System</b>	Integrated PC, Software, and Touch Screen Display	Same
<b>Automated Assay</b>	Yes. Sample preparation, amplification, detection and result interpretation	Same
<b>Sample Type</b>	Direct Nasal Swab	Nasopharyngeal Swab in UTM
<b>Influenza A Viral Target</b>	PB2 gene segment	Matrix Gene
<b>Influenza B Viral Target</b>	PA gene segment	Non-Structural Protein (NSP) Gene
<b>Amplification and Detection Technology</b>	Isothermal nucleic acid amplification and detection of specific amplification products using molecular beacon probes	Real-time RT-PCR using TaqMan probes
<b>Internal Control</b>	Yes	Same

## **K. Standard/Guidance Document Referenced (if applicable):**

None


## **L. Test Principle:**

Individual reactions for influenza A and B and IC, provided as lyophilized pellets contained within the Test Base, specifically amplify unique target regions within the RNA genomes of these two viruses and the IC. Alere™ i Influenza A & B employs fluorescently-labeled molecular beacon probes for real-time detection of amplified products. Alere™ i Influenza A & B amplification technology, Nicking and Extension Amplification Reaction (NEAR) is isothermal and does not require a heat denaturation step to generate a single strand target for amplification. Instead a thermostable, strand-displacing DNA polymerase, a thermostable nicking endonuclease, and two oligonucleotides (typically known as primers, but are referred to as templates in this technology) are utilized. Reverse transcriptase is also present in the reaction. The products of Alere™ i Influenza A & B are two complementary oligonucleotides 30 – 40 nucleotides in length. A product of this size is of sufficient length to be highly unique amongst genomes. The exact length of the products depends on the specific placement of the templates relative to one another and relative to their complementary sequences within the target genome. Alere™ i Influenza A & B assay specificity is achieved through two levels of sequence-specific interrogation, template annealing to complementary target sequence at an elevated temperature to drive target-specific amplification, and amplification product detection as the molecular beacon anneals to Alere™ i Influenza A & B amplification product in a sequence-specific fashion.

### **Alere™ i Influenza A & B Isothermal Amplification Mechanism**


Alere™ i Influenza A & B isothermal amplification, NEAR, is driven by a thermostable strand displacing DNA polymerase, thermostable nicking endonuclease, two oligonucleotide templates, plus a reverse transcriptase. The Alere™ i Influenza A & B oligonucleotide template (primer) consists of three distinct regions. The 5'-end of the template is known as the stabilizing region (SR) which is not involved in target-specific complementary sequence-dependent annealing, but rather in stabilizing the double stranded products generated during amplification. The central region of each template maintains the nicking endonuclease binding site (NEBS) and nicking endonuclease cut site (NECS), which also are not complementary to the target sequence, but provide the sequence required for single strand nicking activity by the nicking endonuclease. At the 3'-end of each template is the recognition region (RR), which provides target-specific, sequence-dependent annealing. Typically, the RR is 11 – 12 nucleotides in length on each template. As shown in Figure 1 below, step 1a, the recognition region (RR) of template T2 (labeled as Y), anneals to a complimentary sequence present in the single strand DNA target (labeled as Y'). In step 2a, a DNA polymerase extends T2 from its 3'-OH group generating sequence X which is complementary to target sequence X'. In this step only the 3' end of T2 anneals to a target sequence in a complementary fashion. The SR, NEBS and NECS regions of the template are not complementary to the target sequence. In step 3a, a second T2 template displaces the original, partially annealed T2. DNA polymerase extension from this second T2 displaces the first synthesized strand. T1 recognizes, anneals to and extends from sequence X of the displaced strand (step 4a). This leaves a product containing a double strand NEBS and NECS (step 5a). A nicking endonuclease will then bind to its recognition site and nick the structure at its cut site. This generates an accessible 3'-OH that can be extended by the DNA polymerase, resulting in the release of a single strand DNA (containing regions X and Y) and generation of a double strand DNA product containing NEBS and NECS on both ends (step 6a), known as the Amplification Duplex. Figure 1 depicts the linear phase of the mechanism, essentially the generation of the NEAR Amplification Duplex.

**Figure 1 – Alere™ i Influenza A & B Mechanism – Linear Phase**



In Figure 2 below, the geometric amplification phase is depicted. The NEAR Amplification Duplex is a double strand DNA structure containing two complete nicking endonuclease binding and cut sites (1b). As shown in 2b, nicking of the top and bottom strands occurs, providing the DNA polymerase two 3'-OH groups for initiating strand displacement extension. This results in the release of X' Y' (specific product 1) and X Y (specific product 2) products, and the generation of two double strand DNA moieties containing single nicking endonuclease binding and cut sites (3b). The two moieties undergo reiterative nicking and extension, generating additional specific products 1 and 2 (3b). Specific products 1 and 2 are substrates for T2 and T1 (4b). Extension of these complexes (4b) generates additional double strand DNA moieties containing single nicking endonuclease binding and cut sites (3b), with the result of this reiterative nicking and extension being geometric amplification.

**Figure 2 – Alere™ i Influenza A & B Amplification Duplex – Geometric**



Phase

## **Alere™ i Influenza A & B Amplification Detection Mechanism**

The intentional skewing of templates leading to asymmetric amplification product formation is necessary for optimal performance of the real-time fluorescence-based molecular beacon detection system coupled to the Alere™ i Influenza A & B amplification technology. Molecular beacons are single-stranded DNA sequences that are able to form hairpin loop structures. The hairpin is formed by complementary 5' and 3' ends that are typically 4 – 7 nucleotides in length. Molecular beacons contain a fluorophore and a quencher, typically at opposite ends of the molecular beacon's linear nucleotide sequence. When the molecular beacon forms a hairpin structure, the fluorophore and quencher are in close proximity, and fluorescence is quenched. Upon opening of the hairpin, which can occur when a complementary DNA or RNA species anneals to the molecular beacon's loop region, the distance between the fluorophore and quencher is lengthened, reducing the ability of the quencher to quench the fluorescence signal. As a result, fluorescence is emitted and can be measured using an optical reader. The amount of fluorescence released is directly proportional to the amount of molecular beacon that has been opened, and concomitantly, the amount of complementary DNA or RNA species present in the sample. The result of the Alere™ i Influenza A & B NEAR amplification reaction is two specific products which are perfect complements to one another. These two specific products form a very stable double strand DNA structure at assay temperature. In order to combine the use of the Alere™ i Influenza A & B NEAR technology and molecular beacons, the concentrations of the two templates present in the Alere™ i Influenza A & B assay are skewed relative to one another to drive asymmetric product formation, i.e., one product is made to a larger extent than the other. The product that is generated in excess is then targeted by the molecular beacon. Annealing of the loop region of the molecular beacon to the Alere™ i Influenza A & B product generated in excess results in opening of the hairpin stem, separation of the fluorophore and quencher, and ultimately fluorescence signal release.

### **M. Performance Characteristics (if/when applicable):**

#### 1. Analytical performance:

##### *a. Precision/Reproducibility:*

##### Reproducibility Study

Reproducibility performance was evaluated for the Alere™ i Influenza A & B at three testing sites. Six Alere™ i instruments (per study site) were utilized at Site 1 and Site 2, and four Alere™ i instruments were utilized at Site 3. Three operators (Operator 1, 2, and 3) participated in the study at each site. Two Alere™ i Influenza A & B reagent lots were used in the Reproducibility Study.

The reproducibility panel consisted of seven samples at virus concentrations near the respective LoD (i.e., Moderate Positive Influenza A, Low Positive Influenza A, High Negative Influenza A, Moderate Positive Influenza B, Low Positive Influenza B, High Negative Influenza B, and True Negative). Each sample was prepared using the influenza A and B strains (described above) spiked into simulated nasal swab matrix, the Universal Transport Medium (UTM). The influenza A strain used in the reproducibility study was Influenza A/Puerto Rico/8/34 (Lot J1215A, inactivated), and the influenza B strain used in the study was Influenza B/Malaysia/2506/2004 (Lot B1012C, inactivated). The targeted concentrations for the Moderate Positive (MP) samples were approximately 2 to 3 X the respective LoD, the targeted concentrations for the Low Positive (LP) samples were approximately 1 X the respective LoD (C<sub>95</sub> concentration), and the targeted concentrations for the High Negative (HN) samples were less than the respective LoD (C<sub>5</sub>

concentration). The True Negative samples contain no influenza virus.

Contrived nasal swab samples were prepared using the influenza A and B strains (described above) dilutions in UTM to the targeted concentrations. The concentrations of the viral stocks (in TCID<sub>50</sub>/mL) were determined by standard virologic method prior to inactivation by the vendors. The targeted concentration for each dilution (in genome equivalents/mL) was also assessed using laboratory developed and validated influenza A and influenza B quantitative real-time PCR assays.

Swabs were prepared by coating 10 µl of the virus dilution onto the swab. Swab coating was done by lightly scratching the swab surface with the pipette tip as the dilution was pipetted to ensure liquid was absorbed into the swab tip. Once all solution was absorbed, swabs were allowed to dry and assembled into blinded panels.

Samples in each of the seven categories were tested twice by each of the three operators per study site, on five distinct days, at three clinical sites, using two lots of reagents (Site-to-Site: 2 tests/day/operator x 3 operators/per site x 5 days x 3 sites = 90 data points for each sample category).

The influenza A and influenza B strains were tested as follows:

- Moderate Positive (MP) influenza A:  $1.47 \times 10^4$  TCID<sub>50</sub>/Swab ( $2.38 \times 10^4$  genome equivalents/Swab); [Note: 10 µl of  $1.47 \times 10^6$  TCID<sub>50</sub>/mL ( $2.38 \times 10^6$  genome equivalents/mL) virus dilution coated onto a swab]
- Low Positive (LP) influenza A:  $4.91 \times 10^3$  TCID<sub>50</sub>/Swab ( $7.95 \times 10^3$  genome equivalents/Swab); [Note: 10 µl of  $4.91 \times 10^5$  TCID<sub>50</sub>/mL ( $7.95 \times 10^5$  genome equivalents/mL) virus dilution coated onto a swab]
- High Negative (HN) influenza A :  $6.80 \times 10^2$  TCID<sub>50</sub>/Swab ( $1.10 \times 10^3$  genome equivalents/Swab); [Note: 10 µl of  $6.80 \times 10^4$  TCID<sub>50</sub>/mL ( $1.10 \times 10^5$  genome equivalents/mL) virus dilution coated onto a swab]
- Moderate Positive (MP) influenza B:  $5.04 \times 10^1$  TCID<sub>50</sub>/Swab ( $6.00 \times 10^2$  genome equivalents/Swab); [Note: 10 µl of  $5.04 \times 10^3$  TCID<sub>50</sub>/mL ( $6.00 \times 10^4$  genome equivalents/mL) virus dilution coated onto a swab]
- Low Positive (LP) influenza B:  $1.68 \times 10^1$  TCID<sub>50</sub>/Swab ( $2.00 \times 10^2$  genome equivalents/Swab); [Note: 10 µl of  $1.68 \times 10^3$  TCID<sub>50</sub>/mL ( $2.00 \times 10^4$  genome equivalents/mL) virus dilution coated onto a swab]
- High Negative (HN) influenza B:  $4.20 \times 10^{-1}$  TCID<sub>50</sub>/Swab ( $5.00 \times 10^0$  genome equivalents/Swab); [Note: 10 µl of  $4.20 \times 10^1$  TCID<sub>50</sub>/mL ( $5.00 \times 10^2$  genome equivalents/mL) virus dilution coated onto a swab]
- True negative (TN): Negative samples (no virus)

All samples were tested using the Alere™ i Influenza A & B on the Alere™ i instruments, according to the product instructions. Each operator tested one external positive and one external negative control swab on each instrument used for the study, on each day of testing prior to performing the study testing.

Operator 2 at Site 2 generated a single invalid result for the external negative control swab tested on Day 1 of the Reproducibility Study. The expected negative result was obtained upon repeat testing per the product instructions. All other external positive and negative control swabs tested in the Reproducibility Study generated the expected results each day of testing at all three sites.

On Day 1 of the Reproducibility Study, Operator 2 at Site 1 generated a “Flu A Invalid/Flu B

Negative” result for a TN sample, and a “Flu A Invalid/Flu B Invalid” result for a HN influenza B and a MP influenza B sample. Upon repeat testing per the product instructions, valid expected results were obtained for all three samples. No other invalid results were obtained. These repeated test results were included in the Reproducibility Study results presented in Table 3 and Table 4 below.

For the Site-to-Site Reproducibility Study, the overall percent agreement with expected result (i.e., negative for TN and HN samples, positive for MP and LP samples) was 100% (90/90, 95% CI, 95.9% - 100%) for the influenza A MP, the influenza A LP, the influenza B MP, and the TN categories; 92.2% (83/90, 95% CI, 84.8% - 96.2%) for the influenza B LP category; and 70.0% (63/90, 95% CI, 59.9% - 78.5%) and 90.0% (81/90, 95% CI, 82.1% - 94.6%) for the influenza A HN and the influenza B HN categories, respectively.

The Reproducibility Study site-to-site qualitative results (agreements with expected results) are presented in Table 3 below:

**Table 3: Reproducibility Study Site-To-Site Qualitative Results**

Sample Category	SITE						Overall Percent Agreement and 95% CI	
	Site 1		Site 2		Site 3			
	Percent Agreement	Count	Percent Agreement	Count	Percent Agreement	Count		
HN <sup>1</sup> Influenza A	66.7%	20/30	80.0%	24/30	63.3%	19/30	70.0% (63/90)	(59.9%, 78.5%)
LP Influenza A	100%	30/30	100%	30/30	100%	30/30	100% (90/90)	(95.9%, 100%)
MP Influenza A	100%	30/30	100%	30/30	100%	30/30	100% (90/90)	(95.9%, 100%)
HN <sup>1</sup> Influenza B	86.7%	26/30	100%	30/30	83.3%	25/30	90.0% (81/90)	(82.1%, 94.6%)
LP Influenza B	93.3%	28/30	86.7%	26/30	96.7%	29/30	92.2% (83/90)	(84.8%, 96.2%)
MP Influenza B	100%	30/30	100%	30/30	100%	30/30	100% (90/90)	(95.9%, 100.0%)
TN	100%	30/30	100%	30/30	100%	30/30	100% (90/90)	(95.9%, 100%)

<sup>1</sup>Percent Agreement correlates to the percent of negative results.

The Reproducibility Study operator-to-operator qualitative results (agreements with expected results) are presented in Table 4 below:

**Table 4: Reproducibility Study Operator-To-Operator Qualitative Results**

Sample Category		HN <sup>1</sup> Influenza A	LP Influenza A	MP Influenza A	HN <sup>1</sup> Influenza B	LP Influenza B	MP Influenza B	TN	Total	
Site	Operator	Count (% Agreement)	Count (% Agreement)	Count (% Agreement)	Count (% Agreement)	Count (% Agreement)	Count (% Agreement)	Count (% Agreement)	Count (% Agreement)	95% CI
Site 1	Operator 1	6/10 (60.0%)	10/10 (100%)	10/10 (100%)	9/10 (90.0%)	9/10 (90.0%)	10/10 (100%)	10/10 (100%)	64/70 (91.4%)	82.5% - 96.0%
	Operator 2	6/10 (60.0%)	10/10 (100%)	10/10 (100%)	9/10 (90.0%)	9/10 (90.0%)	10/10 (100%)	10/10 (100%)	64/70 (91.4%)	82.5% - 96.0%
	Operator 3	8/10 (80.0%)	10/10 (100%)	10/10 (100%)	8/10 (80.0%)	10/10 (100%)	10/10 (100%)	10/10 (100%)	66/70 (94.3%)	86.2% - 97.8%
Site 2	Operator 1	8/10 (80.0%)	10/10 (100%)	10/10 (100%)	10/10 (100%)	9/10 (90%)	10/10 (100%)	10/10 (100%)	67/70 (95.7%)	88.1% - 98.5%
	Operator 2	6/10 (60.0%)	10/10 (100%)	10/10 (100%)	10/10 (100%)	7/10 (70%)	10/10 (100%)	10/10 (100%)	63/70 (90.0%)	80.8% - 95.1%

	Operator 3	10/10 (100%)	70/70 (100%)	94.8% - 100%						
Site 3	Operator 1	5/10 (50%)	10/10 (100%)	10/10 (100%)	8/10 (80.0%)	10/10 (100%)	10/10 (100%)	10/10 (100%)	63/70 (90.0%)	80.8% - 95.1%
	Operator 2	8/10 (80%)	10/10 (100%)	10/10 (100%)	8/10 (80.0%)	9/10 (90%)	10/10 (100%)	10/10 (100%)	65/70 (92.9%)	84.3% - 96.9%
	Operator 3	6/10 (60%)	10/10 (100%)	10/10 (100%)	9/10 (90%)	10/10 (100%)	10/10 (100%)	10/10 (100%)	65/70 (92.9%)	84.3% - 96.9%

<sup>1</sup>Percent Agreement correlates to the percent of negative results.

The Reproducibility Study instrument-to-instrument qualitative results (agreements with expected results) are presented in Table 5 below:

**Table 5: Reproducibility Study Instrument-To-Instrument Qualitative Results**

Sample Category		HN <sup>1</sup> Influenza A	LP Influenza A	MP Influenza A	HN <sup>1</sup> Influenza B	LP Influenza B	MP Influenza B	TN	Total	
Site	Instrument	Count (% Agreement)	Count (% Agreement)	Count (% Agreement)	Count (% Agreement)	Count (% Agreement)	Count (% Agreement)	Count (% Agreement)	Count (% Agreement)	95% CI
Site 1	Instrument 1	N/A	N/A	N/A	N/A	4/4 (100%)	N/A	N/A	4/4 (100%)	51.0% - 100%
	Instrument 2	3/9 (33.3%)	6/6 (100%)	13/13 (100%)	7/8 (87.5%)	4/5 (90%)	9/9 (100%)	8/8 (100%)	50/58 (86.2%)	75.1% - 92/8%
	Instrument 3	N/A	N/A	N/A	N/A	5/5 (100%)	N/A	N/A	5/5 (100%)	56.6% - 100%
	Instrument 4	7/8 (87.5%)	11/11 (100%)	8/8 (100%)	9/10 (90%)	6/6 (100%)	9/9 (100%)	13/13 (100%)	63/65 (96.9%)	89.5% - 99.2%
	Instrument 5	7/10 (70%)	11/11 (100%)	9/9 (100%)	7/9 (77.8%)	3/4 (75%)	11/11 (100%)	5/5 (100%)	53/59 (89.8%)	79.5% - 95.3%
	Instrument 6	3/3 (100%)	2/2 (100%)	N/A	3/3 (100%)	6/6 (100%)	1/1 (100%)	4/4 (100%)	19/19 (100%)	83.2% - 100%
Site 2	Instrument 1	2/3 (66.7%)	4/4 (100%)	4/4 (100%)	3/3 (100%)	2/3 (66.7%)	4/4 (100%)	5/5 (100%)	24/26 (92.3%)	75.9% - 97.9%
	Instrument 2	4/6 (66.7)	6/6 (100%)	6/6 (100%)	7/7 (100%)	6/6 (100%)	10/10 (100%)	4/4 (100%)	44/46 (95.7%)	85.5% - 98.8%
	Instrument 3	2/3 (66.7%)	2/2 (100%)	3/3 (100%)	1/1 (100%)	2/3 (66.7%)	1/1 (100%)	3/3 (100%)	14/16 (87.5%)	64.0% - 96.5%
	Instrument 4	11/13 (84.6%)	5/5 (100%)	8/8 (100%)	6/6 (100%)	6/6 (100%)	5/5 (100%)	6/6 (100%)	47/49 (95.9%)	86.3% - 98.9%
	Instrument 5	2/2 (100%)	4/4 (100%)	2/2 (100%)	4/4 (100%)	4/6 (66.7%)	4/4 (100%)	5/5 (100%)	25/27 (92.6%)	76.6% - 97.9%
	Instrument 6	3/3 (100%)	9/9 (100%)	7/7 (100%)	9/9 (100%)	6/6 (100%)	6/6 (100%)	7/7 (100%)	47/47(100%)	92.4% - 100%
Site 3	Instrument 1	6/8 (75%)	10/10 (100%)	12/12 (100%)	6/8 (75%)	14/14 (100%)	12/12 (100%)	10/10 (100%)	70/74 (94.6%)	86.9% - 97.9%
	Instrument 2	7/12 (58.3%)	10/10 (100%)	8/8 (100%)	10/12 (83.3%)	15/16 (93.8%)	8/8 (100%)	10/10 (100%)	68/76 (89.5%)	80.6% - 94.6%
	Instrument 3	4/6 (66.7%)	5/5 (100%)	4/4 (100%)	6/6 (100%)	N/A	4/4 (100%)	5/5 (100%)	28/30 (93.3%)	78.7% - 98.2%
	Instrument 4	2/4 (50%)	5/5 (100%)	6/6 (100%)	3/4 (75%)	N/A	6/6 (100%)	5/5 (100%)	27/30 (90.0%)	74.4% - 96.5%

<sup>1</sup>Percent Agreement correlates to the percent of negative results.

Due to the fact that evaluating reagent kit lot to lot reproducibility was not addressed by the Reproducibility Study design, although two reagent kit lots were included in the Reproducibility Study, most of the data points were generated using Lot 1. To assess reagent

kit lot to lot reproducibility of the Alere™ i Influenza A & B, the results of the “Near the Cut Off” study were utilized since the same sample categories that were tested in the Reproducibility Study were also tested in the “Near the Cut Off” study. The reproducibility assessment lot-to-lot qualitative results (agreements with expected results) are presented in Table 6 below:

**Table 6: Reproducibility Study Lot-To-Lot Qualitative Results**

Sample Category	Reagent Lot						Overall Count and Percent Agreement	95% CI
	Lot 1		Lot 2		Lot 3			
	Percent Agreement	Count	Percent Agreement	Count	Percent Agreement	Count		
HN <sup>1</sup> Influenza A	70.0%	63/90	N/A	N/A	46.7%	28/60	91/150 (60.7%)	52.7% - 68.1%
LP Influenza A	100%	90/90	N/A	N/A	100%	60/60	150/150 (100%)	97.5% - 100%
MP Influenza A	100%	90/90	N/A	N/A	N/A	N/A	90/90 (100%)	95.9% - 100%
HN <sup>1</sup> Influenza B	90.0%	81/90	N/A	N/A	91.7%	55/60	136/150 (90.7%)	84.9% - 94.4%
LP Influenza B	96.8%	60/62	82.1%	23/28	96.7%	58/60	141/150 (94.0%)	89.0% - 96.8%
MP Influenza B	100%	90/90	N/A	N/A	N/A	N/A	90/90 (100%)	95.9% - 100%
TN	100%	90/90	N/A	N/A	100%	60/60	150/150 (100%)	97.5% - 100%

<sup>1</sup>Percent Agreement correlates to the percent of negative results.

*b. Linearity/assay reportable range:*

Not applicable

*c. Traceability, Stability, Expected values (controls, calibrators, or methods):*

**Controls**

External Positive and Negative Controls

Alere™ i Influenza A & B kits contain Positive and Negative Control Swabs. These swabs monitor the entire assay. The Alere™ i Influenza A & B Instructions for Use requires testing these control swabs once with each new shipment received and once for each new operator. Further controls may be tested in order to conform to local, state and/or federal regulations, accrediting groups, or laboratory’s standard Quality Control procedures.

During the prospective clinical study, control testing using external positive and negative controls was performed prior to sample testing each day and on each instrument Alere™ i Influenza A & B testing was performed, at all eight study sites. A total of 335 external positive controls and 335 external negative controls were tested during the prospective clinical study. Seven external positive controls (7/335, 2.1%, 95% CI, 1.0% - 4.2%) and nine external negative controls (9/335, 2.7%, 95% CI, 1.4% - 5.0%) failed initially, but succeeded on repeat testing per the product instructions.

## Natural Nasal Swab Matrix and Simulated Nasal Swab Matrix Equivalency

Negative natural clinical nasal swab matrix is difficult to obtain in large volumes to be used as diluent in various analytical studies and the Reproducibility Study. In order to compensate for the limited availability of negative natural clinical nasal swab matrix, a simulated nasal swab matrix, the Universal Transport Medium (UTM), was used as diluent in the Reproducibility Study and in a number of analytical studies (i.e., Carryover Contamination Study, Swab Type Validation Study, Reagent Shelf Life Study, Interfering Substances Study, Competitive Interference with Non-Influenza Respiratory Viruses Study, and Competitive Interference with Influenza Viruses Study).

Limit of Detection Studies (UTM vs. natural clinical nasal swab matrix in UTM) were carried out to demonstrate the equivalency between the natural clinical nasal swab matrix and the simulated nasal swab matrix.

The anchoring influenza strains used in the LoD studies were Influenza A/Puerto Rico/8/34 (Lot D1002B - inactivated, and Lot G1112B), and Influenza B/Malaysia/2506/2004 (Lot B1012C - inactivated, and Lot C1320E).

For the LoD study using UTM as the diluent, dilutions of each of the two anchoring influenza strains (described above) were prepared in UTM. 10 µl of each dilution was dispensed onto a swab, using the pipette tip to scratch and rub liquid onto the swab tip. The swab sample was then tested per the product instructions. A total of 20 replicates were tested at the estimated LoD level.

For the LoD study using natural clinical nasal swab matrix in UTM as the diluent, presumed negative natural clinical nasal swabs were eluted in 3 mL of UTM. Swab eluates were combined and mixed thoroughly to create a swab eluate pool. Dilutions of each of the two anchoring influenza strains (described above) were prepared in this swab eluate pool. 10 µl of each dilution was dispensed onto each swab, using the pipette tip to scratch and rub liquid onto the swab tip. The swab sample was then tested per the product instructions. A total of 20 replicates were tested at the estimated LoD level.

The concentrations of the viral stocks (in TCID<sub>50</sub>/mL) were determined by standard virologic method by the vendors and confirmed by the sponsor (except for the inactivated virus lots, where the TCID<sub>50</sub>/mL values were determined by standard virologic method by the vendors prior to inactivation). The concentration for each dilution (in genome equivalents/mL) was also assessed using laboratory developed and validated influenza A and influenza B quantitative real-time PCR assays.

The estimated LoD concentration (defined as the lowest concentration at which ≥ 95% of all replicates tested positive), presented in both calculated TCID<sub>50</sub>/reaction and genome equivalents/reaction, per virus strain tested, is presented in Table 7 below.

**Table 7: Natural Nasal Swab Matrix and Simulated Nasal Swab Matrix Equivalency Demonstration**

A/Puerto Rico/8/34 (H1N1)	Lot	TCID <sub>50</sub> /Reaction	Genome Equivalents/Reaction
LoD Study –Clinical Matrix	G1112B	75	1688
LoD Study –Without Clinical Matrix	D1002B (inactivated)	50	304
B/Malaysia/2506/2004	Lot	TCID <sub>50</sub> /Reaction	Genome Equivalents/Reaction
LoD Study –Clinical Matrix	C1320E	0.8	29
LoD Study –Without Clinical Matrix	B1012C (inactivated)	0.7	8

This analysis demonstrated that the Alere™ i Influenza A & B assay LoDs using simulated nasal swab matrix (i.e., UTM) are similar to those using natural nasal swab matrix. Therefore, the simulated nasal swab matrix is appropriate and adequate to be used in place of natural nasal swab matrix for the analytical studies, including the Reproducibility Study.

### Specimen Stability Study

To provide data supporting the specimen storage recommendations stated in the product package insert, an analytical study was carried out to evaluate specimen stability.

Dilutions of Influenza A/Puerto Rico/8/34 (Lot D1002B, inactivated) and Influenza B/Malaysia/2506/2004 (Lot B1012C, inactivated) were used in this study to prepare contrived samples for testing. The concentrations of the viral stocks (in TCID<sub>50</sub>/mL) were determined by standard virologic method prior to inactivation by the vendors. The targeted concentration for each dilution (in genome equivalents/mL) was also assessed using laboratory developed and validated influenza A and influenza B quantitative real-time PCR assays.

Contrived positive nasal swab samples were prepared using the influenza A and B strains (described above) dilutions in saline with targeted concentrations near the respective LoD levels (approximately 2 to 3 X the respective LoD). Mock negative nasal swab samples were prepared using saline. Swabs were prepared by coating 10 µl of the virus dilution or saline onto the swab. Swab coating was done by lightly scratching the swab surface with the pipette tip as the dilution was pipetted to ensure liquid was absorbed into the swab tip. All prepared swabs were tested in triplicates at each of the following three time points/conditions: Time 0, Hour 3 (after being stored at approximately 18-22°C for 3 hours), and Hour 25 (after being stored at approximately 18-22°C for 3 hours, and then placed at 2-8°C for 22 hours of storage).

The contrived influenza A and influenza B positive swabs were tested at each time point as follows:

- Contrived influenza A swab: 3.30 x 10<sup>3</sup> TCID<sub>50</sub>/Swab (2.00 x 10<sup>4</sup> genome equivalents/Swab); [Note: 10 µl of 3.30 x 10<sup>5</sup> TCID<sub>50</sub>/mL (2.00 x 10<sup>6</sup> genome equivalents/mL) virus dilution coated onto a swab]
- Contrived influenza B swab: 4.20 x 10<sup>1</sup> TCID<sub>50</sub>/Swab (5.00 x 10<sup>2</sup> genome equivalents/Swab); [Note: 10 µl of 4.20 x 10<sup>3</sup> TCID<sub>50</sub>/mL (5.00 x 10<sup>4</sup> genome equivalents/mL) virus dilution coated onto a swab]

Direct nasal swab testing results are presented in Table 8 below.

**Table 8: Direct Nasal Swab Specimen Stability Study Results**

Time Point	Mock Negative Swab		Contrived Flu A Positive Swab		Contrived Flu B Positive Swab	
	Flu A Result	Flu B Result	Flu A Result	Flu B Result	Flu A Result	Flu B Result
Time 0	Negative (3/3)	Negative (3/3)	Positive (3/3)	Negative (3/3)	Negative (3/3)	Positive (3/3)
Hour 3 (approximately 18-22°C)	Negative (3/3)	Negative (3/3)	Positive (3/3)	Negative (3/3)	Negative (3/3)	Positive (3/3)
Hour 25 (22 hrs at 2-8°C)	Negative (3/3)	Negative (3/3)	Positive (3/3)	Negative (3/3)	Negative (3/3)	Positive (3/3)

This analytical study data supports Alere™ i Influenza A & B product insert sample storage recommendations for direct nasal swab samples of up to 3 hours at 18-22°C, or at 2-8°C for up to 25 hours (22 hours at 2-8°C and 3 hours at 18-22°C). This analytical study data also supports Alere™ i Influenza A & B product insert sample storage recommendations for nasal swab of up to 2 hours at 18-22°C, or at 2-8°C for up to 24 hours.

## Product Shelf Life Determination

### Alere™ i Influenza A & B Assay Components Stability

A product stability study was conducted (and currently on-going) to generate real-time stability data to support the expiry dating assigned to Alere™ i Influenza A & B Kit.

Two lots of Alere™ i Influenza A & B were included in this stability study. All test components (Sample Receiver, Test Base, and Transfer Cartridge) were placed at room temperature (28-32°C), 2-8°C and “Customer Storage Condition” (i.e., Test Base at 2-8°C, Sample Receiver and Transfer Cartridge at 28-32°C), and held at these controlled temperatures for the duration of the stability study.

Dilutions of Influenza A/Puerto Rico/8/34 (Lot D1002B, inactivated) and Influenza B/Malaysia/2506/2004 (Lot B1012C, inactivated, and A1012A) were used in this study to prepare contrived positive samples for testing. The concentrations of the viral stocks (in TCID<sub>50</sub>/mL) were determined by standard virologic method prior to inactivation by the vendors. The targeted concentration for each dilution (in genome equivalents/mL) was also assessed using laboratory developed and validated influenza A and influenza B quantitative real-time PCR assays. Elution buffer, Alere™ i Influenza A & B Negative Control Swab (polyester swab hand coated with approximately  $1 \times 10^7$  *Streptococcus dysgalactiae* (ATCC 12388) organisms), and presumed negative natural nasal swabs (collected from asymptomatic donors) were tested as negative samples in this study.

Contrived positive nasal swab samples were prepared using the influenza A and B strains (described above) diluted in UTM with targeted concentrations near the respective LoD levels. 50 µl of each dilution was added to the Sample Receiver and pipetted up and down to mix. The contrived influenza A and influenza B positive samples were tested in triplicates at each specified time point for each stability condition as follows:

- Contrived influenza A #1 (tested at Time 0, Month 4, Month 7, Month 10, and Month 13):  $3.30 \times 10^3$  TCID<sub>50</sub> ( $2.00 \times 10^4$  genome equivalents) in the Sample Receiver; [Note: 50 µl of  $6.59 \times 10^4$  TCID<sub>50</sub>/mL ( $4.00 \times 10^5$  genome equivalents/mL) virus dilution added to the Sample Receiver]
- Contrived influenza A #2 (tested at Time 0, Month 7, and Month 13):  $6.60 \times 10^3$  TCID<sub>50</sub> ( $4.00 \times 10^4$  genome equivalents) in the Sample Receiver; [Note: 50 µl of  $1.32 \times 10^5$  TCID<sub>50</sub>/mL ( $8.01 \times 10^5$  genome equivalents/mL) virus dilution added to the Sample Receiver]
- Contrived influenza B #1A (tested at Time 0 and Month 4):  $2.10 \times 10^2$  TCID<sub>50</sub> ( $2.50 \times 10^3$  genome equivalents) in the Sample Receiver; [Note: 50 µl of  $4.20 \times 10^3$  TCID<sub>50</sub>/mL ( $5.00 \times 10^4$  genome equivalents/mL) virus dilution added to the Sample Receiver]
- Contrived influenza B #2A (tested at Time 0):  $4.20 \times 10^2$  TCID<sub>50</sub> ( $5.00 \times 10^3$  genome equivalents) in the Sample Receiver; [Note: 50 µl of  $8.40 \times 10^3$  TCID<sub>50</sub>/mL ( $1.00 \times 10^5$  genome equivalents/mL) virus dilution added to the Sample Receiver]
- Contrived influenza B #1B (tested at Month 7):  $1.20 \times 10^2$  TCID<sub>50</sub> ( $1.50 \times 10^3$  genome equivalents) in the Sample Receiver; [Note: 50 µl of  $2.39 \times 10^3$  TCID<sub>50</sub>/mL ( $3.00 \times 10^4$  genome equivalents/mL) virus dilution added to the Sample Receiver]
- Contrived influenza B #2B (tested at Month 7):  $2.39 \times 10^2$  TCID<sub>50</sub> ( $3.00 \times 10^3$  genome

equivalents) in the Sample Receiver; [Note: 50  $\mu\text{l}$  of  $4.78 \times 10^3$  TCID<sub>50</sub>/mL ( $6.00 \times 10^4$  genome equivalents/mL) virus dilution added to the Sample Receiver]

- Contrived influenza B #1C (tested at Time 10 and Month 13):  $4.20 \times 10^1$  TCID<sub>50</sub> ( $5.00 \times 10^2$  genome equivalents) in the Sample Receiver; [Note: 50  $\mu\text{l}$  of  $8.40 \times 10^2$  TCID<sub>50</sub>/mL ( $1.00 \times 10^4$  genome equivalents/mL) virus dilution added to the Sample Receiver]
- Contrived influenza B #2C (tested at Month 13):  $8.40 \times 10^1$  TCID<sub>50</sub> ( $1.00 \times 10^3$  genome equivalents) in the Sample Receiver; [Note: 50  $\mu\text{l}$  of  $1.68 \times 10^3$  TCID<sub>50</sub>/mL ( $2.00 \times 10^4$  genome equivalents/mL) virus dilution added to the Sample Receiver]

To perform testing of the Elution Buffer, no sample was added to the Sample Receiver. The Elution Buffer contained in the Sample Receiver was transferred to the Test Base using the Transfer Cartridge when prompted by the instrument according to the product instructions. For testing the Alere™ i Influenza A & B Negative Control Swabs and the presumed negative natural nasal swabs, the swabs were tested according to the product instructions.

The Elution Buffer and the Alere™ i Influenza A & B Negative Control Swab sample were tested in triplicates at Time 0, Month 3, Month 4, Month 7, Month 10, and Month 13, for each stability condition. A total of 20 presumed negative natural nasal swabs were tested at Time 0, Month 8, and Month 13, for each stability condition.

The Alere™ i Influenza A & B assay components were considered “stable” as long as they continued to produce the expected influenza A and B positive results for the positive A & B samples and expected negative results for the Elution Buffer, the Alere™ i Influenza A & B Negative Control Swabs, and the presumed negative natural nasal swabs.

The stability study at 2- 8°C and at the “Customer Storage Condition” (i.e., Test Base at 2-8°C, Sample Receiver and Transfer Cartridge at 28-32°C) is currently ongoing, and available real-time data supports a shelf life of up to 12 months at 2- 8°C and at the “Customer Storage Condition”. The stability study for the entire kit at room temperature (28-32°C) was discontinued after four months due to test results not meeting the acceptance criteria. No stability claim for the entire kit at room temperature (28-32°C) can be supported by real-time data.

#### Alere™ i Influenza A & B External Positive and Negative Control Swabs Stability

The Alere™ i Influenza A & B Positive Control Swab is the same swab control currently used in the FDA cleared BinaxNOW® Influenza A & B test kit. The Alere™ i Influenza A & B Negative Control Swab is the same swab control currently used in the FDA cleared BinaxNOW® Strep A test kit. An analytical study was carried out to confirm the stability claim for the Alere™ i Influenza A & B Positive Control and Negative Control Swab, expiry dating of 2 years when stored at room temperature (approximately 18-22°C), which was established based on real-time stability testing data in the context of the FDA cleared BinaxNOW® Influenza A & B test kit and the BinaxNOW® Strep A test kit.

For the Alere™ i Influenza A & B Positive Control Swab, 20 Positive Control swabs were tested 25 months postdate of manufacture (i.e., 1 month post expiry) according to the product instructions. Expected Influenza A and B dual positive results were obtained for all 20 Positive Control swabs tested.

For the Alere™ i Influenza A & B Positive Control Swab, 20 Negative Control swabs were tested 16 months postdate of manufacture according to the product instructions. Expected Influenza A and B negative results were obtained for all 20 Negative Control swabs tested.

This part of the study is on-going and 20 Negative Control swabs will be tested again at expiry (i.e., December 2014).

## Swab Type Compatibility Study

To assess the compatibility of various swabs to be used with the Alere™ i Influenza A & B, an analytical study was conducted testing the following swabs:

Puritan PurFlock® Ultra flocked swab (standard tip), Puritan PurFlock® Ultra flocked swab (mini tip), Puritan PurFlock® Ultra flocked swab (large tip), HydraFlock® Flocked swab (standard tip), Foam swabs, Polyester swabs, Rayon swabs, Puritan Calcium Alginate swabs, and Copan Regular Flocked swabs.

Dilutions of Influenza A/Puerto Rico/8/34 (Lot D1002B, inactivated) and Influenza B/Malaysia/2506/2004 (Lot B1012C, inactivated) were used in this study to prepare contrived swab samples for testing. The concentrations of the viral stocks (in TCID<sub>50</sub>/mL) were determined by standard virologic method prior to inactivation by the vendors. The targeted concentration for each dilution (in genome equivalents/mL) was also assessed using laboratory developed and validated influenza A and influenza B quantitative real-time PCR assays.

Contrived positive nasal swab samples for each swab type were prepared using the influenza A and B strains (described above) dilutions in UTM with targeted concentrations near the respective LoD levels (approximately 2 to 3 X the respective LoD). Mock negative nasal swab samples for each swab type were prepared using UTM. Swabs were prepared by coating 10 µl of the virus dilution or UTM onto the swab. Swab coating was done by lightly scratching the swab surface with the pipette tip as the dilution was pipetted to ensure liquid was absorbed into the swab tip. All prepared swabs were tested in triplicates according to the product instructions for testing direct swabs.

The influenza A and influenza B positive swab samples with concentration approximately 2 to 3 X times the respective LoD were tested for each swab type as follows:

- Contrived influenza A swab:  $3.30 \times 10^3$  TCID<sub>50</sub>/Swab ( $2.00 \times 10^4$  genome equivalents/Swab); [Note: 10 µl of  $3.30 \times 10^5$  TCID<sub>50</sub>/mL ( $2.00 \times 10^6$  genome equivalents/mL) virus dilution coated onto a swab]
- Contrived influenza B swab:  $4.20 \times 10^1$  TCID<sub>50</sub>/Swab ( $5.00 \times 10^2$  genome equivalents/Swab); [Note: 10 µl of  $4.20 \times 10^3$  TCID<sub>50</sub>/mL ( $5.00 \times 10^4$  genome equivalents/mL) virus dilution coated onto a swab]

A swab type was considered to be compatible for use with the Alere™ i Influenza A & B assay if expected results for positive and negative samples were achieved for all replicates tested.

Of the nine swab types evaluated in this study, all three Puritan PurFlock® Ultra swab types (with standard, mini, or large tips) and the Calcium Alginate swabs were found to be not compatible or appropriate for use with Alere™ i Influenza A & B. The following statement is included in the “Sample Collection and Handling” section of the Product Insert: “*Calcium Alginate and Puritan Purflock® Ultra flocked swabs are not suitable for use in this assay.*”

The other five swab types evaluated in this study, HydraFlock® Flocked swab (standard tip), Foam swabs, Polyester swabs, Rayon swabs, and Copan Regular Flocked swabs were found to be compatible for use with Alere™ i Influenza A & B. The following statements have been added to the “Sample Collection and Handling” section of the Product Insert: “*For optimal performance, use the swabs provided in the test kit. Alternatively, rayon, foam, HydraFlock® Flocked Swab (standard tip), Copan Regular Flocked Swab or polyester nasal swabs can be used to collect nasal swab samples.*”

d. *Detection Limit:*

LoD Study – With Natural Clinical Matrix

An analytical study was carried out to determine the LoD for each Alere™ i Influenza A & B assay targeted influenza viruses in natural nasal swab matrix. Three strains of influenza A virus representing each of the three common currently or recently circulating influenza A subtypes (i.e., A/H1N1, A/H3N2 seasonal, and A/H1N1 pandemic (pdm)) and two strains of influenza B virus representing each of the two influenza B genetic lineages (i.e., Victoria and Yamagata) were included in this study.

Presumed negative nasal swab specimens were eluted in 3 mL of UTM. Swab elutes were combined and mixed thoroughly by gentle vortex to create a clinical matrix pool, and this clinical matrix pool was used as the diluent throughout this study. Vender provided stocks of Influenza A/Puerto Rico/8/34 (A/H1N1, Lot G1112B), A/Perth/16/2009 (A/H3N2, Lot A1121B), A/California/7/2009 (A/H1N1 pdm, Lot E1020A1), Influenza B/Malaysia/2506/2004 (Victoria lineage, Lot C1320E), and B/Bangladesh/3333/2007 (Yamagata lineage, Lot H1012A) were diluted in the natural clinical nasal swab matrix pool as described previously to generate virus dilutions for testing. The vender provided virus strains were re-titered and the concentrations (in TCID<sub>50</sub>/mL) were determined by standard virologic method. The concentration for each dilution (in genome equivalents/mL) was also assessed using laboratory developed and validated influenza A and influenza B quantitative real-time PCR assays.

Contrived swab specimens were prepared by coating 10 µl of virus dilution onto each swab. Swab coating was carried out by lightly scratching the swab surface with the pipette tip as virus dilution was pipetted to ensure liquid was absorbed into the swab tip. The contrived swab samples were tested according to the product instructions for testing direct swabs.

A range finding experiment was conducted initially testing triplicates per dilution per strain. If 100% positive results were obtained with the starting dilution, a series of 2-fold dilutions beginning with the dilution that produced 100% positive results were prepared and tested in triplicates until a negative result was obtained. 20 replicates of the lowest 2-fold dilution that produced 100% (3/3) positive results were then tested for the confirmation of the estimated limit of detection (LoD). If ≥95% positive results were achieved, the level was reported as the confirmed LoD. If <95% positive results were obtained, 20 replicates of the next 2-fold dilution up from the dilution that produced a <95% positive results were tested.

The confirmed LoD in natural nasal swab matrix, defined as the lowest virus level that was detected ≥95% of the time (i.e., virus level at which at least 19 out of 20 replicates tested positive), for each influenza strain tested, were presented in Table 9 below:

**Table 9: Limit of Detection (LoD) Study Results – Natural Nasal Swab Matrix**

Influenza Strain	Influenza A Subtype or Influenza B Genetic Lineage	LoD (TCID <sub>50</sub> /mL)	LoD (TCID <sub>50</sub> /Swab)*	LoD (Genome Equivalents/mL)	LoD (Genome Equivalents/Swab)*
A/Puerto Rico/8/34	A/H1N1	1.88 x 10 <sup>5</sup>	1.88 x 10 <sup>3</sup>	4.22 x 10 <sup>6</sup>	4.22 x 10 <sup>4</sup>
A/Perth/16/2009	A/H3N2	8.60 x 10 <sup>2</sup>	8.60 x 10 <sup>0</sup>	7.91 x 10 <sup>4</sup>	7.91 x 10 <sup>2</sup>
A/California/7/2009	A/H1N1 (pdm)	1.25 x 10 <sup>4</sup>	1.25 x 10 <sup>2</sup>	5.20 x 10 <sup>6</sup>	5.20 x 10 <sup>4</sup>
B/Malaysia/2506/2004	B Victoria lineage	1.90 x 10 <sup>3</sup>	1.90 x 10 <sup>1</sup>	7.24 x 10 <sup>4</sup>	7.24 x 10 <sup>2</sup>
B/Bangladesh/3333/2007	B Yamagata lineage	5.55 x 10 <sup>2</sup>	5.55 x 10 <sup>0</sup>	7.36 x 10 <sup>4</sup>	7.36 x 10 <sup>2</sup>

\*Note: 10 ul of each virus dilution was coated onto a swab

## LoD Study – With Simulated Matrix

An additional analytical study was also carried out to determine the LoD for each Aleré™ i Influenza A & B assay targeted influenza viruses in simulated nasal swab matrix, UTM. Two strains of influenza A virus (A/H1N1) and two strains of influenza B virus were included in this study.

Vender provided stocks of Influenza A/Puerto Rico/8/34 (A/H1N1, Lot D1002B, inactivated), A/Solomon Islands/3/2006 (A/H1N1, Lot G0931A1), Influenza B/Malaysia/2506/2004 (Victoria lineage, Lot B1012C, inactivated), and B/Brisbane/60/2008 (Victoria lineage, Lot D1029C1) were diluted in UTM to generate virus dilutions for testing.

The concentrations of the viral stocks (in TCID<sub>50</sub>/mL) for A/Puerto Rico/8/34 and B/Malaysia/2506/2004 were determined by standard virologic method prior to inactivation by the vender. For B/Brisbane/60/2008, the vender provided virus strain was re-titered and the concentration (in TCID<sub>50</sub>/mL) was determined by standard virologic method. For A/Solomon Islands/3/2006, the vender provided concentration (in TCID<sub>50</sub>/mL) was not confirmed by re-titering due to unavailability of residual viral stock. The concentration for each dilution (in genome equivalents/mL) was also assessed using laboratory developed and validated influenza A and influenza B quantitative real-time PCR assays.

Contrived swab specimens were prepared by coating 10 µl of virus dilution onto each swab. Swab coating was carried out by lightly scratching the swab surface with the pipette tip as virus dilution was pipetted to ensure liquid was absorbed into the swab tip. The contrived swab samples were tested according to the product instructions for testing direct swabs.

A range finding experiment was conducted initially testing triplicates per dilution per strain. If 100% positive results were obtained with the starting dilution, a series of 2-fold dilutions beginning with the dilution that produced 100% positive results were prepared and tested in triplicates until a negative result was obtained. 20 replicates of the lowest 2-fold dilution that produced 100% (3/3) positive results were then tested for the confirmation of the estimated limit of detection (LoD). If ≥95% positive results were achieved, the level was reported as the confirmed LoD. If <95% positive results were obtained, 20 replicates of the next 2-fold dilution up from the dilution that produced a <95% positive results were tested.

The confirmed LoD in simulated nasal swab matrix (UTM), defined as the lowest virus level that was detected ≥95% of the time (i.e., virus level at which at least 19 out of 20 replicates tested positive), for each influenza strain tested, were presented in Table 10 below:

**Table 10: Limit of Detection (LoD) Study Results – Simulated Nasal Swab Matrix**

Influenza Strain	Influenza A Subtype or Influenza B Genetic Lineage	LoD (TCID <sub>50</sub> /mL)	LoD (TCID <sub>50</sub> /Swab)*	LoD (Genome Equivalents/mL)	LoD (Genome Equivalents/Swab)*
A/Puerto Rico/8/34	A/H1N1	1.25 x 10 <sup>5</sup>	1.25 x 10 <sup>3</sup>	7.59 x 10 <sup>5</sup>	7.59 x 10 <sup>3</sup>
A/Solomon Islands/3/2006	A/H1N1	2.34 x 10 <sup>3</sup>	2.34 x 10 <sup>1</sup>	1.12 x 10 <sup>5</sup>	1.12 x 10 <sup>3</sup>
B/Malaysia/2506/2004	B Victoria lineage	1.64 x 10 <sup>3</sup>	1.64 x 10 <sup>1</sup>	1.95 x 10 <sup>4</sup>	1.95 x 10 <sup>2</sup>
B/Brisbane/60/2008	B Victoria lineage	3.98 x 10 <sup>2</sup>	3.98 x 10 <sup>0</sup>	8.08 x 10 <sup>4</sup>	8.08 x 10 <sup>2</sup>

\*Note: 10 ul of each virus dilution was coated onto a swab

*e. Analytical Reactivity:*

## Analytical Reactivity Testing

An analytical reactivity (inclusivity) study was performed to determine whether the Alere™ i Influenza A & B assay is able to detect a variety of influenza A and B strains that represent temporal and geographic diversity.

Vender provided stocks of influenza A and B strains were diluted in UTM to generate virus dilutions for testing. The concentration (in TCID<sub>50</sub>/mL, CEID<sub>50</sub>/mL or EID<sub>50</sub>/mL) for each was determined by standard virologic method. The concentration for each dilution (in genome equivalents/mL) was also assessed using laboratory developed and validated influenza A and influenza B quantitative real-time PCR assays.

Contrived swab specimens were prepared by coating 10 µl of virus dilution onto each swab. Swab coating was carried out by lightly scratching the swab surface with the pipette tip as virus dilution was pipetted to ensure liquid was absorbed into the swab tip. The contrived swab samples were tested according to the product instructions.

The starting dilution concentration selected for testing in this study was higher than the established LoDs in the Limit of Detection studies. Each starting dilution per strain was tested in triplicates initially. If the initial testing concentration tested positive for all three replicates, the strain was further diluted 10-fold and tested in triplicates until at least one out three replicates generated a negative result. When a negative result was obtained, additional 2-fold dilutions were tested, starting from the highest dilution that produced 100% (3/3) positive results. A concentration level was considered “reactive/positive” in this study for all but one strain tested (i.e., B/Texas/06/2011 – see footnote C under Table 11 below) if all three replicates generated a positive result for the expected influenza virus. A concentration level was considered “non-reactive/negative” for all but one strain tested (i.e., B/Texas/06/2011 – see footnote C under Table 11 below) in this study if at least one out of the three replicates generated a negative result for the expected influenza virus.

The Alere™ i Influenza A & B assay was evaluated with 17 Influenza A and nine Influenza B strains. Influenza A strains included four pre-2009 pandemic seasonal A/H1N1 strains, eight seasonal H3N2 strains, three 2009 A/H1N1 pandemic strains, one A/H3N2v strain, and one avian influenza A/H7N9 virus (detected in China in 2013). Influenza B strains included four Victoria lineage strains and five Yamagata lineage strains. The Alere™ i Influenza A & B assay detected all strains tested at the concentrations indicated in Table 11 below:

**Table 11: Analytical Reactivity Study Results**

Influenza Strain	Influenza A Subtype or Influenza B Genetic Lineage	Test Concentration (in TCID <sub>50</sub> or Genome Equivalents , unless indicated otherwise)				Flu A Result (n=3, unless indicate otherwise)	Flu B Result (n=3, unless indicate otherwise)
		TCID <sub>50</sub> /mL	TCID <sub>50</sub> /Swab*	Genome Equivalents /mL	Genome Equivalents/Swab*		
A/New Caledonia/20/1999 <sup>a</sup>	A/H1N1	9.19 x 10 <sup>5</sup>	9.19 x 10 <sup>3</sup>	4.09 x 10 <sup>6</sup>	4.09 x 10 <sup>4</sup>	+	-
A/New Jersey/8/76 <sup>a</sup>	A/H1N1	3.41 x 10 <sup>1</sup>	3.41 x 10 <sup>-1</sup>	1.52 x 10 <sup>5</sup>	1.52 x 10 <sup>3</sup>	+	-
A/Brisbane/59/2007 <sup>a</sup>	A/H1N1	2.11 x 10 <sup>4</sup>	2.11 x 10 <sup>2</sup>	3.39 x 10 <sup>5</sup>	3.39 x 10 <sup>3</sup>	+	-
A/WSN/33 <sup>a</sup>	A/H1N1	2.11 x 10 <sup>2</sup>	2.11 x 10 <sup>0</sup>	2.43 x 10 <sup>5</sup>	2.43 x 10 <sup>3</sup>	+	-
A/Port Chalmers/1/73	A/H3N2	4.22 x 10 <sup>4</sup>	4.22 x 10 <sup>2</sup>	1.31 x 10 <sup>6</sup>	1.31 x 10 <sup>4</sup>	+	-
A/Hong Kong/8/68	A/H3N2	7.03 x 10 <sup>0</sup>	7.03 x 10 <sup>-2</sup>	2.70 x 10 <sup>5</sup>	2.70 x 10 <sup>3</sup>	+	-
A/Aichi/2/68	A/H3N2	2.08 x 10 <sup>5</sup>	2.08 x 10 <sup>3</sup>	7.47 x 10 <sup>6</sup>	7.47 x 10 <sup>4</sup>	+	-
A/Victoria/3/75	A/H3N2	3.68 x 10 <sup>5</sup>	3.68 x 10 <sup>3</sup>	3.39 x 10 <sup>6</sup>	3.39 x 10 <sup>4</sup>	+	-
A/Wisconsin/67/2005	A/H3N2	6.81 x 10 <sup>4</sup>	6.81 x 10 <sup>2</sup>	2.57 x 10 <sup>6</sup>	2.57 x 10 <sup>4</sup>	+	-

A/Brisbane/10/2007	A/H3N2	$3.16 \times 10^2$	$3.16 \times 10^0$	$3.37 \times 10^5$	$3.37 \times 10^3$	+	-
A/Texas/50/2012 <sup>b</sup>	A/H3N2	$2.50 \times 10^0$	$2.50 \times 10^{-2}$	$6.35 \times 10^3$	$6.35 \times 10^1$	+	-
A/Victoria/361/2011	A/H3N2	$1.56 \times 10^1$	$1.56 \times 10^{-1}$	$3.53 \times 10^5$	$3.53 \times 10^3$	+	-
A/California/4/2009	A/H1N1 (pdm)	$1.47 \times 10^4$	$1.47 \times 10^{-2}$	$1.07 \times 10^6$	$1.07 \times 10^4$	+	-
A/Maryland/04/2011	A/H1N1 (pdm)	$7.88 \times 10^4$	$7.88 \times 10^{-2}$	$3.81 \times 10^6$	$3.81 \times 10^4$	+	-
A/New York/18/2009	A/H1N1 (pdm)	$1.25 \times 10^2$	$1.25 \times 10^0$	$9.16 \times 10^5$	$9.16 \times 10^3$	+	-
A/Anhui/1/2013 (Inactivated) <sup>a</sup>	A/H7N9 (Detected in China in 2013)	$4.00 \times 10^6$ EID <sub>50</sub> /mL	$4.00 \times 10^4$ EID <sub>50</sub> /Swab	$1.72 \times 10^6$	$1.72 \times 10^4$	+	-
A/Indiana/10/2011 <sup>a</sup>	A/H3N2v	$2.00 \times 10^8$ EID <sub>50</sub> /mL	$2.00 \times 10^6$ EID <sub>50</sub> /Swab	$2.00 \times 10^6$	$2.00 \times 10^4$	+	-
B/Lee/40	Victoria Lineage	$5.00 \times 10^1$ CEID <sub>50</sub> /mL	$5.00 \times 10^{-1}$ CEID <sub>50</sub> /Swab	$5.40 \times 10^4$	$5.40 \times 10^2$	-	+
B/Victoria/504/2000	Victoria Lineage	$1.19 \times 10^3$	$1.19 \times 10^1$	$8.29 \times 10^4$	$8.29 \times 10^2$	-	+
B/Nevada/03/2011	Victoria Lineage	$1.75 \times 10^3$	$1.75 \times 10^1$	$1.13 \times 10^5$	$1.13 \times 10^3$	-	+
B/Montana/05/2012	Victoria Lineage	$9.00 \times 10^1$	$9.00 \times 10^{-1}$	$2.55 \times 10^4$	$2.55 \times 10^2$	-	+
B/Maryland/1/59	Yamagata Lineage	$8.51 \times 10^2$	$8.51 \times 10^0$	$1.13 \times 10^5$	$1.13 \times 10^3$	-	+
B/Russia/69 <sup>c</sup>	Yamagata Lineage	$4.44 \times 10^1$	$4.44 \times 10^{-1}$	$2.96 \times 10^6$	$2.96 \times 10^4$	-	+
B/Wisconsin/01/2010 <sup>d</sup>	Yamagata Lineage	$3.68 \times 10^4$	$3.68 \times 10^2$	$1.16 \times 10^6$	$1.16 \times 10^4$	-	+
B/Massachusetts/2/2012	Yamagata Lineage	$6.25 \times 10^1$	$6.25 \times 10^{-1}$	$2.28 \times 10^5$	$2.28 \times 10^3$	-	+
B/Texas/06/2011 <sup>d</sup>	Yamagata Lineage	$2.89 \times 10^5$	$6.25 \times 10^3$	$2.00 \times 10^6$	$2.00 \times 10^4$	-	+

\*Note: 10 ul of each virus dilution was coated onto a swab

<sup>a</sup> Although this test has been shown to detect A/H1N1 (pre-2009 pandemic), A/H7N9 (detected in China in 2013) and A/H3N2v viruses cultured from positive human respiratory specimens, the performance characteristics of this device with clinical specimens that are positive for the A/H1N1 (pre-2009 pandemic), A/H7N9 (detected in China in 2013) and A/H3N2v viruses have not been established.

<sup>b</sup> Influenza A/Texas/50/2012 lowest level in which 3/3 replicates were positive is approximately 0.01 to 0.06 x the LoD (as comparing to the Genome Equivalents/Swab values generated in the LoD with simulated matrix study testing A/Puerto Rico/8/34 and A/Solomon Islands/3/2006, respectively). Based on *in silico* analysis, no polymorphism within the target sequence was revealed that can account for the observed increased reactivity. A plausible cause for the observation of increased reactivity is the impact that polymorphisms found outside of the target region can have on amplification and detection. Influenza virus mutates at a very high frequency, and as a result conservation across the entirety of the Influenza genome is very limited. The Influenza genome is RNA, and RNA readily folds into complex secondary and tertiary structures based on the primary sequence of the nucleic acid. Polymorphisms introduced into the primary structure can result in changes in secondary/tertiary structure, potentially making specific regions more or less accessible to nucleic acid binding moieties. The Alere™ i Influenza A & B test requires access to the target nucleic acid sequence in order to initiate amplification and detection, therefore it is plausible that specific changes to the genome's primary structure (polymorphisms) may result in secondary/tertiary structural changes that could potentially increase (more accessible to nucleic acid binding moieties) or decrease (less accessible to nucleic acid binding moieties) the test's analytical reactivity.

<sup>c</sup> Influenza B/Russia/69 lowest level in which 3/3 replicates were positive is approximately 40 to 150 x the LoD (as comparing to the Genome Equivalents/Swab values generated in the LoD with simulated clinical matrix study testing B/Malaysia/2506/2004 and B/Brisbane/60/2008, respectively). A polymorphism within segment PA of the Influenza B genome was identified at a position which is 4 nucleotides from the 3'-end of template 2. This G to A polymorphism results in a G/C (product/template) match to an A/C (product/template) mismatch. An A/C mismatch is determined to be moderately destabilizing, and coupled to its position only 4 nucleotides from the 3'-end of the template 2 recognition region, its impact on annealing is potentially great. The frequency of this G to A polymorphism is determined to be very low. In analyzing the strains present in the NCBI Influenza Virus Resource database from 2/2005 to 3/2014 (N=986), no strains contained this polymorphism, suggesting that it has not been circulating for an extended period of time.

<sup>d</sup> Influenza B/Wisconsin/01/2010 lowest level in which 3/3 replicates were positive is approximately 15 to 60 x the LoD, and Influenza B/Texas/06/2011 lowest level in which at least 1/3 replicates were positive is approximately 25 to 100 x the LoD (as comparing to the Genome Equivalents/Swab values generated in the LoD with simulated clinical matrix study testing B/Malaysia/2506/2004 and B/Brisbane/60/2008, respectively). A single G to A polymorphism within segment PA of the Influenza B genome was identified at a position which is 5 nucleotides from the 3'-end of the molecular beacon annealing region in both strains. The G to A polymorphism results in a C/G match to a C/A mismatch between the molecular beacon and product 1. The C/A mismatch is determined to be moderately destabilizing that can significantly reduce assay sensitivity. An assessment of what impact this polymorphism would have on the melting temperature (T<sub>m</sub>) of the molecular beacon/product 1 annealing was performed and the results showed a T<sub>m</sub> drop from 62.3°C to 55.6°C, just below the assay running temperature. This suggests that annealing would occur, but at a greatly reduced level, with a concomitant loss of assay sensitivity. The frequency of this G to A polymorphism is found at a frequency of approximately 5% within the NCBI Influenza Virus Resource database covering the time frame from 2/2005 through 3/2014.

Supplemental Reactivity Information for Influenza Strains of Human, Swine and Avian Origin (*in silico* analysis):

In addition to laboratory testing, bioinformatics resources were also used to predict reactivity of additional influenza A strains of human, swine, or avian origin with the Alere™ i Influenza A & B.

Predicted reactivity of the Alere™ i Influenza A&B assay with additional influenza A strains of human, avian, and swine origin was generated using a bioinformatics approach. Assay template and probe sequences were aligned with GenBank sequences corresponding to the appropriate gene targets and reactivity were predicated based on the number and location of mismatches in the targeted region. For each subtype, multiple GenBank IDs were evaluated, corresponding to the gene segments targeted by the Alere™ i Influenza A&B assay (PB2 gene for Flu A and PA gene for Flu B). The strains listed in Table 12 below are predicted to react with the Alere™ i Influenza A&B assay as indicated:

**Table 12: Predicted Alere i Influenza A&B Reactivity for Influenza A**

Origin	Subtype	Strain	GenBank ID	Predicted Alere i Influenza A&B Reactivity
Human	H1N1pdm	A/Washington/24/2012(H1N1pdm)	EPI439208	Influenza A (100% sequence homology)
	H1N1	A/Fujian-Gulou/1896/2009(H1N1)	EPI244244	Influenza A (100% sequence homology)
	H3N2v	A/Pennsylvania/09/2011(H3N2v)	JN655532	Influenza A (100% sequence homology)
		A/Pennsylvania/13/2012(H3N2v)	EPI392732	Influenza A (100% sequence homology)
		A/Indiana/123/2012(H3N2v)	EPI397375	Influenza A (100% sequence homology)
		A/Indiana/99/2012(H3N2v)	EPI397401	Influenza A (100% sequence homology)
		A/Indiana/114/2012(H3N2v)	EPI397910	Influenza A (100% sequence homology)
		A/Indiana/21/2013(H3N2v)	EPI467215	Influenza A (100% sequence homology)
		A/Ohio/13/2012(H3N2v)	EPI381854	Influenza A (100% sequence homology)
		A/Ohio/74/2012(H3N2v)	EPI393779	Influenza A (100% sequence homology)
	H2N2	A/Japan/305-MA12/1957(H2N2)	CY087807	Influenza A (100% sequence homology)
		A/Montevideo/2208/1967(H2N2)	CY125829	Influenza A (100% sequence homology)
		A/Ann Arbor/7/1967(H2N2)	CY125845	Influenza A (100% sequence homology)
		A/Ann Arbor/6/60(H2N2)	AY209938	Influenza A (100% sequence homology)
		A/Albany/1/65(H2N2)	AY209942	Influenza A (100% sequence homology)
		A/Berkeley/1/66(H2N2)	AY209946	Influenza A (100% sequence homology)
		A/Albany/3/1958(H2N2)	CY021076	Influenza A (100% sequence homology)
		A/Ann Arbor/6/60(H2N2)	M23970	Influenza A (100% sequence homology)
		A/North Carolina/1/1968(H2N2)	CY031594	Influenza A (100% sequence homology)
	H5N1	A/Cambodia/Q0321176/2006(H5N1)	HQ200422	Influenza A (100% sequence homology)
		A/Shanghai/1/2006(H5N1)	CY098631	Influenza A (100% sequence homology)
A/Anhui/1/2006(H5N1)		CY098665	Influenza A (100% sequence homology)	

			homology)
		A/Guangdong/1/2008(H5N1)	CY098699 Influenza A (100% sequence homology)
		A/Beijing/1/2009(H5N1)	CY098706 Influenza A (100% sequence homology)
		A/Hong Kong/481/97(H5N1)	AF115290 Influenza A (100% sequence homology)
		A/Hong Kong/481/97(H5N1)	AF258837 Influenza A (100% sequence homology)
		A/Cambodia/V0606311/2011(H5N1)	JN588929 Influenza A (100% sequence homology)
	H7N7	A/Netherlands/219/03 (H7N7)	AY342413 Influenza A (100% sequence homology)
		A/England/AV877/1996 (H7N7)	GU053116 Influenza A (100% sequence homology)
	H9N2	A/Hong Kong/1073/99 (H9N2)	AF258835 Influenza A (100% sequence homology)
		A/Hong Kong/1073/99 (H9N2)	AJ404630 Influenza A (100% sequence homology)
		A/Hong Kong/1074/1997 (H9N2)	GU053185 Influenza A (100% sequence homology)
		A/Hong Kong/33982/2009 (H9N2)	KF188313 Influenza A (sequence mismatch – moderate impact on reactivity)
Avian	H2N2	A/duck Hong Kong/319/1978(H2N2)	CY005545 Influenza A (100% sequence homology)
		A/mallard New York/6750/1978(H2N2)	CY130005 Influenza A (sequence mismatch – strong impact on reactivity)
	H5N1	A/Brown-headed Gull Qinghai/3/05(H5N1)	DQ095756 Influenza A (100% sequence homology)
		A/mallard California/2531P/2011(H5N1)	CY133956 Influenza A (100% sequence homology)
	H5N2	A/wild bird Minnesota/460613/2006(H5N2)	GQ923532 Influenza A (100% sequence homology)
		A/mallard Arkansas/473507-9/2006(H5N2)	GQ923348 Influenza A (100% sequence homology)
		A/northern pintail Utah/469648/2006(H5N2)	GQ923300 Influenza A (100% sequence homology)
		A/turkey MB/FAV11/2010(H5N2)	JX570826 Influenza A (100% sequence homology)
		A/mallard MD/865/2002(H5N2)	EU980493 Influenza A (100% sequence homology)
		A/American black-duck Illinois/08OS2688/2008(H5N2)	CY079459 Influenza A (100% sequence homology)
		A/spur-winged goose Nigeria/2/2008(H5N2)	FR771844 Influenza A (100% sequence homology)
		A/tundra swan Shimane/3211A001/2011(H5N2)	AB830625 Influenza A (100% sequence homology)
		A/Muscovy duck New York/62095-1/2006(H5N2)	CY036041 Influenza A (sequence mismatch – strong impact on reactivity)
	H7N2	A/Chicken New York/6777-3/97(H7N2)	EU084916 Influenza A (100% sequence homology)
		A/guinea fowl New York/23165-12/2005(H7N2)	CY031730 Influenza A (100% sequence homology)
		A/chicken NY/8030-2/1997(H7N2)	EU742937 Influenza A (100% sequence homology)
		A/Guinea fowl New York/11646-3/2005(H7N2)	CY035961 Influenza A (100% sequence homology)
		A/chicken New York/88291-14/2005(H7N2)	CY037054 Influenza A (100% sequence homology)
		A/turkey New York/88291-7/2005(H7N2)	CY095131 Influenza A (100% sequence homology)
		A/chicken Hebei/1/2002(H7N2)	AY724251 Influenza A (sequence mismatch – minimal impact on reactivity)
		A/goose NY/8600-3/1998(H7N2)	EU743272 Influenza A (sequence mismatch – strong impact on reactivity)
H7N3	A/turkey Italy/3620/2003 (H7N3)	CY021364 Influenza A (100% sequence	

			homology)
	A/turkey Italy/2685/2003 (H7N3)	CY028675	Influenza A (100% sequence homology)
	A/ruddy turnstone NJ/65/1985 (H7N3)	CY004412	Influenza A (100% sequence homology)
	A/pheasant Minnesota/917/1980 (H7N3)	CY005895	Influenza A (100% sequence homology)
	A/American black duck NB/2538/2007 (H7N3)	EU500849	Influenza A (100% sequence homology)
	A/ruddy turnstone Delaware Bay/121/2007 (H7N3)	CY036790	Influenza A (100% sequence homology)
	A/turkey Minnesota/1138/1980 (H7N3)	CY041838	Influenza A (100% sequence homology)
	A/shorebird Delaware Bay/552/2006 (H7N3)	CY103022	Influenza A (100% sequence homology)
	A/American green-winged teal Mississippi/11OS250/2011 (H7N3)	CY133676	Influenza A (100% sequence homology)
	A/ruddy turnstone New Jersey/175/2006 (H7N3)	GU051674	Influenza A (100% sequence homology)
H7N7	A/northern shoveler Mississippi/09OS643/2009 (H7N7)	CY079419	Influenza A (100% sequence homology)
	A/ruddy turnstone Delaware/892/2006 (H7N7)	GQ257395	Influenza A (100% sequence homology)
	A/wild bird Korea/A72/10 (H7N7)	JN244116	Influenza A (100% sequence homology)
	A/American black duck New Brunswick/00344/2010 (H7N7)	CY138776	Influenza A (100% sequence homology)
	A/common teal Hong Kong/MPL634/2011 (H7N7)	KF260797	Influenza A (100% sequence homology)
	A/chicken Wenzhou/299/2013 (H7N7)	KF260930	Influenza A (sequence mismatch – strong impact on reactivity)
	A/chicken Wenzhou/645/2013 (H7N7)	KF260919	Influenza A (sequence mismatch – strong impact on reactivity)
H9N2	A/avian Saudi Arabia/910135/2006 (H9N2)	GU050294	Influenza A (100% sequence homology)
	A/chicken Guangdong/TS/2004 (H9N2)	JQ639775	Influenza A (100% sequence homology)
	A/bird Guangxi/A1/2006 (H9N2)	EU086301	Influenza A (100% sequence homology)
	A/Bewick's swan Netherlands/5/2007 (H9N2)	CY041281	Influenza A (100% sequence homology)
	A/Rock Pigeon Vietnam/6/2009 (H9N2)	AB753198	Influenza A (100% sequence homology)
	A/chicken Bangladesh/5209/2009 (H9N2)	KC757995	Influenza A (100% sequence homology)
	A/chicken China/Guangxi/14/2000 (H9N2)	DQ485213	Influenza A (sequence mismatch – moderate impact on reactivity)
H10N7	A/mallard Ohio/99/1989 (H10N7)	CY017788	Influenza A (100% sequence homology)
	A/chicken Germany/n/1949 (H10N7)	CY014678	Influenza A (100% sequence homology)
	A/mallard Minnesota/g-00056/2007(H10N7)	CY064080	Influenza A (100% sequence homology)
	A/American green-winged teal Illinois/08OS2713/2008(H10N7)	CY079475	Influenza A (100% sequence homology)
	A/American black duck New Brunswick/00906/2010(H10N7)	CY139392	Influenza A (100% sequence homology)
	A/American black duck New Brunswick/00906/2010(H10N7)	CY139392	Influenza A (100% sequence homology)
	A/mallard Alberta/289/2007(H10N7)	CY137721	Influenza A (100% sequence homology)
	A/mallard California/2954/2012(H10N7)	CY157237	Influenza A (100% sequence homology)
	A/American green-winged teal Mississippi/12OS405/2012(H10N7)	CY167247	Influenza A (100% sequence homology)
	A/American green-winged teal	CY166193	Influenza A (100% sequence

		Wisconsin/11OS3499/2011(H10N7)		homology)
		A/green-winged teal California/10197/2008(H10N7)	CY093893	Influenza A (100% sequence homology)
Swine	H5N1	A/swine Fujian/1/2003(H5N1)	AY747616	Influenza A (100% sequence homology)
Swine	H7N2	A/swine KU/16/2001(H7N2)	CY067683	Influenza A (100% sequence homology)

*f. Analytical Specificity:*

The Alere™ i Influenza A & B was performed on samples containing high levels of non-target organisms using the Alere™ i instrument, to demonstrate the specificity of the assay for detection of influenza A and B. The analytical specificity study testing included 37 bacteria, 15 viruses and 1 yeast common to the nasal cavity/nasopharynx. All bacterial stocks and the yeast stock were obtained from ATCC, with the exception of *Chlamydia pneumoniae* (obtained from the University of Washington). All viral stocks were obtained from ATCC, with the exception of human metapneumovirus (obtained from the Zeptomatrix Corporation).

The majority of the bacteria and yeast were re-grown by Alere Scarborough, Inc. using standard microbiological culture techniques, reconstituted into saline dilution buffer and then quantified by spectrophotometer (cells/ml) versus a saline dilution buffer blank. The *Mycoplasma pneumoniae* received from the vendor was resuspended by Alere Scarborough, Inc. per vendor instructions and tested. Re-growth did not yield enough for harvesting and reconstitution in saline. The *Chlamydia pneumoniae* was supplied by the University of Washington as an irradiated stock and it was, therefore, not re-grown and titered by Alere Scarborough, Inc. prior to testing. The titers of the quantified viruses indicated on certificates of analysis from vendors were not confirmed by Alere Scarborough, Inc. prior to testing.

Swab samples were prepared using the bacterial, viral and yeast strains at the stock concentration indicated in Table 13 below. 75 µl of each organism was dispensed into each of three 1.5 ml tubes. A swab was placed into each tube and stirred to absorb all the liquid into the swab head. The samples were tested in triplicate in Alere™ i Influenza A & B, according to product instructions for testing direct swabs.

**Table 13: Alere™ i Influenza A & B Analytical Specificity Study Results**

Bacteria	Concentration (Cells/mL, Unless Indicated Otherwise)	Concentration (Cells/Reaction, Unless Indicated Otherwise)	Influenza A Result			Influenza B Result		
<i>Acinetobacter calcoaceticus</i>	2.00 x 10 <sup>10</sup>	6.00 x 10 <sup>7</sup>	-	-	-	-	-	-
<i>Bacteroides fragilis</i>	4.00 x 10 <sup>9</sup>	1.20 x 10 <sup>7</sup>	-	-	-	-	-	-
<i>Bordetella pertussis</i>	1.00 X 10 <sup>10</sup>	3.00 x 10 <sup>7</sup>	-	-	-	-	-	-
<i>Chlamydia pneumoniae</i>	9.00 x 10 <sup>8</sup> IFU/ml	2.70 x 10 <sup>6</sup> IFU/ml	-	-	-	-	-	-
<i>Corynebacterium diphtheria</i>	1.00 X 10 <sup>10</sup>	3.00 x 10 <sup>7</sup>	-	-	-	-	-	-
<i>Enterococcus faecalis</i>	1.00 X 10 <sup>10</sup>	3.00 x 10 <sup>7</sup>	-	-	-	-	-	-
<i>Escherichia coli</i>	2.00 X 10 <sup>10</sup>	6.00 x 10 <sup>7</sup>	-	-	-	-	-	-
<i>Gardnerella vaginalis</i>	3.00 X 10 <sup>10</sup>	9.00 x 10 <sup>7</sup>	-	-	-	-	-	-
<i>Haemophilus influenzae</i>	1.00 X 10 <sup>10</sup>	3.00 x 10 <sup>7</sup>	-	-	-	-	-	-
<i>Klebsiella pneumoniae</i>	2.00 X 10 <sup>9</sup>	6.00 x 10 <sup>6</sup>	-	-	-	-	-	-
<i>Lactobacillus casei</i>	7.00 X 10 <sup>9</sup>	2.10 x 10 <sup>7</sup>	-	-	-	-	-	-
<i>Lactobacillus plantarum</i>	5.00 X 10 <sup>9</sup>	1.50 x 10 <sup>7</sup>	-	-	-	-	-	-
<i>Legionella pneumophila</i>	1.00 X 10 <sup>10</sup>	3.00 x 10 <sup>7</sup>	-	-	-	-	-	-

<i>Listeria monocytogenes</i>	7.00 X 10 <sup>9</sup>	2.10 x 10 <sup>7</sup>	-	-	-	-	-	-
<i>Moraxella/Branhamella catarrhalis</i>	2.00 X 10 <sup>9</sup>	6.00 x 10 <sup>6</sup>	-	-	-	-	-	-
<i>Mycobacterium avium</i>	1.00 X 10 <sup>9</sup>	3.00 x 10 <sup>6</sup>	-	-	-	-	-	-
<i>Mycobacterium intracellulare</i>	1.00 X 10 <sup>9</sup>	3.00 x 10 <sup>6</sup>	-	-	-	-	-	-
<i>Mycobacterium tuberculosis</i>	1.00 X 10 <sup>8</sup>	3.00 x 10 <sup>5</sup>	-	-	-	-	-	-
<i>Mycoplasma pneumoniae</i>	1.50 x 10 <sup>3</sup> CFU/ml	4.50 x 10 <sup>0</sup> CFU/ml	-	-	-	-	-	-
<i>Neisseria gonorrhoeae</i>	2.00 X 10 <sup>10</sup>	6.00 x 10 <sup>7</sup>	-	-	-	-	-	-
<i>Neisseria meningitidis</i>	8.00 x 10 <sup>9</sup>	2.40 x 10 <sup>7</sup>	-	-	-	-	-	-
<i>Neisseria sicca</i>	1.00 X 10 <sup>10</sup>	3.00 x 10 <sup>7</sup>	-	-	-	-	-	-
<i>Neisseria subflava</i>	2.00 X 10 <sup>10</sup>	6.00 x 10 <sup>7</sup>	-	-	-	-	-	-
<i>Proteus vulgaris</i>	7.00 X 10 <sup>8</sup>	2.10 x 10 <sup>6</sup>	-	-	-	-	-	-
<i>Pseudomonas aeruginosa</i>	1.00 X 10 <sup>10</sup>	3.00 x 10 <sup>7</sup>	-	-	-	-	-	-
<i>Serratia marcescens</i>	3.00 X 10 <sup>9</sup>	9.00 x 10 <sup>6</sup>	-	-	-	-	-	-
<i>Staphylococcus aureus</i>	3.00 X 10 <sup>10</sup>	9.00 x 10 <sup>7</sup>	-	-	-	-	-	-
<i>Staphylococcus epidermidis</i>	5.00 x 10 <sup>10</sup>	1.50 x 10 <sup>8</sup>	-	-	-	-	-	-
<i>Streptococcus, Group A</i>	3.00 x 10 <sup>9</sup>	9.00 x 10 <sup>6</sup>	-	-	-	-	-	-
<i>Streptococcus, Group B</i>	2.00 x 10 <sup>10</sup>	6.00 x 10 <sup>7</sup>	-	-	-	-	-	-
<i>Streptococcus, Group C</i>	1.00 x 10 <sup>10</sup>	3.00 x 10 <sup>7</sup>	-	-	-	-	-	-
<i>Streptococcus, Group F</i>	3.00 x 10 <sup>9</sup>	9.00 x 10 <sup>6</sup>	-	-	-	-	-	-
<i>Streptococcus, Group G</i>	2.00 x 10 <sup>10</sup>	6.00 x 10 <sup>7</sup>	-	-	-	-	-	-
<i>Streptococcus mutans</i>	1.00 x 10 <sup>10</sup>	3.00 x 10 <sup>7</sup>	-	-	-	-	-	-
<i>Streptococcus pneumoniae</i>	2.00 X 10 <sup>9</sup>	6.00 x 10 <sup>6</sup>	-	-	-	-	-	-
<i>Streptococcus salivaris</i>	6.00 X 10 <sup>9</sup>	1.80 x 10 <sup>7</sup>	-	-	-	-	-	-
<i>Streptococcus sanguis</i>	6.00 X 10 <sup>9</sup>	1.80 x 10 <sup>7</sup>	-	-	-	-	-	-
<b>Yeast</b>	<b>Concentration (Cells/mL, Unless Indicated Otherwise)</b>	<b>Concentration (Cells/Reaction, Unless Indicated Otherwise)</b>	<b>Influenza A Result</b>			<b>Influenza B Result</b>		
<i>Candida albicans</i>	2.00 X 10 <sup>8</sup>	6.00 x 10 <sup>5</sup>	-	-	-	-	-	-
<b>Virus</b>	<b>Concentration (TCID<sub>50</sub>/mL, Unless Indicated Otherwise)</b>	<b>Concentration (TCID<sub>50</sub>/Reaction, Unless Indicated Otherwise)</b>	<b>Influenza A Result</b>			<b>Influenza B Result</b>		
Adenovirus 1	1.58 X 10 <sup>7</sup>	4.47 X 10 <sup>4</sup>	-	-	-	-	-	-
Adenovirus 7	2.81 X 10 <sup>6</sup>	8.43 X 10 <sup>3</sup>	-	-	-	-	-	-
Cytomegalovirus	8.89 x 10 <sup>5</sup>	2.67 X 10 <sup>3</sup>	-	-	-	-	-	-
Coronavirus 229E	2.81 x 10 <sup>4</sup>	8.43 X 10 <sup>1</sup>	-	-	-	-	-	-
Coronavirus OC43	2.81 x 10 <sup>5</sup>	8.43 X 10 <sup>2</sup>	-	-	-	-	-	-
Enterovirus/Coxsackievirus B4	2.81 x 10 <sup>7</sup>	8.43 X 10 <sup>4</sup>	-	-	-	-	-	-
Epstein Barr Virus*	Not Available	Not Available	-	-	-	-	-	-
Human metapneumovirus	2.19 x 10 <sup>6</sup>	6.57 X 10 <sup>3</sup>	-	-	-	-	-	-
Parainfluenza 1	1.58 x 10 <sup>8</sup> CEID <sub>50</sub> /mL	4.74 X 10 <sup>5</sup> CEID <sub>50</sub> /ml	-	-	-	-	-	-
Parainfluenza 2	2.81 x 10 <sup>6</sup>	8.43 X 10 <sup>3</sup>	-	-	-	-	-	-
Parainfluenza 3	2.81 x 10 <sup>7</sup>	8.43 X 10 <sup>4</sup>	-	-	-	-	-	-
Measles (Edmonston)	1.58 x 10 <sup>4</sup>	4.74 X 10 <sup>1</sup>	-	-	-	-	-	-
Mumps (Enders)	1.58 x 10 <sup>6</sup>	4.74 X 10 <sup>3</sup>	-	-	-	-	-	-
Respiratory Syncytial Virus Type B	8.89 x 10 <sup>5</sup>	2.67 X 10 <sup>3</sup>	-	-	-	-	-	-
Rhinovirus 1A	1.58 x 10 <sup>7</sup>	4.74 X 10 <sup>4</sup>	-	-	-	-	-	-

\* Epstein Barr Virus was received from ATCC and the titer was listed as N/A on the Certificate of Analysis (CoA).

None of the bacteria, yeast, or viruses tested cross-reacted in the Alere™ i Influenza A & B at these concentrations.

*g. Potentially Interfering Substances:*

An analytical study was performed to assess the potential interference effects of 25 substances naturally present in respiratory specimens or that may be artificially introduced into the nasal cavity/nasopharynx.

Dilutions of Influenza A/Puerto Rico/8/34 (Lot D1002B, inactivated), Influenza A/Solomon Islands/3/2006 (Lot D1326A), Influenza B/Malaysia/2506/2004 (Lot B1012C, inactivated), and Influenza B/Brisbane/60/2008 (Lot D1326B) were used in this study to prepare contrived positive samples for testing. The concentrations of the viral stocks (in TCID<sub>50</sub>/mL) were determined by standard virologic method. The targeted concentration for each dilution (in genome equivalents/mL) was also assessed using laboratory developed and validated influenza A and influenza B quantitative real-time PCR assays.

Using UTM as the diluent, positive samples were prepared using the influenza A and B strains (described above) diluted to a concentration approximately 2 to 3 times the respective limit of detection (LoD) levels and the various interfering substances at the concentrations specified in Table 14 below. (Note: During preparation of the B/Brisbane/60/2008 dilution a calculation error was made. This strain was inadvertently diluted to a concentration approximately 5 to 10 times the LoD level.) Each influenza strain, combined with each interfering substance, was evaluated separately. Negative samples were prepared using UTM and the various interfering substances at the concentrations specified in Table 10 below.

50 µl of each positive and negative sample with interfering substance was added to the Sample Receiver and pipetted up and down to mix. Each sample was tested in triplicates per the product instructions.

The influenza A and influenza B positive samples were tested as follows:

- Influenza A/Puerto Rico/8/34:  $3.30 \times 10^3$  TCID<sub>50</sub> ( $2.00 \times 10^4$  genome equivalents) in the Sample Receiver; [Note: 50 µl of  $6.59 \times 10^4$  TCID<sub>50</sub>/mL ( $4.00 \times 10^5$  genome equivalents/mL) virus dilution contain an interfering substance added to the Sample Receiver]
- Influenza A/Solomon Islands/3/2006:  $4.10 \times 10^2$  TCID<sub>50</sub> ( $2.80 \times 10^3$  genome equivalents) in the Sample Receiver; [Note: 50 µl of  $8.19 \times 10^3$  TCID<sub>50</sub>/mL ( $5.60 \times 10^4$  genome equivalents/mL) virus dilution contain an interfering substance added to the Sample Receiver]
- Influenza B/Malaysia/2506/2004:  $2.70 \times 10^1$  TCID<sub>50</sub> ( $5.00 \times 10^2$  genome equivalents) in the Sample Receiver; [Note: 50 µl of  $8.40 \times 10^2$  TCID<sub>50</sub>/mL ( $1.00 \times 10^4$  genome equivalents/mL) virus dilution contain an interfering substance added to the Sample Receiver]
- Influenza B/Brisbane/60/2008:  $7.50 \times 10^1$  TCID<sub>50</sub> ( $4.06 \times 10^3$  genome equivalents) in the Sample Receiver; [Note: 50 µl of  $1.50 \times 10^3$  TCID<sub>50</sub>/mL ( $8.12 \times 10^4$  genome equivalents/mL) virus dilution contain an interfering substance added to the Sample Receiver]

The Interfering Substances Study Results are presented in Table 14 below.

**Table 14: Interfering Substances Study Results**

Substance Tested	Concentration Tested	Negative Samples		Influenza Positive Samples			
		Flu A Result	Flu B Result	Flu A Result (A/PR/8/34)	Flu A Result A/Solomon Islands/3/2006)	Flu B Result (B/Malaysia/2506/2004)	Flu B Result (B/Brisbane/6/2008)
Mucin	20 µg /mL	-	-	+	+	+	+
		-	-	+	+	+	+
		-	-	+	+	+	+
Whole Blood	50 µl/mL	-	-	+	+	+	+
		-	-	+	+	+	+
		-	-	+	+	+	+
Sinus Buster Nasal Spray	200 µl/mL	-	-	+	+	+	+
		-	-	+	+	+	+
		-	-	+	+	+	+
NeoSynephrine Cold & Sinus Extra Strength Spray	200 µl/mL	-	-	+	+	+	+
		-	-	+	+	+	+
		-	-	+	+	+	+
Zicam Extreme Congestion Relief	200 µl/mL	-	-	+	+	+	+
		-	-	+	+	+	+
		-	-	+	+	+	+
4-acetamidophenol	200 µg/mL	-	-	+	+	+	+
		-	-	+	+	+	+
		-	-	+	+	+	+
Acetylsalicylic acid	650 µg/mL	-	-	+	+	+	+
		-	-	+	+	+	+
		-	-	+	+	+	+
Chlorpheniramine	145 ng/mL	-	-	+	+	+	+
		-	-	+	+	+	+
		-	-	+	+	+	+
Dexamethasone	0.80 mg/mL	-	-	+	+	+	+
		-	-	+	+	+	+
		-	-	+	+	+	+
Dextromethorphan	1 µg/mL	-	-	+	+	+	+
		-	-	+	+	+	+
		-	-	+	+	+	+
Diphenhydramine	5 µg/mL	-	-	+	+	+	+
		-	-	+	+	+	+
		-	-	+	+	+	+
Doxylamine Succinate	236 ng/mL	-	-	+	+	+	+
		-	-	+	+	+	+
		-	-	+	+	+	+
Ephedrine	273 ng/mL	-	-	+	+	+	+
		-	-	+	+	+	+
		-	-	+	+	+	+
Flunisolide	6.8 ng/mL	-	-	+	+	+	+
		-	-	+	+	+	+
		-	-	+	+	+	+
Guaiacol glycerol ether	3.5 ng/mL	-	-	+	+	+	+
		-	-	+	+	+	+
		-	-	+	+	+	+
Mupirocin	12 mg/mL	-	-	+	+	+	+
		-	-	+	+	+	+
		-	-	+	+	+	+
Oxymetazoline	0.6 mg/mL	-	-	+	+	+	+
		-	-	+	+	+	+
		-	-	+	+	+	+
Phenylephrine	12 mg/mL	-	-	+	+	+	+
		-	-	+	+	+	+
		-	-	+	+	+	+
Rebetol	4.5 µg /mL	-	-	+	+	+	+
		-	-	+	+	+	+
		-	-	+	+	+	+

Relenza	282 ng/mL	-	-	+	+	+	+
		-	-	+	+	+	+
Rimantadine	282 ng/mL	-	-	+	+	+	+
		-	-	+	+	+	+
		-	-	+	+	+	+
Tamiflu	1.1 µg /mL	-	-	+	+	+	+
		-	-	+	+	+	+
		-	-	+	+	+	+
Tobramycin	2.43 mg/ml	-	-	+	+	+	+
		-	-	+	+	+	+
		-	-	+	+	+	+
Triamcinolone	40 µg /mL	-	-	+	+	+	+
		-	-	+	+	+	+
		-	-	+	+	+	+

None of the 25 potential interference substances tested produced false positive or false negative test results in Aleré™ i Influenza A & B at the concentrations specified.

One important substance, FluMist™ was not tested in this study. The following statements have been added to the “Limitations” section of the Product Insert: “*Potential interference effects from FluMist™ have not been evaluated. Individuals who have received nasally administered influenza vaccine reportedly may test positive in commercially available influenza diagnostic tests for up to three days after vaccination.*”

#### *h. Potential Microbial Interference Study:*

##### Low Loads of Influenza A or Influenza B Co-Spiked with Selected Non-influenza Respiratory Viruses

An analytical study was performed to assess any potential interference of high concentrations of selected non-influenza respiratory viruses when co-spiked with low loads of influenza A or influenza B.

The non-influenza respiratory viral stocks used in this study are listed in Table 15 below.

**Table 15: Selected Non-Influenza Respiratory Viral Stocks Tested**

Virus	ATCC Strain	Lot Number	Stock Concentration (TCID <sub>50</sub> /mL)
Adenovirus Type 1	VR-1	59772750	1.58 x 10 <sup>7</sup>
Rhinovirus Type 1A	VR-1559	57685054	1.58 x 10 <sup>7</sup>
Respiratory Syncytial Virus (RSV) Type B, Strain 18537	VR-1580	58214862	8.89 x 10 <sup>5</sup>

Dilutions of Influenza A/Puerto Rico/8/34 (Lot D1002B, inactivated) and Influenza B/Malaysia/2506/2004 (Lot A1012A) were used in this study to prepare contrived positive swab samples for testing. The concentrations of the viral stocks (in TCID<sub>50</sub>/mL) were determined by standard virologic method. The targeted concentration for each dilution (in genome equivalents/mL) was also assessed using laboratory developed and validated influenza A and influenza B quantitative real-time PCR assays.

Using UTM as the diluent, positive samples were prepared using the influenza A and B strains (described above) diluted to a concentration approximately 2 to 3 times the respective limit of detection (LoD) levels. Swabs were prepared by coating 10 µl of the virus dilution onto the swab.

Swab coating was done by lightly scratching the swab surface with the pipette tip as the solution was pipetted to ensure liquid was absorbed into the swab tip.

The contrived influenza A and influenza B positive swabs were prepared as follows:

- Contrived influenza A swab:  $3.30 \times 10^3$  TCID<sub>50</sub>/Swab ( $2.00 \times 10^4$  genome equivalents/Swab); [Note: 10 µl of  $3.30 \times 10^5$  TCID<sub>50</sub>/mL ( $2.00 \times 10^6$  genome equivalents/mL) virus dilution coated onto a swab]
- Contrived influenza B swab:  $5.19 \times 10^1$  TCID<sub>50</sub>/Swab ( $5.00 \times 10^2$  genome equivalents/Swab); [Note: 10 µl of  $5.19 \times 10^3$  TCID<sub>50</sub>/mL ( $5.00 \times 10^4$  genome equivalents/mL) virus dilution coated onto a swab]

For each of the three non-influenza respiratory viruses, 75 µl of each viral stock (described in Table 15 above) was dispensed into each of six 1.5-ml Eppendorf tubes. One contrived influenza A or influenza B swab (prepared as described above) was placed into each tube. Swabs were stirred to absorb all of the liquid into the swab head. A total of three mock co-infection swabs per non-influenza virus were prepared using the influenza A or influenza B contrived positive swabs. The mock co-infection swabs were tested immediately according to the product instructions. Three contrived influenza A positive swabs and three contrived influenza B swabs without co-infection were also tested as positive experimental controls.

The Potential Microbial Interference – Non-Influenza Respiratory Viruses Study Results are presented in Table 16 below.

**Table 16: Microbial Interference Study Results – Selected Non-Influenza Respiratory Viruses**

Non-Influenza Respiratory Virus	Non-Influenza Respiratory Virus Concentration Tested	Contrived Influenza A Positive Swab		Contrived Influenza B Positive Swab	
		Flu A Result	Flu B Result	Flu A Result	Flu B Result
None	None	+	-	-	+
		+	-	-	+
		+	-	-	+
Adenovirus	$1.19 \times 10^6$ TCID <sub>50</sub> /Swab	+	-	-	+
		+	-	-	+
		+	-	-	+
Rhinovirus	$1.19 \times 10^6$ TCID <sub>50</sub> /Swab	+	-	-	+
		+	-	-	+
		+	-	-	+
RSV	$6.67 \times 10^4$ TCID <sub>50</sub> /Swab	+	-	-	+
		+	-	-	+
		+	-	-	+

No competitive inhibitory effect from co-infecting non-influenza respiratory virus was observed for this study.

Low Loads of Influenza A or Influenza B Co-Spiked with Influenza Virus

An analytical study was performed to assess any potential competitive inhibitory effects of high concentrations of influenza virus when co-spiked with low loads of influenza A or influenza B.

Dilutions of Influenza A/Puerto Rico/8/34 (Lot D1002B, inactivated) and Influenza B/Malaysia/2506/2004 (Lot A1012A) were used in this study to prepare contrived positive swab

samples for testing. The concentrations of the viral stocks (in TCID<sub>50</sub>/mL) were determined by standard virologic method. The targeted concentration for each dilution (in genome equivalents/mL) was also assessed using laboratory developed and validated influenza A and influenza B quantitative real-time PCR assays.

Using UTM as the diluent, positive samples were prepared using the influenza A and B strains (described above) diluted to a concentration approximately 2 to 3 times the respective limit of detection (LoD) levels. Swabs were prepared by coating 10 µl of the virus dilution onto the swab. Swab coating was done by lightly scratching the swab surface with the pipette tip as the solution was pipetted to ensure liquid was absorbed into the swab tip.

The contrived influenza A and influenza B positive swabs were prepared as follows:

- Contrived influenza A swab: 3.30 x 10<sup>3</sup> TCID<sub>50</sub>/Swab (2.00 x 10<sup>4</sup> genome equivalents/Swab); [Note: 10 µl of 3.30 x 10<sup>5</sup> TCID<sub>50</sub>/mL (2.00 x 10<sup>6</sup> genome equivalents/mL) virus dilution coated onto a swab]
- Contrived influenza B swab: 5.19 x 10<sup>1</sup> TCID<sub>50</sub>/Swab (5.00 x 10<sup>2</sup> genome equivalents/Swab); [Note: 10 µl of 5.19 x 10<sup>3</sup> TCID<sub>50</sub>/mL (5.00 x 10<sup>4</sup> genome equivalents/mL) virus dilution coated onto a swab]

To prepare the co-infection swabs, 10 µl of influenza A virus dilution (1.65 x 10<sup>6</sup> TCID<sub>50</sub>/mL or 1.00 x 10<sup>7</sup> genome equivalents/mL) was added to the contrived influenza B positive swabs prepared above by lightly scratching the swab surface with the pipette tip as the solution was pipetted to ensure liquid was absorbed into the swab tip (1.65 x 10<sup>4</sup> TCID<sub>50</sub>/swab or 1.00 x 10<sup>5</sup> influenza A genome equivalents/Swab). Likewise, 10 µl of influenza B virus dilution (1.59 x 10<sup>5</sup> TCID<sub>50</sub>/mL or 2.00 x 10<sup>6</sup> genome equivalents/mL) was added to the contrived influenza A positive swabs prepared above using the same method (1.59 x 10<sup>3</sup> TCID<sub>50</sub>/swab or 2.00 x 10<sup>4</sup> influenza B genome equivalents/Swab).

A total of three mock co-infection swabs per influenza virus were prepared using the influenza A or influenza B contrived positive swabs. The mock co-infection swabs were tested immediately according to the product instructions. Three contrived influenza A positive swabs and three contrived influenza B swabs without co-infection were also tested as positive experimental controls.

The Potential Microbial Interference – Influenza Viruses Study Results are presented in Table 17 below.

**Table 17: Microbial Interference Study Results –Influenza Viruses**

Co-Infecting Influenza Virus	Co-Infecting Influenza Virus Concentration Tested	Contrived Influenza A Positive Swab		Contrived Influenza B Positive Swab	
		Flu A Result	Flu B Result	Flu A Result	Flu B Result
None	None	+	-	-	+
		+	-	-	+
		+	-	-	+
Influenza A	1.00 x 10 <sup>5</sup> genome equivalents/Swab	ND	ND	+	+
		ND	ND	+	+
		ND	ND	+	+
Influenza B	2.00 x 10 <sup>4</sup> genome equivalents/Swab	+	+	ND	ND
		+	+	ND	ND
		+	+	ND	ND

Note: ND = Not Done

No competitive inhibitory effect from co-infecting influenza virus was observed for this study.

*i. Carry-Over Study:*

A study was conducted to evaluate the risk of carry-over contamination with the Alere™ i Influenza A & B assay on the Alere™ i instrument.

Dilutions of Influenza A/Puerto Rico/8/34 (Lot D1002B, inactivated) and Influenza B/Malaysia/2506/2004 (Lot B1012C, inactivated) were used in this study to prepare contrived samples for testing. The concentrations of the viral stocks (in TCID<sub>50</sub>/mL) were determined by standard virologic method prior to inactivation by the vendors. The targeted concentration for each dilution (in genome equivalents/mL) was also assessed using laboratory developed and validated influenza A and influenza B quantitative real-time PCR assays.

Contrived positive nasal swab samples were prepared using the influenza A and B strains (described above) dilutions in UTM with targeted concentrations of a minimum of 10<sup>6</sup> TCID<sub>50</sub>/mL. Mock negative nasal swab samples were prepared using UTM. Swabs were prepared by coating 10 µl of the virus dilutions or UTM onto the swab. Swab coating was done by lightly scratching the swab surface with the pipette tip as the dilution was pipetted to ensure liquid was absorbed into the swab tip. Alere™ i Influenza A & B Negative Control swabs were tested as negative swab samples.

The contrived influenza A and influenza B positive swabs were tested as follows:

- Contrived influenza A and B positive swab:
  - For Flu A: 1.00 x 10<sup>5</sup> TCID<sub>50</sub>/Swab (6.07 x 10<sup>5</sup> genome equivalents/Swab); [Note: 10 µl of 1.00 x 10<sup>7</sup> TCID<sub>50</sub>/mL (6.07 x 10<sup>7</sup> genome equivalents/mL) virus dilution coated onto a swab]
  - For Flu B: 1.00 x 10<sup>4</sup> TCID<sub>50</sub>/Swab (1.19 x 10<sup>5</sup> genome equivalents/Swab); [Note: 10 µl of 1.00 x 10<sup>6</sup> TCID<sub>50</sub>/mL (1.19 x 10<sup>7</sup> genome equivalents/mL) virus dilution coated onto a swab]

Positive and negative swab samples were tested per the product instructions. Testing between positive and negative samples alternated for a total of 15 rounds (15 positives and 15 negatives) on one Alere™ i instrument.

The Carry-over Study results are presented in Table 18 below.

**Table 18: Carry-Over Study Results**

Sample	# Tested	Influenza A and B Positive Sample		Influenza A and B Negative Sample	
		Flu A Result	Flu B Result	Flu A Result	Flu B Result
Direct Nasal Swab	15	+	+	-	-

No carry-over/contamination event was observed in this study. This data demonstrated that there is low risk of carryover/contamination when testing positive and negative samples on the Alere™ i instrument when the product insert instructions are followed.

*j. Assay cut-off:*

Alere™ i Influenza A & B consists of three sets of fluorescence response curves for Flu A, Flu B, and IC. The response curves are analyzed within a decision algorithm to form a channel result. There

are nine critical parameters (Table 19 below) associated with each channel in the decision algorithm. Most important among these parameters is the gradient threshold (parameter – *GradThreshold*). For each response curve, the gradient values within the threshold window (parameter – *ThresholdWindow*) are compared to the gradient threshold (parameter – *GradThreshold*) and the results of this comparison are what ultimately determine the channel result, “Asserted” or “Not Asserted”. In addition, if the calculated *NormAve* in a *NormWindow* is not between *NormAveMax* and *NormAveMin*, the channel result is “Indeterminate”.

The values of the nine parameters used within the channel decision algorithm are listed in Table 19 below.

**Table 19: Channel Decision Algorithm Parameter Values**

Parameter Name	Value
<i>NormWindow</i> (s)	40 – 110 (measurement blocks 5 – 12)
<i>ThresholdWindow</i> (s)	120 – End of acquisition (measurement blocks 13 – 58)
<i>GradWindow</i>	5 consecutive data points (50s, 5 measurement blocks)
<i>NormAveMin</i> (mV)	Flu A: 1, Flu B: 1, IC: 1
<i>NormAveMax</i> (mV)	Flu A:1200, Flu B: 1200, IC: 1200
<i>GradThreshold</i> (mV/s)	Flu A: 1.5, Flu B: 1.5, IC: 1.5
<i>GradPlateau</i> (mV/s)	0.2
<i>GradCount</i>	3
<i>AmpMAx</i> (mV)	4000

The gradient thresholds (parameter – *GradThreshold*) (cut-offs) in the channel decision algorithm were derived through analysis of clinical testing data (20 Flu A positives and 105 Flu A negatives for the Flu A *GradThreshold* ; 9 Flu B positives and 117 Flu B negatives for the Flu B *GradThreshold* ) collected during the 2011-2012 respiratory season. The gradients of the true positive and true negative clinical samples (as determined by viral culture and confirmed by PCR when there are discrepancies) were analyzed.

These thresholds were subsequently verified during the prospective clinical study conducted during the 2012-2013 respiratory season.

2. Comparison studies:

a. *Method comparison with predicate device:*

Not applicable. Performance of Alere™ i Influenza A & B was evaluated against the reference method of viral culture in a prospective clinical study.

b. *Matrix comparison:*

Not applicable.

3. Clinical studies:

a. *Prospective Study:*

Clinical performance characteristics of the Alere™ i Influenza A & B were evaluated in a multi-site prospective study during the 2012-2013 flu season in the U.S. A total of eight investigational sites throughout the U.S. participated in the study. To be enrolled in the study, patients had to be

presenting at the participating study centers with flu-like symptoms. Direct nasal swab specimens (without elution in VTM) from patients with flu-like symptoms were collected and tested using the Alere™ i Influenza A & B at the eight study sites. Viral culture performed according to standard virology culture procedures, was utilized as the reference method for this study.

Two nasal swabs were collected from one nostril from each subject using standard collection methods. At all sites, one nasal swab was tested directly on Alere™ i Influenza A & B, according to product instructions, without prior elution in VTM. The other nasal swab was eluted in 3-mL of viral transport media (VTM). Six of the eight sites (Site 1, Site 4, Site 8, Site 10, Site 11, and Site 12) shipped nasal swab samples in VTM to a central testing laboratory for viral culture testing. This central testing laboratory was located at Site 1, which also participated as a sample collection and Alere™ i Influenza A & B testing site. The nasal swab samples in VTM from Site 2 and Site 9 were cultured on site by a local laboratory.

External control testing, using Alere™ i Influenza A & B Positive and Negative Controls, was performed prior to sample testing each day and on each Alere™ i instrument the testing was performed, at all study sites.

All specimens generating discrepant Alere™ i Influenza A & B and viral culture results were investigated by testing using an FDA cleared Influenza RT-PCR assay at a central testing laboratory located at Site 1.

#### Performance - Direct Nasal Swab (Without Elution in VTM)

A total of 612 nasal swab specimens were enrolled in this study. Of those, 27 direct nasal swab specimens did not meet eligibility criteria (i.e., five samples failed to meet inclusion/exclusion criteria, eight samples were not transported to the viral culture testing laboratory per the conditions required by the clinical study protocol, two samples for which site did not perform external QC controls testing before sample testing, nine samples were observed with blood, and three samples had indeterminate viral culture results). A total of 585 direct nasal swab specimens were considered evaluable. Patient age and gender distribution for the 585 specimens is presented in Table 20 below.

**Table 20: Age and Gender Distribution – Direct Nasal Swab Study**

Age Group	Female	Male
<1 year	56	59
1 to 5 years	108	117
6 to 10 years	55	47
11 to 15 years	21	20
16 to 21 years	14	8
>21 to 60 years	51	20
>60 years	5	4
Total	310	275

Of the evaluable 585 specimens, the Alere™ i Influenza A & B generated influenza A invalid results for 14 specimens and influenza B invalid results for 16 specimens, resulting in a total of 571 specimens for influenza A performance analysis and 569 specimens for influenza B performance analysis.

Compared to the viral culture reference method, the performance of Alere™ i Influenza A & B for influenza A and influenza B are presented in Table 21 and Table 22 below, respectively.

**Table 21: Alere™ i Influenza A & B Influenza A Direct Nasal Swab Performance against Viral Culture**

Alere™ i Influenza A & B – Flu A	Culture		
	Positive	Negative	Total
Positive	92	66 <sup>a</sup>	158
Negative	2 <sup>b</sup>	411	413
Total	94	477	571
Sensitivity: 92/94 97.9% (95% CI: 92.6%-99.4%)			
Specificity: 411/477 86.2% (95% CI: 82.8%-89.0%)			

<sup>a</sup> Flu A nucleic acid was detected in 58/66 False Positive specimens using an FDA-cleared molecular test

<sup>b</sup> Flu A nucleic acid was not detected in 1/2 False Negative specimens using an FDA-cleared molecular test

**Table 22: Alere™ i Influenza A & B Influenza B Direct Nasal Swab Performance against Viral Culture**

Alere™ i Influenza A & B – Flu B	Culture		
	Positive	Negative	Total
Positive	74	17 <sup>a</sup>	91
Negative	6 <sup>b</sup>	472	478
Total	80	489	569
Sensitivity: 74/80 92.5% (95% CI: 84.6%-96.5%)			
Specificity: 472/489 96.5% (95% CI: 94.5%-97.8%)			

<sup>a</sup> Flu B nucleic acid was detected in 15/17 False Positive specimens using an FDA-cleared molecular test

<sup>b</sup> Flu B nucleic acid was not detected in 4/6 False Negative specimens using an FDA-cleared molecular test

Site to site performance of Alere™ i Influenza A & B for influenza A and influenza B comparing to viral culture are presented in Table 23 and Table 24 below, respectively.

**Table 23: Site by Site Direct Nasal Swab Performance Obtained for Influenza A with the Alere™ i Influenza A & B in Comparison to Viral Culture**

Clinical Sites	Prevalence <sup>a</sup>	Sensitivity with 95% CI <sup>b</sup>	Specificity with 95% CI <sup>b</sup>
Site 1	0.0% (0/58)	0.0% (0/0) (NA)	89.7% (52/58) <sup>c</sup> (79.2%, 95.2%)
Site 2	15.5% (9/58)	100% (9/9) (70.1%, 100%)	91.8% (45/49) <sup>d</sup> (80.8%, 96.8%)
Site 4	3.2% (1/31)	100% (1/1) (20.7%, 100%)	93.3% (28/30) <sup>e</sup> (78.7%, 98.2%)
Site 8	3.4% (2/58)	100% (2/2) (34.2%, 100%)	67.9% (38/56) <sup>f</sup> (54.8%, 78.6%)
Site 9	27.4% (81/296)	98.8% (80/81) (93.3%, 99.8%)	85.6% (184/215) <sup>g</sup> (80.3%, 89.7%)
Site 10	0.0% (0/33)	0.0% (0/0) (NA)	97.0% (32/33) (84.7%, 99.5%)
Site 11	0.0% (0/15)	0.0% (0/0) (NA)	93.3% (14/15) <sup>h</sup> (70.2%, 98.8%)
Site 12	4.5% (1/22)	0.0% (0/1) <sup>i</sup> (NA)	85.7% (18/21) <sup>1</sup> (65.4%, 95.0%)
Overall	16.5% (94/571)	97.9% (92/94) <sup>1</sup> (92.6%, 99.4%)	86.2% (411/477) <sup>k</sup> (82.8%, 89.0%)

<sup>a</sup> Prevalence based on the reference method only

<sup>b</sup> Confidence interval

<sup>c</sup> Flu A nucleic acid was detected in 6/6 False Positive specimens using an FDA-cleared molecular test

<sup>d</sup> Flu A nucleic acid was detected in 2/4 False Positive specimens using an FDA-cleared molecular test

<sup>e</sup> Flu A nucleic acid was detected in 2/2 False Positive specimens using an FDA-cleared molecular test

<sup>f</sup> Flu A nucleic acid was detected in 16/18 False Positive specimens using an FDA-cleared molecular test

<sup>g</sup> Flu A nucleic acid was detected in 28/31 False Positive specimens using an FDA-cleared molecular test

<sup>h</sup> Flu A nucleic acid was detected in 1/1 False Positive specimen using an FDA-cleared molecular test

<sup>i</sup> Flu A nucleic acid was detected in 3/3 False Positive specimens using an FDA-cleared molecular test

<sup>1</sup> Flu A nucleic acid was not detected in this False Negative specimen using an FDA-cleared molecular test

<sup>j</sup> Flu A nucleic acid was not detected in this False Negative specimen using an FDA-cleared molecular test

<sup>k</sup> Flu A nucleic acid was detected in 58/66 False Positive specimens using an FDA-cleared molecular test

<sup>1</sup>Flu A nucleic acid was not detected in 1/2 False Negative specimens using an FDA-cleared molecular test

**Table 24: Site by Site Direct Nasal Swab Performance Obtained for Influenza B with the Alere™ i Influenza A & B in Comparison to Viral Culture**

Clinical Sites	Prevalence <sup>a</sup>	Sensitivity with 95% CI <sup>b</sup>	Specificity with 95% CI <sup>b</sup>
Site 1	1.7% (1/58)	100% (1/1) (20.7%, 100%)	98.2% (56/57) <sup>c</sup> (90.7%, 99.7%)
Site 2	15.8% (9/57)	100% (9/9) (70.1%, 100%)	95.8% (46/48) <sup>d</sup> (86.0%, 98.8%)
Site 4	19.4% (6/31)	100% (6/6) (61.0%, 100%)	92.0% (23/25) (75.0%, 97.8%)
Site 8	10.3% (6/58)	100% (6/6) (61.0%, 100%)	94.2% (49/52) <sup>e</sup> (84.8%, 98.0%)
Site 9	15.9% (47/295)	87.2% (41/47) <sup>g</sup> (74.8%, 94.0%)	98.4% (244/248) <sup>f</sup> (95.9%, 99.4%)
Site 10	27.3% (9/33)	100% (9/9) (70.1%, 100%)	95.8% (23/24) <sup>h</sup> (79.8%, 99.3%)
Site 11	13.3% (2/15)	100% (2/2) (34.2%, 100%)	92.3% (12/13) <sup>i</sup> (66.7%, 98.6%)
Site 12	0.0% (0/22)	0.0% (0/0) (NA)	86.4% (19/22) <sup>j</sup> (66.7%, 95.3%)
Overall	14.1% (80/569)	92.5% (74/80) <sup>1</sup> (84.6%, 96.5%)	96.5% (472/489) <sup>k</sup> (94.5%, 97.8%)

<sup>a</sup> Prevalence based on the reference method only

<sup>b</sup> Confidence interval

<sup>c</sup> Flu B nucleic acid was detected in 1/1 False Positive specimen using an FDA-cleared molecular test

<sup>d</sup> Flu B nucleic acid was detected in 2/2 False Positive specimens using an FDA-cleared molecular test

<sup>e</sup> Flu B nucleic acid was detected in 3/3 False Positive specimens using an FDA-cleared molecular test

<sup>f</sup> Flu B nucleic acid was detected in 4/4 False Positive specimens using an FDA-cleared molecular test

<sup>g</sup> Flu B nucleic acid was not detected in 4/6 False Negative specimens using an FDA-cleared molecular test

<sup>h</sup> Flu B nucleic acid was detected in 1/1 False Positive specimen using an FDA-cleared molecular test

<sup>i</sup> Flu B nucleic acid was detected in 1/1 False Positive specimen using an FDA-cleared molecular test

<sup>j</sup> Flu B nucleic acid was detected in 3/3 False Positive specimens using an FDA-cleared molecular test

<sup>k</sup> Flu A nucleic acid was detected in 15/17 False Positive specimens using an FDA-cleared molecular test

<sup>1</sup> Flu A nucleic acid was not detected in 4/6 False Negative specimens using an FDA-cleared molecular test

Performance of Alere™ i Influenza A & B for the detection of influenza A and influenza B versus culture is presented in Table 25 below stratified by patient age.

**Table 25: Nasal Swab Performance Obtained for Influenza A and Influenza B with the Alere™ i Influenza A & B in Comparison to Viral Culture – Stratified by Patient Age**

Influenza Type	≤ 5 Years of Age (n = 332)		6 - ≤ 21 Years of Age (n = 162)		≥ 22 Years of Age (n = 77)	
	Sensitivity 95% CI	Specificity 95% CI	Sensitivity 95% CI	Specificity 95% CI	Sensitivity 95% CI	Specificity 95% CI
<b>Flu A</b>	98.3% (58/59) 91.0% - 99.7%	89.0% (243/273) 84.7% - 92.2%	100% (31/31) 89.0% - 100%	85.5% (112/131) 78.5% - 90.5%	75.0% (3/4) 30.1% - 95.4%	76.7% (56/73) 65.8% - 84.9%
<b>Flu B</b>	88.9% (32/36) 74.7% - 95.6%	98.0% (288/294) 95.6% - 99.1%	94.4% (34/36) 81.9% - 98.5%	96.8% (122/126) 92.1% - 98.8%	100% (8/8) 67.6% - 100%	89.9% (62/69) 80.5% - 95.0%

The Alere™ i Influenza A & B detected one mixed influenza A and B infections in the prospective clinical evaluation. This sample tested positive for influenza B only by viral culture, but tested positive for influenza A only by an FDA cleared Influenza RT-PCR assay. Based on this result, the following have been added to the “Result Interpretation” section of the Product Insert: *“Dual infections of Flu A and Flu B are rare. Repeat testing using new test components. Contact Technical Support during normal business hours if multiple samples provide this result.”*

During the prospective clinical study, the initial invalid rate (before repeat testing per the product instructions) was 5.8% (34/585) (95% CI: 4.2% to 8.0%) for Flu A, and 3.6% (21/585) (95% CI: 2.4% to 5.4%) for Flu B. After repeat testing per the product instructions, the invalid rate was 2.4% (14/585) (95% CI: 1.4%, 4.0%) for Flu A, and 2.7% (16/585) (95% CI: 1.7%, 4.4%) for Flu B.

*b. Retrospective Study:*

None.

4. Clinical cut-off:

Not applicable.

5. Expected values:

In the Alere™ i Influenza A & B prospective clinical study (described in the “Clinical Studies” section above), a total of 585 direct nasal swab specimens were determined to be evaluable. The number and percentage of influenza A and influenza B positive cases per specified age group, as determined by the Alere™ i Influenza A & B assay, are presented in Table 26 and Table 27 below:

**Table 26: Influenza A Positives by the Alere™ i Influenza A & B Assay per Age Group**

Age Group	Number of Direct Nasal Swab Specimens	Number of Influenza A Positives	Influenza A Positivity Rate
<1 year	121	25	20.1%
1 to 5 years	219	63	28.8%
6 to 10 years	102	35	34.3%
11 to 15 years	41	12	29.3%
16 to 21 years	22	3	13.6%
>21 to 60 years	71	18	25.4%
>60 years	9	2	22.2%
Total	585	158	27.0%

**Table 27: Influenza B Positives by the Alere™ i Influenza A & B Assay per Age Group**

Age Group	Number of Direct Nasal Swab Specimens	Number of Influenza B Positives	Influenza B Positivity Rate
<1 year	121	12	9.9%
1 to 5 years	219	26	11.9%
6 to 10 years	102	26	25.5%
11 to 15 years	41	7	17.1%
16 to 21 years	22	5	22.7%
>21 to 60 years	71	14	19.7%
>60 years	9	1	11.1%
Total	585	91	15.6%

## N. Instrument Name:

Alere™ i Instrument

## O. System Descriptions:

### 1. Modes of Operation:

The Alere™ i Instrument is a portable bench-top unit designed to perform in laboratory and point-of care environments. Alere™ i Influenza A & B components are color-coded to match the corresponding holder on the instrument and designed to only fit in the corresponding holder in one direction. Once properly placed into the instrument, the Alere™ i Instrument uses image analysis to automatically detect the presence of the Test Base and to check the test type and expiry date on the barcode present on the Test Base. This triggers the test sequence. Heating, mixing, and detection by fluorescence are provided by the instrument. The Alere™ i Instrument is equipped with two optical fluorescence modules, FAM™ and ROX™. Each optical module consists of a LED excitation source and photodiode receiver configured to the fluorescent channel. Fluorescence signals are optically filtered by a confocal lens and mirror arrangement and processed by software running on an on-board computer. The results are calculated automatically and presented on the instrument graphical user interface and available for export in JSON format or printing.

### 2. Software:

FDA has reviewed applicant's Hazard Analysis and software development processes for this line of product types:

Yes  or No

### 3. Specimen Identification:

Specimen ID is entered using on screen keyboard or barcode scanner.

### 4. Specimen Sampling and Handling:

Not applicable. The specimens are manually inserted in the Sample Receiver in the instrument.

### 5. Calibration:

The Alere™ i Instrument is factory calibrated and does not require any further calibration and verification at user site. However, if the instrument was transported or moved, a performance check using Alere™ i Positive and Negative Controls is recommended to ensure proper functionality by the manufacture.

### 6. Quality Control:

Quality control is addressed for each specific FDA-cleared assay to be run on the instrument (separately cleared).

**P. Proposed Labeling:**

The labeling is sufficient and it satisfies the requirements of 21 CFR Part 809.10.

**Q. Conclusion:**

The submitted information in this premarket notification is complete and supports a substantial equivalence decision.